SUMMARY OF PEER-REVIEWED LITERATURE

Fixation Strength Comparison of Rotator Cuff Anchors in Porcine Bone

Overview

Following is a summary of the ultimate load-to-failure results as reported in the Journal of Arthroscopy in May 2013 and July 2011 for select 5.5mm rotator cuff knot-tying anchors.

| Anchor | Manufacturer | Anchor Material | Sutures* | Size |
|----------------|-----------------|-------------------------|-----------------------|-------|
| CrossFT™ (PK) | ConMed Linvatec | PEEK | #2 Hi-Fi [®] | 5.5mm |
| TwinFix (HA) | Smith & Nephew | Biocomposite (HA/PLLA) | #2 Ultrabraid | 5.5mm |
| TwinFix (PK) | Smith & Nephew | PEEK | #2 Ultrabraid | 5.5mm |
| Healicoil (PK) | Smith & Nephew | PEEK | #2 Ultrabraid | 5.5mm |
| Healix (BR) | DePuy Mitek | Biocomposite (TCP/PLGA) | #2 OrthoCord | 5.5mm |
| Healix (PK) | DePuy Mitek | PEEK | #2 OrthoCord | 5.5mm |
| ALLthread (PK) | Biomet | PEEK | #2 MaxBraid | 5.5mm |

Table 1: Rotator cuff anchor properties

Methods

Anchors were implanted in porcine metaphyseal femoral cortex and tensile loads were applied parallel to axis of insertion at 12.5mm/s until failure. Mean failure loads are compared for anchors listed in Table 1.

Results

The CrossFT™ (PK) suture anchor exhibited 569.5 \pm 139.1 N¹ mean ultimate failure load, which was higher than any non-metallic anchor of any size in either study. The comparable size PEEK anchors exhibited the following failure loads: ALLthread (PK) $(476.5 \pm 22.2 \text{ N})^{1}$, TwinFix (PK) $(469.4 \pm 48.7 \text{ N})^2$, Healix (PK) $(404.3 \pm 24.4 \text{ N})^2$, and Healicoil (PK) $(298.7 \pm 37.4 \text{ N})^2$. The comparable size biocomposite anchors exhibited the following failure loads: TwinFix (HA) (382.6 ± 38.2 N)2 and Healix (BR) (312.1 ± 31.0 N)2.




Figure 1: Cortical loads-to-failure (N) in porcine bone for 5.5mm non-metallic anchors.^{1,2}

Clinical Relevance

The CrossFT (PK) anchor had 20% higher pullout than the next highest anchor, ALLthread (PK), and may help reduce the risk of anchor pullout postoperatively.³

References

- 1. Barber FA, et al. Biomechanical Analysis of Pullout Strengths of Rotator Cuff and Glenoid Anchors: 2011 Update. Arthroscopy 2011; 27:895-905.
- 2. Barber FA, et al. Cyclic Loading Biomechanical Analysis of the Pullout Strengths of Rotator Cuff and Glenoid Anchors: 2013 Update. Arthroscopy 2013; 29:832-844.
- 3. Tashjian RZ, et al. Initial Fixation Strength of Massive Rotator Cuff Tears: In Vitro Comparison of Single-Row Suture Anchor and Transosseous Tunnel Constructs. Arthroscopy 2007; 23:710-716