

Getinge Clean

Pre Treatment Foam

A ready to use pre-treatment, pH neutral, waterbased foam to be applied on medical devices immediately after use, to prevent drying of surgical soils. Adheres to medical device surfaces and keeps soils moist for easy removal. Suitable for stainless steel, ceramics, glass, plastic and soft metal surfaces. Water based and easily rinsed off in automated or manual pre-wash applications.



Instructions for use

- Spray on soiled instruments immediately after use, prior to transport
- Shake before use

Cautions

- Check material for compatibility before use
- Do not mix with other chemicals
- Always wear correct PPE

Symbols



Instruments



Pre treatment foam

Technical data

Specific gravity (20 °C):	1.016
pH as supplied (typical):	7
Viscosity (typical):	110 cSt

Ingredients

Surfactants, rheology modifier, preservative and corrosion inhibitors.

Storage information

Store in a cool, dry place protected from frost and away from acids and strong oxidizing agents. Store upright in original containers. Optimum storage temperature 5–25 °C. Shelf life: 18 months in original, sealed containers.

Safety and handling

Refer to Safety Data Sheet.

CE Mark

CE marked in accordance with current European medical devices regulatory requirements.

Service

For technical support in regards to appropriate detergent use, please contact your local Getinge representative.

Environmental awareness

Biodegradable formulation and recyclable packaging.

Ordering number

Size (ml)	Language	Ordering number	Comments
500	en, sv, da, no, fr, de, es, it, pt, tr, ru, zh, ja, pl	6001752131	6-pack
500	en, sv, da, no, fr, de, es, it, pt, tr, ru, zh, ja, pl	6001752231	20-pack

Material compatibility

Always follow the cleaning instructions from the manufacturer of the items and instruments and check material for compatibility before use. The following chart serves only as a guide.

	Aluminium	Anodized aluminium	Titanium	Non ferrous metals	Matt chrome plated	Stainless steel	Glass	Ceramics	Rubber	Plastic
Getinge Clean Universal Detergent	M	M	M	M	✓	✓	✓	✓	✓	M
Getinge Clean Instrument Lubricant	✓	✓	✓	✓	✓	✓	-	-	-	-
Getinge Clean Instrument Lubricant Plus	✓	✓	✓	✓	✓	✓	-	-	-	-
Getinge Clean Pre Treatment Foam	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Getinge Clean Enzymatic Detergent	✓	✓	✓	✓	✓	✓	✓	✓	✓	M
Getinge Clean Manual Plus	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Getinge Clean MIS Detergent	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Getinge Clean Prion Detergent	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Getinge Clean Heavy Soil Detergent	-	-	M	-	✓	✓	✓	✓	✓	M
Getinge Clean Rinse Aid	✓	✓	✓	✓	✓	✓	✓	✓	M	M
Getinge Clean Chamber Foam	-	-	-	-	-	✓	-	-	-	-
Getinge Clean Neutralizer	M	M	✓	✓	✓	✓	✓	✓	M	M
Getinge Clean Neutralizer Plus	-	-	M	M	M	✓	✓	M	M	M

006070

✓ = The product is suitable for the material

- = The product is not suitable for the material

M = Pay special attention to the guidelines from the instruments/items manufacturer