

Getinge Clean

Instrument Lubricant Plus

A water based, silicone and mineral oil free lubricant for application after cleaning to lubricate joints and moving parts of instruments. The product can be applied both manually and in automated washer-disinfector processes. Instrument Lubricant Plus has been tested according to ISO 10993 for biocompatibility and found not to be irritating, sensitising or cytotoxic. Instrument Lubricant Plus will not interfere with steam or ethylene oxide sterilization processes.



Instructions for use

Manual process: In a suitable container, mix 20-40 ml/l in clean potable water, preferably deionized, RO or distilled at <25 °C. Open the instrument joints and place in the basket. Submerge the instruments in the solution for 30 seconds. Remove the instruments from the solution and allow to drain and dry prior to sterilization. The solution must be discarded daily or after each shift or if it becomes visibly contaminated.

Automated process: Dose 1–2 ml/l in the final rinse water, preferably deionized or RO.

Cautions

- Check material for compatibility before use
- Do not mix with other chemicals
- Always wear correct PPE

Symbols



Instrument lubricant



Instruments

Technical data

Specific gravity (20 °C):	1.004 (typical)
Viscosity (concentrate, 20 °C):	4.2 cSt (typical)

Ingredients

Contains nonionic surfactants 5–15 %, anionic surfactants <5 % and preservative.

Storage information

Store in a cool, dry place protected from frost and away from acids and strong oxidizing agents. Store upright in original containers. Optimum storage temperature 5–25 °C. Shelf life: 36 months in original, sealed containers.

Safety and handling

Refer to Safety Data Sheet.

CE Mark

CE marked in accordance with current European medical devices regulatory requirements.

Service

For technical support in regards to appropriate detergent use, please contact your local Getinge representative.

Environmental awareness

Biodegradable formulation and recyclable packaging.

Ordering number

Size (l)	Language	Ordering number	Comments
5	en, da, fr, es, it, tr, ru, zh	6001750331	
25	en, fr, ru, zh	6001750531	
205	en, da, fr, es, it, tr, ru, zh	6001750631	

Material compatibility

Always follow the cleaning instructions from the manufacturer of the items and instruments and check material for compatibility before use. The following chart serves only as a guide.

	Aluminium	Anodized aluminium	Titanium	Non ferrous metals	Matt chrome plated	Stainless steel	Glass	Ceramics	Rubber	Plastic
Getinge Clean Universal Detergent	M	M	M	M	✓	✓	✓	✓	✓	M
Getinge Clean Instrument Lubricant	✓	✓	✓	✓	✓	✓	-	-	-	-
Getinge Clean Instrument Lubricant Plus	✓	✓	✓	✓	✓	✓	-	-	-	-
Getinge Clean Pre Treatment Foam	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Getinge Clean Enzymatic Detergent	✓	✓	✓	✓	✓	✓	✓	✓	✓	M
Getinge Clean Manual Plus	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Getinge Clean MIS Detergent	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Getinge Clean Prion Detergent	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Getinge Clean Heavy Soil Detergent	-	-	M	-	✓	✓	✓	✓	✓	M
Getinge Clean Rinse Aid	✓	✓	✓	✓	✓	✓	✓	✓	M	M
Getinge Clean Chamber Foam	-	-	-	-	-	✓	-	-	-	-
Getinge Clean Neutralizer	M	M	✓	✓	✓	✓	✓	✓	M	M
Getinge Clean Neutralizer Plus	-	-	M	M	M	✓	✓	M	M	M

006070

✓ = The product is suitable for the material

- = The product is not suitable for the material

M = Pay special attention to the guidelines from the instruments/items manufacturer

Performance

This product has been developed, tested and validated on Getinge washers, meeting the appropriate ISO 15883 standards.