




Paradigm™ Retractable Straight Blade and Hook Blade

Featuring surgeon-designed handles, both the Paradigm™ Retractable Straight Blade and Hook Blade are designed for safety. They feature retractable sheaths that allow surgeons to safely probe with the blunt, sheathed tip and then use a lever to reveal the sharp blade when it is time to cut.



Features and Benefits

Retractable Straight Blade

Extremely sharp, piercing blade enables ability to “push-cut” through the hip capsule.

Retractable Hook Blade

Retrograde blade enables effective cutting while providing a soft tissue-protecting blunt tip.

Designed for Safety

Blunt retractable sheath protects the sharp blade when entering the joint, reducing the risk of iatrogenic damage to the articular cartilage and soft tissue.

Ease of Use

The Paradigm™ Retractable Straight Blade and Hook Blade handle can be screwed into CONMED’s DoubleLatch™ Hip Cannula for ease of use.

Experience the Paradigm shift. Take the next step toward simple, safe, and efficient hip arthroscopy. Call **1-866-4CONMED** or visit www.CONMED.com/Paradigm to learn more.



*Paradigm™ Retractable
Straight Blade*



*Paradigm™ Retractable
Hook Blade*

Ordering Information

Description

**Catalog
Number**

Paradigm™ Retractable Blades

Paradigm™ Retractable Blade, Straight

PRBS

Paradigm™ Retractable Blade, Hook

PRBH