

#### High Performance Brilliantly Controlled

You make the choice... the 24k<sup>®</sup> delivers with superior performance



#### HIGH PERFORMANCE DER HIGH PERFORMANCE D

# The 24k<sup>®</sup> delivers high performance with custom features that allow you to get the job done.

- Independent pressure and flow controls allow the 24k<sup>®</sup> Pump to keep the joint clear of debris even at lower pressure settings.
- Independent rotor control allows the 24k Pump to operate as a simple inflow-only pump or a complete arthroscopic fluid management system.
- Four Presets Knee, Shoulder, Small Joint and Hip simplify setting up the pump.
- Custom Wave Mode in the Knee Preset for excellent joint clearing even with a double stopcock scope sheath.
- Simple intuitive tubing set up. Hang the cassettes, spike the bags, prime and go.
- Independent shaver sensor allows the 24k Pump to control suction with any ConMed Linvatec or competitive shaver handpiece.
- Fully autoclavable multifunction remote control gives you complete control of the pump.

## ANDED... DELIVERED.

### functions and flexible



Viove Mode Suction Higher Outlow TIPS Performance's Clearfield Performance

24k® F	PUMP ORDERING INFORMATION
Cat. No.	Description
24k	24k® Pump
24k900	24k Shaver Sensor
24k191	24k Shaver Sensor Extension Cable
C8115	24k Remote
24k100	24k Outflow/Shaver Suction Tube Set
24kTIPS	24k TIPS Tube Set
10k100	10k® Arthroscopy Tube Set (Inflow Only)
10k100D	10k Arthroscopy Tube Set for Dyonics
10k150	10k Arthroscopy Tube Set with Suction Line
10k500	10k Laparoscopic Tubing Set
10k505	10k Laparoscopic Tube Set with Probe
10k600	Eco-Flow <sup>™</sup> Day-use Pump Tubing (Must be used with 10k605)
10k605	Eco-Flow <sup>™</sup> Single-use Patient Tubing (Must be used with 10k600)
C7104	AC Power Cord 115VAC
C7105	AC Power Cord 230 VAC

#### ConMed Linvatec Fluid Management System

11311 Concept Boulevard Largo, Florida 33773-4908

Local 727-392-6464 Toll Free 800-237-0169

Website www.linvatec.com Email customer\_service@linvatec.com