


## Paradigm™ EZ Switch™ Portal Saver System

Say goodbye to the sleds and switching sticks of rigid, motion-restricting cannulas. The innovative EZ Switch™ Portal Saver System is flexible to allow for a wide range of motion while allowing you to easily and efficiently switch portals.



## Features and Benefits

### Designed to Save Time

Easily switch portals without sleds or switching sticks.

### Flexible Design

Allows for a wide range of motion.

### Customizable

Flexible body length (120mm) can be easily trimmed to match patient anatomy.

### Optimal Access

Up to 150mm working length is designed to enable easy access to the surgical site.

### 30mm Threads

Designed to provide strong dermal fixation.

- **PARADIGM™ EZ SWITCH™** Portal Saver System is available in several kit configurations to meet the needs of your technique and your facility.

Experience the Paradigm shift. Take the next step toward simple, safe, and efficient hip arthroscopy. Call **1-866-4CONMED** or visit **CONMED.com/Paradigm** to learn more.



*Paradigm™ EZ Switch™  
Portal Saver System*

## Ordering Information

| Description   | Catalog Number |
|---|----------------|
| <b>NEW</b> Paradigm™ Access                         |                |
| Paradigm™ EZ Switch™ Portal Saver System            | PEZS01         |
| Paradigm™ Access Kit and EZ Switch™ System*         | PEZS02         |
| Paradigm™ Tactile Access Kit and EZ Switch™ System* | PEZS03         |

\*Not for sale or distribution in all markets.