



HandStand™

Traction Tower

The HandStand™ Traction Tower takes all the features you loved about CONMED's TractionTower® and adds modern innovations. With improved maneuverability, sterile strapping, a more intuitive assembly, and a lightweight design, the HandStand™ Traction Tower is designed to ensure an easy day in the OR.

Features and Benefits

Engineered specifically for the small, hard-to-reach anatomic space of the wrist, the HandStand™ device allows for precise distraction, flexion, and extension of the wrist during arthroscopy.

Designed for Optimal Flexion and Extension

The HandStand™ arm can be rotated 180 degrees in 15-degree increments for full access to the forearm and to accommodate the use of fluoroscopy during arthroscopic use. The HandStand™ arm can also bend up to 70 degrees in either direction to accommodate flexion and extension of the wrist.



Intuitive Quick-Connect Assembly

Designed to reduce time and complexity, this system features fewer parts and an intuitive assembly process that makes both setting up the system and adjusting it during the case quick and easy.¹

QuickCool™ Technology

QuickCool™ Technology means the HandStand™ is never too hot to handle. In fact, when compared to our legacy Traction Tower®, the HandStand™ cools almost 10% faster.^{2,3} In situations where a quick turnover is required, place the system in saline post-autoclave and the metal will become cool to the touch within seconds.

Sterile Traps and Straps

Unique arm straps are designed for patient safety by eliminating the need for strapping and unstrapping outside of the sterile field. The updated Digit Trap® design allows for maximum grasp on a variety of finger sizes during arthroscopy. The customizable sizing feature enables secure fixation while the braided hose configuration is designed to distribute compression over the finger, alleviating concerns of digit strangulation.



HandStand™ Digit Trap®

Lightweight Design

Constructed of lightweight aluminum, the entire HandStand™ system weighs less than 11 lbs. for hassle-free sterilization, transportation, assembly, and disassembly.

¹ Memo on File

² Data on File: PDD1435217

³ Data on File: PDD1511206

Need a Hand with your next wrist arthroscopy case? Call **1-866-4CONMED** to schedule an evaluation with your local CONMED Sales Representative or visit www.CONMED.com/HandStand for more information.



525 French Road, Utica, NY 13502

HANDSTAND™
TRACTION TOWER

Ordering Information

Description	Disposable or Reusable	Catalog Number
HandStand™ Traction Tower	Reusable	HS-TT
HandStand™ Disposable Kit (10 Kits Per SKU) Includes: 2 x Digit Trap®, 1 x Forearm Strap, and 1 x Bicep Strap	Disposable	HS-KIT10
HandStand™ Disposable Digit Trap® Includes: 2 x Digit Trap®	Disposable	HS-DT2
HandStand™ Disposable Forearm Strap Includes: 1 x Forearm Strap	Disposable	HS-FS
HandStand™ Disposable Bicep Strap Includes: 1 x Bicep Strap	Disposable	HS-BS
HandStand™ Sterilization Tray	Reusable	HS-ST
HandStand™ Container	Reusable	INST-C352

CONMED Corporation
525 French Road, Utica, NY 13502

Toll Free: 1-866-4CONMED
International: 727-214-3000
www.CONMED.com

customerexperience@conmed.com
internationalorders@conmed.com