


Bexen cardio
Ready for life

Reanibex 800

The unique modular
Monitor Defibrillator
“plug & play”
of the market



Monitor Defibrillator
multiparameter

Manual and AED up to 360J

Universal AED algorithm for
adult and paediatric patients

12 leads acquisition

ECG interpretation

Non invasive pacing

SpO2, SpCO, SpMet, SpHb,
SpOC, PVI

NIBP

EtCO2

IBP

TEMP

CPR assistance

Data transmission
(Bluetooth, GSM)

Reanibex 800.

Modular Monitor Defibrillator



Communications (Reanibex Data Manager, DataLink, DataCloud)



Operates from the mains (AC), from a vehicle battery (DC) and with its battery

12 ECG acquisition key



Modular easily upgradeable. Plug & Play system

50 mm or 106 mm printer

Non invasive pacing

External reusable paddles with contact indicator on the APEX paddle



Status indicator

Large colour display (8.4")

Therapy selector

AED for adult or paediatric

FEATURES

- Biphasic technology.
- Large colour display (8.4").
- Monitoring and diagnostic capabilities.
- 3, 5 or 10 leads patient cable.
- Non invasive pacing (optional).
- 12 leads ECG data transmission from the field to the Emergency Service (optional).
- 12 leads Interpretive Algorithm (University of Glasgow) (Optional).
- Synchronized Cardioversion.
- CPR Metronome with feedback.
- Masimo Rainbow SET technology that monitors SpO2, SpCO and SpMet (optional).
- Capnography from Oridion with Microsteam technology.
- Non-invasive Blood Pressure from SunTech Medical.
- Invasive Blood Pressure (two lines).
- Temperature (two lines).
- Operates from the mains (AC), from a vehicle battery (DC) and with its battery.
- Possibility to operate with two batteries.
- Possibility of setting the screen on high contrast mode.

ACCESSORIES

- Adult and paediatric disposable electrodes for defibrillation.
- Adult and paediatric ECG electrodes.
- Carrying case (1).
- Neoprene upper bag (2)
- Ambulance bracket (EN 1789:2007 certified) (3).
- 3, 5 or 10 leads patient cable.
- Internal paddles (4).
- Rechargeable battery.
- External battery charger.
- Connection cable at 12V (DC).
- "REANIBEX Data Manager, Data Link, Data Cloud" software application to manage and organize all the data gathered.



(1)



(3)



(2)



(4)

Bexen cardio

Osatu, S.Coop.
Edificio Zearrekobueta
Subida de Areitio, 5
48260 Ermua (Spain)

Tel: +34 943 170 220
Fax: +34 943 170 227
info@bexencardio.com
www.bexencardio.com

