

WarwickSASCo
Moulding The Future of Medical Plastics



5 reasons to choose Warwick SASCo

Recognised as the market leader in the design and
manufacture of medical polyware

Wide range of high quality re-usable and single-use polypropylene
products for all areas of healthcare

Unrivalled reputation for providing quality products, delivery and
service to healthcare professionals worldwide

Re-usable products are guaranteed for a minimum of 1000 wash
cycles – delivering excellent value for money

Certification since 1997 demonstrates our continued commitment
to maintaining industry standards. Our quality management system
complies to ISO13485 and ISO9001

MADE IN BRITAIN



Warwick SASCo

For over 30 years Warwick SASCo has been part of evolving supply of
decontamination and nursing products worldwide to hospitals, healthcare and
patient nursing.

Recognised as a market leader, we have a reputation for high quality products
designed for hospital decontamination and sterile services, operating and
surgical applications, patient care and home nursing - where infection
management and good user application is paramount.

Warwick SASCo's low cost reusable and single use product ranges are chosen
by healthcare professionals and end users worldwide.

Our growing family owned business is proud to be represented by leading
distributor partners in over 60 countries.

All our products are manufactured within the UK with Quality Assured
Suppliers. Our customer support and distribution centre is in Warwick, England.

The team at Warwick SASCo is immensely proud of delivering high quality
products and a professional friendly service to you - our customer. You will find
our experienced team are knowledgeable, caring and quick to respond to
your requirements.



Darby Booth,
Managing Director – Warwick SASCo

Products Overview

Our wide range of reusable and disposable medical polyware products are chosen by healthcare professionals for their high quality and low cost.

Reusable Medical Polyware

WARWICK RANGE

The Warwick Range for use in Decontamination Services, Operating Rooms, Surgeries, Clinics and healthcare services. Products are designed for repeated washing, disinfection and sterilization.

The Warwick Range of products can be processed by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H2O2).

Products comply with BS 5452 for Hospital Hollowware where applicable. Material used is polypropylene unless otherwise stated.

Single Use Medical Polyware

STERILITE RANGE

Sterilizable lightweight polyware for single-use and disposable sets. SteriLite Range is designed for Custom Pack Manufacturers and Decontamination Services.

The Sterilite Range of products can be processed by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H2O2).

Material used is polypropylene unless otherwise stated.

Nursing & Care Daily Living Aids & Accessories

SPA RANGE

Spa Range of personal living aids, patient care, toileting and bedside products.

Products comply with BS 2588 for Portable Sanitary Products and BS 3215 for Reusable Urinals where applicable.


Material used is polypropylene unless otherwise stated.

We work with a wide range of industry sectors including:




STERILE SERVICES & DECONTAMINATION

Operating Room, Instrument Packing



PERSONAL USE & HOME CARE

Daily Living Aids, Community Nurse, Incontinent Care




HOSPITALS

Surgery, Clinics, Wards, Veterinary, Dentistry



NURSING HOMES

Patient Care, Rehabilitation



CUSTOM PACK MANUFACTURERS

Single Use Packs, Contract Instrument Processing

WARWICK - reusable range

Compartment Trays	page	6
DIN Style Baskets	page	7-8
Disinfectant & Soaking Containers	page	6
Gallipots	page	9-10
Instrument Trays	page	7-8
Instrument Tray Lids	page	9
Jugs	page	10
Kidney Dishes	page	11
Lotion & Sponge Bowls	page	9
Medical Boxes, Instrument Boxes & Storage Containers	page	12
Scope Containers	page	13
Silicone Protection	page	13
Theatre Bowls	page	11
Tray Tags	page	14
Wash Bowls	page	11

STERILITE - single use range

Compartment Trays	page	16
Gallipots	page	17
Jugs	page	18
Kidney Dishes	page	18
Packing Trays	page	16
Quivers	page	18
Sponge Bowls	page	17
Theatre Bowls	page	18

SPA - care and living aid range

Bedpans	page	20-21
Bedpan Holders	page	21
Bidets	page	21
Commode Pans	page	21
Denture Cups	page	22
Drinking Beakers and Cups	page	22
Drinking Cup Lids	page	22
Jugs	page	24
Jug Sets	page	24
Kidney Dishes	page	23
Medicine Measures	page	24
Patient Trays	page	23
Urinal Bottles	page	25
Urinal Pans	page	25
Vomit Bowls	page	23

How to Order

To place an order, request a quotation or for general enquiries:

CALL us on +44 (0)1926 422 427
 FAX +44 (0)1926 315 405
 EMAIL enquiries@sasco.co.uk



To find out more about our service, product ranges, product decontamination cycle and industry news:

visit www.sasco.co.uk

WARWICK RANGE

Reusable Medical Polyware

The Warwick Range is a wide range of sterilizable long-life plastic polyware products.

Designed as an alternative to stainless steel products, they are ideal for use in Decontamination Services, Operating Rooms, Clinics, Surgeries and other surgical or healthcare services.

Products are designed for repeated washing, disinfection and sterilization. The Warwick Range of products can be processed by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H2O2).

Products comply with BS 5452 for Hospital Holloware where applicable.

Material used is polypropylene unless otherwise stated.

Warwick Benefits

- Low cost.
- Robust and hard wearing - guaranteed for 1000 cycles.
- Easy to clean - blood protein is easily removed from the surface, leaving no residue.
- 50% lighter than stainless steel, assisting users with storage, lifting and handling.
- Can be sterilized by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H2O2).
- Can be used as standalone items or as a support component in surgical sets.
- Widely used in Sterile Services for packing instrument sets and procedure kits for the autoclave.
- Disinfection sets can be used for manual sterilization in Sterile Services, Clinics, etc.
- Used non-sterile in hospital wards and clinics.
- Smooth corners and surfaces, reducing risk of paper rips and the need for re-processing.
- Autoclave sterilization temperature is 134.6°C + 2°C (generally referred to as 135°C).

Compartment Trays, Disinfectant and Soaking Containers



Compartment Tray - 3 compartments

Code: CIT1814

Dimensions: 181x134x54mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 38

R W S E HP LF



Compartment Tray - 3 compartments

Code: CIT2013

Dimensions: 205x130x30mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 120

R W S E HP LF



Compartment Tray - 4 compartments

Code: CIT2718

Dimensions: 270x180x41mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 25

R W S E HP LF



Compartment Tray - 4 compartments

Code: LCT2215

Dimensions: 220x150x24mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 160

R W S E HP LF



Disinfection Set - Soaking container, strainer and hinged lid
Code: DSSW1000

Dimensions: 260x130x65mm
Volume: 1 litre
Material: Polypropylene
Colour: White
Pack Size: 10

R W S E HP LF



Disinfection Set - Soaking container, strainer and lid
Code: DSS2711

Dimensions: 257x180x128mm
Volume: 3litre
Material: Polypropylene
Colour: Blue
Pack Size: 5

R W S E HP LF



Disinfection Set - Soaking container, strainer and lid
Code: DSSW2711

Dimensions: 257x180x128mm
Volume: 3litre
Material: Polypropylene
Colour: White
Pack Size: 5

R W S E HP LF



Disinfection Set - Soaking container, strainer and vented lid
Code: DSSV2711

Dimensions: 257x180x128mm
Volume: 3litre
Material: Polypropylene
Colour: Blue
Pack Size: 5

R W S E HP LF



Disinfection Set - Soaking container, strainer, drip tray and lid
Code: DSS4215

Dimensions: 426x156x138mm
Volume: 4litre
Material: Polypropylene
Colour: Blue
Pack Size: 1

R W S E HP LF

Disinfection Set - Soaking container, strainer, drip tray and lid
Code: DSS8515

Dimensions: 890x200x145mm
Volume: 12litre
Material: Polypropylene
Colour: Blue
Pack Size: 1

R W S E HP LF



Disinfection Set - 2 x soaking containers, 1 x strainer, 2 x lids and 1 x drip tray
Code: DSS8515TP

Dimensions: 890x200x145mm
Volume: 12litre
Material: Polypropylene
Colour: Blue
Pack Size: 1

R W S E HP LF

R Reusable

W Washer Decontamination

S Suitable for Steam Sterilization

E Suitable for Ethylene Oxide Sterilization

HP Suitable for Hydrogen Peroxide Sterilization

LF Latex Free

DIN Style Baskets and Instrument Trays



Instrument Tray -
Solid base

Code: IT1813

Dimensions: 181x134x54mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 40

R W S E HP LF



Instrument Tray -
Solid base

Code: IT1913

Dimensions: 195x135x23mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 115

R W S E HP LF



Instrument Tray -
Solid base

Code: T2015

Dimensions: 200x150x51mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 60

R W S E HP LF



Instrument Tray -
Solid ribbed base

Code: IT2015

Dimensions: 200x150x51mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 60

R W S E HP LF



Instrument Tray -
Solid base

Code: T3025

Dimensions: 300x250x52mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 25

R W S E HP LF



Instrument Tray -
Solid base

Code: IT2715

Dimensions: 270x150x30mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 60

R W S E HP LF

How to order

To place an order,
request a quotation or for
general enquiries:

CALL us on
+44 (0)1926 422 427

FAX
+44 (0)1926 315 405

EMAIL
enquiries@sasco.co.uk



Instrument Tray -
Solid ribbed base

Code: IT3025

Dimensions: 300x250x52mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 25

R W S E HP LF



Instrument Tray -
Solid ribbed base

Code: IT4030

Dimensions: 424x305x75mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 12

R W S E HP LF

DIN Style Baskets and Instrument Trays



Instrument Tray -
Mesh base and sides

Code: MIT1515

Dimensions: 200x180x30mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 50

R W S E HP LF



Instrument Tray -
Mesh base

Code: MIT2015

Dimensions: 200x150x51mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 60

R W S E HP LF



Instrument Tray -
Mesh base and sides

Code: MIT2709

Dimensions: 275x110x30mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 60

R W S E HP LF



Instrument Tray -
Mesh base

Code: MIT2715

Dimensions: 270x150x30mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 60

R W S E HP LF



Instrument Tray -
Mesh base and sides

Code: MIT2817

Dimensions: 280x170x55mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 25

R W S E HP LF



Instrument Tray -
Mesh base and sides

Code: MIT3025

Dimensions: 300x250x52mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 24

R W S E HP LF



Instrument Tray -
Mesh base

Code: MIT4030

Dimensions: 424x305x75mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 12

R W S E HP LF



Instrument Tray -
Perforated base

Code: PIT3025

Dimensions: 300x250x52mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 25

R W S E HP LF



Instrument Tray -
Perforated base

Code: PIT4030

Dimensions: 424x305x75mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 12

R W S E HP LF



DIN Style Tray -
Mesh base and sides

Code: DINT2525

Dimensions: 255x255x50mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 20

R W S E HP LF



DIN Style Tray -
Mesh base and sides

Code: DINT4825

Dimensions: 480x255x50mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 25

R W S E HP LF

R Reusable W Washer Decontamination S Suitable for Steam Sterilization
E Suitable for Ethylene Oxide Sterilization HP Suitable for Hydrogen Peroxide Sterilization LF Latex Free

R Reusable W Washer Decontamination S Suitable for Steam Sterilization
E Suitable for Ethylene Oxide Sterilization HP Suitable for Hydrogen Peroxide Sterilization LF Latex Free

Instrument Tray Lids



Instrument Tray Lid -
For trays IT1813/CIT1814

Code: ITL181

Dimensions: N/A
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 10

R W S E HP LF



Instrument Tray Lid -
For trays IT2015/MIT2015/
T2015

Code: ITL201

Dimensions: N/A
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 10

R W S E HP LF



Instrument Tray Lid -
For trays IT2715/MIT2715

Code: ITL271

Dimensions: N/A
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 10

R W S E HP LF



Instrument Tray Lid -
For tray CIT2718

Code: ITL281

Dimensions: N/A
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 10

R W S E HP LF



Instrument Tray Lid -
For tray MIT2817

Code: ITL287

Dimensions: N/A
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 10

R W S E HP LF



Instrument Tray Lid -
For trays IT3025/MIT3025/
PIT3025/T3025

Code: ITL302

Dimensions: N/A
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 10

R W S E HP LF

How to order

To place an order, request a quotation or for
general enquiries:

CALL us on

+44 (0)1926 422 427

FAX

+44 (0)1926 315 405

EMAIL

enquiries@sasco.co.uk

Gallipots, Lotion and Sponge Bowls



Gallipot

Code: GP40

Dimensions: 40x30mm
Volume: 35ml
Material: Polypropylene
Colour: Blue
Pack Size: 100

R W S E HP LF



Gallipot

Code: GP60

Dimensions: 60x40mm
Volume: 95ml
Material: Polypropylene
Colour: Blue
Pack Size: 300

R W S E HP LF



Gallipot

Code: GP80

Dimensions: 80x60mm
Volume: 280ml
Material: Polypropylene
Colour: Blue
Pack Size: 80

R W S E HP LF



Gallipot -
Flat base

Code: GPLB60

Dimensions: 60x40mm
Volume: 50ml
Material: Polypropylene
Colour: Blue
Pack Size: 400

R W S E HP LF

Gallipots, Lotion and Sponge Bowls



Gallipot -
Flat base

Code: GPLB80

Dimensions: 80x60mm
Volume: 200ml
Material: Polypropylene
Colour: Blue
Pack Size: 200

R W S E HP LF



Lotion Bowl -
Flat base

Code: GPLB100

Dimensions: 100x45mm
Volume: 300ml
Material: Polypropylene
Colour: Blue
Pack Size: 135

R W S E HP LF



Lotion Bowl -
Flat base

Code: GPLB150

Dimensions: 150x70mm
Volume: 500ml
Material: Polypropylene
Colour: Blue
Pack Size: 60

R W S E HP LF



Sponge Bowl -
Flat base

Code: SB140

Dimensions: 140x65mm
Volume: 300ml
Material: Polypropylene
Colour: Blue
Pack Size: 160

R W S E HP LF



Lotion Bowl -
Graduated

Code: LB100

Dimensions: 100x45mm
Volume: 280ml
Material: Polypropylene
Colour: Blue
Pack Size: 200

R W S E HP LF



Lotion Bowl

Code: LB150

Dimensions: 150x70mm
Volume: 900ml
Material: Polypropylene
Colour: Blue
Pack Size: 60

R W S E HP LF



Lotion Bowl

Code: LB200

Dimensions: 200x90mm
Volume: 2000ml
Material: Polypropylene
Colour: Blue
Pack Size: 30

R W S E HP LF



Lotion Bowl

Code: LB250

Dimensions: 250x105mm
Volume: 3000ml
Material: Polypropylene
Colour: Blue
Pack Size: 20

R W S E HP LF

Gallipots, Lotion and Sponge Bowls



Gallipot

Code: GP40

Dimensions: 40x30mm
Volume: 35ml
Material: Polypropylene
Colour: Blue
Pack Size: 100

R W S E HP LF



Gallipot

Code: GP60

Dimensions: 60x40mm
Volume: 95ml
Material: Polypropylene
Colour: Blue
Pack Size: 300

R W S E HP LF



Gallipot

Code: GP80

Dimensions: 80x60mm
Volume: 280ml
Material: Polypropylene
Colour: Blue
Pack Size: 80

R W S E HP LF



Gallipot -
Flat base

Code: GPLB60

Dimensions: 60x40mm
Volume: 50ml
Material: Polypropylene
Colour: Blue
Pack Size: 400

R W S E HP LF



Jug

Code: JG500

Dimensions: N/A
Volume: 500ml
Material: Polypropylene
Colour: Natural
Pack Size: 50

R W S E HP LF



Jug

Code: JG1000

Dimensions: N/A
Volume: 1000ml
Material: Polypropylene
Colour: Natural
Pack Size: 45

R W S E HP LF



Jug

Code: JG2000

Dimensions: N/A
Volume: 2000ml
Material: Polypropylene
Colour: Natural
Pack Size: 20

R W S E HP LF

Kidney Dishes, Theatre Bowls and Wash Bowls



Kidney Dish

Code: KD150

Dimensions: 150x35mm
Volume: 300ml
Material: Polypropylene
Colour: Blue
Pack Size: 210

R W S E HP LF



Kidney Dish

Code: KD200

Dimensions: 200x45mm
Volume: 500ml
Material: Polypropylene
Colour: Blue
Pack Size: 150

R W S E HP LF



Kidney Dish

Code: KD250

Dimensions: 250x55mm
Volume: 750ml
Material: Polypropylene
Colour: Blue
Pack Size: 60

R W S E HP LF



Kidney Dish

Code: KDY250

Dimensions: 250x55mm
Volume: 750ml
Material: Polypropylene
Colour: Yellow
Pack Size: 60

R W S E HP LF



Kidney Dish

Code: KD300

Dimensions: 300x65mm
Volume: 1500ml
Material: Polypropylene
Colour: Blue
Pack Size: 25

R W S E HP LF



Kidney Dish Tray -
Shallow kidney shaped tray

Code: KDT185

Dimensions: 185x92mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 70

R W S E HP LF



Washbowl -
Graduated

Code: WB3112

Dimensions: 310x105mm
Volume: 3000ml
Material: Polypropylene
Colour: Blue
Pack Size: 25

R W S E HP LF



Wash Bowl

Code: WBN3112

Dimensions: 310x105mm
Volume: 3000ml
Material: Polypropylene
Colour: Natural
Pack Size: 25

R W S E HP LF



Theatre Bowl

Code: WB3513

Dimensions: 350x120mm
Volume: 5000ml
Material: Polypropylene
Colour: Blue
Pack Size: 25

R W S E HP LF



Wash Bowl

Code: WBN3513

Dimensions: 350x120mm
Volume: 5000ml
Material: Polypropylene
Colour: Natural
Pack Size: 25

R W S E HP LF

R Reusable W Washer Decontamination S Suitable for Steam Sterilization
E Suitable for Ethylene Oxide Sterilization HP Suitable for Hydrogen Peroxide Sterilization LF Latex Free

Medical Boxes, Instrument Boxes & Storage Containers



Instrument Box -
Storage box with lid

Code: IBX2718

Dimensions: 257x180x128mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 20

R W S E HP LF



Instrument box -
Storage box with vented lid

Code: IBXV2718

Dimensions: 257x180x128mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 20

R W S E HP LF



Instrument Box -
Storage box with lid

Code: IBXN290

Dimensions: 300x224x140mm
Volume: N/A
Material: Polypropylene
Colour: Natural
Pack Size: 12

R W S E HP LF

How to order

To place an order,
request a quotation or for
general enquiries:

CALL us on
+44 (0)1926 422 427

FAX
+44 (0)1926 315 405

EMAIL
enquiries@sasco.co.uk



Transportation Box -
With hinged lid

Code: MB2318G

Dimensions: 230x170x200mm
Volume: N/A
Material: ABS
Colour: Green
Pack Size: 10

R W S E HP LF



Transportation Box -
With hinged lid


Code: MB2318

Dimensions: 230x170x200mm
Volume: N/A
Material: ABS
Colour: Orange
Pack Size: 10

R W S E HP LF

R Reusable W Washer Decontamination S Suitable for Steam Sterilization
E Suitable for Ethylene Oxide Sterilization HP Suitable for Hydrogen Peroxide Sterilization LF Latex Free

Scope Containers



Scope Container -
With silicone dividers and lid

Code: SMT5510

Dimensions: 660x110x66mm
Volume: N/A
Material: Polypropylene
Colour: Lilac
Pack Size: 1

R W S E HP LF



Silicone Plug -
For SMT5510

Code: SG5510

Dimensions: N/A
Volume: N/A
Material: Silicone
Colour: Pink
Pack Size: 10

R W S E HP LF

Silicone Protection



Silicone Pad -
For trays IT/MIT2715 and
small metal DIN tray
Code: SP2313

Dimensions: 230x132mm
Volume: N/A
Material: Silicone
Colour: Pink
Pack Size: 10

R W S E HP LF



Silicone Pad -
For trays IT/MIT/PIT/T3025

Code: SP2823

Dimensions: 280x230mm
Volume: N/A
Material: Silicone
Colour: Pink
Pack Size: 10

R W S E HP LF



Silicone Pad -
Fits 1/2 tray IT/MIT/PIT4030

Code: SP2719

Dimensions: 270x190mm
Volume: N/A
Material: Silicone
Colour: Pink
Pack Size: 10

R W S E HP LF



Silicone Pad -
Fits 1/2 tray DINT4825

Code: SP4622

Dimensions: 230x220mm
Volume: N/A
Material: Silicone
Colour: Pink
Pack Size: 10

R W S E HP LF



Silicone Pad -
For medium metal DIN tray

Code: SP2324

Dimensions: 235x245mm
Volume: N/A
Material: Silicone
Colour: Pink
Pack Size: 10

R W S E HP LF




Silicone Pad -
For large metal DIN tray

Code: SP4724

Dimensions: 470x245mm
Volume: N/A
Material: Silicone
Colour: Pink
Pack Size: 10

R W S E HP LF



Silicone Roll -
1 metre length with 10
divisible sections
Code: SR10

Dimensions: 100x245mm
Volume: N/A
Material: Silicone
Colour: Pink
Pack Size: 1

R W S E HP LF

Tray Tags



Tray Tag

Code: TAG-010

Dimensions: 105x40mm
Volume: N/A
Material: Nylon
Colour: White
Pack Size: 100

R W S E HP LF



Tray Tag

Code: TAG-015

Dimensions: 105x40mm
Volume: N/A
Material: Nylon
Colour: Yellow
Pack Size: 100

R W S E HP LF



Tray Tag

Code: TAG-020

Dimensions: 105x40mm
Volume: N/A
Material: Nylon
Colour: Blue
Pack Size: 100

R W S E HP LF



Tray Tag

Code: TAG-022

Dimensions: 105x40mm
Volume: N/A
Material: Nylon
Colour: Light Green
Pack Size: 100

R W S E HP LF



Tray Tag

Code: TAG-025

Dimensions: 105x40mm
Volume: N/A
Material: Nylon
Colour: Orange
Pack Size: 100

R W S E HP LF



Tray Tag

Code: TAG-030

Dimensions: 105x40mm
Volume: N/A
Material: Nylon
Colour: Red
Pack Size: 100

R W S E HP LF



Tray Tag

Code: TAG-050

Dimensions: 105x40mm
Volume: N/A
Material: Nylon
Colour: Green
Pack Size: 100

R W S E HP LF



Tray Tag

Code: TAG-055

Dimensions: 105x40mm
Volume: N/A
Material: Nylon
Colour: Grey
Pack Size: 100

R W S E HP LF



Tray Tag

Code: TAG-060

Dimensions: 105x40mm
Volume: N/A
Material: Nylon
Colour: Lavender
Pack Size: 100

R W S E HP LF



Tray Tag

Code: TAG-065

Dimensions: 105x40mm
Volume: N/A
Material: Nylon
Colour: Tan
Pack Size: 100

R W S E HP LF



Tray Tag

Code: TAG-070

Dimensions: 105x40mm
Volume: N/A
Material: Nylon
Colour: Black
Pack Size: 100

R W S E HP LF



Stainless Steel Attachment Ring

Code: TR30

Dimensions: 30mm
Volume: N/A
Material: Steel
Colour: Silver
Pack Size: 100

R W S E HP LF

R Reusable

W Washer Decontamination

S Suitable for Steam Sterilization

E Suitable for Ethylene Oxide Sterilization

HP Suitable for Hydrogen Peroxide Sterilization

LF Latex Free

STERILITE RANGE

Single Use Medical Polyware

The SteriLite Range is a wide range of single use plastic polyware products.

Designed to speed up packing and sterilizing procedures. SteriLite is ideal for use in Sterile Pack Manufacturing, Decontamination Services, Operating Rooms and Clinics due to its low cost, product design and performance in hospitals.

The SteriLite Range of products can be processed by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H2O2).

Material used is polypropylene unless otherwise stated.

SteriLite Benefits

- Low cost and disposable.
- No need to invest time and money into retrieval and decontamination after use.
- Clean and ready for sterilization.
- Can be sterilized by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H2O2).
- Can be used as stand-alone items or as a support (supplementary) component in operating rooms and surgical procedures.
- Widely used in Sterile Services for packing instrument sets and procedure kits.
- Used non-sterile in hospital wards and clinics.



Compartment Trays and Packing Trays



Packing Tray -
4 compartments

Code: LCT2215S

Dimensions: 220x150x24mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 160



Packing Tray

Code: ST4330

Dimensions: 430x330x64mm
Volume: N/A
Material: Polypropylene
Colour: Blue
Pack Size: 63



- Single Use
- Washer Decontamination
- Suitable for Steam Sterilization
- Suitable for Ethylene Oxide Sterilization
- Suitable for Hydrogen Peroxide Sterilization
- Latex Free

Gallipots and Sponge Bowls



Gallipot

Code: GPS60

Dimensions: 62x30mm
Volume: 60ml
Material: Polypropylene
Colour: Natural
Pack Size: 2000



Gallipot

Code: GPS120

Dimensions: 85x30mm
Volume: 120ml
Material: Polypropylene
Colour: Natural
Pack Size: 1000



Gallipot

Code: GP100

Dimensions: 100x53mm
Volume: 250ml
Material: Polypropylene
Colour: Natural
Pack Size: 750



Gallipot

Code: GP130

Dimensions: 130x66mm
Volume: 500ml
Material: Polypropylene
Colour: Natural
Pack Size: 735



Sponge Bowl

Code: SB170

Dimensions: 170x84mm
Volume: 1000ml
Material: Polypropylene
Colour: Natural
Pack Size: 340



Sponge Bowl

Code: SB200

Dimensions: 210x90mm
Volume: 2000ml
Material: Polypropylene
Colour: Natural
Pack Size: 220



Jugs and Quivers



Jug

Code: J1000

Dimensions: N/A
Volume: 1000ml
Material: Polypropylene
Colour: Natural
Pack Size: 100



Quiver

Code: QV30

Dimensions: 190x65mm
Volume: N/A
Material: Polypropylene
Colour: Natural
Pack Size: 75



Kidney Dishes and Theatre Bowls



Kidney Dish

Code: KB250

Dimensions: 252x100x48mm
Volume: 750ml
Material: Polypropylene
Colour: Natural
Pack Size: 450



Kidney Dish

Code: KB300

Dimensions: 300x140x57mm
Volume: 1500ml
Material: Polypropylene
Colour: Natural
Pack Size: 135



Theatre Bowl

Code: TB350

Dimensions: 352x116mm
Volume: 5000ml
Material: Polypropylene
Colour: Blue
Pack Size: 50



How to order

To place an order,
request a quotation or for
general enquiries:

CALL us on
+44 (0)1926 422 427

FAX
+44 (0)1926 315 405

EMAIL
enquiries@sasco.co.uk

Single Use
 Washer Decontamination
 Suitable for Steam Sterilization
 Suitable for Ethylene Oxide Sterilization
 Suitable for Hydrogen Peroxide Sterilization
 Latex Free

Single Use
 Washer Decontamination
 Suitable for Steam Sterilization
 Suitable for Ethylene Oxide Sterilization
 Suitable for Hydrogen Peroxide Sterilization
 Latex Free

SPA RANGE

Nursing & Care, Daily Living Aids and Accessories

The Spa range includes reusable products for personal living, patient care, toileting and bedside comfort.

The Spa Range is designed for use in Hospital Wards, Rehabilitation, Nursing Homes, Community Care and for Personal Use in the Home.

Products comply with BS 2588 for Portable Sanitary Products and BS 3215 for Reusable Urinals where applicable.

Material used is polypropylene unless otherwise stated.

Spa Benefits

- Toileting products that preserve patient dignity and independence by empowering patients to toilet themselves with little or no assistance.
- Urinal Bottles specifically designed to the male and female anatomy for a comfortable and safe fit during use.
- Product design reduces risk of cross infection, unpleasant odours and sights.
- Easy to use, clean and disinfect, meaning both patient and carer can be confident that products are safe to re-use.
- Beakers and Drinking Cups are wide, easy to grip, non-scratch, drop-proof and have a nondrip rim.
- Handled beakers are ideal for patients with unsteady or infirm hands.
- Medicine Measures are sturdy, ideal for long-term use and regularly checked by Trading Standards for measurement accuracy and compliance.
- All Spa products can be cleaned in automatic washer disinfectors. Also suitable for hand washing using detergents and disinfectants.
- Urinal Pan SPN and Bedpans FPN, VPN, DPNF, BPN and CMP can be autoclaved.

Bedpans

Bedpan - Adult

Code: BPN

Dimensions: 370x310x90mm
Volume: 2500ml
Material: Polypropylene
Colour: Green
Pack Size: 30

R W S E HP LF



How to order

To place an order, request a quotation or for general enquiries:

CALL us on
+44 (0)1926 422 427
FAX
+44 (0)1926 315 405

EMAIL
enquiries@sasco.co.uk



Slipper Pan - Bedpan with handle and lid

Code: FPN

Dimensions: 449x309x99mm
Volume: 1500ml
Material: Polypropylene
Colour: Green
Pack Size: 12

R W S E HP LF



Slipper Pan - Bedpan with handle and lid

Code: FPNW

Dimensions: 449x309x99mm
Volume: 1500ml
Material: Polypropylene
Colour: White
Pack Size: 12

R W S E HP LF



Bedpan - Hospital style bedpan

Code: VP

Dimensions: 355x308x108mm
Volume: 2500ml
Material: Polypropylene
Colour: Green
Pack Size: 18

R W S E HP LF



Bedpan - Hospital bedpan with lid

Code: VPN

Dimensions: 355x308x108mm
Volume: 2500ml
Material: Polypropylene
Colour: Green
Pack Size: 15

R W S E HP LF



Bedpan - Hospital bedpan with lid

Code: VPNW

Dimensions: 355x308x108mm
Volume: 2500ml
Material: Polypropylene
Colour: White
Pack Size: 15

R W S E HP LF



Bedpan - European bedpan with handle and lid

Code: DPNF

Dimensions: 420x315x100mm
Volume: 1500ml
Material: Polypropylene
Colour: Green
Pack Size: 15

R W S E HP LF

R Reusable W Washer Decontamination S Suitable for Steam Sterilization
E Suitable for Ethylene Oxide Sterilization HP Suitable for Hydrogen Peroxide Sterilization LF Latex Free

Bedpans, Commode Pans, Bidets and Bedpan Holders

Slipper Pan -
Female urinal and
infants bedpan
Code: SPN

Dimensions: 335x245mm
Volume: 1500ml
Material: Polypropylene
Colour: Green
Pack Size: 25

R W S E HP LF



Commode Pan / Bidet

Code: CMP

Dimensions: 359x338x130mm
Volume: 3000ml
Material: Polypropylene
Colour: Green
Pack Size: 20

R W S E HP LF



Commode Pan / Bidet

Code: CMPW

Dimensions: 359x338x130mm
Volume: 3000ml
Material: Polypropylene
Colour: White
Pack Size: 20

R W S E HP LF



Commode Pan -
Commode pan with lid

Code: CBL

Dimensions: 259x259x144mm
Volume: 2500ml
Material: Polypropylene
Colour: White
Pack Size: 25

R W S E HP LF



Bedpan Support -
For pulp disposable bedpan

Code: PBPS

Dimensions: 370x310x90mm
Volume: 2500ml
Material: Polypropylene
Colour: Cream
Pack Size: 30

R W S E HP LF



Bedpan Support -
For pulp disposable midi
slipper pan
Code: PSPS

Dimensions: 470x290x75mm
Volume: 1500ml
Material: Polypropylene
Colour: Cream
Pack Size: 30

R W S E HP LF

How to order

To place an order,
request a quotation or for
general enquiries:

CALL us on
+44 (0)1926 422 427
FAX
+44 (0)1926 315 405

EMAIL
enquiries@sasco.co.uk

Single Use

W Washer Decontamination

S Suitable for Steam Sterilization

E Suitable for Ethylene Oxide Sterilization HP Suitable for Hydrogen Peroxide Sterilization LF Latex Free

Drinking Beakers, Cups, Lids and Denture Cups



Beaker -
Feeder/drinking cup

Code: BFC250

Dimensions: N/A
Volume: 250ml
Material: Polypropylene
Colour: Natural
Pack Size: 300

R W S E HP LF



Beaker -
Feeder/drinking cup with
handles

Code: BFH250

Dimensions: N/A
Volume: 250ml
Material: Polypropylene
Colour: Natural
Pack Size: 140

R W S E HP LF



Lid -
Narrow Spout for BFC250/
BFH250

Code: BNS

Dimensions: 4x4mm spout
Volume: N/A
Material: Polypropylene
Colour: Natural
Pack Size: 300

R W S E HP LF



Lid -
Medium Spout for BFC250/
BFH250

Code: BMS

Dimensions: 7x4mm spout
Volume: N/A
Material: Polypropylene
Colour: Natural
Pack Size: 300

R W S E HP LF



Lid -
Wide spout for BFC250
& BFH250

Code: BWS

Dimensions: 13x8mm spout
Volume: N/A
Material: Polypropylene
Colour: Natural
Pack Size: 300

R W S E HP LF



Lid -
Rim feeder for beakers BFC250
& BFH250

Code: BRF

Dimensions: N/A
Volume: N/A
Material: Polypropylene
Colour: Natural
Pack Size: 200

R W S E HP LF



Drinking Cup -
With handles, spout lid and
feeder lid

Code: GPMC6250

Dimensions: N/A
Volume: 250ml
Material: Polypropylene
Colour: Cream
Pack Size: 90

R W S E HP LF



Denture Cup -
With strainer and lid

Code: DNT

Dimensions: 90x80x60mm
Volume: N/A
Material: Polypropylene
Colour: Green
Pack Size: 75

R W S E HP LF



Beaker Set -
Feeder/drinking cup with han-
dle, spout lid and feeder lid
Code: HDMS

Dimensions: N/A
Volume: 300ml
Material: SAN
Colour: Clear
Pack Size: 40

R W S E HP LF



Teapot Style Beaker

Code: TFB

Dimensions: N/A
Volume: 250ml
Material: Polypropylene
Colour: White
Pack Size: 45

R W S E HP LF



Clear Glass Style Beaker

Code: TBC

Dimensions: N/A
Volume: 280ml
Material: SAN
Colour: Clear
Pack Size: 120

R W S E HP LF

R Reusable

W Washer Decontamination

S Suitable for Steam Sterilization

E Suitable for Ethylene Oxide Sterilization HP Suitable for Hydrogen Peroxide Sterilization LF Latex Free

Patient Trays, Kidney Dishes and Vomit Bowls



Dispensing/Dressing Tray

Code: ITC2715

Dimensions: 270x150x30mm
Volume: N/A
Material: Polypropylene
Colour: Green
Pack Size: 60

R W S E HP LF



Kidney Dish

Code: KDC250

Dimensions: 250x130mm
Volume: 750ml
Material: Polypropylene
Colour: Green
Pack Size: 60

R W S E HP LF



Vomit Bowl - With handle

Code: VB300

Dimensions: 300x144mm
Volume: 1500ml
Material: Polypropylene
Colour: Green
Pack Size: 25

R W S E HP LF

R Reusable W Washer Decontamination S Suitable for Steam Sterilization
E Suitable for Ethylene Oxide Sterilization HP Suitable for Hydrogen Peroxide Sterilization LF Latex Free

Jugs and Jug Sets



Jug

Code: JG500

Dimensions: N/A
Volume: 500ml
Material: Polypropylene
Colour: Natural
Pack Size: 50

R W S E HP LF



Jug

Code: JG1000

Dimensions: N/A
Volume: 1000ml
Material: Polypropylene
Colour: Natural
Pack Size: 45

R W S E HP LF



Jug with lid

Code: JGL1000

Dimensions: N/A
Volume: 1000ml
Material: Polypropylene
Colour: Natural
Pack Size: 40

R W S E HP LF



Jug Set - Includes tray, carafe with lid and beaker
Code: LST

Dimensions: N/A
Volume: N/A
Material: Polypropylene
Colour: N/A
Pack Size: 20

R W S E HP LF

How to order

To place an order, request a quotation or for general enquiries:

CALL us on
+44 (0)1926 422 427
FAX
+44 (0)1926 315 405

EMAIL
enquiries@sasco.co.uk

Medicine Measures



Medicine Measure - Graduated

Code: MM50

Dimensions: N/A
Volume: 50ml
Material: Polypropylene
Colour: Natural
Pack Size: 1000

R W S E HP LF



Medicine Measure - Graduated

Code: MM60

Dimensions: N/A
Volume: 60ml
Material: Polypropylene
Colour: Natural
Pack Size: 1000

R W S E HP LF



Medicine Measure - Graduated black lines

Code: MM60BL

Dimensions: N/A
Volume: 60ml
Material: Polypropylene
Colour: Natural
Pack Size: 1000

R W S E HP LF

Urinal Bottles and Urinal Pans



Slipper Pan - Female urinal

Code: SPN

Dimensions: 335x245mm
Volume: 1500ml
Material: Polypropylene
Colour: Green
Pack Size: 25

R W S E HP LF



Urinal Bottle - Female urinal with handle

Code: UF1000

Dimensions: N/A
Volume: 1000ml
Material: Polyethylene
Colour: Natural
Pack Size: 30

R W S E HP LF



Urinal Pan - Female Urinal

Code: PUF1000

Dimensions: N/A
Volume: 1000ml
Material: Polyethylene
Colour: White
Pack Size: 20

R W S E HP LF



Urinal Bottle - Male urinal with handle

Code: UMA1000

Dimensions: N/A
Volume: 1000ml
Material: Polyethylene
Colour: Natural
Pack Size: 40

R W S E HP LF



Urinal Bottle - Male urinal with handle, hospital style

Code: UBS1000

Dimensions: N/A
Volume: 1200ml
Material: Polyethylene
Colour: Natural
Pack Size: 45

R W S E HP LF



Urinal Bottle - Male urinal with handle

Code: UML1000

Dimensions: N/A
Volume: 1000ml
Material: Polyethylene
Colour: Natural
Pack Size: 50

R W S E HP LF



Urinal Bottle - Male urinal, hospital style

Code: UM1000

Dimensions: N/A
Volume: 1000ml
Material: Polypropylene
Colour: Natural
Pack Size: 25

R W S E HP LF



Urinal Bottle - Male urinal, hospital style

Code: UM650

Dimensions: N/A
Volume: 650ml
Material: Polypropylene
Colour: Natural
Pack Size: 30

R W S E HP LF



Non Leak Valve - For use with UMA1000 & UML1000

Code: UMT

Dimensions: N/A
Volume: N/A
Material: Polypropylene
Colour: White
Pack Size: 25

R W S E HP LF



Urinal Bottle with UMT - Male urinal bottle with non leak valve

Code: USP1200

Dimensions: N/A
Volume: 1000ml
Material: Polypropylene and Polyethylene
Colour: N/A
Pack Size: 10

R W S E HP LF



Urinal Bottle Holder - For UBS1000, UML1000 and UM1000

Code: UBH2002

Dimensions: N/A
Volume: 1200ml
Material: Plastic Coated Wire
Colour: White
Pack Size: 10

R W S E HP LF

How to order

To place an order, request a quotation or for general enquiries:

CALL us on

+44 (0)1926 422 427

FAX

+44 (0)1926 315 405

EMAIL

enquiries@sasco.co.uk

R Reusable W Washer Decontamination S Suitable for Steam Sterilization
E Suitable for Ethylene Oxide Sterilization HP Suitable for Hydrogen Peroxide Sterilization LF Latex Free



Warwick House, Heathcote Way, Warwick, CV34 6TE

T: +00 44 (0) 1926 422 427

F: +00 44 (0) 1926 315 405

E: enquiries@sasco.co.uk

www.sasco.co.uk