



*Stainless Steel*

**MEDICAL - HOSPITAL WARES**





## SURGICAL TRAYS



### SURGICAL TRAYS

|          |                          | SIZE (mm.) |     |     |
|----------|--------------------------|------------|-----|-----|
|          |                          | L          | W   | H   |
| UPL-010  | SURGICAL TRAY with cover | 209        | 158 | 76  |
| UPL-011  | SURGICAL TRAY with cover | 223        | 126 | 45  |
| UPL-011D | SURGICAL TRAY with cover | 223        | 126 | 76  |
| UPL-012  | SURGICAL TRAY with cover | 306        | 196 | 50  |
| UPL-0125 | SURGICAL TRAY with cover | 306        | 196 | 15  |
| UPL-013  | SURGICAL TRAY with cover | 264        | 172 | 47  |
| UPL-013A | SURGICAL TRAY with cover | 270        | 218 | 50  |
| UPL-014  | SURGICAL TRAY with cover | 440        | 320 | 64  |
| UPL-015  | SURGICAL TRAY with cover | 355        | 254 | 50  |
| UPL-016  | SURGICAL TRAY with cover | 380        | 304 | 50  |
| UPL-017  | SURGICAL TRAY with cover | 306        | 196 | 63  |
| UPL-018  | SURGICAL TRAY with cover | 323        | 265 | 63  |
| UPL-018D | SURGICAL TRAY with cover | 323        | 265 | 100 |
| UPL-019  | SURGICAL TRAY with cover | 370        | 320 | 63  |
| UPL-019D | SURGICAL TRAY with cover | 370        | 320 | 100 |

## SURGICAL TRAYS w/o Cover



### SURGICAL TRAYS w/o Cover

|          |                             | SIZE (mm.) |     |     |
|----------|-----------------------------|------------|-----|-----|
|          |                             | L          | W   | H   |
| UPL-020  | SURGICAL TRAY without cover | 209        | 158 | 76  |
| UPL-021  | SURGICAL TRAY without cover | 223        | 126 | 45  |
| UPL-021D | SURGICAL TRAY without cover | 223        | 126 | 76  |
| UPL-022  | SURGICAL TRAY without cover | 306        | 196 | 50  |
| UPL-0225 | SURGICAL TRAY without cover | 306        | 196 | 15  |
| UPL-023  | SURGICAL TRAY without cover | 264        | 172 | 47  |
| UPL-023A | SURGICAL TRAY without cover | 270        | 218 | 50  |
| UPL-024  | SURGICAL TRAY without cover | 440        | 320 | 64  |
| UPL-025  | SURGICAL TRAY without cover | 355        | 254 | 50  |
| UPL-026  | SURGICAL TRAY without cover | 380        | 304 | 50  |
| UPL-027  | SURGICAL TRAY without cover | 306        | 196 | 63  |
| UPL-028  | SURGICAL TRAY without cover | 323        | 265 | 63  |
| UPL-028D | SURGICAL TRAY without cover | 323        | 265 | 100 |
| UPL-029  | SURGICAL TRAY without cover | 370        | 320 | 63  |
| UPL-029D | SURGICAL TRAY without cover | 370        | 320 | 100 |

## INSTRUMENT TRAYS



### INSTRUMENT TRAYS

|         |                               | SIZE (mm.) |     |    |
|---------|-------------------------------|------------|-----|----|
|         |                               | L          | W   | H  |
| UPL-031 | INSTRUMENT TRAY without cover | 210        | 160 | 25 |
| UPL-032 | INSTRUMENT TRAY without cover | 268        | 190 | 30 |
| UPL-033 | INSTRUMENT TRAY without cover | 313        | 218 | 31 |
| UPL-034 | INSTRUMENT TRAY without cover | 365        | 265 | 33 |

\* Note : Surgical, Instrument Trays also available with perforation

## DENTAL TRAYS



|         |             | DENTAL TRAYS |            |    |
|---------|-------------|--------------|------------|----|
|         |             | L            | SIZE (mm.) |    |
|         |             |              | W          | H  |
| UPL-036 | DENTAL TRAY | 208          | 109        | 15 |
| UPL-037 | DENTAL TRAY | 209          | 150        | 15 |

\* Note : Also available with perforation.

## CATHETER TRAYS



|         |                          | CATHETER TRAYS |            |    |
|---------|--------------------------|----------------|------------|----|
|         |                          | L              | SIZE (mm.) |    |
|         |                          |                | W          | H  |
| UPL-041 | CATHETER TRAY with cover | 217            | 78         | 38 |
| UPL-042 | CATHETER TRAY with cover | 318            | 91         | 44 |
| UPL-043 | CATHETER TRAY with cover | 400            | 75         | 49 |
| UPL-044 | CATHETER TRAY with cover | 438            | 113        | 54 |

## MAYO TRAYS



|          |                   | MAYO TRAYS |            |    |
|----------|-------------------|------------|------------|----|
|          |                   | L          | SIZE (mm.) |    |
|          |                   |            | W          | H  |
| UPL-045  | MAYO TRAY shallow | 254        | 165        | 18 |
| UPL-046  | MAYO TRAY shallow | 356        | 252        | 16 |
| UPL-046A | MAYO TRAY shallow | 432        | 295        | 19 |
| UPL-047  | MAYO TRAY shallow | 487        | 323        | 13 |
| UPL-047D | MAYO TRAY shallow | 487        | 323        | 25 |
| UPL-048  | MAYO TRAY shallow | 542        | 418        | 22 |

\* Note : Also available with perforation.

## KIDNEY TRAYS/EMESIS BASIN



|                                  |                                | KIDNEY TRAYS/EMESIS BASIN |            |    |
|----------------------------------|--------------------------------|---------------------------|------------|----|
|                                  |                                | L                         | SIZE (mm.) |    |
|                                  |                                |                           | W          | H  |
| <b>DEEP</b> [with/without cover] |                                |                           |            |    |
| UPL-051                          | KIDNEY TRAY 6"                 | 162                       | 77         | 31 |
| UPL-052                          | KIDNEY TRAY 8"                 | 207                       | 98         | 39 |
| UPL-053                          | KIDNEY TRAY 10"                | 247                       | 122        | 43 |
| UPL-053A                         | KIDNEY TRAY 10" extra deep     | 263                       | 124        | 53 |
| UPL-054                          | KIDNEY TRAY 12"                | 309                       | 149        | 59 |
| <b>SHALLOW</b>                   |                                |                           |            |    |
| UPL-071                          | KIDNEY TRAY shallow            | 207                       | 128        | 33 |
| UPL-072                          | KIDNEY TRAY shallow            | 254                       | 141        | 33 |
| UPL-072F                         | KIDNEY TRAY shallow [flatbase] | 254                       | 141        | 33 |



## GALLIPOT J-TYPE



### GALLIPOT J-TYPE

|         |                 |       | Cap. (ml.) | SIZE (mm.) |    |
|---------|-----------------|-------|------------|------------|----|
|         |                 |       |            | Dia Ø      | H  |
| UPL081  | GALLIPOT J-TYPE | 2 oz  | 60         | 60         | 33 |
| UPL082  | GALLIPOT J-TYPE | 3 oz  | 90         | 69         | 41 |
| UPL083  | GALLIPOT J-TYPE | 4 oz  | 120        | 75         | 42 |
| UPL084  | GALLIPOT J-TYPE | 6 oz  | 180        | 83         | 51 |
| UPL084A | GALLIPOT J-TYPE | 8 oz  | 240        | 91         | 54 |
| UPL084B | GALLIPOT J-TYPE | 10 oz | 300        | 97         | 54 |
| UPL085  | GALLIPOT J-TYPE | 12 oz | 360        | 110        | 56 |
| UPL085A | GALLIPOT J-TYPE | 16 oz | 480        | 122        | 56 |
| UPL085B | GALLIPOT J-TYPE | 18 oz | 540        | 122        | 63 |
| UPL086  | GALLIPOT J-TYPE | 22 oz | 660        | 143        | 58 |
| UPL086A | GALLIPOT J-TYPE | 28 oz | 840        | 143        | 70 |
| UPL086B | GALLIPOT J-TYPE | 30 oz | 900        | 163        | 58 |
| UPL086C | GALLIPOT J-TYPE | 36 oz | 1080       | 163        | 68 |

## GALLIPOT with LIP



### GALLIPOT with LIP

|         |                               |  | Cap. (ml.) | SIZE (mm.) |    |
|---------|-------------------------------|--|------------|------------|----|
|         |                               |  |            | Dia Ø      | H  |
| UPL087  | GALLIPOT LOTION BOWL with lip |  | 30         | 46         | 29 |
| UPL088  | GALLIPOT LOTION BOWL with lip |  | 50         | 56         | 32 |
| UPL089  | GALLIPOT LOTION BOWL with lip |  | 90         | 66         | 36 |
| UPL090  | GALLIPOT LOTION BOWL with lip |  | 120        | 76         | 41 |
| UPL091  | GALLIPOT LOTION BOWL with lip |  | 170        | 88         | 46 |
| UPL091A | GALLIPOT LOTION BOWL with lip |  | 250        | 95         | 51 |
| UPL092  | GALLIPOT LOTION BOWL with lip |  | 325        | 107        | 56 |
| UPL093  | GALLIPOT LOTION BOWL with lip |  | 450        | 128        | 60 |
| UPL093A | GALLIPOT LOTION BOWL with lip |  | 600        | 139        | 65 |
| UPL094  | GALLIPOT LOTION BOWL with lip |  | 940        | 158        | 75 |
| UPL095  | GALLIPOT LOTION BOWL with lip |  | 1390       | 185        | 81 |
| UPL096  | GALLIPOT LOTION BOWL with lip |  | 1940       | 208        | 84 |
| UPL097  | GALLIPOT LOTION BOWL with lip |  | 3200       | 258        | 92 |

## GALLIPOT ROUND / LOTION BOWL



### GALLIPOT ROUND / LOTION BOWL

|          |                              |  | Cap. (ml.) | SIZE (mm.) |    |
|----------|------------------------------|--|------------|------------|----|
|          |                              |  |            | Dia Ø      | H  |
| UPL-100  | GALLIPOT Round / Lotion Bowl |  | 30         | 46         | 29 |
| UPL-100A | GALLIPOT Round / Lotion Bowl |  | 50         | 56         | 32 |
| UPL-100B | GALLIPOT Round / Lotion Bowl |  | 90         | 66         | 36 |
| UPL-100C | GALLIPOT Round / Lotion Bowl |  | 120        | 76         | 41 |
| UPL-101  | GALLIPOT Round / Lotion Bowl |  | 170        | 88         | 46 |
| UPL-101A | GALLIPOT Round / Lotion Bowl |  | 250        | 95         | 51 |
| UPL-102  | GALLIPOT Round / Lotion Bowl |  | 325        | 107        | 56 |
| UPL-103  | GALLIPOT Round / Lotion Bowl |  | 450        | 128        | 60 |
| UPL-103A | GALLIPOT Round / Lotion Bowl |  | 600        | 139        | 65 |
| UPL-104  | GALLIPOT Round / Lotion Bowl |  | 940        | 158        | 75 |
| UPL-105  | GALLIPOT Round / Lotion Bowl |  | 1390       | 185        | 81 |
| UPL-106  | GALLIPOT Round / Lotion Bowl |  | 1940       | 208        | 84 |
| UPL-107  | GALLIPOT Round / Lotion Bowl |  | 3200       | 258        | 92 |

## IODINE CUP



| IODINE CUP |            | Cap. (ml.) | SIZE (mm.) |    |
|------------|------------|------------|------------|----|
|            |            |            | Dia. Ø     | H  |
| UPL-111    | IODINE CUP | 170        | 84         | 50 |
| UPL-111A   | IODINE CUP | 60         | 60         | 30 |
| UPL-111B   | IODINE CUP | 400        | 120        | 64 |

## MEDICINE CUP



| MEDICINE CUP |                               | Cap. (ml.) | SIZE (mm.) |    |
|--------------|-------------------------------|------------|------------|----|
|              |                               |            | Dia. Ø     | H  |
| UPL-134      | MEDICINE CUP, Graduated, 2 oz | 60         | 53         | 53 |

## MIXING BOWL-TAPER



| MIXING BOWL-TAPER |                   | Cap. (ml.) | SIZE (mm.) |    |
|-------------------|-------------------|------------|------------|----|
|                   |                   |            | Dia. Ø     | H  |
| UPL-108           | MIXING BOWL-TAPER | 100        | 80         | 33 |
| UPL-109           | MIXING BOWL-TAPER | 150        | 90         | 38 |
| UPL-109A          | MIXING BOWL-TAPER | 200        | 95         | 53 |
| UPL-109B          | MIXING BOWL-TAPER | 300        | 105        | 53 |
| UPL-110           | MIXING BOWL-TAPER | 510        | 130        | 71 |

## WASH BASIN



| WASH BASIN     |                 | Cap. (lt.) | SIZE (mm.) |     |
|----------------|-----------------|------------|------------|-----|
|                |                 |            | Dia. Ø     | H   |
| <b>SHALLOW</b> |                 |            |            |     |
| UPL-121        | WASH BASIN      | 3.0        | 310        | 72  |
| UPL-122        | WASH BASIN      | 4.3        | 355        | 79  |
| UPL-123        | WASH BASIN      | 7.0        | 405        | 95  |
| <b>DEEP</b>    |                 |            |            |     |
| UPL-124        | WASH BASIN Deep | 4.0        | 318        | 84  |
| UPL-124A       | WASH BASIN Deep | 4.5        | 318        | 90  |
| UPL-124B       | WASH BASIN Deep | 5.0        | 318        | 100 |
| UPL-125        | WASH BASIN Deep | 4.5        | 350        | 88  |
| UPL-125A       | WASH BASIN Deep | 5.0        | 350        | 98  |
| UPL-125B       | WASH BASIN Deep | 6.0        | 350        | 115 |
| UPL-125C       | WASH BASIN Deep | 6.6        | 350        | 125 |

## BUCKET (PAIL)



| BUCKET (PAIL) |                        | Cap. (lt.) | SIZE (mm.) |     |
|---------------|------------------------|------------|------------|-----|
|               |                        |            | Dia. Ø     | H   |
| UPL-116       | BUCKET (PAIL) with lid | 8.0        | 264        | 235 |
| UPL-116A      | BUCKET (PAIL) with lid | 12.0       | 280        | 250 |
| UPL-117       | BUCKET (PAIL) with lid | 15.0       | 315        | 275 |

## IRRIGATOR (DOUCH CAN)



### IRRIGATOR (DOUCH CAN)

|         |                       | Cap. (lit.) | SIZE (mm.) |     |
|---------|-----------------------|-------------|------------|-----|
|         |                       |             | Dia Ø      | H   |
| UPL-146 | IRRIGATOR (DOUCH CAN) | 1.0         | 108        | 141 |
| UPL-147 | IRRIGATOR (DOUCH CAN) | 1.5         | 127        | 137 |
| UPL-148 | IRRIGATOR (DOUCH CAN) | 2.0         | 127        | 169 |

## BEAKER



### BEAKER

|                     |        | Cap. (ml.) | SIZE (mm.) |     |
|---------------------|--------|------------|------------|-----|
| with/without handle |        |            | Dia Ø      | H   |
| UPL-151             | BEAKER | 125        | 54         | 67  |
| UPL-152             | BEAKER | 250        | 70         | 67  |
| UPL-153             | BEAKER | 600        | 86         | 117 |
| UPL-154             | BEAKER | 1000       | 102        | 68  |
| UPL-155             | BEAKER | 1200       | 109        | 140 |
| UPL-156             | BEAKER | 2000       | 128        | 165 |
| UPL-157             | BEAKER | 3000       | 148        | 182 |

## BEAKER-TAPER



### BEAKER-TAPER

|                     |              | Cap. (ml.) | SIZE (mm.) |     |
|---------------------|--------------|------------|------------|-----|
| with/without handle |              |            | Dia Ø      | H   |
| UPL-163             | BEAKER-TAPER | 500        | 94         | 116 |
| UPL-164             | BEAKER-TAPER | 1000       | 120        | 133 |
| UPL-166             | BEAKER-TAPER | 2000       | 140        | 170 |

## TUMBLER



### TUMBLER

|         |               | Cap. (ml.) | SIZE (mm.) |     |
|---------|---------------|------------|------------|-----|
|         |               |            | Dia Ø      | H   |
| UPL-235 | TUMBLER 5 oz  | 150        | 65         | 75  |
| UPL-236 | TUMBLER 7 oz  | 210        | 75         | 78  |
| UPL-237 | TUMBLER 10 oz | 300        | 80         | 98  |
| UPL-238 | TUMBLER 12 oz | 360        | 83         | 110 |

## SPITTEN MUG



### SPITTEN MUG

|         |                      | Cap. (ml.) | SIZE (mm.) |    |
|---------|----------------------|------------|------------|----|
|         |                      |            | Dia Ø      | H  |
| UPL-112 | SPITTEN MUG with lid | 430        | 96         | 65 |



## BED PANS

### BED PANS

| <b>ROUND</b> |                          | SIZE (mm.) |    |
|--------------|--------------------------|------------|----|
|              |                          | Dia. Ø     | H  |
| UPL-141      | BED PAN ROUND with cover | 320        | 85 |



| <b>PERFECTION</b> |                    | SIZE (mm.) |     |     |
|-------------------|--------------------|------------|-----|-----|
|                   |                    | L          | W   | H   |
| UPL-142           | BED PAN PERFECTION | 354        | 289 | 100 |



| <b>FRACTURE</b> |                                | SIZE (mm.) |     |    |
|-----------------|--------------------------------|------------|-----|----|
|                 |                                | L          | W   | H  |
| UPL-143         | BED PAN FRACTURE               | 418        | 292 | 85 |
| UPL-143T        | BED PAN FRACTURE (tube handle) | 418        | 292 | 85 |



## URINALS

### URINALS

| <b>MALE HORIZONTAL</b> |                        | Cap. (lt.) | SIZE (mm.) |     |
|------------------------|------------------------|------------|------------|-----|
|                        |                        |            | Dia. Ø     | H   |
| UPL-144                | URINAL MALE HORIZONTAL | 1.0        | 110        | 229 |



| <b>MALE VERTICAL</b> |                      | Cap. (lt.) | SIZE (mm.) |     |
|----------------------|----------------------|------------|------------|-----|
|                      |                      |            | Dia. Ø     | H   |
| UPL-144A             | URINAL MALE VERTICAL | 1.5        | 114        | 220 |



| <b>FEMALE</b> |               | Cap. (lt.) | SIZE (mm.) |     |     |
|---------------|---------------|------------|------------|-----|-----|
|               |               |            | L          | W   | H   |
| UPL-145       | URINAL FEMALE | 0.6        | 225        | 120 | 140 |

## DRESSING JAR



### DRESSING JAR

|          |                       | Cap. (lit.) | SIZE (mm.) |     |
|----------|-----------------------|-------------|------------|-----|
|          |                       |             | Dia Ø      | H   |
| UPL-130  | DRESSING JAR with lid | 0.5         | 102        | 65  |
| UPL-131  | DRESSING JAR with lid | 1.0         | 103        | 126 |
| UPL-132  | DRESSING JAR with lid | 2.0         | 127        | 161 |
| UPL-133  | DRESSING JAR with lid | 3.0         | 149        | 180 |
| UPL-133A | DRESSING JAR with lid | 4.0         | 172        | 196 |



## FORCEPS JAR

### FORCEPS JAR

|         |             | SIZE (mm.) |     |
|---------|-------------|------------|-----|
|         |             | Dia Ø      | H   |
| UPL-135 | FORCEPS JAR | 55         | 110 |
| UPL-136 | FORCEPS JAR | 55         | 140 |
| UPL-137 | FORCEPS JAR | 55         | 170 |
| UPL-138 | FORCEPS JAR | 55         | 180 |



## THERMOMETER JAR

### THERMOMETER JAR

|         |                 | SIZE (mm.) |    |
|---------|-----------------|------------|----|
|         |                 | Dia Ø      | H  |
| UPL-139 | THERMOMETER JAR | 33         | 80 |



## OINTMENT JAR

### OINTMENT JAR

|         |                         | Cap. (ml.) | SIZE (mm.) |    |
|---------|-------------------------|------------|------------|----|
|         |                         |            | Dia Ø      | H  |
| UPL-113 | OINTMENT JAR with cover | 250        | 76         | 57 |
| UPL-114 | OINTMENT JAR with cover | 500        | 102        | 65 |



## STERILIZING/DRESSING DRUM

### STERILIZING/DRESSING DRUM



|          |               | SIZE (mm.) |     |
|----------|---------------|------------|-----|
|          |               | Dist. Ø    | H   |
| UPL-301  | DRESSING DRUM | 100        | 100 |
| UPL-301A | DRESSING DRUM | 100        | 120 |
| UPL-302  | DRESSING DRUM | 125        | 80  |
| UPL-303  | DRESSING DRUM | 150        | 100 |
| UPL-303A | DRESSING DRUM | 150        | 120 |
| UPL-303B | DRESSING DRUM | 150        | 130 |
| UPL-303C | DRESSING DRUM | 150        | 150 |
| UPL-304  | DRESSING DRUM | 180        | 180 |
| UPL-305  | DRESSING DRUM | 190        | 145 |
| UPL-305A | DRESSING DRUM | 190        | 190 |
| UPL-306  | DRESSING DRUM | 200        | 140 |
| UPL-307  | DRESSING DRUM | 225        | 115 |
| UPL-307A | DRESSING DRUM | 225        | 225 |
| UPL-308  | DRESSING DRUM | 228        | 152 |
| UPL-309  | DRESSING DRUM | 240        | 145 |
| UPL-309A | DRESSING DRUM | 240        | 190 |
| UPL-309B | DRESSING DRUM | 240        | 240 |
| UPL-309C | DRESSING DRUM | 240        | 290 |
| UPL-310  | DRESSING DRUM | 280        | 240 |
| UPL-311  | DRESSING DRUM | 290        | 145 |
| UPL-311A | DRESSING DRUM | 290        | 190 |
| UPL-311B | DRESSING DRUM | 290        | 240 |
| UPL-311C | DRESSING DRUM | 290        | 290 |
| UPL-312  | DRESSING DRUM | 304        | 254 |
| UPL-313  | DRESSING DRUM | 340        | 240 |
| UPL-313A | DRESSING DRUM | 340        | 340 |
| UPL-314  | DRESSING DRUM | 354        | 254 |
| UPL-315  | DRESSING DRUM | 380        | 304 |
| UPL-316  | DRESSING DRUM | 390        | 290 |
| UPL-316A | DRESSING DRUM | 390        | 340 |

## FASTENING SYSTEMS



CHAIN TYPE



ITALIAN TYPE



BUTTON TYPE



CLIP LOCK



## SURGICAL BOXES Stainless Steel



### SURGICAL BOXES Stainless Steel

|          |                             | L   | SIZE (mm) |     |
|----------|-----------------------------|-----|-----------|-----|
|          |                             |     | W         | H   |
| UPL-211  | SURGICAL BOX 5.5 with cover | 180 | 80        | 30  |
| UPL-212  | SURGICAL BOX 5.5 with cover | 180 | 80        | 40  |
| UPL-213  | SURGICAL BOX 5.5 with cover | 180 | 80        | 50  |
| UPL-214  | SURGICAL BOX 5.5 with cover | 200 | 90        | 30  |
| UPL-215  | SURGICAL BOX 5.5 with cover | 200 | 90        | 40  |
| UPL-216  | SURGICAL BOX 5.5 with cover | 200 | 90        | 50  |
| UPL-217  | SURGICAL BOX 5.5 with cover | 250 | 100       | 40  |
| UPL-218  | SURGICAL BOX 5.5 with cover | 250 | 100       | 50  |
| UPL-219  | SURGICAL BOX 5.5 with cover | 250 | 100       | 60  |
| UPL-220  | SURGICAL BOX 5.5 with cover | 300 | 125       | 50  |
| UPL-220A | SURGICAL BOX 5.5 with cover | 300 | 125       | 60  |
| UPL-220B | SURGICAL BOX 5.5 with cover | 300 | 125       | 70  |
| UPL-220C | SURGICAL BOX 5.5 with cover | 300 | 125       | 80  |
| UPL-220D | SURGICAL BOX 5.5 with cover | 300 | 150       | 80  |
| UPL-220E | SURGICAL BOX 5.5 with cover | 300 | 150       | 90  |
| UPL-220F | SURGICAL BOX 5.5 with cover | 300 | 150       | 100 |
| UPL-220G | SURGICAL BOX 5.5 with cover | 360 | 160       | 80  |
| UPL-220H | SURGICAL BOX 5.5 with cover | 360 | 160       | 90  |
| UPL-220I | SURGICAL BOX 5.5 with cover | 360 | 160       | 100 |
| UPL-220J | SURGICAL BOX 5.5 with cover | 400 | 160       | 80  |
| UPL-220K | SURGICAL BOX 5.5 with cover | 400 | 160       | 90  |
| UPL-220L | SURGICAL BOX 5.5 with cover | 400 | 160       | 100 |
| UPL-220M | SURGICAL BOX 5.5 with cover | 500 | 200       | 80  |
| UPL-220N | SURGICAL BOX 5.5 with cover | 500 | 200       | 100 |
| UPL-220O | SURGICAL BOX 5.5 with cover | 500 | 200       | 120 |

## SURGICAL BOXES Aluminium



### SURGICAL BOXES Aluminium

|          |                             | L   | SIZE (mm) |    |
|----------|-----------------------------|-----|-----------|----|
|          |                             |     | W         | H  |
| UPL-201  | 5. BOX ALUMINIUM with cover | 175 | 76        | 20 |
| UPL-202  | 5. BOX ALUMINIUM with cover | 175 | 76        | 30 |
| UPL-202A | 5. BOX ALUMINIUM with cover | 175 | 76        | 40 |
| UPL-203  | 5. BOX ALUMINIUM with cover | 185 | 95        | 20 |
| UPL-204  | 5. BOX ALUMINIUM with cover | 185 | 95        | 30 |
| UPL-205  | 5. BOX ALUMINIUM with cover | 185 | 95        | 40 |
| UPL-205A | 5. BOX ALUMINIUM with cover | 185 | 95        | 50 |
| UPL-206  | 5. BOX ALUMINIUM with cover | 218 | 106       | 30 |
| UPL-207  | 5. BOX ALUMINIUM with cover | 218 | 106       | 40 |
| UPL-208  | 5. BOX ALUMINIUM with cover | 218 | 106       | 50 |

Also available in colours :

\* Height of 5.5 & Aluminium Boxes includes cover.



BLUE



RED



GOLDEN YELLOW



## SCOOP

|         |       | SIZE (mm) |     |
|---------|-------|-----------|-----|
|         |       | Di Ø      | H   |
| UPL-115 | SCOOP | 130       | 200 |

\* Dimensions excluding handle.



## PETRI DISH

|         |            | SIZE (mm) |    |
|---------|------------|-----------|----|
|         |            | Di Ø      | H  |
| UPL-231 | PETRI DISH | 55        | 22 |
| UPL-232 | PETRI DISH | 71        | 22 |



## JAR COVERS

|          |           | SIZE (mm) |    |
|----------|-----------|-----------|----|
|          |           | Di Ø      | H  |
| UPL-221  | JAR COVER | 73        | 23 |
| UPL-221A | JAR COVER | 76        | 23 |
| UPL-222  | JAR COVER | 95        | 26 |
| UPL-222A | JAR COVER | 90        | 26 |
| UPL-222B | JAR COVER | 97        | 26 |
| UPL-223  | JAR COVER | 106       | 22 |
| UPL-224  | JAR COVER | 119       | 24 |
| UPL-225  | JAR COVER | 143       | 28 |



## INSTRUMENT STERILIZER

|         |                       | SIZE (mm) |     |     |
|---------|-----------------------|-----------|-----|-----|
|         |                       | L         | W   | H   |
| UPL-171 | INSTRUMENT STERILIZER | 300       | 150 | 125 |
| UPL-172 | INSTRUMENT STERILIZER | 430       | 200 | 150 |
| UPL-173 | INSTRUMENT STERILIZER | 510       | 200 | 150 |

### ALL PRODUCTS

- Made from Stainless Steel, AISI-304 Grade (18-8).
- Available in both Bright & Matt finish.
- Made in Seamless construction.

ISO 9001 Certified

CE Marked

FDA Registered

**United Poly Engineering Pvt. Ltd.** was incorporated in the year 1980 with the intention of providing clinics and hospitals with high quality Stainless Steel Hollowares. This company is a part of the Grover Group of companies, having a successful and growing presence in the surgical and medical field for the past 45 years.



**UPL's** strength lies in its most updated production facilities including in-house Development & Tooling centre, backed by trained & committed workforce, with the sole intention of providing products of highest quality to meet the exacting requirements of the medical professional. Ensuing its goal of manufacturing high quality product, which are made from Stainless Steel AISI-304 grade (18/8), it has put in place a comprehensive Quality Management System and is an ISO 9001 certified, FDA registered company and its products are CE marked.

The Company with its Customer Oriented Focus has been successful in establishing its presence both in the developing and the developed markets and its products are widely accepted & available in every continent. The performance of the company can be gauged by its achievements, reflected by it receiving the state award for outstanding export performance.

The ideology of the company is to foster growth by earning trust.

[www.unitedpoly.com](http://www.unitedpoly.com)



**UNITED POLY ENGINEERING PVT. LTD.**

D-13/3, Okhla Industrial Area, Phase-II,  
New Delhi - 110 020 (INDIA)

Tel. : +91-11- 26388564, 65, 66

Fax : +91-11- 26389947, 26384454

E-mail : [upl@unitedpoly.com](mailto:upl@unitedpoly.com)

Complete Range made from 18 : 8 non magnetic stainless steel  
All the sizes are approximate (measured from outside without cover)  
We reserve the right to Amend Specification without prior notice.