

Infinity™ Knee System

Modular ACL & PCL Reconstruction Platform

Imagine performing AM or all-inside techniques, with the flexibility of straight, curved, or retro-reaming, and achieving suspensory or aperture fixation with one system. Imagine an adjustable loop button that allows for graft re-positioning after initial femoral fixation. Developed with an international team of world-renowned knee surgeons, this system has been designed to do just that.

Introducing the Infinity™ Knee System – one system with infinite possibilities.



Infinity™ ACL Tibial Footprint Guide with All-Inside Guide Sleeve

Features and Benefits

The World's First Adjustable, Reversible Button

After initial tensioning of the graft in the femoral tunnel, the Infinity™ Femoral Adjustable Loop Button allows surgeons to re-lengthen the graft loop and re-position the graft within the femoral and tibial tunnels.

One Modular System

The versatile Infinity™ Knee System supports:

- All-inside and anteromedial techniques.
- Straight, curved, and retro-reaming preferences.
- Fixation for both soft tissue and soft tissue and bone grafts.

Infinity™ Femoral Adjustable Loop Button

Designed by Innovators


We've partnered with a team of the world's most renowned knee surgeons, to develop products that overcome the challenges surgeons face in the OR every day. The result is a complete system designed to provide versatility, procedural efficiency, and an easier day in the OR.



Infinity™ Anteromedial Guide with Spade Tip Guide Pin

Infinity™ Retro-Reamer

To simplify your next ACL /PCL reconstruction case, call **1-866-4CONMED** to schedule a demonstration or visit CONMED.com/Infinity.



Infinity™ – One system with infinite possibilities.

Ordering Information

	Description	Dimensions	Use	Catalog Number
NEW	Infinity™ Implants			
	Infinity™ Femoral Adjustable Loop Button	–	Disposable	KFB035
	Infinity™ Adjustable Loop Button Cradle	–	Disposable	KFB135
	Infinity™ Standard Tibial Button	–	Disposable	KTB014
	Infinity™ Large Tibial Button	–	Disposable	KTB017
	Infinity™ Adjustable Free Loop	–	Disposable	KFL100
	Infinity™ Guides			
	Infinity™ ACL/PCL Femoral Footprint Guide	–	Reusable	KFA100
	Infinity™ ACL Tibial Footprint Guide	–	Reusable	KTA100
	Infinity™ ACL Tibial Elbow Guide	–	Reusable	KTE100
	Infinity™ ACL Tibial Tip Guide	–	Reusable	KTT100
	Infinity™ PCL Tibial Guide	–	Reusable	KTP100
	Infinity™ Guide Sleeve Straight	2.4mm	Reusable	KTS124
	Infinity™ Guide Sleeve Angled	2.4mm	Reusable	KTS224
	Infinity™ All-Inside Guide Sleeve	3.5mm	Reusable	KTS135
	Infinity™ Guide Body	–	Reusable	KGB100
	Infinity™ Anteromedial Guide Left	7/8mm	Reusable	KBL178
	Infinity™ Anteromedial Guide Left	9/10mm	Reusable	KBL191
	Infinity™ Anteromedial Guide Right	7/8mm	Reusable	KBR178
	Infinity™ Anteromedial Guide Right	9/10mm	Reusable	KBR191
	Infinity™ Accessories			
	Infinity™ Spade Tip Guide Pin	3.5mm	Disposable	KGP035
	Infinity™ Tunnel Dilator	3.5mm	Reusable	KFD035
	Infinity™ Suture Shuttle	–	Disposable	KSP100
	KPS Infinity™ Insert Tray	–	Reusable	INST-ACL04
	Infinity™ Reamers			
	Infinity™ Retro-Reamer	6.0mm	Disposable	KRR060
	Infinity™ Retro-Reamer	6.5mm	Disposable	KRR065
	Infinity™ Retro-Reamer	7.0mm	Disposable	KRR070
	Infinity™ Retro-Reamer	7.5mm	Disposable	KRR075
	Infinity™ Retro-Reamer	8.0mm	Disposable	KRR080
	Infinity™ Retro-Reamer	8.5mm	Disposable	KRR085
	Infinity™ Retro-Reamer	9.0mm	Disposable	KRR090
	Infinity™ Retro-Reamer	9.5mm	Disposable	KRR095
	Infinity™ Retro-Reamer	10.0mm	Disposable	KRR100
	Infinity™ Retro-Reamer	11.0mm	Disposable	KRR110
	Infinity™ Retro-Reamer	12.0mm	Disposable	KRR120