

Tastan-Cakir Sägen / Tastan-Cakir Saws

Zur Durchtrennung von Knochen bei Radix- (konkave Säge)
und transversalen (konvexe Säge) Osteotomien

Entwicklungspartner:

Dr. Eren Taştan und Dr. Barış Çakır, Istanbul

For cutting through bone during radix (concave saw) and
transversal (convex saw) osteotomies

Development partner:


Dr. Eren Taştan and Dr. Barış Çakır, Istanbul



konvex / convex



konkav / concave



||| SURGICAL

medicon[®]

Medicon eG · Gänsäcker 15 · D-78532 Tuttlingen
Tel.: +49 (0) 74 62 / 20 09-0 · Fax: +49 (0) 74 62 / 20 09-50
E-Mail: sales@medicon.de · Internet: www.medicon.de · Germany

Tastan-Cakir Sägen / Tastan-Cakir Saws

Entwicklungspartner / development partner:
Dr. Eren Taştan & Dr. Barış Çakır, Istanbul



67.60.15
Tastan-Cakir Säge konvex /
Tastan-Cakir Saw convex



67.60.16
Tastan-Cakir Säge konkav /
Tastan-Cakir Saw concave

Vorteile:

- » **Speziell angepasst auf die Anatomie** der knöchernen Nasenstrukturen
- » Auf Grund der manuellen Anwendung ist **keine Kühlung notwendig**, was bei elektrisch betriebenen Sägen hingegen unumgänglich ist
- » Präzise Osteotomie, wodurch - im Gegensatz zur Osteotomie mit Meißel - **ungewollte Risse** im Knochen **vermieden** werden

Advantages:

- » **Especially adapted to the anatomy** of the bony structures of the nose
- » Due to the manual application, **no cooling is required**, which is unavoidable when working with electrically operated saws
- » Precise osteotomy, which - in contrast to an osteotomy with a chisel - **avoids unwanted cracks** of the bone

Weitere interessante Produkte für die plastische Chirurgie finden Sie in unserem Prospekt **Moderne Rhinoplastik 451.02.04** und **Advanced Rhinoplasty 451.04.76** oder unter www.medicon.de

Other interesting products for plastic surgery can be found in our brochure **Modern Rhinoplasty 451.02.04** and **Advanced Rhinoplasty 451.04.76** or on www.medicon.de

