


© Copyright 05/2017  Not all products illustrated in this documentation are available for sale in all countries. Please contact your local distributor for more information.



Medicon cannulas for autologous fat transfer

Developed in cooperation with
Dr. med. Henning von Gregory Clinic for Aesthetic-Plastic Surgery and Reconstructive Nose- and Maxillofacial Surgery, Berlin

- Atraumatic treatment of the fat in all steps of the procedure
 - The treatment of scars and the correction of small wrinkles is possible in all areas of the body, fat preparation and fat injection that are possible in the area of the face and the hands.
- Cannula openings:
Smooth on injection cannulas so that fat cells are not destroyed.
 - Atraumatic but anti-friction cannula tips thanks to conical design
- Connectors – handy – ergonomic - direction markings
 - During the injection or aspiration of the fat the cannulas are blocked less often compared to disposable cannulas with a lower product quality.
- The cannulas not only exhibit a very high material quality but also perfect workmanship. They are exclusively produced in highly-specialised factories in Germany.
 - Sterilisable in compliance with DIN 17664 – smooth inner surface which allows effective cleaning.



medicon[®]

Medicon eG · Gänsäcker 15 · D-78532 Tuttlingen
Tel.: +49 (0) 74 62 / 20 09-0 · Fax: +49 (0) 74 62 / 20 09-50
E-Mail: sales@medicon.de · Internet: www.medicon.de · Germany

Medicon cannulas for autologous fat transfer
 Developed in cooperation with
 Dr. med. Henning von Gregory Clinic for Aesthetic-Plastic Surgery
 and Reconstructive Nose- and Maxillofacial Surgery, Berlin

Liposuction Cannula, for lipoinjection in the facial area




a = 40 mm, Ø 0,7 mm	09.89.15	Lipo injection cannula
a = 40 mm, Ø 0,9 mm	09.89.16	Lipo injection cannula
a = 50 mm, Ø 1,2 mm	09.89.17	Lipo injection cannula
a = 60 mm, Ø 1,5 mm	09.89.18	Lipo injection cannula

Liposuction Cannula, to separate scarred connective tissue



sharp



a = 7 cm, Ø 1,4 mm	09.89.22	Toledo lipo injection cannula
--------------------	-----------------	-------------------------------

Liposuction Cannula, for autologous fat transfer




a = 15 cm, Ø 2,5 mm, 4 holes	09.89.24	von Gregory liposuction cannula
------------------------------	-----------------	---------------------------------

Liposuction Cannula, for autologous fat transfer as well as sweat gland removal




sharp



a = 18 cm, Ø 3 mm, 1,3 mm hole diameter	09.89.26	von Gregory liposuction cannula sharp
---	-----------------	---------------------------------------

Infiltration Cannula



a = 15 cm, Ø 2 mm, 16 holes	09.89.28	von Gregory injection cannula
-----------------------------	-----------------	-------------------------------

For further information please see the
 MEDICON catalogue "Plastic Surgery" Art.
 Nr. 450.00.60.

