



# **GUBiSCH**

---

## **Essentials**

**Essential Instrument Set**

---

**medicon<sup>®</sup>**

---

**Set Overview**



**Which instruments do I really need, to have to perform a successful nose surgery?**

For those of you who have asked themselves this question before, I have compiled an "Essential Instrument Set" in cooperation with the company Medicon, that enables ambitious newcomers and experienced surgeons alike, to achieve excellent results. In addition to these essential instruments, you would only need a few standard instruments to be found in every basic set.

Prof. Dr. med. Dr. h. c.  
Wolfgang Gubisch  
Senior Director  
Klinik für Plastische  
Gesichtschirurgie  
Marienhospital Stuttgart

Internet:  
[www.professor-gubisch.de](http://www.professor-gubisch.de)




This clearly structured and beautifully illustrated atlas, which includes numerous high-quality pre-, intra-, and postoperative pictures, documents the techniques employed in order to deal with the problems typically encountered following a previous failed rhinoplasty.

„Mastering Advanced Rhinoplasty“  
Wolfgang Gubisch  
Springer Verlag, Berlin



Item number	Description	Quantity
<b>06.23.12</b>	MICRO-ADSON Tissue forceps · 1 × 2 teeth · 12 cm	2
<b>06.25.12</b>	ADSON-BROWN Tissue forceps · 7 × 7 teeth · 12 cm	1
<b>20.52.10</b>	GUTHRIE Hook · 2 teeth · sharp · 1.5 mm · 13 cm	1
<b>27.10.32</b>	Mallet · made of Novotex plastic · Ø 45 mm · 240 g · 26 cm	1
<b>44.90.70</b>	GUBISCH Pen holder · 13 cm	1
<b>44.90.72</b>	Feather · no 46 · for marking · packing unit 10 pieces	1
<b>60.00.40</b>	CASTROVIEJO Caliper · measuring range 20 mm · 8 cm	1
<b>60.05.12</b>	DESMARRES Lid retractor · figure 2 · 14 mm · 14 cm	1
<b>63.56.20</b>	Cutting burr · barrel style · Ø 5 mm · shank Ø 2.35 mm · 7 cm	1
<b>66.08.59</b>	KILLIAN Nasal speculum · figure 4 · 90 mm · dismantable · 14 cm	1
<b>67.03.52</b>	COTTLE Nasal speculum · 50 mm · 13 cm	1
<b>67.15.98</b>	HARALDSSON Suction elevator · Ø 3 mm · tube connection · 18.5 cm	1
<b>67.16.27</b>	DANIEL-CAKIR Subperichondrial/subperiosteal elevator · 20 cm	1

Set Overview



Item number	Description	Quantity
<b>67.19.76</b>	Cartilage cutting block · made of PPSU · black · 85 × 55 × 10 mm	1
<b>67.19.79</b>	Cartilage cutting block with grooves · made of aluminium · 85 × 55 × 10 mm	1
<b>67.20.07</b>	GUBISCH-KILNER Alar retractor · 10 mm · 8 cm	1
<b>67.21.32</b>	GUBISCH-AUFRICHT Nasal retractor · 30 × 6/10 mm · 17 cm	1
<b>67.30.15</b>	COTTLE-GUBISCH Chisel · graduated · hollow cut · 5 mm · 18 cm	1
<b>67.30.22</b>	COTTLE-GUBISCH Chisel · graduated · hollow cut · 12 mm · 18 cm	1
<b>67.30.53</b>	GUBISCH Osteotome · with aluminium handle · 3 mm · 18.5 cm	1
<b>67.58.42</b>	AIACH-GUBISCH Forceps · 16 cm	1
<b>67.61.01</b>	Handle for nasal rasps · double-ended · 14 cm	1
<b>67.61.03</b>	Wrench · for nasal rasps	1
<b>67.61.11</b>	Rasps insert "HM" · figure 2	1
<b>67.61.14</b>	Rasps insert "HM" · figure 5	1
<b>68.61.05</b>	LINDEMANN Burr · shank-Ø 2.35 mm · Ø 1.4 mm · 44.5 mm	1
<b>68.76.64</b>	Twist drill · STRYKER Shank end · Ø 1.45 mm · 44 mm	1
<b>202.20.61</b>	SC Plus iris scissors straight · 11.5 cm	1
<b>204.55.66</b>	KELLY Dissecting scissors · supercut plus · curved · 16 cm	1

# medicon<sup>®</sup> courses

Our Medicon training courses on plastic surgery provide a platform for an interchange of scientific ideas combined with practical application scenarios of selected surgical cases.

## More Information:

**+49 (0) 74 62 / 20 09-0**  
**www.medicon-courses.de**



## Medicon eG

Gänsäcker 15  
D-78532 Tuttlingen

P. O. Box 44 55  
D-78509 Tuttlingen

Tel.: +49 (0) 74 62 / 20 09-0  
Fax: +49 (0) 74 62 / 20 09-50

E-Mail: sales@medicon.de  
Internet: www.medicon.de

Germany

Not all products illustrated in this documentation are available for sale in all countries.  
Please contact your local distributor for more information.



Modelländerungen vorbehalten  
Patterns are subject to change  
Salvo modificaciones  
Tous droits réservés des changements de modèle  
Ci riserviamo la facoltà di cambiamenti nei modelli

Gedruckt in Deutschland  
Printed in Germany  
Impresso en Alemania  
Imprimé en Allemagne  
Stampato in R.F.G.  
Trykket i Tyskland

© Copyright 04/2017

451 .02.08

