



## Product Description

### Bowie-Dick Test Pack/Paper

M&GPack Bowie-Dick Pack/paper are pre-assembled, single use test packs designed to detection of defects in the sterilization process.

An indicator sheet printed with water based, non-toxic chemical indicator.

The color change from blue to black of indicator after steam sterilization indicates that all air has been removed and steam penetration did take place



### Advantages of M&GPack Bowie-Dick Test Pack/Paper:

Easy-to-read test sheet with detailed instructions for use and interpretations of results

Small size for easy handling and storage

External process indicator for ready identification of processed test packs

Totally reproducible and the test sheet will stay flat after the sterilization cycle

Standards:

M&GPack Bowie-Dick Test Pack/paper comply with ISO11140-1, EN867-4 and ISO11140-4.

### Specification Sheet

Single use test packs designed to detection of defects in the sterilization process

DESCRIPTION	ITEM NO	SIZE
Bowie-Dick Test Packs/Paper	MGBD - 01	115mm*125mm
bowie-Dick Test Paper	MGBD - 02	210mm*300mm

## Sterilization Color Change

Turn to Black from Blue under steam sterilization:



**Before Sterilization**

**After Sterilization**

Material:	Paper
Pack size:	20 test packs per box
Type:	Surgical Supplies
Application:	Hospital, Dental Clinic, etc.
Sterilization:	Steam Bowie-Dick Test
Color Change:	Blue to Black
Latex Free:	Yes
Lead Free Ink & No Toxic:	Yes