



F2 Monobloc



650/690 MM ALSO
AVAILABLE IN SILVER



5-PART MATTRESS

The F2 Monobloc provides a **safe, ergonomic, and reliable patient handling solution**. The independent leg design supports the patient's weight while loading and unloading, reducing strain on medics.

- **Easy to use controls** and **independent leg operation** provide fast, safe loading and unloading from the ambulance
- **Color-coded left-hand and right-hand control levers** are user-friendly and require minimal training
- **Large 200 mm wheels** easily roll over multiple terrain types
- **All Aluminum construction** is lightweight and durable



Features

- **Includes telescoping IV pole,** sidearms and foot support
- **Independent leg** operation
- **Legs unfold automatically** when unloading
- **Easy to use color-coded controls**
- **Two loading heights:** 560 mm and 600 mm
- **Large 200 mm wheels**
- **Control-end swivel wheel locks**
- **Loading-end fixed wheels** and wheel locks
- **Shock frame** with contour feature
- **Four three-position** telescoping carrying handles
- **Adjustable seven-position** backrest
- **Fold-down sidearms**
- **Low maintenance** design
- **Aluminum** construction

Specifications

	MAX	MIN
LENGTH	89 in 2272 mm	76 in 1921 mm
WIDTH	22 in 570 mm	
HEIGHT	32 in 820 mm	11 in 270 mm
LOAD HEIGHT	24 in 600 mm	22 in 560 mm
WEIGHT	95 lbs 43 kg	
LOAD LIMIT	550 lbs 250 kg	

Certifications



EN 1865

EN 1789

Accessories



- **Instrument Tray**
No. 0832158



- **F1 Fastener**
No. 0573081



- **5-Part Mattress**
No. 031-4193