



INSTRAXX[®]

System Components Guide & Specifications







The iNTRAXX Integrated Vehicle Component System **provides crash-ready safety** for providers and patients and **efficiency that reduces operating costs** for agencies.

The iNTRAXX system is designed to secure loose equipment inside the ambulance and allow medics to deliver patient care from a seated, restrained position. This **reduces the risk of death and injury** during crashes, sudden stops, or hard maneuvers.

This modular system features wall-mounted tracks that secure a series of SafeMount Equipment Mounts, SafeCab Storage Cabinets, and SafePak Supply Bags. iNTRAXX modularity supports better workflows and inventory management which **lowers operating costs**.





Delivery of care is now **safer, more convenient, and more efficient** with **iNTRAXX**

CRASH TESTED SAFETY

- **Crash-ready** Dynamically crash-tested to the SAE J3043 requirement that keeps equipment secure during 26g impacts.
- **Seated & Restrained Delivery of Care** Reach supplies and deliver care without standing in a moving vehicle.
- **Soft-Sided Storage** SafePaks are soft-sided, eliminating potential strike hazards during hard maneuvers or crashes.

REDUCED PROCESSING TIME & MAINTENANCE

- **Faster Turn Time** SafePaks and pouches are easy to remove and clean, making restocking fast and easy.
- **Reduced Downtime** Damaged or out-of-service components can be quickly removed and replaced without taking a vehicle out of service.

REDUCED SYSTEM WASTE

- **Improved Inventory Control** SafePaks and pouches make inventory control easier with better visibility of supplies, reducing duplication, waste, and expired supplies.

IMPROVED ERGONOMICS

- **Easy Access** Modularity allows supplies and equipment to be placed next to the provider in the delivery of care position, eliminating the need to get up in a moving vehicle to access supplies.
- **Provider Specific Configuration** Quickly reposition equipment and supplies exactly where you want them.

IMPROVED FLEXIBILITY

- **Interoperability** Supplies and equipment are interchangeable across vehicles and vehicle types.
- **Mission-Specific Customization** Equipment and supplies can be removed or repositioned to adapt a vehicle to different response scenarios (e.g., viral outbreak or mass casualty).
- **Expandable System** Components can be added as your needs change.

IMPROVED INFECTION CONTROL

- **Cleaning and Decontamination** Components are removable and power washable for thorough, easy cleaning.
- **Containment** Individual pouches isolate potential contamination to specific areas, limiting the exposure of expensive supplies.

REDUCE COST OF OWNERSHIP

- **Reduced Space Requirements** Requires less space for storage and equipment, allowing use of smaller, less expensive vehicles.
- **Increased Utilization** Faster turn times so your vehicle is more readily available for revenue-generating activity.

All products come with a **2 Year Warranty**, unless stated otherwise.

CERTIFICATIONS Fully compliant with SAE J3043 testing requirements. Meets KKK-A-1822F, NFPA-1917, and CAAS GVS v1.0 Standards.

Learn more at www.fernoems.com | intraxx



Scalable to Any Size Vehicle

Tracks




Track System | 1016mm



Part Number	0480001
Length	W 7.4 in 189 mm
	D .5 in 12.7 mm
	L 40 in 1016 mm
Weight	7.5 lbs 3.4 kg
Load Capacity	220.4 lbs 100 kg
Crash Tested	  
• Limited Lifetime Warranty	




**Track System | 340mm,
Black End Caps**



Part Number	0480062
Length	W 7.6 in 193 mm
	D .7 in 17 mm
	L 13.4 in 340 mm
Weight	2.6 lbs 1.2 kg
Load Capacity	55.1 lbs 25 kg
Crash Tested	  
• Limited Lifetime Warranty	

**Track | 1016 No Holes,
(Bonding)**



Part Number	0480073
Length	W 7.4 in 189 mm
	D .5 in 12.7 mm
	L 40 in 1016 mm
Weight	7.5 lbs 3.4 kg
Load Capacity	220.4 lbs 100 kg
Crash Tested	  

KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

Tracks

Ceiling Track | Right



Part Number	0480107
Length	W 7.4 in 189mm D .5 in 12.7 mm L 40 in 1016 mm
Weight	7.5 lbs 3.4 kg
Load Capacity	220.4 lbs 100 kg
Crash Tested	

Ceiling Track | Left



Part Number	0480131
Length	W 7.4 in 189mm D .5 in 12.7 mm L 40 in 1016 mm
Weight	7.5 lbs 3.4 kg
Load Capacity	220.4 lbs 100 kg
Crash Tested	

Ceiling Track | Center



Part Number	0480108
Length	W 7.4 in 189mm D .5 in 12.7 mm L 40 in 1016 mm
Weight	7.5 lbs 3.4 kg
Load Capacity	220.4 lbs 100 kg
Crash Tested	

KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

Caps

End Caps | Regular



Part Number Black: 0480002 | Grey: 0480003 | White: 0480004

Length	W 7.6 in 193 mm
	D 1.9 in 48 mm
	L .7 in 17 mm

Weight .605 lbs | 275 g

Crash Tested			
--------------	--	--	--

End Caps | Thin




Part Number Black: 0480049 | Grey: 0480050 | White: 0480051

Length	W 7.4 in 193 mm
	D 1.5 in 48mm
	L .7 in 17 mm

Weight .605 lbs | 275 g

Crash Tested			
--------------	--	--	--

T-Splice Caps



Part Number Black: 0480005 | Grey: 0480006 | White: 0480007

Length	W 7.6 in 193 mm
	D .9 in 24 mm
	L .6 in 16 mm

Weight .14 lbs | 64 g




KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

Track Connectors




SafeMount | Threaded (25kg)



Part Number	0480010
Length	W 5.9 in 151 mm
	D 1.4 in 34.5 mm
	L 6.6 in 167 mm
Weight	1.3 lbs .59 kg
Load Capacity	55.1 lbs 25 kg
Crash Tested	  




SafeMount | Threaded (12kg)



Part Number	0480146
Length	W 5.9 in 151 mm
	D 1.6 in 40 mm
	L 6.6 in 167 mm
Weight (without connector)	1.6 lbs 0.74 kg
Load Capacity	26.5 lbs 12 kg
Crash Tested	12kg   

SafeMount | With Connector



Part Number	0480011
Length	W 7.4 in 187 mm
	D 1.9 in 48.5 mm
	L 6.6 in 167 mm
Weight (without connector)	1.65 lbs .75 kg
Load Capacity	22 lbs 10 kg
Crash Tested	  

KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

Track Connectors

Track Clip



Part Number 0480008

Length	W	2.7 in 69 mm
	D	1.6 in 40.5 mm
	L	1.6 in 41.5 mm

Weight .2 lbs | .08 kg

Load Capacity 5 lbs | 2.27 kg

Crash Tested



KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track




(D) Depth is measured as 'Z' axis in relation to mounting to the track

(L) Length is measured as 'X' axis in relation to mounting to the track

SafePak Mount, SafePaks & SafeCabs




SafePak | Mount



Part Number	0480022
Length	W 43.9 in 1114 mm
	D 2 in 50 mm
	L 19 in 485 mm
Weight	12.79 lbs 5.8 kg
Load Capacity	25 lbs 11.38 kg
Crash Tested	  




SafePak | Small



Part Number	Red: 0480023 Blue: 0480077 White: 0480090
Length	W 12.8 in 325 mm
	D 9.4 in 240 mm
	L 18.1 in 460 mm
Weight	11.1 lbs 5.03 kg
Load Capacity	17 lbs 8 kg
Crash Tested	  

SafePak | Regular



Part Number	Red: 0480024 Blue: 0480059 White: 0480060
Length	W 18.9 in 481 mm
	D 9.4 in 240 mm
	L 18.1 in 460 mm
Weight	16.9 lbs 7.66 kg
Load Capacity	25 lbs 11.34 kg
Crash Tested	  

• Pictured with label sheets

KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

SafePak Mount, SafePaks & SafeCabs

SafeCab | Small



Part Number	0480085
Length	W 13.8 in 350 mm D 13 in 330 mm L 21.5 in 545 mm
Weight	29.8 lbs 13.5 kg
Load Capacity	11 lbs 5 kg
Crash Tested	

SafeCab | Medium



Part Number	0480086
Length	W 13.8 in 350 mm D 13 in 330 mm L 31.5 in 800 mm
Weight	36.4 lbs 16.5 kg
Load Capacity	17.6 lbs 8 kg
Crash Tested	

SafeCab | Large



Part Number	0480087
Length	W 13.8 in 350 mm D 13 in 330 mm L 46.5 in 1180 mm
Weight	51.8 lbs 23.5 kg
Load Capacity	26.5 lbs 12 kg
Crash Tested	

KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

SafePak Pouches

SafePak Pouch | Small



Part Number	0480026
Length	W 4.7 in 120 mm
	D 2.8 in 70 mm
	L 3.9 in 100 mm

Weight .3 lbs | .136 kg

Load Capacity .5 lbs | .226 kg

Crash Tested



- Pack of 2
- Label Sheets for Small Pouch available (2381858)

SafePak Pouch | Medium



Part Number	0480027
Length	W 4.7 in 120 mm
	D 2.8 in 70 mm
	L 7.9 in 200 mm

Weight 0.6 lbs | 0.25 kg

Load Capacity 2.5 lbs | 1.13 kg

Crash Tested



- Pack of 2
- Label Sheets for Medium Pouch available (2381859)

SafePak Pouch | Large Vertical



Part Number	0480028
Length	W 10.2 in 260 mm
	D 2.8 in 70 mm
	L 7.9 in 200 mm

Weight 0.8 lbs | 0.35 kg

Load Capacity 4.5 lbs | 2.04 kg

Crash Tested



- Pack of 2
- Label Sheets for Large Vertical Pouch available (2381859)

KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

SafePak Pouches

SafePak Pouch | Large Horizontal



Part Number	0480029
Length	W 4.7 in 120 mm
	D 2.8 in 70 mm
	L 15.7 in 400 mm

Weight	1.9 lbs 0.85 kg
--------	-------------------

Load Capacity	4.5 lbs 2.04 kg
---------------	-------------------

Crash Tested



- Pack of 2
Label Sheets for Large Horizontal Pouch available (2381860)

SafePak Pouch | Medium Mesh



Part Number	0480092
Length	W 4.7 in 120 mm
	D Expandable to 4 in 101.6 mm
	L 7.9 in 200 mm

Weight	1.5 lbs .7 kg
--------	-----------------

Load Capacity	2.5 lbs 1.13 kg
---------------	-------------------

Crash Tested



- Pack of 2
Pictured with Spiracle disposable BVM

SafePak Pouch Divider



Part Number	0480072
Length	W 3.9 in 99 mm
	D 2.55 in 64.8 mm
	L 1.57 in 40 mm

- Pack of 8

KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

SafePak Pouch Mounts

Pouch Mount Carrier | Single



Part Number	0480030
Length	W 4.6 in 118 mm
	D 0.6 in 15 mm
	L 7.5 in 190 mm
Weight	2.5 lbs 1.13 kg
Crash Tested	

Pouch Track Mount | Single



Part Number	0480031
Length	W 4.9 in 125 mm
	D 2.1 in 53.9 mm
	L 11.4 in 290 mm
Weight	5.3 lbs 2.41 kg
Load Capacity	2.5 lbs 1.13 kg
Crash Tested	

Pouch Track Mount | Double Vertical



Part Number	0480032
Length	W 10.6 in 268 mm
	D 2.1 in 53.9 mm
	L 11.4 in 290 mm
Weight	10.7 lbs 4.85 kg
Load Capacity	4.5 lbs 2.04 kg
Crash Tested	

KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

SafePak Pouch Mounts

Pouch Track Mount | Double Horizontal




Part Number	0480033
Length	W 4.9 in 125 mm
	D 2.1 in 53.9 mm
	L 21.5 in 545 mm
Weight	10.5 lbs 7.78 kg
Load Capacity	4.5 lbs 2.04 kg

Crash Tested



Pouch Mount Plate | 3T2




Part Number	0480155
Length	W 17.8 in 453 mm
	D 1.1 in 28.5 mm
	L 16.4 in 416 mm
Weight	4.2 lbs 1.9 kg
Load Capacity	15 lbs 6.8 kg

Crash Tested



Pouch Mount Plate | 4T2



Part Number	0480160
Length	W 23.1 in 588 mm
	D 1.1 in 28.5 mm
	L 16.4 in 416 mm
Weight	4.2 lbs 1.9 kg
Load Capacity	19.9 lbs 9.04 kg

Crash Tested






KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

SafePak Pouch Mounts




Pouch Mount Plate | 2T2



Part Number	0480161
Length	W 11.7 in 298 mm
	D 1.1 in 28.5 mm
	L 16.4 in 416 mm
Weight	2.8 lbs 1.25 kg
Load Capacity	10 lbs 4.52 kg
Crash Tested	  

Pouch Mount Plate | 1T2



Part Number	0480162
Length	W 7.2 in 183 mm
	D 1.1 in 28.5 mm
	L 16.4 in 416 mm
Weight	1.6 lbs .72 kg
Load Capacity	5.0 lbs 2.26 kg
Crash Tested	  

KEY



(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

Defibrillator Mounts

Defibrillator Wall Mount |

LifePak 15






Part Number	0480013
Length	W 14 in 356 mm
	D 12.6 in 321 mm
	L 8.3 in 210 mm
Weight	14.8 lbs 6.7 kg
Crash Tested	  

- Compatible with LifePak 15 (with or without case)

Defibrillator Swivel Mount |

LifePak 15






Part Number	0480014
Length	W 16.3 in 415 mm
	D 12 in 304 mm
	L 8.3 in 210 mm
Weight	16.3 lbs 7.4 kg
Crash Tested	  

- Compatible with LifePak 15 (with or without case)

Defibrillator Wall + Swivel Mount |

LifePak 15



Part Number	0480015
Length	W 16.3 in 415 mm
	D 13.5 in 344 mm
	L 8.3 in 210 mm
Weight	18.3 lbs 8.3 kg
Crash Tested	  

- Compatible with LifePak 15 (with or without case)

KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

Defibrillator Mounts

Defibrillator Wall Mount |

Zoll X Series



Part Number 0480019

Length	W 11.1 in 282 mm
	D 11.2 in 285 mm
	L 8.3 in 210 mm

Weight 12.3 lbs | 5.6 kg

Crash Tested	EN 1789 COMPLIANT 10g CERTIFIED	SAE J2261 COMPLIANT 22g CERTIFIED	SAE J2265 COMPLIANT 26g CERTIFIED
--------------	--	--	--

- Compatible with Zoll X Soft Case (including rear pouch)

Defibrillator Swivel Mount |

Zoll X Series



Part Number 0480020

Length	W 13.5 in 342 mm
	D 10.4 in 265 mm
	L 8.3 in 210 mm

Weight 17.4 lbs | 7.9 kg

Crash Tested	EN 1789 COMPLIANT 10g CERTIFIED	SAE J2261 COMPLIANT 22g CERTIFIED	SAE J2265 COMPLIANT 26g CERTIFIED
--------------	--	--	--

- Compatible with Zoll X Soft Case (including rear pouch)

Defibrillator Wall + Swivel Mount |

Zoll X Series



Part Number 0480021

Length	W 13.5 in 342 mm
	D 12 in 305 mm
	L 8.3 in 210 mm

Weight 19.4 lbs | 8.8 kg

Crash Tested	EN 1789 COMPLIANT 10g CERTIFIED	SAE J2261 COMPLIANT 22g CERTIFIED	SAE J2265 COMPLIANT 26g CERTIFIED
--------------	--	--	--

- Compatible with Zoll X Soft Case (including rear pouch)




KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

Defibrillator Mounts


Wall Mount | Corpuls 3



Part Number	0480056
Length	W 7.7 in 196 mm
	D 3.9 in 98 mm
	L 10.3 in 261 mm
Weight	3.5 lbs 1.59 kg
Crash Tested	  




Wall Mount | Corpuls Monitor



Part Number	0480093
Length	W 7.8 in 199 mm
	D 3.9 in 98 mm
	L 6.6 in 167 mm
Weight	2.9 lbs 1.3 kg
Crash Tested	

Swivel Mount | Corpuls 3



Part Number	0480094
Length	W 10.3 in 261 mm
	D 9.8 in 250 mm
	L 10.3 in 262 mm
Weight	4.9 lbs 2.24 kg
Crash Tested	  

KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track




(D) Depth is measured as 'Z' axis in relation to mounting to the track

(L) Length is measured as 'X' axis in relation to mounting to the track

Defibrillator Mounts




Wall + Swivel Mount |
Corpuls 3



Part Number	0480095
Length	W 10.3 in 261 mm D 10.5 in 266 mm L 10.3 in 262 mm
Weight	7.6 lbs 3.44 kg
Crash Tested	  




Wall Mount | Physio LP1000



Part Number	0480084
Length	W 7.7 in 195 mm D 2.4 in 62 mm L 10 in 253 mm
Weight	3.1 lbs 1.4 kg
Crash Tested	  

Wall Mount | PD360



Part Number	0480151
Length	W 5.9 in 151 mm D 15.1 in 384 mm L 6.6 in 167 mm
Weight	4 lbs 1.8 kg
Crash Tested	  

• Compatible with D360 ProPac. Sold separately.



KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

Defibrillator Mounts

D360 ProPac/ZOLL Adaptor |
P/N FA868A151



Part Number	0480156
Length	W 1.6 in 41 mm
	D 4.5 in 114 mm
	L 7.1 in 180 mm
Weight	3.7 lbs 1.7 kg
Crash Tested	 

• Compatible with D360 ProPac. Sold separately.




KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

Oxygen

Oxygen Mount | Triple






Part Number	0480035
Length	W 28.7 in 730 mm
	D 9.1 in 230 mm
	L 27 in 686 mm
Weight	50 lbs 22.68 kg
Crash Tested	  

- Weight does not include oxygen cylinders
- Compatible with Spiracle EE-Lite tanks (081-9589) combined with Model 804 Regulator (081-9567). Both sold separately.

Oxygen Mount | Double






Part Number	0480157 Configured for EE-Lite 3K O2 Cylinder 0480145 Configured for EE-Lite 3K O2 Cylinder w/pres. Transducers 0480126 Configured for M60 Aluminum O2 Cylinder 0480158 Configured for M60 Aluminum O2 Cylinder w/pres. Transducers
Length	W 34.5 in 876 mm
	D 8.9 in 226 mm
	L 17.8 in 451 mm
Weight	34 lbs 15.4 kg
Crash Tested	  

- Weight does not include oxygen cylinders or transducers
- Compatible with EE-Lite 3K O2, or M60 Aluminum O2 Cylinders. Sold separately.

Oxygen Mount | Single



Part Number	0480034
Length	W 18.9 in 480 mm
	D 8.3 in 210 mm
	L 6.6 in 167 mm
Weight	4.9 lbs 2.2 kg
Crash Tested	  

- Weight does not include oxygen cylinder.
- Compatible with 'D' or Jumbo 'D' tanks. Sold separately.




KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

Ventilator Mounts




Vent Mount | Impact 731



Part Number	0480038
Length	W 7.8 in 199 mm
	D 2.5 in 64 mm
	L 6.6 in 167 mm
Weight	2 lbs 0.9 kg
Crash Tested	  




Swivel Mount | HT70 Vent



Part Number	0480052
Length	W 13.5 in 344 mm
	D 2.6 in 65 mm
	L 8.3 in 211 mm
Weight	5.1 lbs 2.3 kg
Crash Tested	  

Wall Mount | HT70 Vent



Part Number	0480058
Length	W 6.1 in 155 mm
	D 1.8 in 45 mm
	L 6.6 in 167 mm
Weight	3.3 lbs 1.5 kg
Crash Tested	  

KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

Ventilator Mounts

Mounting Plate | Pneupac



Part Number 0480074

Length W 7.8 in | 199 mm
D 2.4 in | 60 mm
L 6.6 in | 167 mm

Weight 2.4 lbs | 1.08 kg

Crash Tested   

Vent Mount | MeduMat



Part Number 0480079

Length W 5.9 in | 151 mm
D 2.2 in | 57 mm
L 6.6 in | 167 mm

Weight 2.7 lbs | 1.2 kg

Crash Tested   

Wall Vent Mount | Draeger



Part Number 0480109

Length W 11.2 in | 285 mm
D 10.2 in | 259 mm
L 9.5 in | 242 mm

Weight 6.9 lbs | 3.14 kg

Crash Tested   




KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
(D) Depth is measured as 'Z' axis in relation to mounting to the track
(L) Length is measured as 'X' axis in relation to mounting to the track

Ventilator Mounts

Swivel Vent Mount | Draeger



Part Number	0480110
Length	W 13.8 in 351 mm
	D 9.9 in 251 mm
	L 9.6 in 243 mm
Weight	8.7 lbs 3.93 kg
Crash Tested	
	
	




KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

Suction Mounts

LSU Mount | Laerdal






Part Number	0480042
Length	W 6.8 in 173 mm
	D 2.2 in 57 mm
	L 7.5 in 190 mm
Weight	6.6 lbs 3 kg
Crash Tested	  

• Requires LSU Wall Bracket (with or without power cord).




Suction Mount | SSCOR



Part Number	0480043
Length	W 5.9 in 151 mm
	D 2.4 in 60 mm
	L 7.9 in 200 mm
Weight	2.6 lbs 1.2 kg
Crash Tested	  

Suction Mount | Weinmann



Part Number	0480078
Length	W 5.9 in 151 mm
	D 2.2 in 57 mm
	L 6.6 in 167 mm
Weight	2.6 lbs 1.17 kg
Crash Tested	  




KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

Suction Mounts




HERSILL Mount |
V7 Plus B Emergency



Part Number	0480124
Length	W 5.9 in 151 mm
	D 1.9 in 48 mm
	L 6.6 in 167 mm
Weight	2.3 lbs 1.05 kg
Crash Tested	  

Mount | Weinmann Base 1NG



Part Number	0480088
Length	W 6.4 in 163 mm
	D 2.2 in 57 mm
	L 6.6 in 167 mm
Weight	2.8 lbs 1.25 kg
Crash Tested	  




KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

Pole & Rail Mounts




Pole Mount | Single



Part Number	0480080
Length	W 8.1 in 206 mm
	D 5.5 in 140 mm
	L 6.6 in 167 mm
Weight	3.6 lbs 1.62 kg
Crash Tested	  



Pole Mount | Double



Part Number	0480081
Length	W 12.6 in 320 mm
	D 4.1 in 103 mm
	L 6.6 in 167 mm
Weight	3.7 lbs 1.7 kg
Crash Tested	  

Pole Mount | Vertical



Part Number	0480112
Length	W 5.9 in 151 mm
	D 10 in 255 mm
	L 6.6 in 167 mm
Weight	2.8 lbs 1.29 kg
Crash Tested	  

KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

Pole & Rail Mounts

Medrail Mount | 25mm



Part Number	0480082
Length	W 5.9 in 151 mm
	D 3.2 in 81 mm
	L 12 in 305 mm
Weight	3.1 lbs 1.4 kg
Crash Tested	

Medrail Mount | 30mm



Part Number	0480091
Length	W 5.9 in 151 mm
	D 3.2 in 81 mm
	L 12 in 305 mm
Weight	3.2 lbs 1.44 kg
Crash Tested	

KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track

(D) Depth is measured as 'Z' axis in relation to mounting to the track

(L) Length is measured as 'X' axis in relation to mounting to the track

Accessory Mounts

Litter Arm | Rotating



Part Number	0480036
Length	W 7.2 in 182 mm
	D 28.7 in 730 mm
	L 6.6 in 167 mm
Weight	13.2 lbs 6 kg
Load Capacity	352.7 lbs 160 kg (per pair)

Litter Arm | Loading Wedge



Part Number	0480153
Length	W 1.8 in 45 mm
	D 7.7 in 195 mm
	L 2 in 52 mm
Weight	.8 lbs .37 kg

Ram Adapter | 2.5 SMC



Part Number	0480047
Length	W 2.5 in 64 mm
	D 0.3 in 7 mm
	L 2.5 in 64 mm
Weight	0.1 lbs 0.06 kg
Crash Tested	




KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

Accessory Mounts




Work Tray | With Arm



Part Number	0480039
Length	W 7.8 in 199 mm D 12.7 in 323 mm L 12.7 in 323 mm
Weight	3 lbs 1.37 kg
Crash Tested	  




Work Tray | Hinged 2.0

PHOTO NOT AVAILABLE

Part Number	0480123
Length	W 4.6 in 117 mm D 14.3 in 362 mm L 12.7 in 323 mm
Weight	2.4 lbs 1.11 kg
Crash Tested	  

Tablet Mount



Part Number	0480041
Length	W 7.7 in 196 mm D 10.8 in 275 mm L 10.8 in 275 mm
Weight	3.75 lbs 1.7 kg
Load Capacity	5 lbs 2.75 kg
Crash Tested	  
<ul style="list-style-type: none"> • Compatible with iPad Pro 9.7" 	

KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

Accessory Mounts

Trash Bin Mount



Part Number 0480044

Length	W	10.3 in 262 mm
	D	9.2 in 234 mm
	L	10.6 in 268 mm

Weight 4.4 lbs | 2 kg

Load Capacity	2.2 lbs 1 kg
---------------	----------------

- Complete with Trash Bin (0480075)
- 90-Day Warranty on Trash Bin

Sharps Bin Mount



Part Number 0480045

Length	W	13.8 in 350 mm
	D	6.8 in 173 mm
	L	13 in 330 mm

Weight 4.9 lbs | 2.2 kg

Load Capacity	2.2 lbs 1 kg
---------------	----------------

Crash Tested



- Complete with Sharps Bin (0480076)
- 90-Day Warranty on Sharps Bin

Lock Knob | Rotating



Part Number 0480071

Length	W	1.5 in 37.6 mm
	D	2.1 in 53 mm
	L	1.5 in 37.6 mm

Weight 0.5 lbs | 0.23 kg

Crash Tested	SAE J2265 COMPLIANT 10g CRASH TESTED	SAE J2265 COMPLIANT 22g CRASH TESTED	SAE J2265 COMPLIANT 26g CRASH TESTED
--------------	--------------------------------------	--------------------------------------	--------------------------------------

- Pack of 2




KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

Accessory Mounts




T Slot | Foot



Part Number	0480089
Length	W 0.7 in 19 mm
	D 0.3 in 7 mm
	L 1.2 in 31 mm
Weight	0.02 lbs 0.01 kg
Crash Tested	  


Toughbook Mount | SAE



Part Number	0480040 SAE
Length	W 10 in 255 mm
	D 12.2 in 311 mm
	L 12.2 in 311 mm
Weight	9.15 lbs 4.15 kg
Load Capacity	15 lbs 6.75 kg
Crash Tested	  
<ul style="list-style-type: none"> • Compatible with Panasonic model CF-19 	

Toughbook Mount | CEN



Part Number	0480070 CEN
Length	W 10 in 255 mm
	D 12.2 in 311 mm
	L 12.2 in 311 mm
Weight	5.78 lbs 2.62 kg
Load Capacity	10 lbs 4.5 kg
Crash Tested	
<ul style="list-style-type: none"> • Compatible with Panasonic model CF-19 	




KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

Spare Parts

SafeMount | Connector



Part Number	0480012
Length	W 1.8 in 45.7 mm
	D .6 in 14.9 mm
	L —
Weight	.1 lbs .045 kg
Load Capacity	22 lbs 10kg
Crash Tested	  

Spacer Plate | 0.5 in



Part Number	0480057
Length	W 3.3 in 85 mm
	D .05 in 12.7 mm
	L 3.3 in 85 mm
Weight	0.3 lbs 0.12 kg

Splicing Key



Part Number	0480054
Length	W 0.3 in 8 mm
	D 0.2 in 6 mm
	L 2.4 in 60 mm
Weight	.019 lbs 8.6 g

KEY

(W) Width is measured as 'Y' axis in relation to mounting to the track
 (D) Depth is measured as 'Z' axis in relation to mounting to the track
 (L) Length is measured as 'X' axis in relation to mounting to the track

Spare Parts

V Block | Tall



Part Number	0480065
Length	W 4.4 in 110 mm
	D 0.6 in 15 mm
	L 2.8 in 72 mm

V Block | Short



Part Number	0480066
Length	W 3.4 in 87 mm
	D 0.6 in 15 mm
	L 2.8 in 72 mm

KEY

- (W) Width is measured as 'Y' axis in relation to mounting to the track
- (D) Depth is measured as 'Z' axis in relation to mounting to the track
- (L) Length is measured as 'X' axis in relation to mounting to the track

Product Spec Index

ITEM NUMBER	PRODUCT	WIDTH	DEPTH	LENGTH	WEIGHT	MATERIAL	FINISH	CERTIFICATION STANDARD	LOAD CAPACITY/ SWL	SOLD IN USA
0480001	Track System 1016mm	7.4 in 189 mm	.5 in 12.7 mm	40 in 1016 mm	7.5 lbs 3.4 kg	Aluminum	Clear Natural Anodized	10G 22G 26G	220.4 lbs 100 kg	Yes
0480062	Track System 340mm, Black End Caps	7.6 in 193 mm	.7 in 17 mm	13.4 in 340 mm	2.6 lbs 1.2 kg	6005A Aluminum	Clear Natural Anodized	10G 22G 26G	55.1 lbs 25 kg	Yes
0480073	Track 1016 No Holes, (Bonding)	7.4 in 189 mm	.5 in 12.7 mm	40 in 1016 mm	7.5 lbs 3.4 kg	Aluminum	Clear Natural Anodized	10G 22G 26G	220.4 lbs 100 kg	Yes
0480107	Ceiling Track Right	7.4 in 189 mm	.5 in 12.7 mm	40 in 1016 mm	7.5 lbs 3.4 kg	Aluminum	Clear Natural Anodized	10G 22G 26G	220.4 lbs 100 kg	Yes
0480131	Ceiling Track Left	7.4 in 189 mm	.5 in 12.7 mm	40 in 1016 mm	7.5 lbs 3.4 kg	Aluminum	Clear Natural Anodized	10G 22G 26G	220.4 lbs 100 kg	Yes
0480108	Ceiling Track Center	7.4 in 189 mm	.5 in 12.7 mm	40 in 1016 mm	7.5 lbs 3.4 kg	Aluminum	Clear Natural Anodized	10G 22G 26G	220.4 lbs 100 kg	Yes
0480002 0480003 0480004	Black End Caps Regular Grey End Caps Regular White End Caps Regular	7.6 in 193 mm	1.9 in 48 mm	.7 in 17 mm	.605 lbs 275 g	Zytel	Black, Grey, or White	10G 22G 26G	—	Yes
0480049 0480050 0480051	Black End Caps Thin Grey End Caps Thin White End Caps Thin	7.4 in 193 mm	1.5 in 48mm	.7 in 17 mm	.605 lbs 275 g	Zytel	Black, Grey, or White	10G 22G 26G	—	Yes
0480005 0480006 0480007	Black T-Splice Caps Grey T-Splice Caps White T-Splice Caps	7.6 in 193 mm	.9 in 24 mm	.6 in 16 mm	.14 lbs 64 g	Zytel	Black, Grey, or White	10G 22G 26G	—	Yes
0480010	SafeMount Threaded (25kg)	5.9 in 151 mm	1.4 in 34.5 mm	6.6 in 167 mm	1.3 lbs .59 kg	Aluminum 7075-T6, Stainless Steel, Plastic	Clear Natural Anodized, Powder Coated White	10G 22G 26G	25 kg	Yes
0480146	SafeMount Threaded (12kg)	5.9 in 151 mm	1.6 in 40 mm	6.6 in 167 mm	1.6 lbs .74 kg	Aluminum 7075-T6, Stainless Steel, Plastic	Clear Natural Anodized, Powder Coated White	10G 22G 26G	26.5 lbs 12 kg	Yes
0480011	SafeMount With Connector	7.4 in 187 mm	1.9 in 48.5 mm	6.6 in 167 mm	1.65 lbs .75 kg	Aluminum 7075-T6, Stainless Steel, Plastic	Clear Natural Anodized, Powder Coated White	10G 22G 26G	22 lbs 10 kg	Yes
0480012	SafeMount Connector	1.8 in 45.7 mm	.6 in 14.9 mm	—	.1 lbs .045 kg	Aluminum 7075-T6	Clear Natural Anodized	10G 22G 26G	22 lbs 10 kg	Yes
0480008	Track Clip	2.7 in 69 mm	1.6 in 40.5 mm	1.6 in 41.5 mm	2 lbs .08 kg	Aluminum 7075-T6, Stainless Steel, Plastic	Clear Natural Anodized, Powder Coated Red & White	10G 22G 26G	5 lbs 2.27 kg	Yes
0480022	SafePak Mount	43.9 in 1114 mm	2 in 50 mm	19 in 485 mm	12.79 lbs 5.8 kg	Plastic, Aluminum	Matte, Powder Coated White & Red	10G 22G 26G	25 lbs 11.38 kg	Yes
0480023 0480077 0480090	Red SafePak Small Blue SafePak Small White SafePak Small	12.8 in 325 mm	9.4 in 240 mm	18.1 in 460 mm	11.1 lbs 5.03 kg	Plastic	Red, Blue, or White	10G 22G 26G	17 lbs 7.71 kg	Yes
0480024 0480059 0480060	Red SafePak Regular Blue SafePak Regular White SafePak Regular	18.9 in 481 mm	9.4 in 240 mm	18.1 in 460 mm	16.9 lbs 7.66 kg	Plastic	Red, Blue, or White	10G 22G 26G	25 lbs 11.34 kg	Yes
0480085	SafeCab Small	13.8 in 350 mm	13 in 330 mm	21.5 in 545 mm	29.8 lbs 13.5 kg	CPVC	White with Red Trim	10G 22G 26G	11 lbs 5 kg	Yes
0480086	SafeCab Medium	13.8 in 350 mm	13 in 330 mm	31.5 in 800 mm	36.4 lbs 16.5 kg	CPVC	White with Red Trim	10G 22G 26G	17.6 lbs 8 kg	Yes
0480087	SafeCab Large	13.8 in 350 mm	13 in 330 mm	46.5 in 1180 mm	51.8 lbs 23.5 kg	CPVC	White with Red Trim	10G 22G 26G	26.5 lbs 12 kg	Yes
0480026	SafePak Pouch Small	4.7 in 120 mm	2.8 in 70 mm	3.9 in 100 mm	.3 lbs .136 kg	Plastic	Glossy Grey Shell, Matte Grey Plastics	10G 22G 26G	.5 lbs .226 kg	Yes
0480027	SafePak Pouch Medium	4.7 in 120 mm	2.8 in 70 mm	7.9 in 200 mm	0.6 lbs 0.25 kg	Plastic	Glossy Grey Shell, Matte Grey Plastics	10G 22G 26G	2.5 lbs 1.13 kg	Yes
0480028	SafePak Pouch Large Vertical	10.2 in 260 mm	2.8 in 70 mm	7.9 in 200 mm	0.8 lbs 0.35 kg	Plastic	Glossy Grey Shell, Matte Grey Plastics	10G 22G 26G	4.5 lbs 2.04 kg	Yes
0480029	SafePak Pouch Large Horizontal	4.7 in 120 mm	2.8 in 70 mm	15.7 in 400 mm	1.9 lbs 0.85 kg	Plastic	Glossy Grey Shell, Matte Grey Plastics	10G 22G 26G	4.5 lbs 2.04 kg	Yes
0480092	SafePak Pouch Medium Mesh	4.7 in 120 mm	Expandable to 4 in 101.6 mm	7.9 in 200 mm	1.5 lbs .7 kg	Flexible Mesh	Mesh Fabric, Matte Grey Plastics	10G 22G 26G	2.5 lbs 1.13 kg	Yes
0480072	SafePak Pouch Divider	3.9 in 99 mm	2.55 in 64.8 mm	1.57 in 40 mm	—	Plastic	Grey	—	—	Yes
0480030	Pouch Mount Carrier Single	4.6 in 118 mm	0.6 in 15 mm	7.5 in 190 mm	0.1 lbs 0.04 kg	Plastic	Grey, Red	10G 22G 26G	2.5 lbs 1.13 kg	Yes
0480031	Pouch Track Mount Single	4.9 in 125 mm	2.1 in 53.9 mm	11.4 in 290 mm	5.3 lbs 2.41 kg	Aluminum, ABS Plastic	Powder Coated White, Grey Plastics	10G 22G 26G	2.5 lbs 1.13 kg	Yes
0480032	Pouch Track Mount Double Vertical	10.6 in 268 mm	2.1 in 53.9 mm	11.4 in 290 mm	10.7 lbs 4.85 kg	Aluminum, ABS Plastic	Powder Coated White, Grey Plastics	10G 22G 26G	4.5 lbs 2.04 kg	Yes
0480033	Pouch Track Mount Double Horizontal	4.9 in 125 mm	2.1 in 53.9 mm	21.5 in 545 mm	10.5 lbs 7.78 kg	Aluminum, ABS Plastic	Powder Coated White, Grey Plastics	10G 22G 26G	4.5 lbs 2.04 kg	Yes
0480155	Pouch Mount Plate 3T2	17.8 in 453 mm	1.1 in 28.5 mm	16.4 in 416 mm	4.2 lbs 1.9 kg	Aluminum, ABS Plastic	Powder Coated White	10G 22G 26G	15 lbs 6.8 kg	Yes

Product Spec Index

ITEM NUMBER	PRODUCT	WIDTH	DEPTH	LENGTH	WEIGHT	MATERIAL	FINISH	CERTIFICATION STANDARD	LOAD CAPACITY/ SWL	SOLD IN USA
0480160	Pouch Mount Plate 4T2	23.1 in 588 mm	1.1 in 28.5 mm	16.4 in 416 mm	4.2 lbs 1.9 kg	Aluminum, ABS Plastic	Powder Coated White	10G 22G 26G	19.9 lbs 9.04 kg	Yes
0480161	Pouch Mount Plate 2T2	11.7 in 298 mm	1.1 in 28.5 mm	16.4 in 416 mm	2.8 lbs 1.25 kg	Aluminum, ABS Plastic	Powder Coated White	10G 22G 26G	10 lbs 4.52 kg	Yes
0480162	Pouch Mount Plate 1T2	7.2 in 183 mm	1.1 in 28.5 mm	16.4 in 416 mm	1.6 lbs .72 kg	Aluminum, ABS Plastic	Powder Coated White	10G 22G 26G	5.0 lbs 2.26 kg	Yes
0480013	Defibrillator Wall Mount LifePak 15	14 in 356 mm	12.6 in 321 mm	8.3 in 210 mm	14.8 lbs 6.7 kg	Aluminum 7075-T6, Stainless Steel, Plastic	Matte, Powder Coated White & Red	10G 22G 26G	—	Yes
0480014	Defibrillator Swivel Mount LifePak 15	16.3 in 415 mm	12 in 304 mm	8.3 in 210 mm	16.3 lbs 7.4 kg	Aluminum 7075-T6, Stainless Steel, Plastic	Matte, Powder Coated White & Red	10G 22G 26G	—	Yes
0480015	Defibrillator Wall+Swivel Mount LifePak 15	16.3 in 415 mm	13.5 in 344 mm	8.3 in 210 mm	18.3 lbs 8.3 kg	Aluminum 7075-T6, Stainless Steel, Plastic	Matte, Powder Coated White & Red	10G 22G 26G	—	Yes
0480019	Defibrillator Wall Mount Zoll X Series	11.1 in 282 mm	11.2 in 285 mm	8.3 in 210 mm	12.3 lbs 5.6 kg	Plastic, Aluminum	Matte, Powder Coated White & Red	10G 22G 26G	—	Yes
0480020	Defibrillator Swivel Mount Zoll X Series	13.5 in 342 mm	10.4 in 265 mm	8.3 in 210 mm	17.4 lbs 7.9 kg	Plastic, Aluminum	Matte, Powder Coated White & Red	10G 22G 26G	—	Yes
0480021	Defibrillator Wall+Swivel Mount Zoll X Series	13.5 in 342 mm	12 in 305 mm	8.3 in 210 mm	19.4 lbs 8.8 kg	Plastic, Aluminum	Matte, Powder Coated White & Red	10G 22G 26G	—	Yes
0480056	Wall Mount Corplus 3	7.7 in 196 mm	3.9 in 98 mm	10.3 in 261 mm	3.5 lbs 1.59 kg	Aluminum	Powder Coated White	10G 22G 26G	—	Yes
0480093	Wall Mount Corplus Monitor	7.8 in 199 mm	3.9 in 98 mm	6.6 in 167 mm	2.9 lbs 1.3 kg	Aluminum	Powder Coated White	10G	—	Yes
0480094	Swivel Mount Corplus 3	10.3 in 261 mm	9.8 in 250 mm	10.3 in 262 mm	4.9 lbs 2.24 kg	Aluminum	Powder Coated White	10G 22G 26G	—	Yes
0480095	Wall + Swivel Mount Corplus 3	10.3 in 261 mm	10.5 in 266 mm	10.3 in 262 mm	7.6 lbs 3.44 kg	Aluminum	Powder Coated White	10G 22G 26G	—	Yes
0480084	Wall Mount Physio LP1000	7.7 in 195 mm	2.4 in 62 mm	10 in 253 mm	3.1 lbs 1.4 kg	Aluminum	Powder Coated White	10G 22G 26G	—	Yes
0480151	Wall Mount D360	5.9 in 151 mm	15.1 in 384 mm	6.6 in 167 mm	4 lbs 1.8 kg	Aluminum	Powder Coated White	10G 22G 26G	—	—
0480156	D360 ProPac/Zoll Adaptor P/N FA868A151	1.6 in 41 mm	4.5 in 114 mm	7.1 in 180 mm	3.7 lbs 1.7 kg	Aluminum	Powder Coated White	10G 20G 26G	—	Yes
0480035	Oxygen Mount Triple	28.7 in 730 mm	9.1 in 230 mm	27 in 686 mm	50 lbs 22.68 kg	Plastic, Aluminum	Powder Coated White	10G 22G 26G	—	Yes
0480126 0480145 0480157 0480158	Oxygen Mount Double	34.5 in 876 mm	8.9 in 226 mm	17.8 in 451 mm	34 lbs 15.4 kg	Plastic, Aluminum	Powder Coated White	10G 22G 26G	—	Yes
0480034	Oxygen Mount Single	18.9 in 480 mm	8.3 in 210 mm	6.6 in 167 mm	4.9 lbs 2.2 kg	Plastic, Aluminum	Powder Coated White	10G 22G 26G	25 kg	Yes
0480038	Vent Mount Impact 731	7.8 in 199 mm	2.5 in 64 mm	6.6 in 167 mm	2 lbs 0.9 kg	Plastic, Aluminum	Clear Natural Anodized, Powder Coated White	10G 22G 26G	—	Yes
0480052	Swivel Mount HT70 Vent	13.5 in 344 mm	2.6 in 65 mm	8.3 in 211 mm	5.1 lbs 2.3 kg	Aluminum	Powder Coated White	10G 22G 26G	—	Yes
0480058	Wall Mount HT70 Vent	6.1 in 155 mm	1.8 in 45 mm	6.6 in 167 mm	3.3 lbs 1.5 kg	Aluminum	Powder Coated White	10G 22G 26G	—	Yes
0480074	Mounting Plate PneuPac	7.8 in 199 mm	2.4 in 60 mm	6.6 in 167 mm	2.4 lbs 1.08 kg	Aluminum	Powder Coated White	10G 22G 26G	—	Yes
0480079	Vent Mount MeduMat	5.9 in 151 mm	2.2 in 57 mm	6.6 in 167 mm	2.7 lbs 1.2 kg	Aluminum	Clear Natural Anodized	10G 22G 26G	—	Yes
0480109	Wall Vent Mount Draeger	11.2 in 285 mm	10.2 in 259 mm	9.5 in 242 mm	6.9 lbs 3.14 kg	Aluminum	Powder Coated White	10G 22G 26G	—	Yes
0480110	Swivel Vent Mount Draeger	13.8 in 351 mm	9.9 in 251 mm	9.6 in 243 mm	8.7 lbs 3.93 kg	Aluminum	Powder Coated White	10G 22G 26G	—	Yes
0480042	LSU Mount Laerdal	6.8 in 173 mm	2.2 in 57 mm	7.5 in 190 mm	6.6 lbs 3 kg	Aluminum	Matte Anodized	10G 22G 26G	—	Yes
0480043	Suction Mount SSCOR	5.9 in 151 mm	2.4 in 60 mm	7.9 in 200 mm	2.6 lbs 1.2 kg	Aluminum	Matte Anodized	10G 22G 26G	—	Yes
0480078	Suction Mount Weinmann	5.9 in 151 mm	2.2 in 57 mm	6.6 in 167 mm	2.6 lbs 1.17 kg	Aluminum	Clear Natural Anodized	10G 22G 26G	—	Yes
0480124	HERSILL Mount V7 Plus B Emergency	5.9 in 151 mm	1.9 in 48 mm	6.6 in 167 mm	2.3 lbs 1.05 kg	Aluminum	Powder Coated White	10G 22G 26G	—	Yes
0480088	Mount Weinmann Base 1NG	6.4 in 163 mm	2.2 in 57 mm	6.6 in 167 mm	2.8 lbs 1.25 kg	Aluminum	Powder Coated White	10G 22G 26G	—	Yes
0480080	Pole Mount Single	8.1 in 206 mm	5.5 in 140 mm	6.6 in 167 mm	3.6 lbs 1.62 kg	Aluminum	Clear Natural Anodized	10G 22G 26G	—	Yes
0480081	Pole Mount Double	12.6 in 320 mm	4.1 in 103 mm	6.6 in 167 mm	3.7 lbs 1.7 kg	Aluminum	Clear Natural Anodized	10G 22G 26G	—	Yes
0480112	Pole Mount Vertical	5.9 in 151 mm	10 in 255 mm	6.6 in 167 mm	2.8 lbs 1.29 kg	Aluminum	Clear Natural Anodized	10G 22G 26G	—	Yes
0480082	Medrail Mount 25mm	5.9 in 151 mm	3.2 in 81 mm	12 in 305 mm	3.1 lbs 1.4 kg	Aluminum	Clear Natural Anodized	10G 22G 26G	—	Yes
0480091	Medrail Mount 30mm	5.9 in 151 mm	3.2 in 81 mm	12 in 305 mm	3.2 lbs 1.44 kg	Aluminum	Clear Natural Anodized	10G 22G 26G	—	Yes
0480036	Litter Arm Rotating	7.2 in 182 mm	28.7 in 730 mm	6.6 in 167 mm	13.2 lbs 6 kg	Aluminum	Powder Coated White	—	352.7 lbs 160 kg per pair	Yes
0480153	Litter Arm Loading Wedge	1.8 in 45 mm	7.7 in 195 mm	2 in 52 mm	.8 lbs .37 kg	Plastic	Black	—	—	Yes

Product Spec Index

ITEM NUMBER	PRODUCT	WIDTH	DEPTH	LENGTH	WEIGHT	MATERIAL	FINISH	CERTIFICATION STANDARD	LOAD CAPACITY/ SWL	SOLD IN USA
0480039	Work Tray With Arm	7.8 in 199 mm	12.7 in 323 mm	12.7 in 323 mm	3 lbs 1.37 kg	Aluminum, Plastic	White	10G 22G 26G	—	Yes
0480123	Work Tray Hinged 2.0	4.6 in 117 mm	14.3 in 362 mm	12.7 in 323 mm	2.4 lbs 1.11 kg	Aluminum, Plastic	White	10G 22G 26G	—	Yes
0480041	Tablet Mount	7.7 in 196 mm	10.8 in 275 mm	10.8 in 275 mm	3.75 lbs 1.7 kg	Plastic	Matte Black	10G 22G 26G	5 lbs 2.75 kg	Yes
0480044 0480075	Trash Bin Mount Trash Bin Container	10.3 in 262 mm	9.2 in 234 mm	10.6 in 268 mm	4.4 lbs 2 kg	Plastic, Aluminum	Powder Coated Black	10G 22G 26G	2.2 lbs 1 kg	Yes
0480045 0480076	Sharps Bin Mount Sharps Bin Container	13.8 in 350 mm	6.8 in 173 mm	13 in 330 mm	4.9 lbs 2.2 kg	Plastic, Aluminum	Powder Coated Black	10G 22G 26G	2.2 lbs 1 kg	Yes
0480047	Ram Adapter 2.5 SMC	2.5 in 64 mm	0.3 in 7 mm	2.5 in 64 mm	0.1 lbs 0.06 kg	Aluminum	Black	10G 22G 26G	—	Yes
0480071	Lock Knob Rotating	1.5 in 37.6 mm	2.1 in 53 mm	1.5 in 37.6 mm	0.5 lbs 0.23 kg	Stainless Steel & Plastic	White	10G 22G 26G	—	Yes
0480089	T Slot Foot	0.7 in 19 mm	0.3 in 7 mm	1.2 in 31 mm	0.02 lbs 0.01 kg	Aluminum	Clear Natural Anodized	10G 22G 26G	—	Yes
0480054	Splicing Key	—	—	2.4 in 60 mm	.019 lbs 8.6 g	Aluminum 7075-T6	None	10G 22G 26G	—	Yes
0480040	Toughbook Mount SAE	10 in 255 mm	12.2 in 311 mm	12.2 in 311 mm	9.15 lbs 4.15 kg	Aluminum	Matte Anodized Black	10G 22G 26G	15 lbs 6.75 kg	Yes
0480070	Toughbook Mount CEN	10 in 255 mm	12.2 in 311 mm	12.2 in 311 mm	5.78 lbs 2.62 kg	Aluminum	Matte Anodized Black	10G 22G	10 lbs 4.5 kg	Yes
0480057	Spacer Plate 0.5 in	3.3 in 85 mm	0.5 in 12.7 mm	3.3 in 85 mm	0.3 lbs 0.12 kg	Plastic	White	—	—	Yes
0480065	V Block Tall	4.4 in 110 mm	0.6 in 15 mm	2.8 in 72 mm	—	Plastic	White	—	—	Yes
0480066	V Block Short	3.4 in 87 mm	0.6 in 15 mm	2.8 in 72 mm	—	Plastic	White	—	—	Yes



70 Weil Way
Wilmington, OH 45177

Ambulance Systems

Aviation

Emergency Response System

Military

Mortuary

Rescue
