



## Disposable Laparoscopic Specimen Bags

Self-opening, high-strength bags ensure safe specimen capture with easy, one-handed operation.



**UNIVERSAL PLUS™**  
SUCTION/IRRIGATION



**DETACHATIP® III**  
SYSTEM



**VCARE®**  
UTERINE MANIPULATOR



To learn more about these and other innovative products, call 800-448-6506 or visit [CONMED.com](http://CONMED.com)



525 French Road, Utica, NY 13502

**DISPOSABLE  
LAPAROSCOPIC  
SPECIMEN BAGS**

Sterile, disposable device designed for use as a receptacle for the collection and extraction of tissue, organs, and calculi during laparoscopic surgical procedures. The Specimen Bag System consists of a flexible memory wire specimen bag and introducer. The Specimen Bag System can be used through a 5, 10, or 12mm cannula.

**Unique Design**


The unique patented design ensures ease of use, and can be deployed and manipulated with one hand. The self-opening memory wire, high-strength large bag offers easy and safe specimen capture while resisting tearing, reducing the potential for cell migration.

**Easy to Use**

Indicators on the thumb ring position the opening (top) of the bag for correct placement. Bag can be detached within abdominal cavity for later retrieval.

**Benefits**

- Large bag volume to capture large specimens
- Minimum insertion force enables smooth deployment
- Ergonomic design provides easy one-handed operation
- Insert indicator ensures bag will be in proper position when deployed



**Ordering Information**

Description	Cannula	Diameter	Length	Volume	Qty Per Box	Catalog Number
Attached Memory Wire Bag, Polymer	5mm	3in (7.62cm)	4in (10.2cm)	120ml	10/case 1/pkg	SB534
Independent Memory Wire Bag, Polymer	10mm	3in (7.62cm)	6in (15.24cm)	255ml	10/case 1/pkg	SB836
Independent Memory Wire Bag, Polymer	10mm	5in (12.7cm)	7in (17.78cm)	680ml	10/case 1/pkg	SB857
Attached Memory Wire Bag, Polymer	10mm	3in (7.62cm)	6in (15.24cm)	255ml	5/case 1/pkg	SB936
Attached Memory Wire Bag, Polymer	10mm	5in (12.7cm)	7in (17.78cm)	680ml	5/case 1/pkg	SB957
Attached Memory Wire Bag, Polymer	12mm	7.5in (19.1cm)	9in (22.9cm)	1500ml	5/case 1/pkg	SB979
Attached Memory Wire Bag, Nylon	10mm	3in (7.62cm)	6in (15.24cm)	255ml	5/case 1/pkg	SB1036
Attached Memory Wire Bag, Nylon	10mm	5in (12.7cm)	7in (17.78cm)	680ml	5/case 1/pkg	SB1057
Attached Memory Wire Bag, Nylon	12mm	7.5in (19.1cm)	9in (22.9cm)	1500ml	5/case 1/pkg	SB1079

Customer Service: 1-800-448-6506  
 International Sales: 1-315-797-8375  
 Fax: 1-800-438-3051

CONMED.com  
 info@mail.CONMED.com