




GoldLine® Electrosurgical Pencils

Key Safety Features

- Hex hub fittings accommodate industry standard 3/32" electrodes and eliminate unintended rotation during use
- Fully sealed activation buttons ensure no fluid ingress
- Slip-proof handling provides surgeons with utmost control and minimizes unintended activation

Versatility

CONMED's wide variety of pencils appeal to a wide variety of surgeon preferences. Electrosurgical pencils are available in Push Button, Rocker Switch, and Foot-Control configurations, as well as with stainless steel and UltraClean electrode options.




SYSTEM 5000™
ELECTROSURGICAL UNIT



SUREFIT™
DISPERSIVE ELECTRODES



ULTRACLEAN®
COATED ELECTRODES



To learn more about these and other innovative products, call 800-448-6506 or visit CONMED.com



525 French Road, Utica, NY 13502

GOLDLINE®
ELECTROSURGICAL
PENCILS

Ordering Information

Description	Catalog Number
Hand Switch Pencils, Push Button	
UltraClean® Blade, Holster, 10' Cable, 40/case	131309A
UltraClean® Blade, Holster, 15' Cable, 40/case	131326A
Standard Blade, Holster, 10' Cable, 40/case	130309A
Standard Blade, Holster, 15' Cable, 40/case	130326A
Standard Blade, Holster, Tip Cleaner, 10' Cable, 40/case	130302A
Needle Electrode, Holster, 10' Cable, 40/case	130302A
Extended Insulated Blade, Holster, 10' Cable, 40/case	138114A
Hand Switch Pencils, Rocker Switch	
UltraClean® Blade, Holster, 10' Cable, 40/case	131307A
UltraClean® Blade, Holster, 15' Cable, 40/case	131327A
Standard Blade, Holster, 10' Cable, 40/case	130307A
Standard Blade, Holster, 15' Cable, 40/case	130327A
Standard Blade, Holster, Tip Cleaner, 10' Cable, 40/case	130304A
Needle Electrode, Holster, 10' Cable, 40/case	130301A
Foot Control Pencils	
UltraClean® Blade, Holster, 10' Cable, 40/case	138650A
Needle Electrode, Holster, 10' Cable, 40/case	138652A
Blade Electrode, Holster, 10' Cable, 40/case	138659A
Tip Cleaner	
Tip Cleaner, Single Use, Sterile, 40/case	138029

Customer Service: 1-800-448-6506
International Sales: 1-315-797-8375
Fax: 1-800-438-3051

CONMED.com
info@mail.CONMED.com