




Reflex® ELC530

Disposable Laparoscopic Clip Applier



AIRSEAL® iFS
INSUFFLATION
MANAGEMENT SYSTEM



**UNIVERSAL
PLUS™**
SUCTION/IRRIGATION



**MICROLAP®
INSTRUMENTS**



To learn more about these and
other innovative products,
call 800-448-6506 or
visit CONMED.com



525 French Road, Utica, NY 13502


REFLEX® ELC530

Maximum Visibility

By comparison, the Reflex® ELC jaw design has a 20% greater angle allowing for maximum visibility for precise clip placement.

Maximum Control

The Reflex® ELC has an automatic, direct drive clip feeding system with an audible feedback mechanism indicating complete clip closure. Our superior design provides maximum control, feel and ease-of-use.



Tissue Lock™ Clip

The Tissue Lock™ Clip is designed with opposing male/female interlocking features. This design provides superior on-vessel clip retention and improved hemostasis.

Ordering Information

Description

Qty

Catalog Number

Reflex® ELC530, 10mm, 20 M/L Clips

1pkg, 6/cs

530

Customer Service: 1-800-448-6506
International Sales: 1-315-797-8375
Fax: 1-800-438-3051

CONMED.com
info@mail.CONMED.com