




Ask for the 'purple one'™

NEW

Anchor* Tissue Retrieval System™ by CONMED



AIRSEAL® IFS
INSUFFLATION MANAGEMENT
SYSTEM



**MICROLAP®
INSTRUMENTS**
LESS INVASIVE APPROACH



VCARE®
UTERINE MANIPULATOR



To learn more about these and
other innovative products,
call 800-448-6506 or
visit CONMED.com

Anchor* Tissue Retrieval System™ by CONMED

SIMPLE

- Intuitive design provides ease of use
- Opens automatically when deployed

INNOVATIVE

- Unique deployment mechanism offers safety & security
- Proprietary bag material to reduce breaks & leaks

STRONG

- Rip-stop nylon for exceptional strength & quality
- Provides strong & secure seal for safe tissue retrieval

TRS-ROBO-8
125 ml

TRS100SB2
235 ml


TRS-ROBO-12
300 ml

TRS-VATS-15
1550 ml

TRS175SB2
1550 ml

TRS190SB2
1850 ml

TRS-TV-25
4000 ml





525 French Road, Utica, NY 13502

Ordering Information

Anchor* Tissue Retrieval System™

Introducer/ Trocar Size mm	Mouth Diameter cm / in	Bag Length cm / in	Approx Capacity ml	Qty per Box	Catalog Number
8	5.1 / 2.0	14.0 / 5.5	125	5	TRS-ROB0-8
10	5.6 / 2.2	14.0 / 5.5	235	5	TRS100SB2
12	7.6 / 3.0	14.0 / 5.5	300	5	TRS-ROB0-12
15	14.1 / 5.5	25.4 / 9.9	1550	3	TRS-VATS-15
15	14.1 / 5.5	25.4 / 9.9	1550	3	TRS175SB2
15	15.2 / 6.0	26.5 / 10.4	1850	3	TRS190SB2
25	16.5 / 6.5	38.1 / 15.0	4000	3	TRS-TV-25

Customer Service: 1-800-448-6506 CONMED.com
International Sales: 1-315-797-8375 info@mail.CONMED.com
Fax: 1-800-438-3051

*Anchor is a trademark of the Anchor Products Company, Addison, Illinois.