

PenAdapt® Surgical Smoke Evacuation Adapter

PenAdapt® is an ergonomically designed smoke evacuation adapter for most electrocautery pencils that captures surgical smoke at the point of creation.



PenAdapt® shown with ESU pencil

PenAdapt® shown with ESU (reusable) pencil

SYSTEM 5000®
GENERATOR



GOLDVAC® SLIM
PENCILS



VIROSAFE® VISICLEAR®
FILTER



To learn more about these and other innovative products, call **800-448-6506** or visit **CONMED.com**.

Features & Benefits

Universal Design

- Pliable silicone material allows for proper fit with most electro-surgical pencil brands

Angled Back End Clasp

- Reduces tension placed on surgeon's wrist providing comfort

Efficient Workflow

- Eliminates need for OR personnel to hold tubing

Low Profile Front Nozzle

- Improves visibility of the surgical site



To learn more about this and other innovative products, call **1-800-448-6506** or visit CONMED.com.

Ordering Information

NEW

Description	Qty	Catalog Number
PenAdapt® Ideal with wall suction applications. Includes 3/8 in (9.5 mm) x 1 ft (30.5 cm) tubing with an adapter to 3/8 in (9.5 mm) tube.	25/Case	PA1025
PenAdapt® (Cut-Back) Ideal with wall suction applications. Cut-back Adapter provides increased visibility. Include 3/8 in (9.5 mm) x 1 ft (30.5 cm) tubing with an adapter to 3/8 in (9.5 mm) tube.	25/Case	PA1025CB
PenAdapt® Ideal with wall surgical smoke evacuator applications. Includes visibility. Includes 3/8 in (9.5 mm) x 4 ft (1.22 m) tubing with an adapter to 7/8 in (22 mm) x 6 ft (1.8 m) tube	10/Case	PA2010
PenAdapt® Ideal with wall surgical smoke evacuator applications. Cut-back Adapter provides increased visibility. Includes 3/8 in (9.5 mm) x 4 ft (1.22 m) tubing with an adapter to 7/8 in (22 mm) x 6 ft (1.8 m) tube	10/Case	PA2010CB
PenAdapt® Ideal for attachment to reusable electro-surgical pencil applications. Includes 3/8 in (9.5 mm) x 4 ft (1.22 m) tubing with an adapter to 7/8 in (22 mm) x 6 ft (1.8 m) tube	10/Case	PA2010B
PenAdapt® Ideal for attachment to reusable electro-surgical pencil applications. Includes 3/8 in (9.5 mm) x 4 ft (1.22 m) tubing with an adapter to 7/8 in (22 mm) x 6 ft (1.8 m) tube	10/Case	PA1025BCB

* ESU pencils sold separately

CONMED Corporation
525 French Road, Utica, NY 13502

Customer Service: 1-800-448-6506
International Sales: 1-315-797-8375
Fax: 1-800-438-3051

CONMED.com
info@mail.CONMED.com