



Traction devices

for Getinge OR Tables

This document is intended to provide information to an international audience outside of the US.

GETINGE 



Traction devices

– the ideal addition

Essential accessories to easily transform universal tables into trauma and orthopaedic tables.

Getinge is built on a genuine compassion for people's health, safety, and well-being. Founded in 1904 with roots dating back to 1838, Getinge has grown organically and through acquisitions to become a global market leader. Our portfolio offers solutions and support throughout the clinical pathway, and features well-known and dependable product brands, including

Maquet OR Table Systems that have been trusted for more than 50 years.

Good access and optimal patient positioning are crucial for the outcome of procedures. With the traction devices for Getinge OR Tables, surgical teams can easily position the table in the ideal configuration for both surgeon and patient.



Traction device 1419.01

Tried and tested ergonomic functionality

Traction devices enable easy treatment of extremities under tension. The functional principle of the bar and insertion system has been known for its simple handling for many years. The clip-on device is used in conjunction with a universal OR table.

The simple and safe mounting of the extension device is carried out with the help of the practical trolley. This provides safe storage for all accessories and is part of the scope of delivery.

Including basic accessories

- 1 pair of adapters
- 1 trolley with basket
- 1 telescopic bar, short, with foot-end side rail
- 1 telescopic bar, long, with foot-end side rail
- 1 screw tension device
- 1 foot plate mount
- 2 support bars
- 2 radial setting clamps
- 1 countertraction post, including padding roll
- 1 side rail extension

All benefits at a glance

- Excellent usage options for the C-arm during interventions
- Rust-free and disinfectant-proof version in CrNi steel
- Simple mounting of accessories
- Simple and safe mounting to the OR table with the help of the trolley
- Trolley with basket
- The positioning of the patient may be carried out by a single person

Versions

- **1419.01B0** suitable for
 - Maquet Alphastar 1132 and Maquet Betastar 1131.12
 - Meera 7200.01, Meera ST 7100.01 and Meera CL 7000.01
- **1419.01C0** suitable for
 - Maquet Betaclassic 1118.01 and Maquet Alphaclassic Pro 1118.06
 - Maquet Lyra 1118.10
- **1419.01HC** suitable for Maquet Alphamaxx 1133
- **1419.01JC** suitable for
 - Maquet Alphastar Pro 1132.21/22 and Maquet Betastar 1131.12 on the "normal back" mounting point
 - Meera 7200.01, Meera ST 7100.01 and Meera CL 7000.01



Trolley with basket and accessories



Mounting to the leg plate mounting point

Positioning examples

Treatment of femur and tibia fractures



Femoral neck treatment in the supine position



Femoral treatment in the lateral position



Tibia treatment with countertraction post for tibia and fibula and joint supporting arm




Femur treatment with traction stirrup clamp with rotation

Additional traction accessories

Overview

Accessories for traction devices



1001.65A0	Leg holder with ball joint adjustment, with securing strap, SFC padding
1003.23C0	Radial setting clamp to accommodate 1001.65A0
1003.34A0	Rotating and tilting clamp to accommodate foot plates
1003.35A0	Traction stirrup clamp for traction with wire
1003.75A0* 1003.76A0* 1003.77A0* 1003.78A0*	Traction boot Disposable inlay incl. mounting brackets Re-usable inlay Disposable mounting brackets
1001.97A0	Foot plate , adjustable width, with Velcro® strap and heel pad
1150.64B0 1150.64E0	Pair of leg plates , with SFC padding Pair of leg plates , with SFC padding suitable for Maquet Yuno II and Maquet Otesus
1003.50C0	Countertraction post for tibia and fibula , height adjustable and can be swiveled, including PUR padding roll
1004.93B0	Condyle fixation , for 1003.50, individually adjustable
1007.51A0/B0	Vertical adjustment , long and short
1007.43A0	Screw tension device
1003.59V0	Telescopic bar (steel), reinforced
1007.45A0	Support plate
1007.48A0	Frame for sterile drape
1007.49A0	Traction frame for prone and side position
1004.89B0	Perineal support
1004.8800	Side rail clamp
1003.54A0	Bar elongation
1003.55A0	Traction bar , carbon fiber version
1003.49A0/B0	Foot plate mount

* Legal manufacturer
Condor MedTec GmbH

Femur treatment



1004.86B0

Universal support, can be swiveled, height can be adjusted, leg rest on ball joint, for positioning the calf for transcondyle wire traction

1004.90B0

Countertraction post for locking nailing with the patient in the supine position, with PUR padding

1004.85B0

Femur countertraction post with support plate for lateral position, height adjustable, incl. padding rolls and seat plate, with PUR padding

Storage



1007.15A0

Accessories trolley

1007.47A0

Accessories trolley

1160.57A0

Trolley for 1160.59XC and 1433.41XC

Attachments for Maquet Yuno II and Otesus

Specifically for Yuno II and Otesus an own table top has been developed. This table top enables the use of existing traction accessories.




Maquet Yuno II

- **1433.41AC** Adapter plate for traction without additional accessories
- **1433.41BC** Adapter plate for traction including traction accessories

1433.41XC Adapter plate for traction

- Suitable for Maquet Yuno II
- Can be mounted on top interface
- Quick conversion of the OR table
- Good C-arm access
- For use of existing traction accessories



Maquet Otesus

- **1160.59AC** Support plate for extension without additional accessories
- **1160.59BC** Support plate for extension including traction accessories

1160.59XC Support plate for extension

- Suitable for the table tops 1160
- Can be mounted on top interface
- Quick conversion to an extension table
- Good C-arm access



With a firm belief that every person and community should have access to the best possible care, Getinge provides hospitals and life science institutions with products and solutions aiming to improve clinical results and optimize workflows. The offering includes products and solutions for intensive care, cardiovascular procedures, operating rooms, sterile reprocessing and life science. Getinge employs over 10,000 people worldwide and the products are sold in more than 135 countries.

This document is intended to provide information to an international audience outside of the US.

Maquet GmbH · Kehler Str. 31 · 76437 Rastatt · Germany · +49 7222 932-0

www.getinge.com