

TRIMANO FORTIS

Improving ergonomics

GETINGE GROUP



For years, the TRIMANO 3D support arm has offered a helping hand in the operating theatre. Maquet continues this success story with the TRIMANO FORTIS.

TRIMANO FORTIS brings strength, flexibility and ergonomic advantages to the operating theatre. It is a further development of the TRIMANO 3D support arm. TRIMANO FORTIS is more robust than its predecessor, with significant improvements to load capacity¹ and handling. TRIMANO FORTIS builds upon the popular features of the widely used TRIMANO 3D support arm. It also offers ergonomic benefits for clinicians. The force required to release the handle has been reduced, making repositioning easier.

Getinge Group. Passion for life.

Benefits at a glance:

- **Easy to use:** TRIMANO FORTIS requires no power supply or additional connection
- **Work relief:** thanks to its versatile applications the TRIMANO improves working conditions the surgical team
- **Extremely movable:** the “Click and Move” system enables free movement in all directions while securing a stable fixation when the handle is released
- **Multidisciplinary:** due to a large range of accessories, the TRIMANO is usable in many surgical disciplines
- **Compatibility:** can be used with common OR tables with a side rail
- **Maximize your investment:** backwards compatible with predecessor TRIMANO (with exceptions), meaning existing accessories can be used

¹ max. proportional weight of patient arm of 250 kg

An extra-hand in the operating room

Holding and lifting body extremities during orthopedic interventions is a huge effort for surgical assistants. This task is performed by TRIMANO FORTIS easily and with more stability.

The TRIMANO FORTIS can be connected directly to the OR table accessory rail and is ready for immediate use. No need for any other connections. It can be moved in any direction and securely fixed at any position.

TRIMANO's features – such as multidisciplinary use, flexible positioning options, increased load capacity, and comfortable handling – make it an essential tool in the operating theatre. And as TRIMANO FORTIS is mostly compatible with the previous version, you can continue using your existing equipment.

Ordering numbers and compatible accessories

1002.30A0	TRIMANO FORTIS
1002.31A0	Adapter Fortis
1002.32A0	Armboard Fortis, SFC
1002.43A0	Arm rest, lateral position
1002.43B0	Arm rest, beach chair
1002.52A0	Arm rest, disposable, sterile
1002.45A0	Counter traction post
1002.48A0	Finger fixation
1002.47A0	Wound retractor support



Shoulder treatment in beach chair position



Wrist treatment in supine position



Plastic surgery in lateral position



Elbow treatment in prone position

GETINGE GROUP

MAQUET GmbH
Kehler Str. 31, 76437 Rastatt, Germany

MAQUET
GETINGE GROUP

ARJOHUNTLEIGH
GETINGE GROUP

GETINGE
GETINGE GROUP

www.getingegroup.com