



More than just comfort

Patient warming is about more than providing comfort. The clinical benefits of intra-operative and post-operative patient warming – to prevent hypothermia – have been well documented:

- Reduces surgical site infections by as much as 67%
- Decreases blood loss
- Speeds wound healing and reduces pain
- Reduces incidence of myocardial ischemia and post-operative angina
- Increases cardiac output
- Shortens hospital stays
- Improves patient comfort and satisfaction

Getinge's flexible warming technologies help you maintain patient normothermia from admission to discharge.

Available in a range of sizes and types, Getinge has a warming solution to fit most any budget.



GET250L*

*One basket is standard



GET250 over GET350L

GETINGE *

Getinge is a registered trademark of Getinge AB, its subsidiaries, or affiliates in the United States or other countries · Getinge blanket and fluid warming solutions are manufactured for Getinge by Enthermics · Copyright 2019 Getinge AB or its subsidiaries or affiliates · All rights reserved.

Sales Office · Getinge
1777 East Henrietta Road
Rochester, NY 14623
1-800-475-9040

www.getinge.com

MCV00089748 REVA
EXP 03/21



Getinge blanket and fluid warming cabinets

A better brand of warmth

GETINGE *

Blanket warmers

GET250

- Capacity: 4-6 blankets
- Exterior (H x W x D): 22" (559 mm) x 18.5" (470 mm) x 23.5" (597 mm)
- Interior (H x W x D): 15.25" (387 mm) x 15.5" (394 mm) x 18.75" (476 mm)
- Heavy-duty caster and mobile cart options
- Stacks with GET250, GET250L, GET350 and GET350L warmers
- Net weight: 67 lbs (30 kg)



GET350

- Capacity: 10-12 blankets
- Exterior (H x W x D): 28" (711 mm) x 18.5" (470 mm) x 23.5" (597 mm)
- Interior (H x W x D): 21.25" (540 mm) x 15.5" (394 mm) x 18.75" (476 mm)
- Heavy-duty caster and mobile cart options
- Stacks with GET250, GET250L, GET350 and GET350L warmers
- Net weight: 91 lbs (41 kg)



GET750

- Capacity: 20-25 blankets
- Exterior (H x W x D): 35.5" (902 mm) x 24" (610 mm) x 25.75" (654 mm)
- Interior (H x W x D): 28.25" (718 mm) x 21" (533 mm) x 21.75" (552 mm)
- Optional pass-through door makes it easy to stock and access blankets from front and back
- Stacks with GET750 warmers (not available in CA)
- Net weight: 134 lbs (61 kg)



GET2060

- Capacity: 60-70 blankets
- Exterior (H x W x D): 71" (1803 mm) x 32" (813 mm) x 29.25" (743 mm)
- Interior (H x W x D): 56.75" (1441 mm) x 26.5" (673 mm) x 23.75" (603 mm)
- Stand-alone, floor model
- Heavy-duty casters are standard
- Net weight: 347 lbs (157 kg)



Fluid warmers

GET250L

- Capacity: 12 one-liter bottles or 16 one-liter bags per basket*
- Exterior (H x W x D): 22" (559 mm) x 18.5" (470 mm) x 23.5" (597 mm)
- Interior (H x W x D): 15.25" (387 mm) x 15.5" (394 mm) x 18.75" (476 mm)
- Heavy-duty caster, mobile cart, and glide-out basket options
- Net weight: 67 lbs (30 kg)



*One basket is standard

GET350L

- Capacity: 24 one-liter bottles or 24 one-liter bags
- Exterior (H x W x D): 28" (711 mm) x 18.5" (470 mm) x 23.5" (597 mm)
- Interior (H x W x D): 21.25" (540 mm) x 15.5" (394 mm) x 18.75" (476 mm)
- Heavy-duty caster and mobile cart options
- Stacks with GET250, GET250L, GET350 and GET350L warmers
- Net weight: 91 lbs (41 kg)



Combination warmers

GET1850BL

- Blanket chamber capacity: 30-40 blankets
- Fluid chamber capacity: 30 one-liter bottles or 24 one-liter bags
- Exterior (H x W x D): 75" (1905 mm) x 30" (762 mm) x 28.5" (724 mm)
- Interior (H x W x D): 51" (1295 mm) x 26.5" (673 mm) x 23.75" (603 mm)
- Net weight: 457 lbs (207 kg)



Blankets, irrigation fluids and injection fluids need to be warmed to different temperatures

Match the right temperature to the right product in the right place at the right time.

Getinge offers a variety of fluid and blanket warmers to stack to create the warming configuration to fit your needs.