

Getinge Assured

Record Labels for Steam, EO and VHP

The Getinge Assured Record Labels are double self-adhesive record labels with process indicators for Steam, Ethylene Oxide and Plasma sterilization.

The labels can be printed with the Getinge Assured Label Gun which features three printing lines. Following sterilization, the user can remove the labels and place on the patient's medical record.



Performance Characteristics

The double adhesive label allows you to print:

- Code number of the operator
- Number or code of the autoclave and cycle
- Lot number
- Sterilization date
- Expiration date

Benefits

- Easy printing of identical labels with minimal time and effort
- Eliminates mistakes in batch numbers
- Easy interpretation of the results due to precise color change
- No need for additional process indicators

Environmental Limitations

Storage: 5° - 30°C (41° - 86°F) and 20-70% Relative Humidity.

Shelf Life

24 months. The expiry date is printed on the product packaging.

Standards Classification

EN ISO 11140-1, Type 1

Ordering Information

Getinge Assured	Order Number
Record Label Steam	6005500408
Record Label EO	6005500409
Record Label VHP	6005500410
Record Label Gun	6005500411