

Getinge Assured

Protein Test Instrument Lumen

The Getinge Assured Protein Test Lumen is a quick and easy test used to detect residual protein in hollow lumen surgical instruments.



Performance Characteristics

The protein test can detect protein residuals to 1 µg sensitivity within 10 seconds. When protein is present, the solution will turn from light brown to blue. This test provides a semiquantitative measure of surface cleanliness; the higher the level of contamination of protein (bioburden), the darker and faster the color changes to blue.

Physical Characteristics

The Getinge Assured Protein Test Lumen contains 25 vials of protein detection solution and 25 swabs.

Environmental Limitations

Keep away from light and heat exposure.

Shelf Life

24 months when refrigerated and 6 months when not refrigerated.

Standards Classification

Consistent with ISO 15883, routine monitoring of washer disinfectors.

Ordering Information

Getinge Assured

Protein Test Instrument Lumen

Order Number

503911200