

# Getinge Assured

## Multi-Incubator

The Getinge Assured Multi-Incubator is designed for incubation of the Getinge Assured EZTest® Self-Contained Biological Indicators (SCBI's), 37°C (Bacillus Atrophaeus), 57°C or 60°C (Geobacillus Stearothermophilus).

The Incubator features an integrated thermometer for easy temperature reference. In addition to the 13 incubation wells, it has activation well for crushing the Getinge Assured SCBI's.



### Performance Characteristics

#### Operating Range:

- **Incubator:** 37°C / 57°C / 60°C.
- **Thermometer:** 37°-60°C.

#### Electrical Rating:

- **Incubator:** 50-60 Hz. 12V.
- **Power supply:** 120-240V 50-60Hz.

### Physical Characteristics

The incubator has 13 wells for the placement of the Getinge Assured EZTest SCBI's. Three temperature settings: 37°C for EO Gas, 57°C for steam and 60°C for Getinge Assured Smart-Read SCBI's (10 hour) for steam.

### Material

Incubator and Quick Start Guide

### Incubator Dimensions

- **Incubator:** W 9cm x H 5.08cm x D 10.2cm
- **Well Dimensions:** Diameter 94cm x Depth 2.54cm

### Environmental Limitations

Storage 10°-38°C and 20-80% Relative Humidity.

### Ordering Information

**Getinge Assured**  
Multi-Incubator

**Order Number**  
504054700