

WARMING CABINETS



BETTER “LOOK-IN” DESIGN

GETINGE Warming Cabinets allow you to check your heated blankets and solutions at a glance, without ever having to open a door. Innovative “look-in” glass doors are an industry first, allowing inspection of supplies without heat loss. And tempered glass makes them completely safe to touch.

Getinge USA, Inc.
1777 East Henrietta Road
Rochester, New York 14623-3133
U.S.A.
Phone: (800) 475-9040
Fax: (585) 272-5033
info@getingeusa.com
www.getingeusa.com

MEMBER OF THE GETINGE GROUP

ALL THE FEATURES YOU'RE LOOKING FOR, WITH INDUSTRY-LEADING GETINGE QUALITY

Designed for warm storage of the solutions and blankets used for patient care in surgery, recovery, OB/GYN, ICU and ER, these warming cabinets offer independently controlled, convection-heated chambers. They're available in a variety of sizes and mounting options to fit space requirements. And every Getinge product is backed by an industry-leading team of service professionals.

Clearly, these are the best warming cabinets you've ever looked into.

- **Individual Compartments:** Upper compartment for solution bottles and lower compartment for blankets.
- **Individual "Look-In" Glass Doors:** Innovative glass doors allow you to check inventory visually, without heat loss from opening the door. The tempered glass doors are completely safe to touch.
- **Set-Point Temperature "Lock-In":** Control temperatures in each compartment. Set-point temperatures can be locked in to prevent unauthorized changes.
- **Easy To Install:** Backed by Getinge service professionals.
- **90°F to 160°F Set-Point Range:** Choose the right temperature for your supplies.
- **Control Memory:** Temperature set-point values are retained in the event of a power loss.
- **Forced-Air Convection Heating:** Warms evenly, eliminating hot spots and scorched blankets.
- **Overheat Protection With Alarm:** If the temperature exceeds 10°F above the set-point, an alarm alerts the operator and the heaters are turned off.



GETINGE 5618 Warming Cabinet

This unique, 18-inch deep, recessed cabinet allows in-line installation with recessed stainless steel OR cabinets. Your warm blankets and solutions are conveniently available in the operating room.



GETINGE 5524 Warming Cabinet

When space is at a premium, look to the compact size of the 5524. This model can sit on a table top, hang on a wall, or be recessed into a wall.



GETINGE 5624 Warming Cabinet

A full 24-inches deep, this model provides maximum storage capacity. Available as a freestanding cabinet with a pedestal base, recessed (with or without pedestal base) or wall mounted.