

Getinge Assured

Wash Monitor Ultrasonic

The Getinge Assured Ultrasonic Wash Monitor is a single use synthetic soil test for routine monitoring of an ultrasonic wash process.

The test is inserted in a Wash Monitor Holder and then placed in the ultrasonic washer.



Performance Characteristics

The synthetic test soil (blue diamond) simulates human blood and tissue. The removal of all test soil indicates effective cleaning.

Physical Characteristics

The Wash Monitor Ultrasonic consists of a metal strip with synthetic test soil.

- **Dimensions:** 89 mm x 25 mm (3.5" X 1")
- **Substrate Material:** Metal

Environmental Limitations

Storage: 4° - 30°C (40° - 86°F) and 10-60% Relative Humidity. Keep away from disinfectants and chemical fumes.

Shelf Life

24 months. The expiry date is printed on the product packaging.

Standards Classification

Consistent with ISO 15883, routine monitoring of washer disinfectors.

Ordering Information

Getinge Assured	Order Number
Wash Monitor Ultrasonic 50 pcs	503865200
Wash Monitor Holder	503864800