

Getinge Assured

Plasma Indicator

The Getinge Assured Plasma Indicator gives visible confirmation that the goods have been exposed to hydrogen peroxide.



Performance Characteristics

When exposed to H₂O₂ sterilization process, the Plasma Indicator will change from blue to a shade of pink.

Physical Characteristics

- **Dimensions:** L 101mm x W 19mm
- **Packaging:** 250 indicators in a zip lock moisture barrier bag

Environmental Limitations

Storage: In provided sealed packaging at room temperature 10°-38°C and 10-60% Relative Humidity away from oxidizing agents.

Shelf Life

24 months. The expiry date is printed on the product packaging.

Ordering Information

Getinge Assured

Plasma Indicator

Order Number

504052000