

Customer: [REDACTED]

Reference: [REDACTED]

TAPE DISPENSER

PRODUCT DESCRIPTION


Dispenser for sterilizer tape that holds two (2) reels of tape. The ability to hold two (2) rolls of tape gives the packaging operator a flexible support. The dispenser is designed as a free-standing table model with a built-in cutting knife for easy one-hand operation of indicator tape. A heavy-duty bottom plate, fitted with anti-slip rubber, prevents the dispenser from slipping when tape is torn off. The dispenser is made of high quality coated steel that assures a long lifetime.

APPLICATION

Specially designed for use in reprocessing units such as a Central Sterile Processing/CSSD and point-of-use Sterile Processing/TSSU. The tape dispenser is equally suitable for use in smaller medical clinics and dental clinics.

KEY FEATURES

- Practical
- Solid
- Durable



PRINCIPLE OF OPERATION

Operation is simple. Depending on the size of the goods to be packed, tear off tape either from one or both of the rolls.

TESTING

All tape dispensers have been tested according to Getinge demands.

TECHNICAL DATA

Tick your selection.

Art. No.	Product name	Size (L x W x H), mm
561 57 07-01 <input type="checkbox"/>	Tape dispenser	260 x 60 x 120