


# PlumePen<sup>®</sup> Pro

## Surgical Smoke Evacuation Pencil

The PlumePen<sup>®</sup> Pro features an easy to remove x-ray detectable capture port that allows for effortless electrode replacement to accommodate surgeries where longer electrodes are desired. PlumePen<sup>®</sup> Pro is smaller than competitive surgical smoke evacuation pencils and similar in size to standard electro-surgical pencils.



**SYSTEM 5000<sup>™</sup>**  
GENERATOR



**PLUMEPEN<sup>®</sup> ELITE**  
PENCILS



**PLUMEPEN<sup>®</sup> ULTRA**  
PENCILS



**VISICLEAR<sup>®</sup>**  
SURGICAL SMOKE EVACUATOR



To learn more about these and other innovative products, call **866-4CONMED** or visit **[www.CONMED.com](http://www.CONMED.com)**.

## Features & Benefits

### Sleek, Slim Profile

- Smaller than competitive surgical smoke evacuation pencils; similar in size to standard ESU pencils

### Removable Smoke Capture Port

- Allows for easy electrode and capture port replacement

### Compatible with Standard Electrodes

- Designed for surgeries where longer electrodes are desired

### Full Blade Tip Reveal

- Provides an unobstructed view of surgical site

### 360 Degree Swivel

- Allows for fluid hand movement

### Textured Profile

- Improves surgeon's grip



To learn more about this and other innovative products, call [866-4CONMED](tel:866-4CONMED) or visit [www.CONMED.com](http://www.CONMED.com).

## Ordering Information

Description	Qty	Catalog Number
PlumePen® Pro (Non-Stick ESU Blade) 10 ft (3 m) tubing	20/Case	PLPR04020
PlumePen® Pro (Non-Stick ESU Blade) 15 ft (4.6 m) tubing	20/Case	PLPR04520
Capture Ports (X-ray detectable) Accommodates 4 in (102 mm) electrodes	20/Case	PROTIP420
Capture Ports (X-ray detectable) Accommodates 5 in (127 mm) electrodes	20/Case	PROTIP520
Capture Ports (X-ray detectable) Accommodates 6 in (152 mm) electrodes	20/Case	PROTIP620

CONMED Corporation  
525 French Road, Utica, NY 13502

866-4CONMED  
International: 727-214-3000

[www.CONMED.com](http://www.CONMED.com)  
[Customerexperience@CONMED.com](mailto:Customerexperience@CONMED.com)