




Looking Glass™ Cable Connector Interface (CCI)

Designed for flexibility and convenience during both installation and clinical use, the CCI replaces traditional bulky camera boxes with a compact console where the Looking Glass™ camera head and light cables connect.



Features and Benefits

Easy Connection

Easily connect and disconnect camera heads and light cables.


Integration

Integrated USB port connects to Looking Glass image capture to simplify transferring images and video to external memory devices.

LED Indicator Lights

LED indicators show when the CCI is powered (green) and when the camera head and light cables are plugged in and initialized (blue).

To learn more about the Looking Glass™ System,
visit www.CONMED.com



Ordering Information

Description	Catalog Number
Looking Glass™ Cable Connector Interface (CCI)	AV41CON

CONMED Corporation
11311 Concept Blvd.,
Largo, FL 33773

Toll Free: 1-866-4CONMED
International: 727-214-3000
www.CONMED.com

customerexperience@conmed.com
internationalorders@conmed.com