



Looking Glass™ Camera Heads

Customizable, autoclavable, and adjustable from the sterile field—CONMED's Looking Glass™ camera heads are designed with innovative technology for maximum performance and efficiency.



Features and Benefits

Crisp 4K Technology

4K delivers precise, detailed imaging for multi-specialty visualization.

Innovative 3-CMOS Color Imaging

3-CMOS offers a dedicated chip for each color spectrum (red, green, blue), while a single CMOS chip works to blend all colors.

Customizable 5-Button Design

Unique 5-button design allows you to customize the system to your preferences so you can easily make adjustments from the OR sterile field.

Unique ShockFlex™ Durability

Resilient ShockFlex™ prism mounting prevents damage from accidental drops.

Autoclavable to Keep Costs Down

Autoclavable camera heads eliminate the need for expensive sterilization methods. This improves efficiency while maintaining maximum performance.

7-Year Warranty* Against Autoclave Failure


Provides peace of mind that CONMED stands behind its products.

Learn how you can invest in best-in-class Looking Glass™ innovation and efficiency by scheduling a no-obligation evaluation with your local CONMED Sales Representative at **1-866-4CONMED** or **www.CONMED.com**

*Warranty information can be found at www.conmed.com/LookingGlassWarranty



11311 Concept Blvd., Largo, FL 33773



Ordering Information

Description	Catalog Number
Looking Glass™ Camera Head - 4K 3MOS Eyecup	AV41CHE

CONMED Corporation
11311 Concept Blvd.
Largo, FL 33773

Toll Free: 1-866-4CONMED
International: 727-214-3000
www.CONMED.com

customerexperience@conmed.com
internationalorders@conmed.com