

ProTwin

Professional spray system – Trolley wash unit

APPLICATION

Designed for reprocessing Hospital S/S equipment/furniture e.g. Stainless steel trolleys, close transport cabinets etc.

Use in Central Sterile Processing/CSSD and the OR/point-of-use Sterile Processing/TSSU located within a central surgery department.

The trolley wash unit works equally well at healthcare centers or smaller clinics, where an effective and hygienic trolley wash unit is required.

Time-saving, hygienic manual trolley washer with wall-mounted spray-gun unit, including holder for detergents.

ProTwin is a cleaning and sanitizing station that automatically mixes and dispenses concentrated chemicals. With a clean, professional appearance it comes with a built-in hose storage.



Performance

- ✚ Pressure: Min 1 bar Max6 bar (90 psi)
- ✚ Flow-rate (Rinse) 2.7 gpm (Prod) 3.5 gpm
- ✚ ProTwin foam kit
- ✚ Utilizes the ProTwin's single selector for detergent and sanitizing functions
- ✚ ProTwin's venturi maintain a proper and consistent chemical action/consumption
- ✚ Compressed medical air facilitate drying of items (to be installed)
- ✚ Multi-layer hose, improves flexibility and chemical resistant,

Adjustable spray pattern & thick foam all in one



Features

- ✚ Adjustable spray pattern
- ✚ Separate selectable chemical position eliminate risk of chemical crossover
- ✚ Adjustable accessories – easy to use
- ✚ Venturi ensure accurate dosing, easy to clean and be replaced
- ✚ Intuitive ergonomic design for user friendly operation
- ✚ 3 functions spray/foam and rinse up to 60°C (140°F) tolerance
- ✚ Multi-layer hose, improves flexibility and chemical resistant