



Medhold Redefining Healthcare Through
Technology That Keeps Moving With You.

CLINICAL ACCESSORIES PORTFOLIO



**CAPITAL EQUIPMENT
SUPPLIES & ACCESSORIES**

INDEX

GE Healthcare offers you a simple, reliable one-stop solution for high-quality GE Healthcare-compatible clinical accessories. We provide an expanded inventory of products that have been verified by our engineers.

[CLICK THE PAGE YOU REQUIRE](#)

ECG 1

SpO2

TECHNOLOGY GE TRUSIGNAL 5

TECHNOLOGY NELLCOR™ 7

TECHNOLOGY MASIMO™ 9

NIBP 10

TEMPERATURE PROBES & CABLES 14

INVASIVE BLOOD PRESSURE 16

CARESCAPE RESPIRATORY

MODULES AND COMPACT

AIRWAY GAS MODULE

GAS MONITORING 17

PATIENT SPIROMETRY ACCESSORIES 18

SPIROMETRY TUBES & KITS 19

RESPIRATORY MONITORING

SAMPLE LINE 20

GAS MONITORING

AIRWAY ADAPTERS 21

GAS MONITORING ACCESSORIES 22

MULTI ABSORBERS 24

ENTROPY 25

NEUROMUSCULAR TRANSMISSION - NMT 26

EEG NEUROLOGY 27

CARDIAC OUTPUT 28

PATIENT CIRCUITS 31

PATIENT BACTERIAL & VIRAL FILTERS 33

DEFIBRILLATOR SUPPLIES AND

ACCESSORIES REANIBEX 34

DIAGNOSTIC CARDIOLOGY

AND STRESS ECG SYSTEMS

MAC™ CASE™ AND CARDIOSOFT™ 42

EASYONE CS AND TRUEFLOW

SPIROMETRY 45

SEER® HOLTER MONITORS 46

TONOPORT VI 49

TANGO M2 49

MATERNAL INFANT CARE

GIRAFFE* OMNIBED*, GIRAFFE

INCUBATOR, GIRAFFE WARMER

AND PANDA* WARMER 50

MATERNAL FETAL CARE

COROMETRICS™ & MINI TELEMETRY 51

INTERNAL MONITORING 52

INFANT CARE PHOTOTHERAPY

ACCESSORIES RESUSCITATION

RESUSCITATION 53

LULLABY™ WARMER 53

VENTILATION

CARESCAPE R860 VENTILATOR

FILTERS & PATIENT CIRCUITS 54

GAS MODULE ACCESSORIES 55

AEROGEN NEBULIZERS & ACCESSORIES 56

CABLES, EXHALATION VALVE HEATER

AND SENSORS 57



ECG


GE Healthcare ECG trunk cables and lead wires help you reliably make the correct connection to GE monitoring devices. These high-quality reusable accessories complement your monitoring equipment, enabling accurate measurements.

ECG Accessories with IEC standard


PATIENT MONITORING ACCESSORIES, IEC FOR GE'S MONITORS

ECG CABLES & LEADWIRES

MULTI-LINK® RECTANGULAR CONNECTION - 3/5 LEAD


| Part number | Length | Description | Cable | Quantity |
|-------------|--------|---|--|----------|
| 2017004-003 | 3.6 m | Multi-Link 3-lead ECG care cable neonatal DIN |  | 1 |

ECG TRUNK CABLE ACCESSORIES

| Description | Part number | Quantity |
|---|-------------|-------------|
|  DIN Converter | 3-Lead | 2106403-003 |
| | 5-Lead | 2106403-001 |

ECG PRE-WIRED RADIOTRANSLUCENT ELECTRODES WITH DIN CONNECTOR


| Population | Conductive Medium | Cloth | Dimensions | Part Number | Quantity |
|------------|-------------------|-------|--------------------|------------------|---------------------|
| PICU, NICU | Adhesive Hydrogel | X | 22.225 mm/0.875 in | 2009101-406* IEC | 600 (200 bags of 3) |




PATIENT MONITORING ACCESSORIES, IEC FOR GE'S MONITORS

ECG CABLES & LEADWIRES



MULTI-LINK SERIES - RECTANGULAR CONNECTOR - 3/5-LEAD

| Part number | Length | Description | Trunk cable | Quantity |
|-------------|--------|------------------------------------|--|----------|
| 2017003-003 | 3.6 m | Multi-Link 3/5-lead ECG care cable |  | 1 |

ECG TRUNK CABLES - IEC



| | Description | Length | Part number | Quantity |
|--|-----------------------|-------------|-------------|----------|
|  | 3/5-lead w/ESU filter | 3.6 m/12 ft | 2106310-002 | 1 |

MULTI-LINK SERIES - RECTANGULAR CONNECTOR - 3/5-LEAD




| Part number | Length | Description | Leadwire | Quantity |
|--------------------------|-----------------|--|--|----------|
| 412682-003 412682-004 | 74 cm 130 cm | Multi-Link 3-leadwire set, grouped, grabber |  | 1 |
| 414556-003 414556-004 | 74 cm 130 cm | Multi-Link 5-leadwire set, individually replaceable, grabber |  | 1 |

PATIENT MONITORING ACCESSORIES, IEC FOR GE'S MONITORS



MULTI-LINK SERIES - RECTANGULAR CONNECTOR - 6-LEAD

| Part number | Length | Description | Trunk cable | Quantity |
|---------------------------------------|--------|-----------------------------------|--|----------|
| 2017005-003 | 3.6 m | Multi-Link 6-lead ECG care cable |  | 1 |
| Part number | Length | Description | Leadwire | Quantity |
| 421933-002 | 130 cm | Individually replaceable, grabber |  | 1 |
| Also accepts Multi-Link 3/5-leadwires | | | | |

MULTI-LINK SERIES - RECTANGULAR CONNECTOR - 12-LEAD






| Part number | Length | Description | Trunk cable | Quantity |
|--------------------------|-----------------|---|--|----------|
| 2017006-003 | 3.6 m | Multi-Link 12-lead ECG care cable |  | 1 |
| Part number | Length | Description | Leadwire | Quantity |
| 414556-003 414556-004 | 74 cm 130 cm | Multi-Link 5-leadwire set, individually, replaceable, grabber |  | 1 |
| 416467-004 416467-002 | 74 cm 130 cm | Multi-Link C-leadwire set, individually replaceable, grabber, C2-C6 |  | 1 |

CARESCAPE ONE - 12-LEAD

| Part number | Length | Description | Leadwire | Quantity |
|---------------|--------|---|--|----------|
| 2066468 - 008 | 130 cm | Extension from 6 Leadwires to 10, grabber IEC |  | 1 |
| 2106398 -002 | 130 cm | ECG Leadwire set, 6 - Lead, grouped, grabber, IEC |  | 1 |

PATIENT MONITORING ACCESSORIES, IEC FOR GE'S MONITORS

MULTI-LINK SERIES - RADIOLUCENT LEADWIRES - 5-LEAD, 12-LEAD

| Part number | Length | Description | Trunk cable | Quantity |
|-------------|--------|---|--|----------|
| 2017003-003 | 3.6 m | Multi-Link 3/5-lead ECG care cable |  | 1 |
| 2017006-003 | 3.6 m | Multi-Link 12-lead ECG care cable |  | 1 |
| Part number | Length | Description | Leadwire | Quantity |
| 403751-026 | 150 cm | Multi-Link radiolucent leadwire set, grabber, 10-lead |  | 1 |
| 403751-025 | 150 cm | Multi-Link radiolucent leadwire set, grabber, 5-lead |  | 1 |
| 403751-033 | 150 cm | Multi-Link radiolucent leadwire set, grabber, C2-C6 |  | 1 |






SpO² SENSORS


Through the legacy of Datex-Ohmeda, GE SpO² technology is a story of innovation that has been going on for more than 20 years. TruSignal™ sensors allow continuous non-invasive monitoring of SpO² and heart rate. We offer a wide range of SpO² sensors, designed and tested for their performance and always serving the comfort of your patients. For your Monitors with GE TruSignal Technology.

TECHNOLOGY GE TRUSIGNAL

REUSABLE INTEGRATED SENSORS






| | Sensors | Patient | Patient Weight | Length | Part Number | Quantity |
|--|------------------------------------|------------|----------------|--------|-------------|----------|
|  | TruSignal Integrated Finger sensor | Adult, | ≥ 20 kg | 2 m | TS-F2-GE | 1 |
| | | Paediatric | ≥ 20 kg | 4 m | TS-F4-GE | 1 |
|  | TruSignal Integrated FingerTip | Adult | > 30 kg | 4 m | TS-SA4-GE | 1 |
| | | Paediatric | 15 - 30 kg | 3 m | TS-SP3-GE | 1 |
|  | TruSignal Integrated Ear sensor | Adult, | ≥ 10 kg | 2 m | TS-E2-GE | 1 |
| | | Paediatric | ≥ 10 kg | 4 m | TS-E4-GE | 1 |

INTERCONNECT CABLES





| | Product | Length | Part Number | Quantity |
|---|--|--------|-------------|----------|
|  | TruSignal interconnect cable with GE connector | 3 m | TS-G3 | 1 |

TECHNOLOGY GE TRUSIGNAL

REUSABLE SENSORS (REQUIRES COMPATIBLE INTERCONNECT CABLE)

| | Sensors | Patient | Patient Weight | Length | Part Number | Quantity |
|---|---------------------------------------|---------------------------------------|---------------------|------------|--------------------|----------|
|  | TruSignal Finger sensor | Adult, Paediatric | ≥ 20 kg | 1 m | TS-F-D | 1 |
|  | TruSignal FingerTip sensor | Paediatric Adult | 15-30 kg > 30 kg | 1 m 1 m | TS-SP-D TS-SA-D | 1 |
|  | TruSignal Ear sensor | Adult, Paediatric | ≥ 10 kg | 1 m | TS-E-D | 1 |
|  | TruSignal Wrap sensor | Adult, Paediatric | ≥ 3 kg | 1 m | TS-W-D | 1 |
|  | TruSignal Sensitive Skin sensor | Adult, Paediatric, Infant, Neonate | No limit | 1 m | TS-SE-3 | 3 |

DISPOSABLE SENSORS (REQUIRES COMPATIBLE INTERCONNECT CABLE)

| | Sensors | Patient | Patient Weight | Length | Part Number | Quantity |
|---|--------------------------------------|---------------------------------------|----------------|--------|------------------------|----------|
|  | TruSignal sensor | Paediatric, Infant | ≥ 20 kg | 0,3 m | TS-AP-10 TS-AP-25 | 10 25 |
|  | TruSignal sensor | Adult, Paediatric, Infant, Neonate | No limit | 0,5 m | TS-AF-10 TS-AF-25 | 10 25 |
|  | TruSignal Adhesive Wrap sensor | Adult | > 20 kg | 0,3 m | TS-AAW-10 TS-AAW-25 | 10 25 |
|  | TruSignal Adhesive Wrap sensor | Paediatric, Infant | 3 - 20 kg | 0,3 m | TS-PAW-10 TS-PAW-25 | 10 25 |

TECHNOLOGY GE TRUSIGNAL

GE SPO2 ACCESSORIES

| Product | Size | Part Number | Quantity |
|-----------------------|-------|-------------|----------|
| Foam wrap replacement | Large | TS-RWL | 24 |

TECHNOLOGY NELLCOR™

For your monitors with Nellcor SpO2 technology.



INTERCONNECT CABLES


| Product | Length | Part Number | Quantity |
|--|--------|-------------|----------|
| Interconnect cable with GE connector and Nellcor OxiMax™ connector | 3 m | 2021406-001 | 1 |
| Nellcor SpO ₂ Interface cable | 3 m | DOC-10 | 1 |

NELLCOR™ SpO₂ REUSABLE SENSORS

PRODUCT AND FEATURES

ADULT, PAEDIATRIC

REUSABLE SpO₂ CLIP SENSORS



- Adult reusable finger and ear clip sensor
- Paediatric reusable two-piece clip sensor
- Suitable for short-term monitoring and spot-checks

- Rubber boot sensor
- Limited continuous
- Monitoring for virtually all areas of care for up to six hours
- Large and small sizes

NELLCOR™ SPO₂ REUSABLE SENSORS

| Sensors | Patient Weight | Accuracy 70–100% | Part Number | Quantity |
|---|----------------|------------------------------------|-------------|----------|
| Nellcor SpO ₂ Finger probe | Adult >40 Kg | ± 3 | DS100A | 1 |
| Nellcor SpO ₂ Wrap sensor, ear clip for D-YS | >30 Kg | ± 3.5 | D-YS,D-YSE | 1 |
| Nellcor SpO ₂ Wrap Sensor | >1 Kg | Infant-Adult ± 3.5 Neonatal ± 4 | D-YS | 1 |
| Nellcor SpO ₂ Adult Soft Sensor | Large >20 Kg | ± 2.5 | FLEXMAX | 1 |
| Nellcor SpO ₂ Paediatric Soft Sensor | Small >20 Kg | ± 2.5 | FLEXMAX-P | 1 |

TECHNOLOGY MASIMO™

For your monitors with MASIMO SpO² technology.

RD INTERCONNECT CABLE

| Product | Length | Part Number | Quantity |
|---------------------------------|--------|----------------|----------|
| Masimo RD SET GE cable for V100 | 3,6 m | MASIMO-RD-4085 | 1 |

APPLICATION



RD SET Adt
finger application



RD SET Pdt
finger application



RD SET Inf
thumb application



RD SET Neo
foot application



RD SET NeoPt /
NeoPt-500
foot application

RD SET SENSOR SERIES - TECHNOLOGY MASIMO™

| Sensors | Part Number | Quantity |
|--|-------------|----------|
| Masimo RD set DC-I, Adult Reusable sensor | 4050 | 1 |
| Masimo RD set DC-I, Paediatric Reusable sensor | 4051 | 1 |
| Masimo RD set DC-I, Reusable Tip Clip sensor | 4053 | 1 |
| Masimo RD set DC-I, Multi-site Reusable sensor | 4054 | 1 |
| Masimo RD Disposable Infant sensor for 3-20Kg | 4002 | 1 |



NIBP NON INVASIVE BLOOD PRESSURE

A patient's blood pressure (BP) is often the first step in determining the best clinical course of action. That's why quality data acquisition through your BP cuff is so essential. GE CRITIKON® BP cuffs from GE Healthcare are designed and tested for performance accuracy – from the first touch to the final measurement. By starting with accurate information, you can make well-informed decisions that may ultimately enhance patient outcomes.

NIBP HOSES

HOSES

| CRITIKON AIR HOSE | PATIENT | LENGTH | PART NUMBER | Quantity |
|---|-------------------|--------|-------------|----------|
| Assy Care Hose, 2-Tube DINACLICK™ Connector | Adult, Paediatric | 3.6m | 2058203-002 | 1 |
| Assy Care Hose, 2-Tube NEO-SNAP™ Connector | Neonate | 3.6m | 2058204-003 | 1 |
| Assy Critikon Hose, 2-Tube DINACLICK™ Connector | Adult | 3.6m | 2058205-001 | 1 |
| Assy Critikon Hose, 2-Tube NEO-SNAP™ Connector | Neonate | 3.6m | 2058206-001 | 1 |
| CARESCAPE ONE Hose, 2-Tube DINACLICK™ Connector | Adult, Paediatric | 3.6m | 2087389-002 | 1 |
| CARESCAPE ONE Hose, 2-Tube NEO-SNAP™ Connector | Neonate | 3.6m | 2089791-003 | 1 |

CUFFS DURA-CUF

CRITIKON™ DURA-CUF™ is gentle and designed for durability. They are constructed with materials not made with PVC, making them safer for the environment, and natural rubber latex, making them safer for your patient.



Easy-to-use, fast **DINACLICK connector** enables nurses to link the cuff to the NIBP hose quicker than a pair of traditional tube connectors. The audible click lets users know there is a proper connection.

2-TUBE DINACLICK CUFFS

| Size | Colour | Limb Circumference | 2- Tube DINACLICK | Quantity |
|------------------|------------|--------------------|-------------------|----------|
| Infant | Orange | 8-13 cm | DUR-P1-2A | 1 |
| Child | Green | 12-19 cm | DUR-P2-2A | 1 |
| Child Long | Green | 12-19 cm | DUR-P2-2A-L | 1 |
| Small Adult | Royal Blue | 17-25 cm | DUR-A1-2A | 1 |
| Small Adult Long | Royal Blue | 17-25 cm | DUR-A1-2A-L | 1 |
| Adult w/ Hanger | Navy | 23-33 cm | DUR-A2-2A | 1 |
| Adult Long | Navy | 23-33 cm | DUR-A2-2A-L | 1 |
| Large Adult | Wine | 31-40 cm | DUR-A3-2A | 1 |
| Large Adult Long | Wine | 31-40 cm | DUR-A3-2A-L | 1 |
| Thigh | Brown | 38-50 cm | DUR-T1-2A | 1 |


PRECISE MEASURES FOR INFORMED CHOICES

GE Healthcare offers the **first BP forearm cuff** clinically validated¹ for use on **obese or overly muscular adult patients**.

“By adapting to the anatomy of obese patients, the RADIAL-CUF cuff optimizes the accuracy of the **blood pressure** measurements of this particularly risky patient population. **77% Of patients** also consider it **more comfortable** than traditional blood pressure cuffs*.”

2-TUBE DINACLICK CUFF



| Size | Colour | Limb Circumference | 2-Tube Dinaclick | Quantity |
|-------------|--|--------------------|------------------|----------|
| Obese Adult |  Purple | 26-36 cm (Forearm) | SFT-F1-2A | 1 |













CUFFS SOFT-CUF

Our most popular cuff, the CRITIKON™ SOFT-CUF™ is designed with soft, absorbent material with rounded corners for maximizing patient comfort during extended stays at the hospital. These cuffs are ideal for patients with longer term stays or for protocols where the cuff travels with the patient through their total length of time in the hospital. It can be cleaned and reused on multiple patients, but is priced economically so it can be disposed of after a single patient use, helping you affordably prevent infections.



2-TUBE DINAClick CUFFS



| Size | | Colour | Limb Circumference | 2-Tube Dinaclick | Quantity |
|------------------|---|------------|--------------------|------------------|----------|
| Infant |  | Orange | 8-13 cm | SFT-P1-2A | 1 |
| Child |  | Green | 12-19 cm | SFT-P2-2A | 1 |
| Child Long |  | Green | 12-19 cm | SFT-P2-2A-L | 1 |
| Small Adult |  | Royal Blue | 17-25 cm | SFT-A1-2A | 1 |
| Small Adult Long |  | Royal Blue | 17-25 cm | SFT-A1-2A-L | 1 |
| Adult |  | Navy | 23-33 cm | SFT-A2-2A | 1 |
| Adult Long |  | Navy | 23-33 cm | SFT-A2-2A-L | 1 |
| Large Adult |  | Pink | 31-40 cm | SFT-A3-2A | 1 |
| Large Adult Long |  | Pink | 31-40 cm | SFT-A3-2A-L | 1 |
| Thigh |  | Brown | 38-50 cm | SFT-T1-2A | 1 |






CUFFS NEONATAL

Our neonatal BP cuffs are designed with the performance accuracy you trust and the gentle comfort you expect.


Each of our 5 sizes is offered in a extra soft absorbent material (SOFT-CUF™ range) or easy-to-clean vinyl (CLASSIC-CUF™ range). Economical enough for pairing with a single patient throughout hospital stay.

2-TUBE NEO-SNAP CUFFS



| Size | Colour | Limb Circumference | Soft-Cuf | Quantity |
|------|---|--------------------|-----------|----------|
| 1 |  Orange | 3-6 cm | SFT-N1-2B | 1 |
| 2 |  LT. Blue | 4-8 cm | SFT-N2-2B | 1 |
| 3 |  Green | 6-11 cm | SFT-N3-2B | 1 |
| 4 |  Navy | 7-13 cm | SFT-N4-2B | 1 |
| 5 |  Pink | 8-15 cm | SFT-N5-2B | 1 |

ASSORTMENTS NEO-SNAP CONNECTOR

| Product | 2-Tube Neo-Snap |
|---|---|
| 20 cuffs SOFT-CUF : 2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo # 5 |  2059305-002 |

“Playful numeric and colour-coded animal designs to help you find the right size at a glance.”





TEMPERATURE PROBES & CABLES




Make the right connection

- Disposable probes are not made with natural rubber latex.
- Single-patient-use temperature probes are well suited for infection control environments.
- Temperature probes and cables help you effectively utilise GE Healthcare patient monitoring technologies.
- Select from a product line that provides nasopharyngeal, oral, oesophageal, rectal, bladder, auditory canal, and axillary application options.

REUSABLE TEMPERATURE PROBES

Available in general purpose and skin options, with sizes for both paediatrics and adults, our reusable 400 series probes are designed for accuracy and flexibility.


TEMPERATURE PROBES

| | Application area | Response time | Sensor diameter | Cable length | Quantity |
|--|--|---------------|-----------------|--------------|----------|
| M1024247-General purpose probe, adult (1 per pkg.) | | | | | |
|  | Esophagus, rectum | 30 sec. | 4 mm/0.16" | 3 m/10' | 1 |
| M1024251-General purpose probe, paediatric (1 per pkg.) | | | | | |
|  | Esophagus, rectum | 20 sec. | 3 mm/0.12" | 3 m/10' | 1 |
| M1024254-Skin purpose probe, adult/paediatric (1 per pkg.) | | | | | |
|  | Finger, toe, axillary (armpit) using tape or posey wrap | 15 sec. | 10 mm/0.39" | 3 m/10' | 1 |


TEMPERATURE CARE CABLES

A wide breadth of cable options support disposable and reusable temperature probe use on GE monitoring technologies.


TEMPERATURE CABLES

| | Description and REF code | Cable length | Compatibility | Quantity |
|---|---|--------------|---|----------|
| Temperature care cables, for use with GE Healthcare patient monitors (1 per pkg.) | | | | |
|  | Dual Temperature Cable 2016998-001 | 0.5 m/1.7' | Connects directly to reusable probes or to interconnect cable (165640 and 165641) for disposable probes | 1 |
| | Single Temperature Cable 2021700-001 | 0.5 m/1.7' | | 1 |

Direct connect cable for disposable probes, for use with GE Healthcare patient monitors (1 per pkg.)

| | | | | |
|--|---|------------|--|---|
|  | Direct Connect Cable, Single 2021701-001 | 3.6 m/12' | Connects disposable probes directly to GE healthcare patient monitors (no interconnect cable required) | 1 |
| | 2021702-002 | 1.2 m/3.9' | | |

Interconnect cable for disposable probes (1 per pkg.)

| | | | | |
|--|--------|------|--------------------|---|
|  | 165641 | 2.8m | Interconnect Cable | 1 |
|--|--------|------|--------------------|---|



INVASIVE BLOOD PRESSURE

Featuring durable materials for extended cable life, our invasive pressure cables help provide accurate monitoring and measurement. These products have been verified for use with GE Healthcare equipment.

INVASIVE BLOOD PRESSURE ADAPTER



| Product | Length | Part Number | Quantity |
|---------------------------------------|--------|-------------|----------|
| Invasive Blood Pressure adapter cable | 30 cm | 2103966-001 | 1 |



CARESCAPE RESPIRATORY MODULES AND COMPACT AIRWAY GAS MODULE

GAS MONITORING

GE Healthcare's gas sampling lines are validated and rigorously tested to ensure optimal performance on your GE Healthcare equipment. Just like the other water traps in the D-Fend™ range, the Mini D-Fend protects your monitors from moisture, contaminations and secretions.

D-FEND, D-FEND+, MINI D-FEND, D-FEND PRO AND D-FEND PRO+ WATER TRAPS ARE SPECIFIED FOR USE WITH THE FOLLOWING MODULE

| | E-CO | E-COV | E-COVX | E-CAiO | E-CAiOV | E-CAiOVX | E-MINIC | | |
|-------------|-------|--------|---------|---------|----------|-----------|----------|-----------|--------|
| D-fend | | | | X | X | X | | | |
| D-fend+ | X | X | X | | | | | | |
| Mini D-fend | | | | | | | X | | |
| | E-sCO | E-sCOV | E-sCOVX | E-sCAiO | E-sCAiOV | E-sCAiOVX | E-sCAiOE | E-sCAiOVE | N-CAiO |
| D-fend Pro | | | | X | X | X | X | X | X |
| D-fend Pro+ | X | X | X | | | | | | |

WATER TRAPS



| Product | Speciality | Colour | Part Number | Quantity |
|-------------|-------------------------|------------|-------------|----------|
| D-fend | Anaesthesia | Black | 876446-HEL | 10 |
| D-fend | Pro Anaesthesia | Steel Blue | M1182629 | 10 |
| Mini D-Fend | Anaesthesia/Reanimation | Light Grey | 8002174 | 10 |
| D-fend Pro+ | Reanimation | Green | M1200227 | 10 |
| D-fend+ | Reanimation | Green | 881319-HEL | 10 |

In order to allow optimal operation of your gas module, the D-Fend water trap should be changed every month in the anaesthesia department and every 24 hours in the intensive care unit.

PATIENT SPIROMETRY ACCESSORIES

D-LITE FLOW SENSORS FOR PATIENT SPIROMETRY



733910-HEL



73393

Ventilator independent, the flow sensor obtains measurement at the patients airway. Only one connector is needed for lung mechanics, airway gases and metabolism. Lightweight sensors can be used anywhere GE Healthcare respiratory monitoring technology is applied, including in anesthesia delivery and ICU ventilation.

FOR HUMID CONDITIONS:

The D-lite+ and Pedi-Lite+ flow sensors' hydrophilic surface prevents the condensed water from disturbing the measurement.

SELECTION

Selection of the flow sensor is based on the tidal volume of the patient:

- D-lite(+): 150 – 2000 ml
- Pedi-lite(+): 5 – 300 ml ¹

CLEANING

The disposable D-lite, D-lite+ and Pedi-lite+ sensors are for single patient use only.

The reusable D-lite and Pedi-lite sensors can be:

- hand washed
- machine washed
- disinfected with suitable disinfectants
- steam autoclaved (max. temp. 132° C /270° F)

POSITIONING

Place all ports upwards to avoid condensed water from entering the tubings.

| | | |
|------------|----------------------------|---|
| 733910-HEL | D-Lite sensor, reusable | 1 |
| 73393 | Pedi-Lite sensor, reusable | 1 |


Flow sensors are specified for use with the following modules:

E-COV, E-COVX, E-CAiOV, E-CAiOVX, E-sCOV, E-sCAiOV

TECHNICAL SPECIFICATIONS OF FLOW SENSORS:

| | 733910-HEL D-lite Reusable | 73393 Pedi-lite Reusable |
|-------------------------------------|-------------------------------|--|
| Measurement range (Tidal volume) | 150 – 2000 ml | 5 – 300 ml 15 – 300 ml ¹ |
| Accuracy | ±6% or 30 ml | ±6% or 4 ml |
| Dead space | 9.5 ml | 2.5 ml |

¹ Measurement range 15 – 300 with E-COV, E-COVX, E-CAiOV, E-CAiOVX modules



SPIROMETRY

For measuring the respiratory gases and lung mechanics at the patient's airway, the products in the GE Healthcare patient spirometry range guarantee optimal performance in humid conditions, including in case of low flow anesthesia, thanks to their internal lining and their precise tolerances. Only GE Healthcare patient spirometry supplies are approved for use with your GE Healthcare equipment.

PATIENT SPIROMETRY TUBES

Spirometry tubes are for single patient use only

Note:

Use 2 m sampling lines with the gas exchange measurement

| P/N | Description | Quantity |
|------------|--|----------|
| 891191-HEL | Spirometry Tube, 6 m (20 ft) Note: Use for MRI monitoring | 1 |

Spirometry tubes are specified for use with the following modules:

E-COV, E-COVX, E-CAiOV, E-CAiOVX, E-sCOV, E-sCAiOV

PATIENT SPIROMETRY KITS



Pre-connected kits meet all application areas and are designed to help make set-up quicker and easier. Specific kits designed for humid conditions. Selection between D-lite+ and Pedi-lite+ sensor is based on the tidal volume of the patient.

Patient Spirometry kits are for single patient use only.

Kit Contains:

- Gas sampling line PVC/PE
- D-lite, D-lite+ or Pedi-lite+
- Spirometry tube

| P/N | Description | Quantity |
|----------|---|----------|
| 8004381 | Patient Spirometry kit, humid conditions, D-lite+, 2 m (7 ft) | 1 |
| 8004382 | Patient Spirometry kit, humid conditions, D-lite+, 3 m (10 ft) | 1 |
| 8002718 | Patient Spirometry kit, humid conditions, Pedi-lite+, 2 m (7 ft) | 1 |
| M1032634 | Patient Spirometry kit, humid conditions, Pedi-lite+, 3 m (10 ft) | 1 |

Spirometry tubes are specified for use with the following modules:

E-COV, E-COVX, E-CAiOV, E-CAiOVX, E-sCOV, E-sCAiOV

RESPIRATORY MONITORING SAMPLING LINES

ANESTHESIA GAS SAMPLING LINES, PVC/PE



73318-HEL,
73319-HEL,
73306

| P/N | Description | Quantity |
|-----------|---|----------|
| 73318-HEL | Gas sampling line, 2 m (7 ft) | 10 |
| 73319-HEL | Gas sampling line, 3 m (10 ft) | 10 |
| 73306 | Gas sampling line, 6 m (20 ft) | 5 |
| | Note: Use for MRI monitoring | |
| 8004383 | Anesthesia gas sampling line, 2 m (7 ft) with Elbow Connector | 10 |

Note: Use with D-lite sensors

GE Healthcare's gas sampling lines help ensure accurate measurements and problem-free use. They are tested and validated with GE Healthcare equipment. Gas sampling lines with an elbow connector are well-suited for use with flow sensors.

Anesthesia gas sampling lines are for single patient use only.

Use PVC/PE gas sampling lines when anesthetic agents are present.

Use 2 m gas sampling lines with the gas exchange measurement.

Gas sampling lines are specified for use with the following modules:

E-CO, E-COV, E-COVX, E-CAiO, E-CAiOV, E-CAiOVX, E-sCO, E-sCOV, E-sCAiO, E-sCAiOV, E-miniC

CO2 GAS SAMPLING LINES, PVC



733162-HEL,
733163

| P/N | Description | Quantity |
|------------|------------------------------------|----------|
| 733162-HEL | CO2 gas sampling line, 2 m (7 ft) | 10 |
| 733163 | CO2 gas sampling line, 3 m (10 ft) | 10 |

CO2 gas sampling lines are for single patient use only. Green identification line in CO₂ gas sampling lines. Use PVC/PE sampling lines in Anesthesia.

CO2 gas sampling lines are specified for use with the following modules:

E-CO, E-COV, E-COVX, E-sCO, E-sCOV, E-miniC

GAS MONITORING AIRWAY ADAPTERS

AIRWAY ADAPTERS



Airway adapters are for single patient use only.

Dead space:

Straight T-adapter: 7.5 ml

Elbow adapter: 10 ml

| P/N | Description | Quantity |
|-----------|--------------------|----------|
| 73385-HEL | Straight T-adapter | 10 |
| 73386-HEL | Elbow adapter | 10 |

Airway adapters are specified for use with the following modules:

E-CO, E-COV, E-COVX, E-CAiO, E-CAiOV, E-CAiOVX, E-sCO, E-sCOV, E-sCAiO, E-sCAiOV, E-miniC

LOW DEAD SPACE AIRWAY ADAPTERS



8001542

Low dead space airway adapters are well suited for patients with low tidal volume (TV).

Deadspace:

ET tubes with ID 2.5 mm: 0.1 ml

ET tubes with ID 3.0 mm: 0.15 ml

ET tubes with ID 3.5 mm: 0.2 ml

ET tubes with ID 4.0 mm: 1 ml

ET tubes with ID 4.5 mm: 1 ml

ET tubes with ID 5.0 mm: 1 ml

| P/N | Description | Quantity |
|---------|--|----------|
| 8001540 | Low dead space adapter for ET tubes with ID | 10 |
| 8001541 | Low dead space adapter for ET tubes with ID 3.0 mm | 10 |
| 8001542 | Low dead space adapter for ET tubes with ID 3.5 mm | 10 |
| 8001543 | Low dead space adapter for ET tubes with ID 4.0 mm | 10 |
| 8000062 | Low dead space adapter for ET tubes with ID 4.5 mm | 10 |
| 8000063 | Low dead space adapter for ET tubes with ID 5.0 mm | 10 |

Low dead space airway adapters are specified for use with the following modules:

E-CO, E-COV, E-COVX, E-CAiO, E-CAiOV, E-CAiOVX, E-sCO, E-sCOV, E-sCAiO, E-sCAiOV, E-miniC

GAS MONITORING ACCESSORIES

GAS SCAVENGING



M1031274



8004462, 8004463

| P/N | Description |
|--|---------------------------------|
| M1031274 | Exhaust line for gas scavenging |
| Exhaust lines are specified for use with the following modules: E-CO, E-COV, E-COVX, E-CAiO, E-CAiOV, E-CAiOVX, E-sCO, E-sCOV, E-sCAiO, E-sCAiOV | |

| P/N | Description |
|---------|--|
| 8004462 | Exhaust line with Coulter fitting, 41 in (1 m) Note: Use with S/5 Aespire* Carestion* and S/5 Avance* Carestation |
| 8004463 | Exhaust line with Coulter fitting, 7 in (18 cm) Note: Use with S/5 Avance Carestation |



M1031273

| P/N | Description |
|----------|------------------------|
| M1031273 | Adapter for gas return |

Gas return adapter is specified for use with the following modules:

E-CO, E-COV, E-COVX, E-CAiO, E-CAiOV, E-CAiOV

GAS MONITORING ACCESSORIES

CALIBRATION GASES & REGULATORS

Use the calibration gas recommended for your module and follow the recommended calibration intervals to ensure that the measurement accuracy remains within specifications.

| | ANESTHESIA | CRITICAL CARE | CO ₂ , E-MINIC | GAS EXCHANGE IN CRITICAL CARE |
|----------|---|---|---|--|
| | 755534-HEL Calibration Gas Regulator + 755583-HEL Calibration Gas, blue 5% CO ₂ , 33% N ₂ O, 2% Desfl., 55% O ₂ , Balance N ₂ | 755534-HEL Calibration Gas Regulator + 755581-HEL Calibration Gas, green 5% CO ₂ , 40% N ₂ O, 55% O ₂ | 755534-HEL Calibration Gas Regulator + 755580 Calibration Gas, grey 5% CO ₂ , Balance air | M1006864 Calibration Gas Regulator + 755587 Calibration Gas 5% CO ₂ , Balance O ₂ |
| E-sCO | | X | | |
| E-sCOV | | X | | |
| E-sCAiO | X | | | |
| E-sCAiOV | X | | | |
| E-CO | | X | | |
| E-COV | | X | | |
| E-COVX | | | | X |
| E-CAiO | X | | | |
| E-CAiOV | X | | | |
| E-CAiOVX | X | | | |
| E-miniC | | | X | |



MULTI ABSORBERS

**MEDICAL SODA LIME
SAFE AND EFFECTIVE CO₂ ABSORPTION***

The use of low flow anesthesia is a growing clinical trend throughout the world, driven by the need to reduce cost and improve patient care. The unique, D-profile granule of GE Medisorb® soda lime offers you cost-effectiveness with optimal performance.

MULTI ABSORBER

| Part Number | Description | Quantity |
|-------------|---|----------|
| 8003138 | Multi Absorber, disposable, 0.8 kg, white – violet change | 1 |

Pink to White colour change products not available for sale in the USA. Check with your local GE Healthcare representative for availability in your region.

DISPOSABLE MEDISORB® EX AND MEDISORB EF EX SODA LIME CO₂ ABSORBENT FOR GE CARESTATION™ 600 SERIES

Medisorb® EX and EF EX absorbent in disposable cartridges

Medisorb EX and low alkaline Medisorb EF EX disposable cartridges are designed for use with the GE Carestation™ 600 series to ensure optimal performance.

The EX cartridges have been fully tested and validated for use with the GE Carestation 600 series. The pre-filled, sealed cartridges can be installed quickly and easily to ensure CO₂ is removed from low-flow, closed loop anesthesia breathing circuits.



MEDISORB®

| Part Number | Description | Quantity |
|-------------|----------------|----------|
| 2079796-001 | Medisorb EX | 1 |
| 2079797-001 | Medisorb EF EX | 1 |





ENTROPY

Adequacy of Anaesthesia (AoA) is a concept made up of various parameters to help you assess patients' individual responses to the delivery of inhaled and intravenous hypnotics, opioids, and NMBA's during general anaesthesia.

The GE Healthcare Entropy module, E-ENTROPY, and accessories are indicated for adult and paediatric patients older than 2 years within a hospital for monitoring the state of the brain by data acquisition of electroencephalograph (EEG) and frontal electromyograph (FEMG) signals during general anaesthesia. The use of Entropy parameters may be associated with a reduction of anaesthetic use and faster emergence from anaesthesia¹.

SENSORS



| Product | Use | Part Number | Quantity |
|----------------|------------|-------------|----------|
| Entropy Sensor | Disposable | M1038681 | 1 |



| | | | |
|------------------------|------------|----------|---|
| Entropy EasyFit Sensor | Disposable | M1174413 | 1 |
|------------------------|------------|----------|---|

CABLE



| Product | Use | Part Number | Quantity |
|---------------|-------|-------------|----------|
| Entropy Cable | 3,5 m | M1050784 | 1 |



NEUROMUSCULAR TRANSMISSION - NMT

Quantitative NMT monitoring gives a clear picture of the individual dosage needs and facilitates optimal administration of NMBAs and antagonists. Studies have shown that the implementation of quantitative EMG neuromuscular monitoring resulted in a significant reduction in the incidence of incompletely reversed patients in the PACU¹.

CABLE

| Product | Length | Part Number | Quantity |
|------------------|--------|-------------|----------|
| NMT Sensor Cable | 3 m | 888414 | 1 |

SENSORS



| Product | Use | Patient | Part Number | Quantity |
|-------------------|------------|-------------------|-------------|----------|
| NMT Electrodes | Disposable | Adult, Paediatric | 57268-HEL | 30 |
| NMT Mechanosensor | Reusable | Adult | 888418 | 1 |
| | | Paediatric | 897439 | |
| NMT Electrosensor | Reusable | Adult, Paediatric | 888416 | 1 |

The life of a Mechanosensor in the context of average use is 2 years. Remember to renew them to guarantee the quality of the measurements.



Electromyography (EMG)
Requires 5 electrodes



Kinemyography (KMG)
Requires 2 electrodes

“EMG is the process of recording the specific electrical activity of muscle fibres in response to ulnar nerve stimulation.”



EEG NEUROLOGY


Electroencephalogram(EEG), Frontal Electromyogram (FEMG) and Audio Evoked Potential (AEP) measurements require good electrode contacts for high-quality signal. GE EEG electrodes and accessories are designed to ensure positive contact with the skin and enable reliable connections.

For use with EEG module, E-EEGX module and N-EEGX headbox with compatible CARESCAPE™ monitors.

REUSABLE PRECONFIGURED EEG LEADSETS

| Product | Length | Part Number | Quantity |
|--|--------|-------------|----------|
| Pre-configured EEG Lead Set 1, Basic | 0,6 m | 898050 | 1 |
| Pre-configured EEG Lead Set 2, General | 0,6 m | 898051 | 1 |

REUSABLE PRECONFIGURED EEG LEADSETS



| Product | Part Number | Quantity |
|---------------------------------|-------------|------------------|
| EEG Clip electrodes Ag/AgCl | M1240794 | 40 pouches of 25 |
| EEG Cup electrodes Ag/AgCl | 572686 | 500 |
| EEG Prewired cup electrodes | 2083319-001 | 46 pouches of 10 |
| EEG Prewired wet gel electrodes | 2083320-001 | 35 pouches of 12 |

All electrodes within a set-up should be similar



CARDIAC OUTPUT

Cardiac output cables from GE Healthcare are tested and verified to function with your equipment. Made with premium quality materials that minimize cable curls and tangles, the cables have a connector for an injectate probe and also a standard three-pin catheter connector that fits most Swan-Ganz™ catheters. A molded connector assembly helps keep cables flexible and prevents ingress of fluids.

CARDIAC OUTPUT TRUNK CABLES

| Cable | Length | Part Number | Quantity |
|----------------------|--------|-------------|----------|
| Cardiac Output Cable | 2,4 m | 2104180-002 | 1 |
| | 3,6 m | 2104180-001 | |

For use with the Patient Data Module (PDM)



CARDIAC OUTPUT PROBES

| Description | Length | For use with | Part Number | Quantity |
|-------------------------|------------|----------------------|-------------|----------|
| Injectate Bath Probe | 1.8 m/6 ft | CARESCAPE™ PDM | 9446-909 | 1 |
| In-Line Injectate Probe | .9 m/3 ft | Edwards™ CO Set, PDM | 9446-910 | 1 |
| | .9 m/3 ft | BD CritiKit, PDM | 9446-911 | 1 |
| | 2.4 m/8 ft | Thermoset, PDM | 9446-913 | 1 |



CARDIAC OUTPUT

ARGON™ ACCESSORIES

| Description | For use with | Part Number | Quantity |
|-------------------------------------|----------------|-------------|----------|
| Spectramed CO Set Temperature Probe | E-COP, E-COPSv | 16593-HEL | 1 |

EDWARDS™ ACCESSORIES

| Description | For use with | Part Number | Quantity |
|--|---------------------------------------|-------------|----------|
| Injectate Bath Temperature Probe | E-COP, E-COPSv | 16591 | 1 |
| Flow-through Injectate Temperature Probe | Edwards™ catheters, E-COP, E-COPSv | 16592VS | 1 |

OTHER ACCESSORIES

| Description | For use with | Part Number | Quantity |
|---|----------------|-------------|----------|
| CO Set for Room Temperature Injectate, disposable | E-COP, E-COPSv | 16574VS | 1 |
| Catheter Connecting Cable | E-COP, E-COPSv | 16590-HEL | 1 |
| Thermodilution Catheter, disposable | E-COP, E-COPSv | M1238258 | 5 |

PICCO™ CARDIAC OUTPUT CABLES


| Description | Length | For use with | Part Number | Quantity |
|---------------------------------|--------|----------------|-------------|----------|
| Continuous Cardiac Output Cable | 4.1m | E-PiCCO module | M1154585 | 1 |
| Injectate Sensor Cable | 0.9m | | M1182448 | 1 |



CARDIAC OUTPUT

PICCO™ KITS

| Patient | Description | Includes | PULSIOCATH Arterial Thermodilution | Part Number | Quantity |
|-----------------|-----------------|--|---------------------------------------|-------------|----------|
| Adult | Femoral artery | Five Catheters, one pressure transducer with flush device and one PV4046 PiCCO™ -Technology injectate temperature sensor housing | PV2015L20 Standard | M1219758 | 1 |
| | Brachial artery | | PV2014L22 Brachial | M1219759 | 1 |
| Children/Infant | Femoral artery | | PV2013L07 Femoral | M1219760 | 1 |






PATIENT CIRCUITS

There are many new products which offer the opportunity to choose from an expanded portfolio and to increase your opportunities to purchase from Medhold.


DISPOSABLE PATIENT CIRCUITS

| Part Number | Description | Quantity |
|-------------|---|----------|
| M1019517 | Manual Ventilation Kit for ABS Bag Arm, Adult, Disposable Including: 1 Manual ventilation tube, 0.9 m 1 Latex-free manual bag, 2 L 1 ABS Bag arm connector | 1 |




DISPOSABLE PATIENT CIRCUITS

| Part Number | Description | Quantity |
|-------------|---|----------|
| M1019527 | Patient Circuit Kit, Adult, Disposable Including: 1 Patient circuit, 1.8 m 1 Gas sampling elbow 1 Manual ventilation tube, 1.8 m 1 Latex-free manual bag, 2 L 1 Anaesthesia gas sampling line, 3.0 m | 1 |



DISPOSABLE PATIENT CIRCUITS


| Part Number | Description | Quantity |
|-------------|---|----------|
| M1019507 | ADU Patient Circuit, Adult, Disposable Including: 1 Patient circuit, 3.0 m 1 Standard elbow 1 Manual ventilation tube with quick connector, 3.0 m 1 Ventilation tube, 1.0 m 1 Latex-free manual bag, 2 L | 1 |



PATIENT CIRCUITS


DISPOSABLE PATIENT CIRCUITS

| Part Number | Description | Quantity |
|-------------|---|----------|
| M1019541 | Patient Circuit Kit, Adult, Disposable Including: 1 Patient circuit, 1.5 m 1 Standard elbow | 1 |



DISPOSABLE PATIENT CIRCUITS

| Part Number | Description | Quantity |
|-------------|--|----------|
| M1019544 | Patient Circuit, Paediatric, Disposable Including: 1 Patient circuit, 1.5 m | 1 |



PATIENT FILTERS - BACTERIAL AND VIRAL FILTERS, SINGLE PATIENT USE

Filter media is made of permanently charged hydrophobic fibers that help enhance particle capture efficiency and have a low and stable breathing resistance, even when wet.

Various types are available, with or without sampling port.



| Part Number | Description | Bacterial Efficiency | Viral Efficiency | Weight | Deadspace | Box Quantity |
|-------------|--------------|----------------------|------------------|--------|-----------|--------------|
| 557022500 | Uni filter/S | 99.9999% | 99.999% | 16g | 35ml | 1 |

COMBINED HUMIDIFICATION AND FILTRATION









| Part Number | Product | Filtration efficiency | Weight | Dead Space | Quantity |
|-------------|--------------------|-----------------------|--------|------------|----------|
| 557070100 | HMEF 1000/S | 99.9999 | 24g | 77ml | 1 |
| M1010534 | Angled HMEF 1000/S | 99.9999 | 24g | 77ml | 1 |
| M1004132 | HMEF 750/S | 99.9999 | 17g | 34ml | 1 |
| M1010538 | Angled HMEF 750/S | 99.9999 | 17g | 34ml | 1 |
| 557070500 | HMEF 500/S | 99.999 | 15g | 30ml | 1 |
| 8004231 | HMEF Mini/S | 99.999 | 14g | 21ml | 1 |










DEFIBRILLATOR SUPPLIES AND ACCESSORIES REANIBEX



ELECTRODE PADS FOR DEFIBRILLATION & PACING

| Part Number | Product |
|---|---|
|  TKL 0311 D | Adult Disposable cable-electrodes for defibrillation |
|  TKL 0313 A | Adult Preconnected disposable cable-electrodes. |
|  TKL 0312 B | Paediatric Disposable cable-electrodes |
|  KSA 0501 D | Paediatric Disposable cable-electrodes. |
|  KSA 0108 | Reusable longer cable ONLY for ADULT disposable cable-electrodes for defibrillation (1,5 m long). Valid only for model REANIBEX-Serie 200 and 700. |
|  KSA 0104 A | Reusable longer cable for ADULT and PAEDIATRIC disposable cable-electrodes for defibrillation (1,5 m long). Valid only for models REANIBEX-300,500 and 800. |
|  KSA 0502 | Set of ADULT & PAEDIATRIC disposable cable-electrodes preconnected for defibrillation. |

ECG CABLES

| | |
|---|---|
|  KPA 0400 A | 3 lead patient cable (only for defibrillator R300 with Manual override option) |
|  KPA 0401 A | 4 I lead patient cable (visualisation of Lead I, II and III) |
|  TKC 0616 A | 6 lead patient cable to be attached to the 4 lead patient cable (visualisation of Lead I, II, III, aVF, aVL, aVR and from V1 to V6) |
|  KPA 0311 | 3 leads patient cable |
|  KPA 0309 | 5 lead patient cable |
|  KPA 0310 | 10 lead patient cable |
|  KPA 0306 A | 3 leads patient cable |
|  KPA 0307 A | 5 lead patient cable |
|  KPA 0308 A | 10 lead patient cable |

| REANIBEX MODEL | | | | | |
|----------------|-----|-----|-----|-----|-----|
| 100 | 200 | 300 | 500 | 700 | 800 |
| | X | X | X | X | X |
| | X | X | X | X | X |
| | X | | | | |
| | | X | X | X | X |
| | X | | | X | |
| | | X | X | | X |
| X | | | | | |
| | | X | X | | |
| | | X | X | | |
| | | X | X | | |
| | | | | X | |
| | | | | X | |
| | | | | X | |
| | | | | | X |
| | | | | | X |
| | | | | | X |

DEFIBRILLATOR SUPPLIES AND ACCESSORIES REANIBEX



MASIMO SPO2 ACCESSORIES

REANIBEX MODEL

| | Part Number | Product |
|--|-------------|--|
| | TKC 0201 A | Masimo Rainbow SpO2 transducer longer cable (12') |
| | TKC 0201 B | Masimo Rainbow SpO2 transducer longer cable (4') |
| | TKC 0203 A | Masimo Pulsioxymeter transducer longer cable (4 feet) |
| | TKC 0213 A | Masimo Reusable Adult's finger SpO2 transducer (> 30 kg) |
| | TKC 0214 A | Masimo Multisite Reusable SpO2 transducer (> 1kg) |
| | TKC 0216 A | Masimo Reusable Paediatric SpO2 transducer (10-50 kg) |
| | TKC 0221 A | Rainbow DCI, Adult Reusable Sensor, 3 ft - SpCO, SpMet, SpO2 |
| | TKC 0226 A | Rainbow DCIP, Paediatric Reusable Sensor, 3 ft - SpCO, SpMet, SpO2 |
| | TKC 0222 A | Rainbow DCI, Adult Reusable Sensor, 3ft - SpHb, SpMet, SpO2 |

| 100 | 200 | 300 | 500 | 700 | 800 |
|-----|-----|-----|-----|-----|-----|
| | | | | | X |
| | | | | | X |
| | | | X | X | |
| | | | X | X | X |
| | | | | X | X |
| | | | X | X | X |
| | | | | | X |
| | | | | | X |
| | | | | | X |
| | | | | X | X |
| | | | | X | X |
| | | | | X | X |
| | | | | X | X |
| | | | | X | X |
| | | | | X | X |








SUNTECH NIBP ACCESSORIES

| | | |
|--|------------|---|
| | TKC 0710 A | 3 m NIBP air hose |
| | TKC 0712 A | infant reusable cuff (8 to 13 cm) |
| | TKC 0713 A | child reusable cuff (12 to 19 cm) |
| | TKC 0714 A | small adult reusable cuff (17 to 25 cm) |
| | TKC 0711 A | adult reusable cuff (23 to 33 cm) |
| | TKC 0715 A | large adult reusable cuff (31 to 40 cm) |
| | TKC 0716 A | Thigh reusable cuff (38 to 50 cm) |

DEFIBRILLATOR SUPPLIES AND ACCESSORIES REANIBEX





MASIMO ETCO2 ACCESSORIES & SPARES




| | Part Number | Product |
|---|-------------|--|
|  | TES 0602 A | Module ISA CO2 side stream analyzer |
|  | KIA 0501 A | Communication Cable EtCO2 to R500 EMS |
|  | TKC 0840 B | Adaptor probe EtCO2 Nomoline Multi-patient to R500 EMS |
|  | TKC 0801 A | Adaptor T intubated EtCO2 to R500 EMS |
|  | TKC 0802 A | NOMO 2 m extension to R500 EMS |
|  | TKC 0841 A | Probe ORAL/NASAL Adult EtCO2 to R500 EMS |
|  | TKC 0842 A | Probe ORAL O2/NASAL Adult EtCO2 ' to R500 EMS |

ORIDION ETCO2 ACCESSORIES

ORIDION NON-INTUBATED PATIENTS (NASAL - ORAL)

| | | |
|---|------------|---|
|  | TKC 0810 A | Smart Capnoline Plus Adult /Intermediate O2 (ORAL-NASAL short term duration, adult size, with O2 tubing) pack of 25 units ref: 009822 (TKC 0810) |
|  | TKC 0811 A | Smart Capnoline Paediatric O2 (ORAL-NASAL short term duration, paediatric size, with O2 tubing) pack of 25 units ref: 007269 (TKC 0811) |

ORIDION INTUBATED PATIENTS

| | | |
|---|------------|--|
|  | TKC 0814 A | Filterline Set Adult/Paediatric (short term duration, adult/paediatric size) pack of 25 units ref: XS 04620 (TKC 0814) |
|  | TKC 0815 A | Filterline H Set Adult/Paediatric (long term duration, adult/paediatric size) pack of 25 units ref: XS04624 (TKC 0815) |
|  | TKC 0816 A | Filterline H Set Infant/Neonates (long term duration, infant/neonates size) pack of 25 units ref: 006324 (TKC 0816) |

| REANIBEX MODEL | | | | | |
|----------------|-----|-----|-----|-----|-----|
| 100 | 200 | 300 | 500 | 700 | 800 |
| | | | X | | |
| | | | X | | |
| | | | X | | |
| | | | X | | |
| | | | X | | |
| | | | X | | |
| | | | X | | |
| | | | | | X |
| | | | | | X |
| | | | | | X |
| | | | | | X |



REANIBEX MODEL

37



DEFIBRILLATOR SUPPLIES AND ACCESSORIES REANIBEX




ECG ELECTRODES

| Part Number | Product |
|--|--|
|  TKL 102 A | Adult disposable electrodes for the patient cable (bag of 50 units) |
|  TKL 212 B | Paediatric disposable electrodes for the patient cable (bag of 50 units) |

PAPER

| | |
|---|--|
|  XEG 001 A | Thermosensible paper roll, 50 mm and 28 m for the printer |
|  XEG 0004 A | Thermosensible paper roll, 106 mm and 25 m for the printer |

GEL

| | |
|---|-------------------|
|  XEH 001 A | Gel tube of 250ml |
|---|-------------------|

TEST LOAD

| | |
|-----------|-------------------|
| XEH 001 A | Test Load 50 Ohms |
|-----------|-------------------|

CPR KIT

AED CPR Kit. Contains: 1 disposable CPR mask, 1 hygienic protective disposable mask, scissors, 1 pair nitrile gloves, gauze and razor.










REANIBEX MODEL

| 100 | 200 | 300 | 500 | 700 | 800 |
|-----|-----|-----|-----|-----|-----|
| | | X | X | X | X |
| | | X | X | X | X |
| | | | | X | X |
| | | | | | X |
| | | X | X | X | X |
| | | | | X | X |
| X | X | X | X | | |

DEFIBRILLATOR SUPPLIES AND ACCESSORIES REANIBEX



BATTERIES



| Part Number | Product |
|--|--|
|  DFS 4100 A | Non rechargeable LiSO2 battery of 7,5 A/h. With a shelf life of 5 years, minimum 300 discharges or more than 12 hours of monitoring. |
|  DIS 3100 A | Non rechargeable LiSO2 battery. With a shelf life of 5 years, minimum 270 discharges at 360J. |
|  DIS 3000 A | Rechargeable Li-Ion battery. Minimum 200 discharges at 360 J. |
|  DGS 5000 | Rechargeable NiMH battery. Autonomy of more than 130 discharges at 200J. Charging time: 3 hours |
|  DHS 3000 A | Rechargeable NiMH battery. Autonomy of more than 150 discharges at 360 J. Charging time: 3 hours |
|  DJS 3000 | Non rechargeable battery LiMnO2 for defibrillator model Reanibex-100. Capacity of more than 350 shocks at 200 Joules and more than 36 hours of ECG monitoring. |
|  KAA 109 A | 12VDC battery connection cable R700 |
|  KAA 0207 A | 12VDC battery connection cable R800 |
|  KAA 0501 A | 220V adapter allows the R800 to work with 2 batteries and 220V |

CARRY CASE

200

| | |
|---|----------------------|
|  DFX 0120 A | Simple transport bag |
|---|----------------------|

300 & 500

| | |
|---|--|
|  DIX 8040 A | BASIC Carrying case with shoulder strap |
|  DIS 8010 A | BASIC Thermoconformed Carrying case with shoulder strap. |




REANIBEX MODEL

| 100 | 200 | 300 | 500 | 700 | 800 |
|-----|-----|-----|-----|-----|-----|
| | X | | | | |
| | | X | X | | |
| | | X | X | | |
| | | | | X | |
| | | | | | X |
| X | | | | | |
| | | | | X | |
| | | | | | X |
| | | | | | X |
| | X | | | | |
| | | X | X | | |
| | | X | X | | |


DEFIBRILLATOR SUPPLIES AND ACCESSORIES REANIBEX







CARRY CASE 300 & 500 - CONTINUED

| | Part Number | Product |
|---|-------------|--|
|  | DIS 8000 A | FULL Thermoconformed Carrying case with shoulder strap and side pocket. |
| | DIS 8001 A | FULL Thermoconformed Carrying case with shoulder strap and XL side pocket. |
|  | DIS 8020 A | Side pocket for Basic Thermoconformed Basic Carrying case |
|  | DIS 8030 A | Side pocket XL for Basic Thermoconformed Basic Carrying case |






700

| | | |
|--|-----------|---|
|  | DGX 705 A | Transport bag with side pockets and shoulder strap. |
|--|-----------|---|

800

| | | |
|---|------------|--|
|  | DHS 8100 A | Transport bag with 2 side pockets. Includes shoulder strap, right pouch and left pouch |
|  | DHS 8130 A | Replacement Black shoulder strap for the transport bag |
|  | DHS 8140 A | Backpack strap for the transport bag |
|  | DHX 8020 A | Neoprene upper bag (only for defibrillator without external reusable paddles) |

WALL BRACKETS



| | | |
|---|--------------|---|
|  | VYC 0201 | AED signal (A4) |
|  | DFS 0800 A | Wall bracket for Ambulance for REANIBEX-200 AED. 10 G certified (EN 1789) |
|  | VYC 0210 A/B | REANIBEX-AED Wallcase. Optional Alarm system when opened. Other models available with Communication systems (GSM, GPRS), heater, videovigilance |
|  | DIS 8100 A | Wall bracket for Ambulance for REANIBEX-300 and REANIBEX-500. 10 G certified (EN 1789). Includes plastic backside part for hanging AED. |
|  | DGS 0850 | Wall bracket for ambulance. 10 G certified (EN 1789) |

| REANIBEX MODEL | | | | | |
|----------------|-----|-----|-----|-----|-----|
| 100 | 200 | 300 | 500 | 700 | 800 |
| | | X | X | | |
| | | | X | | |
| | | X | X | | |
| | | | X | | |
| | | | | X | |
| | | | | | X |
| | | | | | X |
| | | | | | X |
| | | | | | X |
| X | X | X | X | | |
| X | X | | | | |
| X | X | X | X | | |
| | | X | X | | |
| | | | | X | |

DEFIBRILLATOR SUPPLIES AND ACCESSORIES REANIBEX





WALL BRACKETS - CONTINUED

| Part Number | Product |
|---|--|
|  DHB 2048 A | Basic Wall bracket support for ambulance. 10 G certified (EN 1789) Wall bracket support for ambulance with 12 VDC charger (without plug connector). 10 G certified (EN 1789) Wall bracket support for ambulance with 12 VDC charger (with plug connector). 10 G certified (EN 1789) Wall bracket support for ambulance with AC charger. 10 G certified (EN 1789) |
|  DHS 8530 A | Wide metal BED RAIL HOOK for hospital use |


INTERNAL DEFIBRILLATION ID SYSTEM

| | |
|-----------|--|
| 13461/202 | ID Smart cable (reusable) |
| 1100 | ID Electrode common for adult and paediatric patients (disposable) |
| 1200 | ID Paddle for adult patient (disposable) |
| 1300 | ID Paddle for paediatric patient (disposable) |

TRAINING MATERIAL

| | |
|---|---|
|  QBS 500 A | ADULT Reusable cable for the REAN IBEX-200& 300 TRAINER without the electrodes. |
|  QBY 030 A (APEX) / QBY 031 A (STERNUN) | Pair of ADULT electrodes for the reusable cable. |
| DIS 8350 A | PAEDIATRIC Reusable cable for the REAN IBEX-300 TRAINER without the electrodes. |
| DIY 8350 A (APEX) / DIY 8351 (STERNUN) | Pair of PAEDIATRIC electrodes for the reusable cable. |

| REANIBEX MODEL | | | | | |
|----------------|-----|-----|-----|-----|-----|
| 100 | 200 | 300 | 500 | 700 | 800 |
| | | | | | X |
| | | | | | X |
| | | | | | X |
| | | | | | X |
| | | | | | X |
| | | | | | X |
| | | | | X | X |
| | | | | X | X |
| | | | | X | X |
| | | | | X | X |
| | X | X | X | | |
| | X | X | X | | |
| | X | X | X | | |
| | X | X | X | | |



DIAGNOSTIC CARDIOLOGY AND STRESS ECG SYSTEMS

MAC™, CASE™
AND CARDIOSOFT™



CABLES AND LEADWIRES – IEC

ECG TRUNK CABLES – IEC



| Description | Length | Part Number | MAC™ 2000 | MAC™ 600 | Quantity |
|-------------------------|--------|-------------|--------------|-------------|----------|
| Trunk Cable 10-lead IEC | 210 cm | 2104722-001 | X | X | 1 |

INTEGRATED CABLE AND LEADWIRE SET – IEC



| Description | Length | Part Number | MAC™ 2000 | MAC™ 600 | Quantity |
|---|---|--------------------------|--------------|-------------|----------|
| ECG Value Cable 10-lead Cable / Leadwire IEC | Chest leads: 80 cm / 32 in Limb leads: 125 cm / 49 in Trunk cable: 210cm | 2029890-001V 38401804 | X | X | 1 |

CABLES AND LEADWIRES – IEC

ECG LEADWIRE SET – GRABBER – IEC



| Description | Length | Part Number | MAC™ 2000 | CASE™ | CardioSoft™ |
|--------------|---------------------|-------------|--------------|-------|-------------|
| Base 10, IEC | Limb lead: 70 cm | 38401808 | X | X | X |

ACCESSORIES

ACCESSORIES

| Description | Part Number | MAC™ 2000 | MAC™ 600 | CASE™ | CardioSoft™ | Quantity |
|-------------------------------|-------------|--------------|-------------|-------|-------------|----------|
| Adapter Clip, Wide Mouth Clip | 38401804 | X | X | X | X | 14 |





ECG RECORDING PAPER

ACCESSORIES

| Folding | White patient data area | Paper size | Red grid wide | Queue | MAC™ 2000 | MAC™ CASE™ 600 | Part Number | Quantity | |
|---------|-------------------------|--------------------------------------|-------------------|-------|-----------|----------------|-------------|-------------|--------------------|
| z-fold | No | A4 * | 195 mm/ 7.6 in | Hole | | | X | 2009828-021 | 300 sheets/ 1 pack |
| | | | 205 mm/ 8 in | Block | X | | | 22616605 | 150 sheets/ 1 pack |
| | | 80 mm x 90 mm/ 3.1 in x 3.5 cm | 75 mm/ 2.9 in | Block | | X | | 2030887-001 | 280 sheets/ 1 pack |



22616605



2009828-021

ECG ELECTRODE

GE Healthcare reusable bulb electrodes with silver-silver chloride (Ag/AgCl) plating are designed to ensure positive contact with the skin and provide maximum signal quality.

REUSABLE BULB ELECTRODE

| Description | Part Number | Quantity |
|-------------------------|-------------|----------|
| Reusable Bulb Electrode | 2104783-001 | 6 |



BULB ELECTRODE ACCESSORIES

| Description | Compatibility | Part Number | Quantity |
|----------------------|--|-------------|----------|
| Clamp Limb IEC Value | MAC™ 600, MAC™ 800, MAC™ 2000 ECG System | 2104785-001 | 4/set |



EASYONE CS AND TRUEFLOW SPIROMETRY

| Description | Part Number | Quantity |
|----------------------|-------------|-----------|
| Disposable Spirettes | 2030431-007 | Box of 50 |





SEER® HOLTER MONITORS

SEER® 12

SEER® 12 – IEC

| Description | Length | Type | Part Number | Quantity |
|-----------------------------|--------------|-------|-------------|----------|
| 12-Channel Holter ECG Cable | 130 cm/51 in | Round | 2104808-001 | 1 |
| | 105 cm/41 in | Flat | 2104809-001 | |



SEER® 12 HOLTER MONITOR ACCESSORIES

| Description | Part Number | Quantity |
|---|-------------|----------|
| SEER® 12 lead portable pouch, waist belt, shoulder strap, 1 battery compartment | 2104803-001 | 1 |
| SEER® 12 Carrying pouch, 2 battery compartment | 2014606-297 | 1 |



SEER® 1000

SEER® 1000 LEADWIRE SET – IEC

| Description | Length | Part Number | Quantity |
|------------------|---------------|-------------|----------|
| 5 leadwire, 2 CH | 75 cm/29.5 in | 2104810-001 | 1 |
| 7 leadwire, 3 CH | | 2104811-001 | 1 |
| 5 leadwire, 2 CH | 105 cm/41 in | 2104812-001 | 1 |
| 7 leadwire, 3 CH | | 2104813-001 | 1 |
| 3 leadwire, 3 CH | 75 cm/29.5 in | 2104818-001 | 1 |
| | 105 cm/41 in | 2104819-001 | 1 |



Product compatibility


SEER® 1000 ECG cables are compatible with SEER® 1000 24-hour Holter monitor, SEER® 1000 48-hour Holter monitor and SEER® 1000 7-day Holter monitors.



SEER® 1000

2067634-073

SEER® 1000 HOLTER MONITOR ACCESSORIES

| Description | Part Number | Quantity | |
|--|-------------|----------|---|
| Pouch with waist belt and shoulder strap | 2067634-073 | 1 |  |
| Pouch Belt Clip | 2067634-166 | 1 | |

2067634-166

SEER® LIGHT KITS

SEER® LIGHT LEADWIRE SET KITS - IEC

| Description | Length | Part Number | Quantity |
|---|--------------|-------------|----------|
| Carry Pouch and 2 Channel Patient cable | 140 cm/55 in | 2042686-003 | 1 |
| Carry Pouch and 3 Channel Patient cable | | 2042686-004 | 1 |
| Carry Pouch and 2 Channel Patient cable | 60 cm/24 in | 2042686-007 | 1 |
| Carry Pouch and 3 Channel Patient cable | | 2042686-008 | 1 |

Product compatibility

SEER® Light and SEER® Light Extend

SEER® LIGHT ACCESSORIES

SEER® LIGHT ACCESSORIES

| Description | Part Number | Quantity | |
|----------------|-------------|----------|--|
| Compact Flash | 2008750-003 | 1 | Product compatibility SEER® Light and SEER® Light Extend |
| Connect Kit | 2020044-292 | 1 | |
| Controller Kit | 2042685-007 | 1 | |

SEER® LIGHT ACCESSORIES

SEER® LIGHT HOLTER MONITOR ACCESSORIES

| Description | Part Number | Quantity |
|--|--------------|----------|
| SEER® Light, shoulder straps, waist belt, portable pouch | 2104827-001* | 1 |



Product compatibility

2104827-001 is compatible with SEER® Light and SEER® Light Extend


TONOPORT VI

ACCESSORIES

| Description | Part Number | Quantity |
|-------------------------------|-------------|----------|
| Supply BP Bag Belt TONOPORT V | 2104824-001 | 1 |




TONOPORT AMBULATORY BLOOD PRESSURE - TONOPORT CUFFS

| Description | Part Number | Quantity |
|--|-------------|----------|
|  TONOPORT V BP Cuff Adult Small | 2001589-211 | 1 |
| TONOPORT BP Cuff Adult Standard | 2001589-212 | 1 |
| TONOPORT BP Cuff Adult Large | 2001589-213 | 1 |
| TONOPORT BP Cuff Adult Extra Large | 2001589-214 | 1 |

TANGO M2

TANGO M2 CUFFS

| Description | Part Number | Quantity |
|--|-------------|----------|
|  SPLY TANGO M2 CUFF ORBIT-K SMALL ADULT | 2071179-006 | 1 |
| SPLY TANGO M2 CUFF ORBIT-K ADULT | 2071179-003 | 1 |
| SPLY TANGO M2 CUFF ORBIT-K ADULT PLUS | 2071179-004 | 1 |
| SPLY TANGO M2 CUFF ORBIT-K LARGE ADULT | 2071179-005 | 1 |
| SPLY TANGO M2 ORBIT-K MICROPHONE | 2071179-012 | 1 |

MATERNAL[®] INFANT CARE

GIRAFFE FAMILY OPTIONS

GIRAFFE* OMNIBED*, GIRAFFE INCUBATOR, GIRAFFE WARMER, AND PANDA* WARMER



PHOTOTHERAPY OPTIONS (NOT SHOWN)

| Description | Part Number | Quantity |
|--|-------------|----------|
| Pad covers, disposable, small BiliSoft | M1093120 | 1 |
| Pad covers, disposable, large BiliSoft | M1093121 | 1 |
| BiliSoft Nest, disposable pad, small | M1093122 | 1 |
| BiliSoft Nest, disposable pad, large | M1093123 | 1 |


GIRAFFE FAMILY OPTIONS (NOT SHOWN)

| Description | Part Number | Quantity |
|--------------------------------------|---------------|----------|
| Disposable patient probe | 2074817-001 | 1 |
| Reusable patient probe | 2075796-001 | 1 |
| Heat reflecting patch | 0203-1980-300 | 50 |
| IRIS cuff | 6600-0499-800 | 1 |
| IRIS cuff | 6600-0500-800 | 1 |
| T-piece, disposable | M1091335 VS | 1 |
| T-piece, disposable with mask size 0 | M1091316 VS | 1 |
| T-piece, disposable with mask size 1 | M1091365 VS | 1 |



MATERNAL FETAL CARE COROMETRICS™ & MINI TELEMETRY

MIC FETAL MONITORING RECORDING THERMAL PAPER



| Scaling | Corometrics™ 170 series | Corometrics™ | Part Number | Quantity |
|--------------------|----------------------------|--------------|-------------|----------|
| 30-240 bpm scaling | X | X | 2009828CAO | 1 |
| 50-210 bpm scaling | X | X | 2009828DAO | 1 |

COROMETRICS™

MECG ACCESSORIES

| Description | | Part Number | Quantity |
|----------------------|---------------------------|-------------|----------|
| Gray MEGC Cable | 3 Lead Set, Model 118 | 1554AAO | 1 |
| | 3 Lead Set IEC, Model 118 | 1554BAO | 1 |
| 120 Green MEGC Cable | 3 Lead Set | 1553AAO | 1 |
| | 3 Lead Set IEC | 1553BAO | 1 |



TRANSDUCER BELTS

| Description | | Part Number | Quantity |
|--------------------|--------------------|-------------|----------|
| Single patient use | Button-style, blue | 2015827-001 | 10 |



2015827-001

ULTRASOUND TRANSDUCERS

| Description | Part Number | Quantity |
|---------------------------------------|-------------|----------|
| Ultrasound Transducer Button and Loop | 5700 HAX | 1 |
| TOCO Transducer Button and Loop | 2264 HAX | 1 |



INTERNAL MONITORING

INTRAUTERINE PRESSURE CATHETER

| Description | Part Number | Quantity |
|--------------------------------------|-------------|----------|
| Intrauterine Pressure Catheter | 2108982-001 | 10 |
| Intrauterine Pressure Catheter Cable | 2108998-001 | 2 |

FETAL SPIRAL ELECTRODE ACCESSORIES

| Description | Type | Part Number | Quantity |
|------------------|----------------------------------|-------------|----------|
| Cable | Rectangluar Connector | 2109005-001 | 1 |
| | Round Connector | 2109004-001 | 1 |
| ATTACHMENT PAD | Large, Single Patient Use | 2109006-001 | 50 |
| SPIRAL ELECTRODE | Single Helix, Single Patient Use | 2109002-001 | 50 |



INFANT CARE PHOTOTHERAPY ACCESSORIES RESUSCITATION



RESUSCITATION

INFANT RESUSCITATOR T-PIECE CIRCUIT KIT



| Description | Mask Size | Part Number | Quantity |
|------------------------------|-----------|-------------|----------|
| Disposable T-Piece with Mask | 0 | 2109060-001 | 1 |
| | 1 | 2109061-001 | 1 |

T-Piece Circuit: Disposable resuscitator designed for neonates, provides accurate reliable ventilation and pulmonary resuscitation.

Product compatibility

With a universal connector, it fits with all commonly used manual resuscitator units and does not require an adapter. Specifications available on request



Panda™ Warmer



StAr

LULLABY™ WARMER

TEMPERATURE PROBES

| Description | Part Number | Quantity |
|---|---------------|----------|
| Disposable Temperature Probe | 6600-0208-700 | 1 |
| Reusable Temperature probe – Lullaby™ Warmer | 6600-0628-700 | 1 |

TEMPERATURE PROBE PATCHES

| Description | Part Number | Quantity |
|----------------------------|---------------|----------|
| Patch Heat Reflecting ROHS | 0203-1980-300 | 50 |

VENTILATION




CARESCAPE R860 VENTILATOR



FILTERS




| FILTERS | | | | ISG FILTER | | | |
|--|--|-------------|-----|--|--------------------------|-------------|-----|
| | Description | Part Number | Qty | | Description | Part Number | Qty |
|  | Uni-Filter | 557021200 | 1 |  | Inspiratory Safety Guard | 5514015 | 1 |
|  | HMEF 500, (Paediatric only) | 557070500 | 1 | | | | |
|  | HMEF 1000/S, with sample port (Adult only) | 557070100 | 1 | | | | |

PATIENT CIRCUITS

| PATIENT CIRCUITS | | | | PATIENT CIRCUITS | | | |
|--|--|-------------|-----|--|---|-------------|-----|
| | Description | Part Number | Qty | | Description | Part Number | Qty |
|  | Patient tube, Hytrel, adult | 587040700 | 1 |  | Patient circuit kit, paediatric, disposable, 1.5 m/60 in | M1019544 | 1 |
| | Patient tube Hytrel, paed | 427010700 | 1 | | | | |
|  | Patient circuit kit, adult, disposable, 1.5 m/60 in, | M1019541 | 1 | | | | |

GAS MODULES ACCESSORIES

GAS MODULES AND ACCESSORIES

| | Description | Part Number |
|---|------------------------------------|-------------|
|  | D-fend+ water trap, single use | 881319-HEL |
|  | D-fend Pro+ water trap, single use | M1200227 |
|  | Mini D-fend water trap | 8002174 |

GAS MODULES AND ACCESSORIES









| | Description | Part Number |
|--|---|-------------|
|  | Spirometry tube, disposable, 2 m/7ft | 8004381 |
|  | CO2 sampling line 2 m/7 ft | 73318 |
|  | Adult reusable D-lite sensor, 1/pk | 733910-HEL |
|  | Adult disposable D-lite+ sensor, for humid conditions | 896952 |
|  | Paediatric reusable Pedi-lite sensor, 1/pk | 73393 |
|  | Paediatric disposable Pedi-lite+ sensor, for humid conditions | 8001948 |

AEROGEN NEBULIZERS AND ACCESSORIES

AEROGEN NEBULIZERS AND ACCESSORIES



| | Description | Part Number |
|--|---|---------------|
|  | Aerogen Pro Nebulizer kit with adult and paed tee and cable | 1505-3846-000 |
|  | Nebulizer head with filler cap | AG-AP1000 |
|  | Nebulizer cable | 1505-5602-000 |
|  | Filler cap for Nebulizer | AG-AP1030 |
|  | Neonatal tee with silicone plug 12 mm/12 mm | AG-AP1035 |
|  | Nebulizer replacement kit, adult | AG-AP1100 |
|  | Nebulizer replacement kit, paediatric | AG-AP1200 |
|  | Conn Luer | AG-AS3400 |
|  | Silicone Plug (1/pk) | AG-AS3450 |

AEROGEN NEBULIZERS AND ACCESSORIES



| | Description | Part Number |
|--|---|-------------|
|  | Plug, Nebulizer tee | AG-AP1005 |
|  | Adult tee with plug | AG-AP1010 |
|  | Paediatric tee with plug | AG-AP1020 |
|  | Aerogen Solo Neb | AG-AS3100-S |
|  | Aerogen Solo Neb | AG-AS3200-S |
|  | Adult Nebulizer tee | AG-AS3010 |
|  | Paediatric Nebulizer tee | AG-AS3020 |
|  | Neonatal Nebulizer tee, 12 mm M/12 mm F | AG-AS3035 |

CABLES

CABLES


| | Description | Part Number | Qty |
|--|--|-------------|-----|
|  | Cable for Networking (Required to Connect ventilator to Networking cable) | 2070435-001 | 1 |
|  | Cable for Nurse Call (Required to Connect ventilator to Nurse Call cable) | 2070436-001 | 1 |

CABLES


| | Description | Part Number | Qty |
|--|-----------------------------|-------------|-----|
|  | Display Unit cable (2 m) | M1236496 | 1 |
|  | Display Unit cable (5 m) | M1236498 | 1 |

EXHALATION VALVE HEATER

EXHALATION VALVE HEATER

| | Description | Part Number | Qty |
|--|--|-------------|-----|
|  | Exhalation valve heater w/cable (includes M1188723) | M1200693 | 1 |

EXHALATION VALVE HEATER



| | Description | Part Number | Qty |
|--|-----------------------------------|-------------|-----|
|  | Cable for Exhalation valve heater | M1188723 | 1 |

SENSORS

SENSORS

| | Description | Part Number | Qty |
|--|---------------------------------------|---------------|-----|
|  | Exhalation valve kit with flow sensor | 1505-3848-000 | 1 |
|  | Exhalation valve kit | 1505-8568-000 | 1 |
|  | Expiratory flow sensor | 1505-3231-000 | 1 |
|  | Neonatal flow sensor | 1505-3272-000 | 1 |

SENSORS

| | Description | Part Number | Qty |
|--|---|---------------|-----|
|  | Neonatal flow sensor cable | 1505-5604-000 | 1 |
|  | Intratracheal pressure sensor, 10/pkg (For use with Spirodynamics) | M1045564 | 1 |