

Ready for anywhere

AMX Traveller is a commercial configuration of the Traveller.

gehealthcare.com

A compact system

111111

with a full range of motion

AMX Traveller is a highly functional, ultra-mobile X-ray system with all-around positioning flexibility.

This compact, mobile X-ray system fits an entire radiographic room in a small, 240 kg unit. Set it in motion with a light push. Then easily maneuver it through tight spots with its oversized rear wheels and swiveling front wheels. Need to make a quick stop? Release the handle to cautiously bring it to a standstill.

Once in place, the X-ray tube rotates on both the X-axis and the Z-axis. The entire column rotates as well. That's three levels of rotation for easy positioning in any situation.

AMX Traveller is small, yet highly functional. It comes with an easy-to-use APR preloaded with 32 different techniques. And with a 32 kW generator, it has the perfect amount of performance to do everything you need it to do.



Moves with just 4-6 kg of force



Three levels of rotation for easy positioning





Short, compact design fits in tight places



GE Healthcare is a leading global medical technology and digital solutions innovator, GE Healthcare enables clinicians to make faster, more informed decisions through intelligent devices, data analytics, applications and services, supported by its Edison intelligence platform. With over 100 years of healthcare industry experience and around 50,000 employees globally, the company operates at the center of an ecosystem working toward precision health, digitizing healthcare, helping drive productivity and improve outcomes for patients, providers, health systems and researchers around the world. Follow us on Facebook, LinkedIn, Twitter and Insights, or visit our website www.gehealthcare.com for more information.

©2020 General Electric Company - All rights reserved.

AMX Traveller is a commercial configuration of the Traveller.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. Get, the GE Monogram, and AMX are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

May 2020 JB79362XX