





Solving Flatfoot

Offering Straightforward Implant and Orthobiologics Solutions for Flatfoot Reconstruction



ExtremiLock™ Foot Plating System*

All-inclusive trauma and reconstruction plating system offering anatomical contoured implants, double-lead screws, and one of the largest variable-angle locking ranges on the market.



InstaFix™ Shape Memory Fixation System*

Body-temperature activated shape memory nitinol implant designed to provide continuous and bicortical compression across the surgical site.



Acutrak® 2 Headless Compression Screw System


First and only fully threaded continuously variable thread pitch screw allowing gradual, consistent compression along the entire length of the screw. Available in 6 different families.



ExtremiFix™ Cannulated Screw System | Mini & Small*

Double-lead partially threaded self-drilling and self-tapping screws available in headed and headless options. Offered in four diameters: 2.0 mm, 2.4 mm, 3.0 mm, and 4.0 mm, including various lengths ranging from 6 mm to 52 mm.

*OsteoMed products



ExtremiFix™ Cannulated Screw System | Midsize & Large*

Double-lead headed and headless partially threaded self-drilling and self-tapping screws available in three diameters: 4.5 mm, 5.5 mm, and 6.5 mm. Offered in lengths ranging from 20 mm to 120 mm in short and long threads.





Cannulated Screw System

Headed partially and fully threaded self-drilling and self-tapping screws, available in three diameters: 4.0 mm, 6.5 mm, and 7.3 mm. Offered in lengths ranging from 30 mm to 150 mm in short and long threads.



Talar-Fit Subtalar Arthroereisis Implant*

Anatomical cannulated cone-shaped implant with blunt, deep threads designed to provide superior blockage of the talus and minimize implant extrusion. Available in multiple sizes to accommodate various patient anatomies.



DBM Putty*

100% bone derived, no inert carrier. Multiple sizes in putty and crunch forms.

PrimaGraft® Evans/Cotton Wedges*

Bicortical wedges designed specifically to address the most common angles for both Evans and Cotton osteotomies.

*OsteoMed products