



RT2 TENDER 2021 ELECTROSURGICAL UNITS

01 DECEMBER 2021
30 NOVEMBER 2024
AWARDS TO MEDHOLD MEDICAL

CONTACTS

- | | |
|--|-------------------------------------|
| • HEAD OFFICE JOHANNESBURG: 011 966 0600 | • GQEBERHA OFFICE: 041 364 0168 |
| • CAPE TOWN OFFICE: 021 976 7499 | • EAST LONDON OFFICE: 043 748 4986 |
| • DURBAN OFFICE: 031 764 3638 | • BLOEMFONTEIN OFFICE: 051 406 2400 |

TENDER - RT2 – 06 – 02

TENDER - RT2-06-02 ELECTROSURGICAL UNIT BASIC

ELECTROSURGICAL UNIT BASIC - WITH STARTER PACK

RT2-2021 ITEM CODE	MEDHOLD CODE	DESCRIPTION
RT2-06-02 SAP NO _ 151185	602450230	ELECTROSURGERY UNIT, GENERAL SURGERY
	OBS6090	FOOTSWITCH, DOUBLE FOR ELECTROSURGERY UNIT, GENERAL
	R1031A	ELECTRODE, REUSABLE STAINLESS STEEL OR CONDUCTIVE RUBBER BIPOLAR PATIENT ELECTRODE SURGERY
	OBS6090	CABLE, PATIENT ELECTRODE SURGERY
	0048R	ELECTRODE HANDLES, MONOPOLAR REUSABLE, FOOT SWITCH CONTROLLED, ELECTRODE WITH CABLES, X2
	0049R	ELECTRODE HANDLES, MONOPOLAR REUSABLE, HANDSWITCH SWITCH CONTROLLED, ELECTRODE WITH CABLES, X2
	OBSNBBL	SET MONOPOLAR REUSABLE ELECTRODES X1
	0550-08	NEEDLE, BLADE, BALL, LOOP
	0550-08	BIPOLAR FORCEPS, REUSABLE, FOOTSWITCH, WITHN CONNECTING CABLE X1
	60-715-2001	DEDICATED TROLLEY



ELECTROSURGICAL UNIT BASIC




SMOKE EVACUATION UNIT

OTHER CONSUMABLES OR ACCESSORIES FOR THE MODEL CONMED 2450 ON OFFER

RT2-2021 ITEM CODE	MEDHOLD CODE	DESCRIPTION
RT2-06-02N	VC220	SMOKE EVACUATION UNIT
RT2-06-0AT	PLP2020	PLUMEPEN ULTRA SURGICAL SMOKE EVACUATION PENCIL QTY1
RT2-06-02AW	VS1352	SMOKE EVACUATION FILTER - VISICLEAR
RT2-06-02AL	60-8084-001	FILTER ASSEMBLY (2 PER CASE) – AER DEFENCE
RT2-06-02BF	OBS2031C	DUAL PLATE WITH CABLE, ADULT AND PEAD
RT2-06-02Z	410-2200	ELECTROSURGERY UNIT, GENERAL SURGERY C18:C39
RT2-06-02AA	425-2508	MACROLYTE DUAL PAED ELECTRODE WITH CABLE
RT2-06-02AB	440-2400	PREEMIE DISPERSIVE ELECTRODE, EACH
RT2-06-02BG	HS102B	HAND SWITCHING PENCIL AND BLADE
RT2-06-02BH	HS102BT	HAND SWITCHING PENCIL AND BLADE (TEFLON COATED)
RT2-06-02X	131309A	HANDSWITCHING PENCIL ULTRACLEAN (131309)
RT2-06-02BH	LB102	ELECTRODE BLADE, STRAIGHT, LONG
RT2-06-02BI	LN102	ELECTRODE NEEDLE, STRAIGHT, LONG
RT2-06-02N	MLBA102	ELECTRODE LONG BALL,
RT2-06-0AT	OBS0301	ELECTRODE BALL, STRAIGHT, 2.5MM DIAMETER
RT2-06-02AW	OBS0302	ELECTRODE BALL, STRAIGHT, 4MM DIAMETER
RT2-06-02AL	OBS0320	ELECTRODE, MICROPOINT, 70MM
RT2-06-02BF	OBS0326C	BLADE ELECTRODE NON STICK 100MM DISP
RT2-06-02Z	OBS0327C	NEEDLE ELECTRODE NON STICK 100MM DISP
RT2-06-02AA	OBS0341	LOOP ELECTRODE, ROUND, 10X5X130MM
RT2-06-02AB	OBS0343	LOOP ELECTRODE, ROUND, 15X15X130MM
RT2-06-02BG	OBS0353	LOOP ELECTRODE, ROUND, 30X12X130MM
RT2-06-02BH	OBS0360	LOOP ELECTRODE, ROUND, 10X5X75MM
RT2-06-02X	OBS0362	LOOP, ROUND, 15X15X86MM
RT2-06-02BH	OBS0364	LOOP ELECTRODE, ROUND, 30X12X88MM
RT2-06-02BI	SCRATCHPAD	SCRATCH PAD
RT2-06-02BI	SN102	ELECTRODE NEEDLE, SHORT
RT2-06-02BI	SW111	FOOT CONTROL ESU PENCIL WITH BLADE



ELECTROSURGICAL PLATE



SMOKE EVACUATION PENCILS


 ELECTROSURGICAL
HANDSWITCHING PENCILS

 ELECTROSURGICAL
HANDSWITCHING ELECTRODES

TENDER - RT2 - 06 - 04

TENDER - RT2-06-04 ELECTROSURGICAL UNIT ADVANCED

ELECTROSURGICAL UNIT BASIC - WITH STARTER PACK

RT2-2021 ITEM CODE	MEDHOLD CODE	DESCRIPTION
RT2-06-04 SAP NO _ 182825	608015SYS	ELECTROSURGERY UNIT, GENERAL SURGERY
	60-510-3001	FOOTSWITCH, DOUBLE FOR ELECTROSURGERY UNIT, GENERAL SURGERY
	R1031A	ELECTRODE, REUSABLE STAINLESS STEEL OR CONDUCTIVE RUBBER BIPOLAR PATIENT ELECTRODE
	OBS6090	CABLE, PATIENT ELECTRODE
	SW111	ELECTRODE HANDLES, MONOPOLAR REUSABLE, FOOT SWITCH CONTROLLED, ELECTRODE WITH CABLES, X2
	0048R	ELECTRODE HANDELS WITH CABLE, REUSABLE AND STERILIZABLE
	0049R	BLADE KNIFE ELECTRODES STRAIGHT REUSABLE X2
	LN102	NEEDLE ELECTRODE THICK X2
	OBS0320	NEEDLE ELECTRODE THIN X2
	OBS0301	BALL ELECTRODE 4MM DIAMETER X2
	OBS0302	BALL ELECTRODE 2.5MM DIAMETER X1
	0550-08	BIPOLAR ELECTRODE REUSABLE STERILIZABLE X2



ELECTROSURGICAL UNIT ADVANCED




SMOKE EVACUATION UNIT

OTHER CONSUMABLES OR ACCESSORIES FOR THE MODEL CONMED 5000 ON OFFER

RT2-2021 ITEM CODE	MEDHOLD CODE	DESCRIPTION
RT2-06-04N	VC220	SMOKE EVACUATION UNIT _ VISICLEAR
RT2-06-04AP	PLP2020	PLUMEPEN ULTRA SURGICAL SMOKE EVACUATION PENCIL QTY1
RT2-06-04AV	VS1352	SMOKE EVACUATION FILTER - VISICLEAR
RT2-06-05AQ	60-8084-001	FILTER ASSEMBLY (2 PER CASE) – AER DEFENCE
RT2-06-04BF	OBS2031C	DUAL PLATE WITH CABLE, ADULT AND PEAD
RT2-06-04Y	410-2200	SUREFIT DUAL DISPERSIVE PLATE ADULT WITH CABLE
RT2-06-04Z	425-2508	MACROLYTE DUAL PAED ELECTRODE WITH CABLE
RT2-06-04AA	440-2400	PREEMIE DISPERSIVE ELECTRODE, EACH
RT2-06-04BG	HS102B	HAND SWITCHING PENCIL AND BLADE
RT2-06-04BH	HS102BT	HAND SWITCHING PENCIL AND BLADE (TEFLON COATED)
RT2-06-04BI	131309A	HANDSWITCHING PENCIL ULTRACLEAN (131309)
RT2-06-04BI	LB102	ELECTRODE BLADE , STRAIGHT, LONG
RT2-06-04BJ	LN102	ELECTRODE NEEDLE , STRAIGHT, LONG
RT2-06-04BK	MLBA102	ELECTRODE LONG BALL,
RT2-06-04BL	OBS0301	ELECTRODE BALL, STRAIGHT, 2.5MM DIAMETER
RT2-06-04BM	OBS0302	ELECTRODE BALL, STRAIGHT, 4MM DIAMETER
RT2-06-04BN	OBS0320	ELECTRODE, MICROPOINT, 70MM
RT2-06-04BO	OBS0326C	BLADE ELECTRODE NON STICK 100MM DISP
RT2-06-04BP	OBS0327C	NEEDLE ELECTRODE NON STICK 100MM DISP
RT2-06-04BQ	OBS0341	LOOP ELECTRODE, ROUND, 10X5X130MM
RT2-06-04BR	OBS0343	LOOP ELECTRODE, ROUND, 15X15X130MM
RT2-06-04BS	OBS0353	LOOP ELECTRODE, ROUND, 30X12X130MM
RT2-06-04BT	OBS0360	LOOP ELECTRODE, ROUND, 10X5X75MM
RT2-06-04BU	OBS0362	LOOP ELECTRODE, ROUND, 15X15X86MM
RT2-06-04BV	OBS0364	LOOP ELECTRODE, ROUND, 30X12X88MM
RT2-06-04BX	SCRATCHPAD	SCRATCH PAD
RT2-06-04BY	SN102	ELECTRODE NEEDLE , SHORT
RT2-06-04BZ	SW111	FOOT CONTROL ESU PENCIL WITH BLADE



ELECTROSURGICAL PLATE



SMOKE EVACUATION PENCILS


 ELECTROSURGICAL
HANDSWITCHING PENCILS

 ELECTROSURGICAL
HANDSWITCHING ELECTRODES

TENDER - RT2 – 06 – 05

TENDER - RT2-06-05 ELECTROSURGICAL UNIT ARGON

ELECTROSURGICAL UNIT BASIC -WITH STARTER PACK ARGON

RT2-2021 ITEM CODE	MEDHOLD CODE	DESCRIPTION
RT2-06-05 HELIX AR_220454	608800SYS	ELECTROSURGERY UNIT, GENERAL SURGERY ARGON
	60-510-3001	FOOTSWITCH, DOUBLE FOR ELECTROSURGERY UNIT, GENERAL SURGERY
	R1031A	ELECTRODE, REUSABLE STAINLESS STEEL OR CONDUCTIVE RUBBER BIPOLAR PATIENT ELECTRODE
	OBS6090	CABLE, PATIENT ELECTRODE
	SW111	ELECTRODE HANDLES, MONOPOLAR REUSABLE, FOOT SWITCH CONTROLLED, ELECTRODE WITH CABLES, X2
	0048R	ELECTRODE HANDELS WITH CABLE, REUSABLE AND STERILIZABLE
	0049R	BLADE KNIFE ELECTRODES STRAIGHT REUSABLE X2
	LN102	NEEDLE ELECTRODE THICK X2
	OBS0320	NEEDLE ELECTRODE THIN X2
	OBS0301	BALL ELECTRODE 4MM DIAMETER X2
	OBS0302	BALL ELECTRODE 2.5MM DIAMETER X1
	0550-08	BIPOLAR ELECTRODE REUSABLE STERILIZABLE X2
	134006	FLEXIBLE ARGON PROBE
	HU100	ARGON ELECTRODE CONNECTION



ELECTROSURGICAL UNIT ARGON



SMOKE EVACUATION UNIT

OTHER CONSUMABLES OR ACCESSORIES FOR THE MODEL HELIX AR ON OFFER


RT2-2021 ITEM CODE	MEDHOLD CODE	DESCRIPTION
RT2-06-05CS	130321	ABC TRIPLE HANDCONTROL HANDPIECE, MONOPOLAR
RT2-06-05CT	133023	ABC FLEXIBLE PROBE 220CM
RT2-06-05CU	133270	ABC FLEXIBLE PROBE 270CM
RT2-06-05CV	134006	ABC BEND-A-BEAM HANDPIECE, HANDC, MONOP
RT2-06-05CW	160655	ABC LAP PROBE 5MM, FOOTCONTROL
RT2-06-05CX	160656	ABC LAPARAROSCOPIC PROBE 5MM, HANDSWITCH
RT2-06-05R	VC220	SMOKE EVACUATION UNIT _ VISICLEAR
RT2-06-05AT	PLP2020	PLUMEPEN ULTRA SURGICAL SMOKE EVACUATION PENCIL QTY1
RT2-06-05AZ	VS1352	SMOKE EVACUATION FILTER - VISICLEAR
RT2-06-05AQ	60-8084-001	FILTER ASSEMBLY (2 PER CASE) – AER DEFENCE
RT2-06-05BH	OBS2031C	DUAL PLATE WITH CABLE, ADULT AND PEAD
RT2-06-05AC	410-2200	SUREFIT DUAL DISPERSIVE PLATE ADULT WITH CABLE
RT2-06-05AD	425-2508	MACROLYTE DUAL PAED ELECTRODE WITH CABLE
RT2-06-05AE	440-2400	PREEMIE DISPERSIVE ELECTRODE, EACH
RT2-06-05AF	HS102B	HAND SWITCHING PENCIL AND BLADE
RT2-06-05BI	HS102BT	HAND SWITCHING PENCIL AND BLADE (TEFLON COATED)
RT2-06-05BJ	131309A	HANDSWITCHING PENCIL ULTRACLEAN (131309)
RT2-06-05BK	LB102	ELECTRODE BLADE , STRAIGHT, LONG
RT2-06-05BL	LN102	ELECTRODE NEEDLE , STRAIGHT, LONG
RT2-06-05BM	MLBA102	ELECTRODE LONG BALL,
RT2-06-05BN	OBS0301	ELECTRODE BALL, STRAIGHT, 2.5MM DIAMETER
RT2-06-05BO	OBS0302	ELECTRODE BALL, STRAIGHT, 4MM DIAMETER
RT2-06-05BP	OBS0320	ELECTRODE, MICROPOINT, 70MM
RT2-06-05BN	OBS0326C	BLADE ELECTRODE NON STICK 100MM DISP
RT2-06-05BQ	OBS0327C	NEEDLE ELECTRODE NON STICK 100MM DISP
RT2-06-05BR	OBS0341	LOOP ELECTRODE, ROUND, 10X5X130MM
RT2-06-05BS	OBS0343	LOOP ELECTRODE, ROUND, 15X15X130MM
RT2-06-05BU	OBS0353	LOOP ELECTRODE, ROUND, 30X12X130MM
RT2-06-05BV	OBS0360	LOOP ELECTRODE, ROUND, 10X5X75MM
RT2-06-05BW	OBS0362	LOOP ELECTRODE, ROUND, 15X15X86MM
RT2-06-05BX	OBS0364	LOOP ELECTRODE, ROUND, 30X12X88MM
RT2-06-05BY	SCRATCHPAD	SCRATCH PAD
RT2-06-05BZ	SN102	ELECTRODE NEEDLE , SHORT
RT2-06-05CA	SW111	FOOT CONTROL ESU PENCIL WITH BLADE

TENDER - RT2 - 06 - 05


TENDER - RT2-06-05 ELECTROSURGICAL UNIT ARGON

OTHER CONSUMABLES OR ACCESSORIES FOR THE MODEL HELIX AR ON OFFER

ELECTROSURGICAL PLATE



SMOKE EVACUATION PENCILS



Argon Beam
Coagulation (ABC)[®]



ELECTROSURGICAL ARGON BEAM ACCESSORIES



ELECTROSURGICAL HANDSWITCHING PENCILS

SYSTEM 2450™ ELECTROSURGICAL GENERATOR

The System 2450™ provides state-of-the-art innovation in an easy-to-use format. The System 2450 is simple on the outside, yet sophisticated and powerful on the inside.

FEATURES:

- Dynamic response technology.
- Continuous ARM safety monitoring.
- Automatic programming restores the System 2450 to the last settings used.
- One hand switch/foot switch combination jack and a separate hand switch jack allows for multiple accessory connections.
- Serial interface port enables advanced communications between the System 2450 and compatible robotic surgical systems or O.R. integration systems.
- Integrated interface for activation of a smoke evacuator.



SYSTEM 5000™ ELECTROSURGICAL GENERATOR

The System 5000™ electrosurgery generator (ESU) unites technology with versatility providing the opportunity to fine tune the electrosurgical generator to specific procedure and surgeon preferences. Unique modes on the System 5000™ include the fluids mode, for an additional boost of energy in procedures taking place in a fluid environment, and lap mode, which limits open circuit voltages to below 2700 volts in all modes and reduces the potential risk of inadvertent burns from capacitive coupling.



HELIXAR™ ELECTROSURGICAL GENERATOR WITH ARGON BEAM COAGULATION (ABC®)



The HelixAR™ combines the advanced specialty modes of a premium ESU with the benefits of CONMED's latest ABC® technology. The HelixAR is intended for both open and laparoscopic procedures, and includes:

NEW MONOPOLAR SPECIALTY MODES

- Laparoscopic: limits voltage to 2700v in lap mode, helping reduce the risks of inadvertent burns from capacitive coupling.
- Fluids: rapid energy initiation for immediate clinical effect in high load/low impedance conditions.
- Pulse coag: pulsing bursts of coagulation for controlled hemostasis with reduced tissue carbonization.

NEW ABC® MODES

- Argon dissect: works with the conmed dissecting electrode to provide argon enhanced cutting with four varying levels of hemostasis.
- Pulse argon: provides controlled hemostasis for small surface bleeders with four pulse timings.

NEW ABC® MODES

- Argon dissect: works with the conmed dissecting electrode to provide argon enhanced cutting with four varying levels of hemostasis.
- Pulse argon: provides controlled hemostasis for small surface bleeders with four pulse timings.

VISICLEAR® SURGICAL SMOKE EVACUATOR



Visiclear® is our most advanced acute care offering for surgical smoke evacuation. Combining Safeport Technology™, an occlusion warning, filter life tracking, and one-touch programmable procedure modes for open tubing, electrosurgical pencil, and laparoscopic procedures, visiclear effectively captures and filters surgical smoke, removes odours, particulates and other potentially hazardous by-products. Ezlink®, an automatic activation device, can be used to remotely turn the evacuator on and off, saving filter life. Virosafe® fluid trap is designed to capture fluids to help maintain the integrity of the filtration system and assist with maintaining filter life.

PROCEDURE OPERATING MODES

- One of three operating modes (laparoscopic, open tubing, and electrosurgical pencil mode) can easily be selected, automatically customizing filter life based on surgical procedure.

ADVANCED SOUND CONTROL DESIGN

- Whisper Technology™ ensures ultra-quiet operation with unprecedented sound quality.

SLEEK SHELF DESIGN

- Ergonomically designed to accommodate stacking.

OCCLUSION WARNING

- Safety feature that alerts medical personnel when an occlusion has occurred.

INTUITIVE COLOUR TOUCH SCREEN

- Easy to use, requires minimal training, and can be operated with gloved hand.

UNIQUE SLIDING PORT CLOSURE

- Safeport Technology™ Ensures reduced risk of staff coming in contact with a contaminated port

Clear the Air[®] SMOKE MANAGEMENT

START YOUR SMOKE FREE JOURNEY TODAY



CONMED

**SMOKE EVACUATION
SYSTEMS**

PENCILS



LAPAROSCOPIC



Why is Surgical Smoke so Hazardous?



CHEMICAL

150 CHEMICALS
AND 40
CARCINOGENIC
CHEMICALS (3)



BIOLOGICAL

VIABLE
BACTERIA
AND VIRUS (4)



PHYSICAL

ULTRA-FINE
PARTICLES
SIMILAR TO
COAL DUST (5)

Whats the impact?

ONE DAY OF SURGICAL SMOKE PRODUCED
IN THE OR IS EQUIVALENT TO

27-30 CIGARETTES



**THERE ARE
DANGERS IN SURGICAL
SMOKE PLUME**

"TOGETHER, WE CAN HELP
ELIMINATE THEM."



MEDHOLD

People. Products. Partnerships.



RT2 TENDER 2021 ANAESTHESIA AND VENTILATION ACCESSORIES

01 DECEMBER 2021
30 NOVEMBER 2024
AWARDS TO MEDHOLD MEDICAL

ANAESTHESIA EQUIPMENT ACCESSORIES

OTHER CONSUMABLES OR ACCESSORIES FOR THE ANAESTHESIA DELIVERY SYSTEM		
ITEM CODE		DESCRIPTION
RT2-01-01A	M1210695	ANAESTHESIA CIRCUIT, DISPOSABLE. 1.8M, COMPLETE WITH Y-PIECE, ELBOW, CATHETER MOUNT, REBREATHING BAG, EACH
RT2-01-01C	2105488-001	ANAESTHESIA CIRCUIT NEONATAL, SINGLE/DOUBLE LIMB DESIGN) DISPOSABLE. 1.8M COMPLETE WITH Y-PIECE, ELBOW, CATHETER MOUNT, REBREATHING BAG, BOX OF 20
RT2-01-01D	M1019527	RE-BREATHING CIRCUIT COMPLETE, ADULT, BOX OF 20
RT2-01-01D	2105488-021	RE-BREATHING CIRCUIT COMPLETE, PAEDIATRIC, EACH
RT2-01-01E	M1047747S	NON-RE-BREATHING CIRCUIT COMPLETE, EACH
RT2-01-01I	4800VS	FACE MASKS: SIZES 0, EACH
RT2-01-01I	4810VS	FACE MASKS: SIZES 1, EACH
RT2-01-01I	4820VS	FACE MASKS: SIZES 2, EACH
RT2-01-01I	4830VS	FACE MASKS: SIZES 3, EACH
RT2-01-01I	4840VS	FACE MASKS: SIZES 4, EACH
RT2-01-01I	4851VS	FACE MASKS: SIZES 5, EACH
RT2-01-01I	4861VS	FACE MASKS: SIZES 6, EACH
RT2-01-01J	557070100	BACTERIAL FILTER, EACH
RT2-01-01J	557022500	VIRAL FILTER, EACH
RT2-01-03A	2096536-003	ANAESTHESIA CIRCUIT, DISPOSABLE. 1.8M, COMPLETE WITH Y-PIECE, ELBOW, CATHETER MOUNT, REBREATHING BAG, PK OF 2
RT2-01-03C	2105488-001	ANAESTHESIA CIRCUIT NEONATAL, SINGLE/DOUBLE LIMB DESIGN) DISPOSABLE. 1.8M COMPLETE WITH Y-PIECE, ELBOW, CATHETER MOUNT, REBREATHING BAG, BOX OF 20
RT2-01-03D	M1019527	RE-BREATHING CIRCUIT COMPLETE, ADULT, BOX OF 20
RT2-01-03D	2105488-021	RE-BREATHING CIRCUIT COMPLETE, PAEDICATRIC, EACH
RT2-01-03E	M1047747S	NON-RE-BREATHING CIRCUIT COMPLETE, EACH
RT2-01-01I	4800VS	FACE MASKS: SIZES 0, EACH
RT2-01-01I	4810VS	FACE MASKS: SIZES 1, EACH
RT2-01-01I	4820VS	FACE MASKS: SIZES 2, EACH
RT2-01-01I	4830VS	FACE MASKS: SIZES 3, EACH
RT2-01-01I	4840VS	FACE MASKS: SIZES 4, EACH
RT2-01-01I	4851VS	FACE MASKS: SIZES 5, EACH
RT2-01-01I	4861VS	FACE MASKS: SIZES 6, EACH
RT2-01-03J	557070100	BACTERIAL FILTER, EACH
RT2-01-03J	557022500	VIRAL FILTER, EACH

ANAESTHESIA EQUIPMENT ACCESSORIES

OTHER CONSUMABLES OR ACCESSORIES FOR THE ANAESTHESIA DELIVERY SYSTEM		
ITEM CODE		DESCRIPTION
RT2-01-05A OFFER 1	2096536-003	ANAESTHESIA CIRCUIT, DISPOSABLE. 1.8M ,COMPLETE WITH Y-PIECE, ELBOW, CATHETER MOUNT, REBREATHING BAG, PK OF 2
RT2-01-05A OFFER 2	2096536-002 (ABS BAG ARM)	ANAESTHESIA CIRCUIT, DISPOSABLE. 1.8M ,COMPLETE WITH Y-PIECE, ELBOW, CATHETER MOUNT, REBREATHING BAG, PACK OF 2
RT2-01-05A OFFER 3	2096536-002 (ABS BAG ARM)	ANAESTHESIA CIRCUIT, DISPOSABLE. 1.8M ,COMPLETE WITH Y-PIECE, ELBOW, CATHETER MOUNT, REBREATHING BAG, PK OF 2
RT2-01-05C OFFER 1,2, 3	2105488-001	ANAESTHESIA CIRCUIT NEONATAL, SINGLE/DOUBLE LIMB DESIGN) DISPOSABLE. 1.8M COMPLETE WITH Y-PIECE, ELBOW, CATHETER MOUNT, REBREATHING BAG, BOX OF 20
RT2-01-05D OFFER 1,2,3	M1019527	RE-BREATHING CIRCUIT COMPLETE, ADULT, BOX OF 20
RT2-01-05D OFFER 1,2,3	2105488-021	RE-BREATHING CIRCUIT COMPLETE, PAEDICATRIC, EACH
RT2-01-05E OFFER 1,2,3	M1047747S	NON-RE-BREATHING CIRCUIT COMPLETE, EACH
RT2-01-05I OFFER 1,2,3	4800VS	FACE MASKS: SIZES 0, EACH
RT2-01-05I OFFER 1,2,3	4810VS	FACE MASKS: SIZES 1, EACH
RT2-01-05I OFFER 1,2,3	4820VS	FACE MASKS: SIZES 2, EACH
RT2-01-05I OFFER 1,2,3	4830VS	FACE MASKS: SIZES 3, EACH
RT2-01-05I OFFER 1,2,3	4840VS	FACE MASKS: SIZES 4, EACH
RT2-01-05I OFFER 1,2,3	4851VS	FACE MASKS: SIZES 5, EACH
RT2-01-05I OFFER 1,2,3	4861VS	FACE MASKS: SIZES 6, EACH
RT2-01-05J OFFER 1,2,3	557070100	BACTERIAL FILTER, EACH
RT2-01-05J OFFER 1,2,3	557022500	VIRAL FILTER, EACH
RT2-01-05L OFFER 2,3	M1019517	MANUAL VENTILATION HOSE, EACH
RT2-02-02G	557022500	BACTERIA FILTER, PK OF 5
RT2-02-02I	M1019541	ADULT DISPOSABLE PATIENT CIRCUIT, PK OF 10
RT2-02-02J	M1019544	PAEDS DISPOSABLE PATIENT CIRCUITS, PK OF 10



CONTACTS



WWW.MEDHOLD.CO.ZA

INFO@MEDHOLD.CO.ZA

FOR MORE INFORMATION CONTACT:

National Sales Manager: Yolandi Berkhout
082 854 7615

Product Manager – Airway: Karen Coles
076 430 1957

Product Manager – ESU: Shelley Olivier
071 702 2251

HEAD OFFICE GAUTENG

Telephone: +27 (11) 966 0600
General Fax: +27 (11) 975 3870
MSI Business Park,
68 Rigger Rd, Spartan,
Kempton Park,
1619

EAST LONDON

Branch Manager: Francois Fourie
Telephone: +27 (043) 748-4986
Fax: +27 (043) 748-4981
Suite No. 18 Beacon Park
89 Beacons Road
East London
5252

GQEBERHA

Branch Manager: Francois Fourie
Telephone: +27 (041) 364-0168
Fax: +27 (041) 365-2591
27 Shirley Street
Newton Park
Port Elizabeth
6025

FREE STATE

Branch Manager: Dawid Venter
Telephone: +27 (051) 406-2400
Fax: +27 (051) 448-9319
Palm Park Suite 4&5
94 Kellner Street
Westdene
Bloemfontein
9301

KWA-ZULU NATAL

Branch Manager: Wayne Gimingham
Telephone: +27 (031) 764-3638
Fax: +27 (031) 764-4891
78 Old Main Road,
Kloof
Durban
3610

WESTERN CAPE

Branch Manager: Hayley Bloomer
Telephone: +27 (021) 976-7499
Fax2mail: 087 231 4668
8 Langenhoven Street
Durbanville
7550