

RT2-2021

**SUPPLY, DELIVERY,
INSTALLATION, COMMISSIONING
AND MAINTENANCE OF
MEDICAL EQUIPMENT.**

**TO THE STATE FOR THE PERIOD
01 DECEMBER 2021
TO 01 DECEMBER 2024.**

CONTACT DETAILS

- HEAD OFFICE JOHANNESBURG: 011 966 0600
- CAPE TOWN OFFICE: 021 976 7499
- DURBAN OFFICE: 031 764 3638
- GQEBERHA OFFICE: 041 364 0168
- EAST LONDON OFFICE: 043 748 4986
- BLOEMFONTEIN OFFICE: 051 406 2400

INDEX

TENDERS

CONTACT: Liza Du Preez l.dupreez@medhold.co.za
Head Office: 011 966 0633

SECTION A - ANAESTHESIA

CONTACT: Filipe Caldeira f.caldeira@medhold.co.za
Head Office: 011 966 0600
Or contact your local branch – details on the back page

SECTION B - CARDIOLOGY

CONTACT: Dave Evans d.evans@medhold.co.za
Head Office: 011 966 0600
Or contact your local branch – details on the back page

SECTION C - DIGITAL SOLUTIONS

CONTACT: Dee-Jay Duvenage d.duvenage@medhold.co.za
Head Office: 011 966 0600
Or contact your local branch – details on the back page

SECTION D - INFECTION CONTROL

CONTACT: Ina Buitenbos i.buitenbos@medhold.co.za
Head Office: 011 966 0600
Or contact your local branch – details on the back page

SECTION E – MATERNAL INFANT CARE & VENTILATION

CONTACT: Natasha Kraamwinkel n.kraamwinkel@medhold.co.za
Head Office: 011 966 0600
Or contact your local branch – details on the back page

SECTION F – PATIENT MONITORING

CONTACT: Johan Roodt j.roodt@medhold.co.za
Head Office: 011 966 0600
Or contact your local branch – details on the back page

SECTION G – THEATRE TABLES & LIGHTS

CONTACT: Mkhusele Mabindla m.mabindla@medhold.co.za
Head Office: 011 966 0600
Or contact your local branch – details on the back page

RT2-2021

**SUPPLY, DELIVERY, INSTALLATION,
COMMISSIONING AND MAINTENANCE OF
MEDICAL EQUIPMENT**

**TO THE STATE FOR THE PERIOD
01 DECEMBER 2021
TO 01 DECEMBER 2024.**

SECTION A

ANAESTHESIA


CONTACT DETAILS

Sibongile Leeu: s.leeu@medhold.co.za
Head Office: 011 966 0600

Or contact your local branch – details on the back page

ANAESTHESIA

ITEM NUMBER: RT2-01-01 OFFER 1
ITEM NAME: ANAESTHETIC MACHINE , BASIC WITHOUT PHYSIOLOGICAL MONITOR

Make/Brand	GE Healthcare	
Model	9100C NXT	
Medhold Part Number	M9100175	
DOH SAP number	202777	
Country of Origin	China	
UNIT INCLUDES: <ul style="list-style-type: none">• Main Anaesthesia workstation• Guided Checklist and calibration• 3 Gasses O2, N2O , AIR• 1 x EZchange CO2 bypass system enables absorbent canister replacement during prolonged surgeries without introducing leaks into the breathing system• condenser unit		
VENTILATION OPERATING MODES: <ul style="list-style-type: none">• VCV, PCV, SIMV, PSVPro™ with apnea backup• Treat neonates to adult patients• Tidal Volume monitoring as low as 5ml in pressure control ventilation and delivery as low as 20ml in volume control ventilation		
EXTRA SOFTWARE OPTIONS INCLUDED: <ul style="list-style-type: none">• Spirometry loops - ability to save and compare reference loops for continuous monitoring of ventilatory status• SIMV & PSVPro ventilation modes• Inspiratory O2 Monitoring on the anaesthetic unit.		
UNIT PRICE INCLUSIVE OF VAT AS PER RT2 CIRCULAR AND ATTACHED PRICELIST See attached pricing sheet Offer price Includes: Start-up accessories & consumables as per tender spec, Plus Yr 1 - 2 Warranty and Preventative Maintenance.		
		VAPOURISERS INCLUDED: <ul style="list-style-type: none">• 1 x Isoflurane Vapouriser• 1 x Sevoflurane Vapouriser
		CONSUMABLES AND ACCESSORIES INCLUDED: <ul style="list-style-type: none">• Dual flow sensors• 1 x Circuit startup kit• 1 x Reusable Absorber canister• 1 x Pack Foam Dust Filters• 4 x Castor Guards / Cable Guard• Manual Breathing Bags 1 to 3L• Paediatric circuit• Neonatal circuit• Jackson Rees circuit• Face masks size 0 to 5• Adult & Paediatric Resuscitator 550ml /1500ml• Filter

OPTIONAL ITEMS
Extended Warranty and Maintenance Contract - 5 years

 Optional Year 3 to 7 Warranty & Preventative Maintenance Including VAT **includes flow sensors & O₂ cells

WEB SITE LINK 9100C NXT
[HTTPS://TINYURL.COM/2P8MDM4M](https://tinyurl.com/2P8MDM4M)

9100C NXT



ANAESTHESIA

ITEM NUMBER: RT2-01-01 OFFER 1 - 9100C NXT

ITEM NAME: ANAESTHETIC MACHINE , BASIC WITHOUT PHYSIOLOGICAL MONITOR

DOH Item Number	Replacement Accessories & Consumables List	Medhold Part Number	Quantity
RT2-01-01a	Anaesthesia circuit (Adult, single/double limb design) Disposable. 1.8m (Complete with Y-piece, Elbow, Catheter mount, Rebreathing bag)	M1210695	Each
RT2-01-01c	Anaesthesia circuit (neonatal, single/double limb design) Disposable. 1.8m (Complete with Y-piece, Elbow, Catheter mount, Rebreathing bag)	2105488-001	Box of 20
RT2-01-01d	Re-breathing circuit complete - Adult - (20 pcs)	M1019527	Box of 20
RT2-01-01d	Re-breathing circuit complete - Neonatal	2105488-021	Each
RT2-01-01e	Non-re-breathing circuit complete	M1047747S	Each
RT2-01-01f	Flow sensor (Pack of 2)	2089610-001-S	Pack of 2
RT2-01-01g	Spirometry line	891191-HEL	Pack 5
RT2-01-01h	Sample line (Pack of 10)	73319-HEL	Pack 10
RT2-01-01i	Face masks: sizes (00,0,1,2,3,4, 5)- End-user to specify size	Various	Each
RT2-01-01j	Bacterial/Viral filter - HMEF	557070100	Each
RT2-01-01j	Viral filter	557022500	Each
RT2-01-01k	Y connector - Adult	8570144	Each
RT2-01-01k	Y connector - Paediatric	8570145	Each
RT2-01-01l	Manual ventilation hose	M1019517	Each
RT2-01-01m	Fresh gas tubing	GAS001	Each
RT2-01-01n	Seal for back-up cylinder connection on the unit	OO5739 LOCAL	Each
RT2-01-01o	Scavenger tubing (as per institution's scavenging outlet)	GAS001	Each
RT2-01-01p	Inspiratory and expiratory unidirectional valves	M1167404-S	Each
RT2-01-01q	Side stream CO2 kit (adult/paed), Box of 10	M1182629 D-fend	Box of 10
RT2-01-01r	Disposable Adult Spirometry Set complete with spirometry tube, sample line and flow sensor, 3 m, Box of 20 pieces	8004382	Box of 20
RT2-01-01s	Disposable Paediatric Spirometry Set complete with spirometry tube, sample line and flow sensor, 3 m, Box of 20 pieces	M1032634	Box of 20
RT2-01-01t	Gas sample tube with luer lock connection, 3m, disposable (10 pcs)	73319-HEL	Pack of 10
RT2-01-01u	Suction Tubing Set inc connector - Serresa	SERRESA	Each
RT2-01-01u	Suction Tubing Set inc connector - Gabler	FSA100191111	Each
RT2-01-01v	Resuscitator Adult, 1500 ml, autoclavable, silicone	BE2101IA	Each
RT2-01-01w	Resuscitator Paediatric, 550 ml, autoclavable, silicone	BE2201IA	Each
RT2-01-01x	Tubing set Paediatric for circle systems	M1019544	Each
RT2-01-01y	Sevoflurane vaporizer	1177-9850-SEV	Each
RT2-01-01z	Isoflurane vaporizer	1177-9850-ISO	Each
RT2-01-01zz	Desflurane vaporizer	1107-9616-000	Each

ANAESTHESIA

ITEM NUMBER: RT2-01-02 OFFER 1

ITEM NAME: ANAESTHETIC MACHINE , BASIC WITH INTEGRATED PHYSIOLOGICAL MONITOR

Make/Brand	GE Healthcare
Model	9100C NXT B105 Datex Ohmeda / Trusignal
Medhold Part Number	M9100175
DOH SAP number	202777
Country of Origin	China

UNIT INCLUDES:

- Main Anaesthesia workstation
- Guided Checklist and calibration
- 3 Gasses O₂, N₂O, AIR
- 1 x EZchange CO₂ bypass system enables absorbent canister replacement during prolonged surgeries without introducing leaks into the breathing system condenser unit

VENTILATION OPERATING MODES:

- VVCV, PCV, SIMV, PSVPro™ with apnea backup
- Treat neonates to adult patients
- Tidal Volume monitoring as low as 5ml in pressure control ventilation and delivery as low as 20ml in volume control ventilation

EXTRA SOFTWARE OPTIONS INCLUDED:

- Spirometry loops - ability to save and compare reference loops for continuous monitoring of ventilatory status
- SIMV & PSVPro ventilation modes
- Inspiratory O₂ Monitoring on the anaesthetic unit.

B105 MONITOR: ADVANCED CAPABILITIES

- B1x5M range of monitors can be deployed seamlessly across a variety of care settings:
- EK-Pro v14 ECG 4-lead simultaneous arrhythmia analysis
- DINAMAP™ SuperSTAT non-invasive blood pressure measurement
- Choice of SpO₂ technologies: GE TruSignal™,
- GE EtCO₂ sidestream measurement, anesthetic agents and cardiac output
- Entropy™ monitoring* for monitoring the state of the brain
- NMT for neuromuscular transmission blockade and reversal monitoring
- Connectivity to GE CARESCAPE™ networks
- Flexibility to share parameter modules and accessories with
- CARESCAPE monitors



HAEMODYNAMIC ACCESSORIES RESPIRATORY MODULE:

- E-sCAIO Respiratory module included
- D-Fends pck
- Sample lines pck

Vapourisers included:

- 1 x Isoflurane Vapouriser
- 1 x Sevoflurane Vapouriser

Consumables and accessories included:

- Dual flow sensors
- 1 x Circuit startup kit
- 1 x Reusable Absorber canister
- 1 x Pack Foam Dust Filters
- 4 x Castor Guards / Cable Guard
- Manual Breathing Bags 1 to 3L
- Paediatric circuit
- Neonatal circuit
- Jackson Rees circuit
- Face masks size 0 to 5
- Adult & Paediatric Resuscitator 550ml / 1500ml
- Filter

UNIT PRICE INCLUSIVE OF VAT AS PER RT2 CIRCULAR AND ATTACHED PRICELIST

See attached pricing sheet

Offer price Includes: Start-up accessories & consumables as per tender spec, Plus Yr 1 - 2 Warranty and Preventative Maintenance.

OPTIONAL ITEMS

Extended Warranty and Maintenance Contract - 5 years
 Optional Year 3 to 7 Warranty & Preventative Maintenance
 Including VAT

**includes Flow sensors and Oxygen cells

WEB SITE LINK

9100C NXT

[HTTPS://TINYURL.COM/2P8MDM4M](https://tinyurl.com/2P8MDM4M)

9100C NXT

Scan me



ANAESTHESIA

ITEM NUMBER: RT2-01-02 (OFFER 1) 9100C NXT WITH B105 OHMEDA SpO2

DOH Item Number	Replacement Accessories & Consumables List	Medhold Part Number	Quantity
RT2-01-01a	Anaesthesia circuit (Adult, single/double limb design) Disposable. 1.8m (Complete with Y-piece, Elbow, Catheter mount, Rebreathing bag)	M1210695	Each
RT2-01-01c	Anaesthesia circuit (neonatal, single/double limb design) Disposable. 1.8m (Complete with Y-piece, Elbow, Catheter mount, Rebreathing bag)	2105488-001	Box of 20
RT2-01-01d	Re-breathing circuit complete - Adult - (20 pcs)	M1019527	Box of 20
RT2-01-01d	Re-breathing circuit complete - Neonatal	2105488-021	Each
RT2-01-01e	Non-re-breathing circuit complete	M1047747S	Each
RT2-01-01f	Flow sensor (Pack of 2)	2089610-001-S	Pack of 2
RT2-01-01g	Spirometry line	891191-HEL	Pack 5
RT2-01-01h	Sample line (Pack of 10)	73319-HEL	Pack 10
RT2-01-01i	Face masks: sizes (00,0,1,2,3,4, 5)- End-user to specify size	Various	Each
RT2-01-01j	Bacterial/Viral filter - HMEF	557070100	Each
RT2-01-01j	Viral filter	557022500	Each
RT2-01-01k	Y connector - Adult	8570144	Each
RT2-01-01k	Y connector - Paediatric	8570145	Each
RT2-01-01l	Manual ventilation hose	M1019517	Each
RT2-01-01m	Fresh gas tubing	GAS001	Each
RT2-01-01n	Seal for back-up cylinder connection on the unit	OO5739 LOCAL	Each
RT2-01-01o	Scavenger tubing (as per institution's scavenging outlet)	GAS001	Each
RT2-01-01p	Inspiratory and expiratory unidirectional valves	M1167404-S	Each
RT2-01-01q	Side stream CO2 kit (adult/paed), Box of 10	M1182629 D-fend	Box of 10
RT2-01-01r	Disposable Adult Spirometry Set complete with spirometry tube, sample line and flow sensor, 3 m, Box of 20 pieces	8004382	Box of 20
RT2-01-01s	Disposable Paediatric Spirometry Set complete with spirometry tube, sample line and flow sensor, 3 m, Box of 20 pieces	M1032634	Box of 20
RT2-01-01t	Gas sample tube with luer lock connection, 3m, disposable (10 pcs)	73319-HEL	Pack of 10
RT2-01-01u	Suction Tubing Set inc connector - Serresa	SERRESA	Each
RT2-01-01u	Suction Tubing Set inc connector - Gabler	FSA100191111	Each
RT2-01-01v	Resuscitator Adult, 1500 ml, autoclavable, silicone	BE2101IA	Each
RT2-01-01w	Resuscitator Paediatric, 550 ml, autoclavable, silicone	BE2201IA	Each
RT2-01-01x	Tubing set Paediatric for circle systems	M1019544	Each
RT2-01-01y	Sevoflurane vaporizer	1177-9850-SEV	Each
RT2-01-01z	Isoflurane vaporizer	1177-9850-ISO	Each
RT2-01-01zz	Desflurane vaporizer	1107-9616-000	Each

ANAESTHESIA

ITEM NUMBER: RT2-01-02 OFFER 2
ITEM NAME: ANAESTHETIC MACHINE , BASIC WITH PHYSIOLOGICAL MONITOR

Make/Brand	GE Healthcare
Model	9100C NXT B105 Nellcor
Medhold Part Number	M9100175
DOH SAP number	202777
Country of Origin	China

UNIT INCLUDES:

- Main Anaesthesia workstation
- Guided Checklist and calibration
- 3 Gasses O₂, N₂O , AIR
- 1 x EZchange CO₂ bypass system enables absorbent canister replacement during prolonged surgeries without introducing leaks into the breathing system
- condenser unit

VENTILATION OPERATING MODES:

- VCV, PCV, SIMV, PSVPro™ with apnea backup
- Treat neonates to adult patients
- Tidal Volume monitoring as low as 5ml in pressure control ventilation and delivery as low as 20ml in volume control ventilation

EXTRA SOFTWARE OPTIONS INCLUDED:

- Spirometry loops - ability to save and compare reference loops for continuous monitoring of ventilatory status
- SIMV & PSVPro ventilation modes
- Inspiratory O₂ Monitoring on the anaesthetic unit
- B105 Monitor:

UNIT PRICE INCLUSIVE OF VAT
AS per RT2 Circular and attached pricelist

Offer price Includes: Start-up accessories & consumables as per tender spec, Plus Yr 1 - 2 Warranty and Preventative Maintenance.

OPTIONAL ITEMS
Extended Warranty and Maintenance Contract - 5 years

Optional Year 3 to 7 Warranty & Preventative Maintenance Including VAT

 **includes flow sensors & O₂ cells

ADVANCED CAPABILITIES

- B1x5M range of monitors can be deployed seamlessly across a variety of care settings:
- EK-Pro v14 ECG 4-lead simultaneous arrhythmia analysis
- DINAMAP™ SuperSTAT non-invasive blood pressure measurement
- SpO₂ technology: Nelcor
- GE EtCO₂ sidestream measurement, anesthetic agents and cardiac output
- Entropy™ monitoring* for monitoring the state of the brain
- NMT for neuromuscular transmission blockade and reversal monitoring
- Connectivity to GE CARESCAPE™ networks
- Flexibility to share parameter modules and accessories with CARESCAPE monitors

HAEMODYNAMIC ACCESSORIES
RESPIRATORY MODULE:

- E-sCAIO Respiratory module included
- D-Fends pck
- Sample lines pck

VAPOURISERS INCLUDED:

- 1 x Isoflurane Vapouriser
- 1 x Sevoflurane Vapouriser

CONSUMABLES AND ACCESSORIES INCLUDED:

- Dual flow sensors
- 1 x Circuit startup kit
- 1 x Reusable Absorber canister
- 1 x Pack Foam Dust Filters
- 4 x Castor Guards / Cable Guard
- Manual Breathing Bags 1 to 3L
- Paediatric circuit
- Neonatal circuit
- Jackson Rees circuit
- Face masks size 0 to 5
- Adult & Paediatric Resuscitator 550ml /1500ml
- Filter

WEB SITE LINK 9100 NXT
[HTTPS://TINYURL.COM/2P8MDM4M](https://tinyurl.com/2P8MDM4M)

9100C NXT



ANAESTHESIA

ITEM NUMBER: RT2-01-02 (OFFER 2) 9100C NXT WITH B105 NELLCOR SpO2

DOH Item Number	Replacement Accessories & Consumables List	Medhold Part Number	Quantity
RT2-01-02a	A full set paediatric 3/5 lead ECG cable	2017003-003	Each
RT2-01-02b	A full set adult 3/5 lead ECG cable	2017003-003	Each
RT2-01-02c	Reusable paediatric cuff	DUR-P2-2A	Each
RT2-01-02d	Reusable adult cuff	DUR-A2-2A	Each
RT2-01-02e	Reusable neonate cuff	DUR-P1-2A	Each
RT2-01-02f	Reusable obese cuff	DUR-A3-2A	Each
RT2-01-02g	Air hose	2058203-002	Each
RT2-01-02h	Finger sensor for paed	D-YS	Each
RT2-01-02i	Finger sensor for adults	Flexmax	Each
RT2-01-02j	Earlobe sensor for paed	D-YS	Each
RT2-01-02k	Earlobe sensor for adults	D-YSE	Each
RT2-01-02l	Finger wrap-around probe for neonates	D-YS	Each
RT2-01-02m	Sensor lead	2021406-001	Each
RT2-01-02n	Reusable temperature probe	M1024247	Each

ANAESTHESIA

ITEM NUMBER: RT2-01-04 OFFER 8

ITEM NAME: ANAESTHESIA DELIVERY SYSTEM: STANDARD WITH INTEGRATED PHYSIOLOGICAL MONITOR

Make/Brand	GE Healthcare
Model	Carestation 650 with B650 Monitor and PDM / Nellcor
Medhold Part Number	1012-9650-000
DOH SAP number	220405
Country of Origin	China



UNIT INCLUDES:

- Main Anaesthesia workstation
- Guided Checklist and calibration
- 3 Gasses O₂, N₂O, AIR
- integrated EZchange CO₂ bypass system enables absorbent canister replacement during prolonged surgeries without introducing leaks into the breathing system
- Condensor unit integrated into Canister and compact breathing system

VENTILATION OPERATING MODES:

- VCV, PCV, SIMV PCV, SIMV VCV, PCV-VG, SIMV-PCV-VG, PSVPro™ with apnea backup
- Treat neonates to adult patients
- Tidal Volume monitoring as low as 5ml in pressure control ventilation and delivery as low as 20ml in volume control ventilation

SOFTWARE OPTIONS INCLUDED:

- Spirometry loops - ability to save and compare reference loops for continuous monitoring of ventilatory status
- SPONTANEOUS BREATHING SUITE
- ADVANCED PRESSURE AND SPONTANEOUS SUITE
- PAUSE GAS SOFTWARE OPTION
- Lung Recruitment :
 - 1) Single step - Vital Capacity
 - 2) Multi step - Lung Cycling

B650 WITH PDM MONITOR:

Providing clinical excellence, dependable data continuity and integration

- Combines the best of Marquette® Electronics and Datex-Ohmeda® legacies into one platform
- Wide set of clinical measurements are available with advanced parameter modules
- CARESCAPE Patient Data Module provides consistent hemodynamic measurement during intra-hospital transport and also transfers trend data
- Innovative algorithms aid in accurate diagnoses, including GE EK-Pro, GE DINAMAP® non-invasive blood pressure, and 12 SL™ diagnostic ECG with direct 2-way MUSE® communication
- Advanced airway modules provide for complete respiratory monitoring with gas exchange measurement

UNITY ID - COMMUNICATION BETWEEN CS650 & B650

HAEMODYNAMIC ACCESSORIES RESPIRATORY MODULE:

- E-sCAIO Respiratory module included
- D-Fends pck
- Sample lines pck

VAPOURISERS INCLUDED:

- 1 x Isoflurane Vapouriser
- 1 x Sevoflurane Vapouriser

CONSUMABLES AND ACCESSORIES INCLUDED:

UNIT PRICE INCLUSIVE OF VAT

- Dual flow sensors
- 1 x Circuit startup kit
- 1 x Reusable Absorber canister
- 1 x Pack Foam Dust Filters
- 4 x Castor Guards / Cable Guard
- Manual Breathing Bags 1 to 3L
- Paediatric circuit
- Neonatal circuit
- Jackson Rees circuit
- Face masks size 0 to 5
- Adult & Paediatric Resuscitator 550ml / 1500ml
- Filter

UNIT PRICE INCLUSIVE OF VAT

AS PER RT2 CIRCULAR AND ATTACHED PRICELIST

See attached pricing sheet

Offer price Includes: Start-up accessories & consumables as per tender spec, Plus Yr 1 - 2 Warranty and Preventative Maintenance.

OPTIONAL ITEMS

Extended Warranty and Maintenance Contract - 5 years

Optional Year 3 to 7 Warranty & Preventative Maintenance Including VAT

**includes Flow sensors and Oxygen cells

WEB SITE LINK

CARESTATION 650

[HTTPS://TINYURL.COM/3X45CHK5](https://tinyurl.com/3X45CHK5)

CARESTATION 650

Scan me



ANAESTHESIA

ITEM NUMBER: RT2-01-04 (OFFER 8) CARESTATION 650 B650 WITH PDM WITH NELLCOR SpO2

DOH Item Number	Replacement Accessories & Consumables List	Medhold Part Number	Quantity
RT2-01-04a	A full set paediatric 3/5 lead ECG cable	2017003-003	Each
RT2-01-04b	A full set adult 3/5 lead ECG cable	2017003-003	Each
RT2-01-04c	Reusable paediatric cuff	DUR-P2-2A	Each
RT2-01-04d	Reusable adult cuff	DUR-A2-2A	Each
RT2-01-04e	Reusable neonate cuff	DUR-P1-2A	Each
RT2-01-04f	Reusable obese cuff	DUR-A3-2A	Each
RT2-01-04g	Air hose	2058203-002	Each
RT2-01-04h	Finger sensor for paed	D-YS	Each
RT2-01-04i	Finger sensor for adults	Flexmax	Each
RT2-01-04j	Earlobe sensor for paed	D-YS	Each
RT2-01-04k	Earlobe sensor for adults	D-YSE	Each
RT2-01-04l	Finger wrap-around probe for neonates	D-YS	Each
RT2-01-04m	Sensor lead	2021406-001	Each
RT2-01-04n	Reusable temperature probe	M1024247	Each

UPGRADES: CS600

RT2-01-04r	Pressure control mode suites	M1809980	Each
RT2-01-04r	Spontaneous breathing suites	M1809981	Each
RT2-01-04r	CPAP-PSV Ventilation	M1809983	Each
RT2-01-04r	VCV Cardiac bypass	M1809984	Each
RT2-01-04r	Ecoflow	M1809985	Each
RT2-01-04r	Active Scav 30mm	M1809986	Each
RT2-01-04r	Pause Gas	M1809987	Each
RT2-01-04r	Lung Recruitment - VITAL CAPACITY AND CYCLING SOFTWARE OPTION	M1809988	Each

ANAESTHESIA

ITEM NUMBER: RT2-01-06 OFFER 1

ITEM NAME: ANAESTHESIA DELIVERY SYSTEM: ADVANCED WITH INTEGRATED PHYSIOLOGICAL MONITOR

Make/Brand	GE Healthcare
Model	Carestation 650 with B650 Monitor and PDM / Nellcor
Medhold Part Number	1012-9650-000
DOH SAP number	220401
Country of Origin	China

UNIT INCLUDES:

- Main Anaesthesia workstation
- Guided Checklist and calibration
- 3 Gasses O₂, N₂O, AIR
- integrated EZchange CO₂ bypass system enables absorbent canister replacement during prolonged surgeries without introducing leaks into the breathing system
- Condensor unit integrated into Canister and compact breathing system

VENTILATION OPERATING MODES:

- VCV, PCV, SIMV PCV, SIMV VCV, PCV-VG, SIMV-PCV-VG, PSVPro™ with apnea backup
- Treat neonates to adult patients
- Tidal Volume monitoring as low as 5ml in pressure control ventilation and delivery as low as 20ml in volume control ventilation

SOFTWARE OPTIONS INCLUDED:

- Spirometry loops - ability to save and compare reference loops for continuous monitoring of ventilatory status
- SPONTANEOUS BREATHING SUITE
- ADVANCED PRESSURE AND SPONTANEOUS SUITE
- PAUSE GAS SOFTWARE OPTION
- Lung Recruitment :

B650 WITH PDM MONITOR:

Providing clinical excellence, dependable data continuity and integration

- Combines the best of Marquette® Electronics and Datex-Ohmeda® legacies into one platform
- Wide set of clinical measurements are available with advanced parameter modules
- CARESCAPE Patient Data Module provides consistent hemodynamic measurement during intra-hospital transport and also transfers trend data
- Innovative algorithms aid in accurate diagnoses, including GE EK-Pro, GE DINAMAP® non-invasive blood pressure, and 12 SL™ diagnostic ECG with direct 2-way MUSE® communication
- Advanced airway modules provide for complete respiratory monitoring with gas exchange measurement



UNITY ID - communication between CS650 & B650

HAEMODYNAMIC ACCESSORIES RESPIRATORY MODULE:

- E-sCAIOV Respiratory module included
- D-Fends pck
- Sample lines pck

VAPOURISERS INCLUDED:

- 1 x Isoflurane Vapouriser
- 1 x Sevoflurane Vapouriser

CONSUMABLES AND ACCESSORIES INCLUDED:

UNIT PRICE INCLUSIVE OF VAT

- Dual flow sensors
- 1 x Circuit startup kit
- 1 x Reusable Absorber canister
- 1 x Pack Foam Dust Filters
- 4 x Castor Guards / Cable Guard
- Manual Breathing Bags 1 to 3L
- Paediatric circuit
- Neonatal circuit
- Jackson Rees circuit
- Face masks size 0 to 5
- Adult & Paediatric Resuscitator 550ml /1500ml
- Filter

UNIT PRICE INCLUSIVE OF VAT

AS per RT2 Circular and attached pricelist

Offer price Includes: Start-up accessories & consumables as per tender spec, Plus Yr 1 - 2 Warranty and Preventative Maintenance.

OPTIONAL ITEMS

Extended Warranty and Maintenance Contract - 5 years

Optional Year 3 to 7 Warranty & Preventative Maintenance Including VAT

**includes Flow sensors and Oxygen cells

WEB SITE LINK

CARESTATION 650

[HTTPS://TINYURL.COM/3X45CHK5](https://tinyurl.com/3X45CHK5)

CARESTATION 650

Scan me



ANAESTHESIA

ITEM NUMBER: RT2-01-06 (OFFER 1) CARESTATION 650 WITH B650 WITH PDM NELLCOR SpO2

DOH Item Number	Replacement Accessories & Consumables List	Medhold Part Number	Quantity
RT2-01-06a	A full set paediatric 3/5 lead ECG cable	2017003-003	Each
RT2-01-06b	A full set adult 3/5 lead ECG cable	2017003-003	Each
RT2-01-06c	Reusable paediatric cuff	DUR-P2-2A	Each
RT2-01-06d	Reusable adult cuff	DUR-A2-2A	Each
RT2-01-06e	Reusable neonate cuff	DUR-P1-2A	Each
RT2-01-06f	Reusable obese cuff	DUR-A3-2A	Each
RT2-01-06g	Air hose	2058203-002	Each
RT2-01-06h	Finger sensor for paedts	D-YS	Each
RT2-01-06i	Finger sensor for adults	Flexmax	Each
RT2-01-06j	Earlobe sensor for paedts	D-YS	Each
RT2-01-06k	Earlobe sensor for adults	D-YSE	Each
RT2-01-06l	Finger wrap-around probe for neonates	D-YS	Each
RT2-01-06m	Sensor lead	2021406-001	Each
RT2-01-06n	Reusable temperature probe	M1024247	Each

ANAESTHESIA

ITEM NUMBER: RT2-01-06 OFFER 4

ITEM NAME: ANAESTHESIA DELIVERY SYSTEM: ADVANCED WITH INTEGRATED PHYSIOLOGICAL MONITOR

Make/Brand	GE Healthcare
Model	Avance CS2 with B650 PDM Nellcor SpO2
Medhold Part Number	21012-9650-000
DOH SAP number	220403
Country of Origin	China

UNIT INCLUDES:

- Main Anaesthesia workstation
- Guided Checklist and calibration
- 3 Gasses O₂, N₂O, AIR
- integrated EZchange CO₂ bypass system enables absorbent canister replacement during prolonged surgeries without introducing leaks into the breathing system
- Condensor unit integrated into Canister and compact breathing system

VENTILATION OPERATING MODES:

- VCV, PCV, SIMV PCV, SIMV VCV, PCV-VG, SIMV-PCV-VG, PSVPro™ with apnea backup
- Treat neonates to adult patients
- Tidal Volume monitoring as low as 5ml in pressure control ventilation and delivery as low as 20ml in volume control ventilation
- Electronic mixer

SOFTWARE OPTIONS INCLUDED:

- Spirometry loops - ability to save and compare reference loops for continuous monitoring of ventilatory status
- SPONTANEOUS BREATHING SUITE
- ADVANCED PRESSURE AND SPONTANEOUS SUITE
- PAUSE GAS SOFTWARE OPTION
- Lung Recruitment :

B650 WITH PDM MONITOR:

Providing clinical excellence, dependable data continuity and integration

- Combines the best of Marquette® Electronics and Datex-Ohmeda® legacies into one platform
- Wide set of clinical measurements are available with advanced parameter modules
- CARESCAPE Patient Data Module provides consistent hemodynamic measurement during intra-hospital transport and also transfers trend data
- Innovative algorithms aid in accurate diagnoses, including GE EK-Pro, GE DINAMAP® non-invasive blood pressure, and 12 SL™ diagnostic ECG with direct 2-way MUSE® communication
- Advanced airway modules provide for complete respiratory monitoring with gas exchange measurement

OPTIONAL ITEMS

Extended Warranty and Maintenance Contract - 5 years
 Optional Year 3 to 7 Warranty & Preventative Maintenance Including VAT

**includes Flow sensors and Oxygen cells



UNITY ID - communication between CS650 & B650

HAEMODYNAMIC ACCESSORIES RESPIRATORY MODULE:

- E-sCAIOV Respiratory module included
- D-Fends pck
- Sample lines pck

VAPOURISERS INCLUDED:

- 1 x Isoflurane Vapouriser
- 1 x Sevoflurane Vapouriser

CONSUMABLES AND ACCESSORIES INCLUDED:

UNIT PRICE INCLUSIVE OF VAT

- Dual flow sensors
- 1 x Circuit startup kit
- 1 x Reusable Absorber canister
- 1 x Pack Foam Dust Filters
- 4 x Castor Guards / Cable Guard
- Manual Breathing Bags 1 to 3L
- Paediatric circuit
- Neonatal circuit
- Jackson Rees circuit
- Face masks size 0 to 5
- Adult & Paediatric Resuscitator 550ml /1500ml
- Filter

UNIT PRICE INCLUSIVE OF VAT

AS per RT2 Circular and attached pricelist

Offer price Includes: Start-up accessories & consumables as per tender spec, Plus Yr 1 - 2 Warranty and Preventative Maintenance.

WEB SITE LINK

AVANCE CS2

[HTTPS://TINYURL.COM/2Z-EE9AZW](https://tinyurl.com/2Z-EE9AZW)

AVANCE CS2

Scan me



ANAESTHESIA

ITEM NUMBER: RT2-01-06 (OFFER 4) AVANCE CS2 WITH B650 WITH PDM NELLCOR SpO2

DOH Item Number	Replacement Accessories & Consumables List	Medhold Part Number	Quantity
RT2-01-06a	A full set paediatric 3/5 lead ECG cable	2017003-003	Each
RT2-01-06b	A full set adult 3/5 lead ECG cable	2017003-003	Each
RT2-01-06c	Reusable paediatric cuff	DUR-P2-2A	Each
RT2-01-06d	Reusable adult cuff	DUR-A2-2A	Each
RT2-01-06e	Reusable neonate cuff	DUR-P1-2A	Each
RT2-01-06f	Reusable obese cuff	DUR-A3-2A	Each
RT2-01-06g	Air hose	2058203-002	Each
RT2-01-06h	Finger sensor for paedts	D-YS	Each
RT2-01-06i	Finger sensor for adults	Flexmax	Each
RT2-01-06j	Earlobe sensor for paedts	D-YS	Each
RT2-01-06k	Earlobe sensor for adults	D-YSE	Each
RT2-01-06l	Finger wrap-around probe for neonates	D-YS	Each
RT2-01-06m	Sensor lead	2021406-001	Each
RT2-01-06n	Reusable temperature probe	M1024247	Each

ANAESTHESIA

ITEM NUMBER: RT2-01-06 OFFER 7

ITEM NAME: ANAESTHESIA DELIVERY SYSTEM: ADVANCED WITH INTEGRATED PHYSIOLOGICAL MONITOR

Make/Brand	GE Healthcare
Model	Aisys CS2 with B650 Monitor and PDM / Nellcor
Medhold Part Number	1011-9050-000
DOH SAP number	220407
Country of Origin	USA

UNIT INCLUDES:

- Main Anaesthesia workstation
- Guided Checklist and calibration
- 3 Gasses O₂, N₂O, AIR
- integrated EZchange CO₂ bypass system enables absorbent canister replacement during prolonged surgeries without introducing leaks into the breathing system
- Condensor unit integrated into Canister and compact breathing system

VENTILATION OPERATING MODES:

- VCV, PCV, SIMV PCV, SIMV VCV, PCV-VG, SIMV-PCV-VG, PSVProTM with apnea backup
- Treat neonates to adult patients
- Tidal Volume monitoring as low as 5ml in pressure control ventilation and delivery as low as 20ml in volume control ventilation

ET CONTROL

Digitally enabled target control

- Et Control** automatically adjusts fresh gas concentrations to quickly and efficiently achieve and maintain end tidal oxygen and end tidal agent targets.
- Estimated MAC display helps you establish end tidal agent targets.

SOFTWARE OPTIONS INCLUDED:

- Spirometry loops - ability to save and compare reference loops for continuous monitoring of ventilatory status
- SPONTANEOUS BREATHING SUITE
- ADVANCED PRESSURE AND SPONTANEOUS SUITE
- PAUSE GAS SOFTWARE OPTION
- Lung Recruitment :

B650 WITH PDM MONITOR:

Providing clinical excellence, dependable data continuity and integration

- Combines the best of Marquette® Electronics and Datex-Ohmeda® legacies into one platform
- Wide set of clinical measurements are available with advanced parameter modules
- CARESCAPE Patient Data Module provides consistent hemodynamic measurement during intra-hospital transport and also transfers trend data
- Innovative algorithms aid in accurate diagnoses, including GE EK-Pro, GE DINAMAP® non-invasive blood pressure, and 12 SL™ diagnostic ECG with direct 2-way MUSE® communication
- Advanced airway modules provide for complete respiratory monitoring with gas exchange measurement



UNITY ID - communication between CS650 & B650

HAEMODYNAMIC ACCESSORIES

RESPIRATORY MODULE:

- E-sCAIOVe Respiratory module included
- D-Fends pck
- Sample lines pck

VAPOURISERS INCLUDED:

- 1 x Isoflurane Vapouriser
- 1 x Sevoflurane Vapouriser

CONSUMABLES AND ACCESSORIES INCLUDED:

UNIT PRICE INCLUSIVE OF VAT

- Dual flow sensors
- 1 x Circuit startup kit
- 1 x Reusable Absorber canister
- 1 x Pack Foam Dust Filters
- 4 x Castor Guards / Cable Guard
- Manual Breathing Bags 1 to 3L
- Paediatric circuit
- Neonatal circuit
- Jackson Rees circuit
- Face masks size 0 to 5
- Adult & Paediatric Resuscitator 550ml /1500ml
- Filter

UNIT PRICE INCLUSIVE OF VAT

AS per RT2 Circular and attached pricelist

Offer price Includes: Start-up accessories & consumables as per tender spec, Plus Yr 1 - 2 Warranty and Preventative Maintenance.

OPTIONAL ITEMS

Extended Warranty and Maintenance Contract - 5 years
 Optional Year 3 to 7 Warranty & Preventative Maintenance
 Including VAT **includes Flow sensors and Oxygen cells

WEB SITE LINK

AISYS CS2

[HTTPS://TINYURL.COM/2Z-FE9AZW](https://tinyurl.com/2Z-FE9AZW)

AISYS CS2

Scan me



ANAESTHESIA

ITEM NUMBER: RT2-01-06 (OFFER 7) AYSIS CS2 WITH B650 WITH PDM NELLCOR SpO2

DOH Item Number	Replacement Accessories & Consumables List	Medhold Part Number	Quantity
RT2-01-06a	A full set paediatric 3/5 lead ECG cable	2017003-003	Each
RT2-01-06b	A full set adult 3/5 lead ECG cable	2017003-003	Each
RT2-01-06c	Reusable paediatric cuff	DUR-P2-2A	Each
RT2-01-06d	Reusable adult cuff	DUR-A2-2A	Each
RT2-01-06e	Reusable neonate cuff	DUR-P1-2A	Each
RT2-01-06f	Reusable obese cuff	DUR-A3-2A	Each
RT2-01-06g	Air hose	2058203-002	Each
RT2-01-06h	Finger sensor for paed	D-YS	Each
RT2-01-06i	Finger sensor for adults	Flexmax	Each
RT2-01-06j	Earlobe sensor for paed	D-YS	Each
RT2-01-06k	Earlobe sensor for adults	D-YSE	Each
RT2-01-06l	Finger wrap-around probe for neonates	D-YS	Each
RT2-01-06m	Sensor lead	2021406-001	Each
RT2-01-06n	Reusable temperature probe	M1024247	Each

ANAESTHESIA

RT2-01-05 -ADV-VAPORIZERS

RT2-2021 Item	Replacement Accessories & Consumables List	SAP Number	Medhold Part Number	Quantity
RT2-01-05y	Sevoflurane vaporizer TEC 850 - CS600/Avance CS2/9100C NXT	150058/150059	1177-9850-SEV	Each
RT2-01-05z	Isoflurane vaporizer TEC 850 - CS600/Avance CS2/9100C NXT	150058/150060	1177-9850-ISO	Each
RT2-01-05aa	Desflurane vaporizer 6+ - CS600/Avance CS2/9100C NXT	150058/150061	1107-9616-000	Each

9100C NXT

Carestation 650

Avance CS2



RT2-2021 Item	Replacement Accessories & Consumables List	SAP Number	Medhold Part Number	Quantity
RT2-01-05y	Sevoflurane vaporizer Aladin 2 Cassette - Aisys CS2	149961	1100-9030-000	Each
RT2-01-05z	Isoflurane vaporizer Aladin 2 Cassette - Aisys CS2	149959	1100-9029-000	Each
RT2-01-05zz	Desflurane vaporizer Aladin 2 Cassette - Aisys CS2		1100-9025-000	Each

AGENT SPECIFIC ALADIN2™ CASSETTES

ELECTRONIC CONTROL – A UNIQUE TECHNOLOGY FOR VAPORIZING ANAESTHETIC AGENT

Our enhanced design provides:

- Accurate overflow protection.
- On screen, electronic level sensing and agent identification
- A prismatic site glass for clear indication of liquid agent level
- Electronic Control allows for automatic record keeping, gas usage calculation, fresh gas flow compensation, and temperature/pressure compensation
- Each cassette is magnetically coded for its specific agent allowing the Aisys Carestation to automatically identify the anaesthetic cassette being used
- Lightweight cassette allows for easier handling
- Virtually maintenance free
- Vaporiser self-check and diagnostics



Aisys CS2

ANAESTHESIA

RT2-01-05 -ADV-VAPORIZERS

TEC 6 PLUS VAPORIZER

- Electronic vaporizer which delivers the anaesthetic agent.
- Heats Suprane to maintain temperature and vapor pressure for consistent output.
- For added convenience and safety, the Tec® 6 Plus LED display indicates vaporizer status, including:
 - Operational
 - Warm-up
 - No output
 - Alarm battery low
 - Low agent
- Attaches to Datex-Ohmeda's interlocking Selecta Tec® manifold which allows you to change vaporizer
- arrangement from case to case or move vaporizers from suite to suite without the need for tools.



TEC 850 VAPORIZER

- The agent specific Tec™ 850 Vaporizer from GE Healthcare
- delivers isoflurane and sevoflurane effectively and efficiently.
- Each Tec 850 presents an improved ergonomic design, with consistent agent delivery and reliability while reducing total cost of ownership.



WEB SITE LINKS

ALADIN

[HTTPS://TINYURL.COM/3J2THNX2](https://tinyurl.com/3J2THNX2)

TEC 6

[HTTPS://TINYURL.COM/4HJXHKHJ](https://tinyurl.com/4HJXHKHJ)

TEC 850

[HTTPS://TINYURL.COM/VMM4W8BM](https://tinyurl.com/VMM4W8BM)

ALADIN

Scan me



TEC 6

Scan me



TEC 850

Scan me



ANAESTHESIA

MODULES

RT2-2021 Item	Replacement Accessories & Consumables List	SAP Number	Medhold Part Number	Quantity
RT2-01-06q	Neuromuscular transmission monitoring - E-NMT	223065	M1127542	Each
RT2-01-06r	Metabolic monitoring - E-sCAIOVX		5697038	Each
RT2-01-06s	Depth of anaesthesia monitoring system - Entropy	3011432	M1235680	Each
RT2-01-06t	Invasive cardiac output monitoring - E-COP		M1206807	Each
RT2-01-06t	Invasive cardiac output monitoring - PICCO		M1187062	Each


THE E-ENTROPY MODULE

is a single-width, plug-in module with the unique Entropy* algorithm that monitors the state of the brain. It is indicated for use within the hospital for adult and **pediatric patients older than two years.**

Benefits of Entropy measurement (1,2)

In adults, Entropy measurement

- May be used as an aid in monitoring the effects of certain anesthetic agents
- May help the user titrate anesthetic drugs according to the individual needs
- May be associated with a reduction of anesthetic use and faster emergence from anesthesia Measurement technology
- Utilizes the intuitive and published Entropy algorithm, a Datex-Ohmeda innovation (3)
- Based on acquisition and processing of EEG and diagnostics



FEMG SIGNALS

- Features two Entropy parameters

Response Entropy (RE) - a fast reacting parameter for detecting activation of facial muscles

State Entropy (SE) - a steady and robust parameter for assessing the effect of anesthetic drugs in the brain in adults

Clinical use of Entropy




THE NEUROMUSCULAR TRANSMISSION MODULE, E-NMT, INTEGRATES

the measurement of the level of neuromuscular block into a single-width, plug-in module.

FEATURES

- Provides a quantitative, automatic measurement of muscle response to an electrical stimulus
- Offers all common stimulation modes: train of four (TOF), single twitch (ST), double burst stimulation (DBS) and tetanic stimulation
- Covers the whole range of neuromuscular blockade
- Provides two measurement sensor options: the easy-to-apply MechanoSensor (KMG) for routine clinical use, and ElectroSensor (EMG) for research use



ANAESTHESIA

MODULES

CONTINUOUS CARDIAC OUTPUT MODULE, E-PICCO

For less invasive continuous cardiac output monitoring. The E-PiCCO module provides continuous cardiac output (CCO) measurement based on pulse contour calculation, transpulmonary thermodilution cardiac output (C.O.) and blood pressure (P). When combined with a CARESCAPE® modular monitor, the E-module enables a graphical view for quick hemodynamic status assessment.

FEATURES

- Direct key on the module for zeroing invasive pressure channels
- Easy insertion/removal of module without interrupting other monitoring
- Uses Pulsion's PiCCO catheters




Figure 6. PICCO graphical view on a CARESCAPE monitor




Figure 1. Transpulmonary thermodilution

WEB SITE LINK

E-NMT

[HTTPS://TINYURL.COM/4XJJN8K7](https://tinyurl.com/4XJJN8K7)

ENTROPY

[HTTPS://TINYURL.COM/YHSY72YH](https://tinyurl.com/YHSY72YH)

ENTROPY

Scan me



NMT

Scan me



RT2-2021

**SUPPLY, DELIVERY, INSTALLATION,
COMMISSIONING AND MAINTENANCE OF
MEDICAL EQUIPMENT**

**TO THE STATE FOR THE PERIOD
01 DECEMBER 2021
TO 01 DECEMBER 2024.**

SECTION B

CARDIOLOGY

CONTACT DETAILS

Sibongile Leeu: s.leeu@medhold.co.za
Head Office: 011 966 0600

Or contact your local branch – details on the back page

DEFIBRILLATOR, ECG & PACING

ITEM NUMBER: RT2-05-04
ITEM NAME: DEFIBRILLATOR, ECG & PACING

Make/Brand	Bexen Cardio
Model	Reanibex 700C
Medhold Part Number	12437/2020 700C with ECG, AED, Pacing
DOH SAP number	220485
Country of Origin	Spain



The REANIBEX 700 is a monitor defibrillator with AED, (Automated External Defibrillator) portable and lightweight ideal for rapid provision of the first treatment to patients suffering a sudden cardiac arrest.

Made with the latest technologies in the field of defibrillation (biphasic waveform up to 200J), incorporates a screen that allows not only the visualization of the user guidance messages but also of the ECG signal in real time.

Operates from the mains (AC), from a vehicle battery (DC) or from it's own battery

FEATURING:

- Biphasic up to 200J
- Control of the defibrillator and the printer from the paddles
- Integrated paediatric paddles
- Autotest when switched on and during operation
- Simple and intuitive operation
- AED with ERC/AHA guidelines 2015 and its review of 2017
- Non Invasive pacing fixed and demand mode (optional)
- Automatic daily selftests ensures the device continual readiness
- Optional 4, 5 or 10 Leads ECG
- Synchronised cardioversion
- Printer
- Optional "Reanibex Data Manager" software to manage and organise all the data gathered

UNIT INCLUDES:

Equipment supplied: Reanibex 700C with ECG, AED, and Pacing P/n RT2-05-04
 Manual and AED (adult and pediatric patients) biphasic defibrillator monitor REANIBEX Serie-700C with Non invasive PACEMAKER, with 50 mm 2 channel printer and with 5.7" TFT color screen 3 waves, ECG up to 12 leads, rechargeable battery, internal battery charger, external adult and pediatric paddles, 5 leads ECG patient cable, compact flash memory card, a set of adult disposable cable-electrodes, a set of pediatric disposable cable-electrodes and standard accessories: main power and grounding cables, fuses, disposable ECG adult electrodes for the patient cable (1 bag of 50 units), 10 rolls of paper (30 m), 1 tube of gel (250 c.c.) and 1 Instruction Manual, & Test Load 50 Ohms commence after formal Commissioning and handover of the equipment.

WARRANTY (2-YEAR): YEAR 1 AND YEAR 2

A 2-year warranty against poor workmanship and latent defects and parts. This must be all inclusive and include, BUT NOT LIMITED TO, amongst others, ALL PARTS (including Batteries), Labour, Traveling and Accommodation. The 2-year warranty must also include all quality check, quality assurance requirements, Preventative Maintenance/Calibrations, Software updates and upgrades to be included. This 2-year warranty will commence after formal Commissioning and handover of the equipment.

UNIT PRICE INCLUSIVE OF VAT
(STANDARD ACCESSORIES & CONSUMABLES, AND 2-YEAR WARRANTY)

See attached pricing sheet

BEXEN CARDIO REANIBEX 700

WEBSITE PRODUCT LINK:
[HTTPS://TINYURL.COM/HD5CY8J3](https://tinyurl.com/HD5CY8J3)

Scan me



DEFIBRILLATOR, ECG & PACING

ITEM NUMBER: RT2-05-04
ITEM NAME: DEFIBRILLATOR, ECG & PACING
ACCESSORIES & OPTIONS:

RT2-05-04	Defibrillator Accessories and Consumables	Medhold Part Number	Quantity
RT2-05-04a	Reusable external adult/paediatric defibrillator paddles	DGS 600 A	1
RT2-05-04b	3-Lead ECG patient cable (x1)	KPA 137 C	1
RT2-05-04c	ECG electrodes (x25)	TKL 102 A	2 x 25
RT2-05-04d	AED pads (x2 pairs)	TKL 0311 A	2 of
RT2-05-04e	Cable for AED pads (x1)	KSA 0103 A	1
RT2-05-04f	Pacing electrodes	TKL 0311 A	1
RT2-05-04g	Pacing cable	KSA 0103 A	1
RT2-05-04h	Electrode gel (x1: 250ml)	XEH 001 A	1
RT2-05-04i	Printer paper (x 10 rolls)	XEG 001 A	10 rolls
RT2-05-04j	Dedicated trolley	70126202DC	1
RT2-05-04k	Internal Defibrillation iD System		
	ID Smart cable (Reusable)	KSS 0730	1
	ID Electrode common for adult and pediatric patients (disposable)	TKC 0523	1
	ID Paddle for adult patient (disposable)	TKC 0520	1
	ID Paddle for pediatric patient (disposable)	1300	1
	Total for one set		1

ACCESSORIES & OPTIONS:

DGX 705 A	Transport bag with side pockets and shoulder strap.	DGX 705 A	1
R700+BAG - R700ONLY - DGS 0830 B	Wall bracket for ambulance. DGS 0850 B / 10 G certified (EN 1789)	R700+BAG - DGS 0850 B / R700ONLY - DGS 0830 B	1
DGS 5600 A	Rechargeable NiMH battery. Autonomy of more than 130 discharges at 200J. Charging time:3 hours	DGS 5600 A	1
DGS 0800 A	External battery charger	DGS 0800 A	1
DII 0020 A	PC application "REANIBEX- Data Manager" software : for transferring the ECG data and audio voices from the Defibrillator to the computer through a compact flash reader.	DII 0020 A	1
70126202DC	Trolley UniCart II with Roll 125mm (local) with basket HU30	70126202DC	1

DISPOSABLES AND SPARE PARTS

TKL 0311 A	Adult Disposable cable-electrodes, for defibrillation or for pacemaker option use.	TKL 0311 A	1
TKL 0313 A (TYCO)	Adult Preconnected disposable cable-electrodes.	TKL 0313 A (TYCO)	1
TKL 0312 B	Pediatric Disposable cable-electrodes	TKL 0312 B	1
KSA 0103 A	Reusable longer cable for ADULT disposable cable-electrodes for defibrillation (1,5 mts long).	KSA 0103 A	1
KAA 109 A	12VDC battery connection cable	KAA 109 A	1
KPA 139 C	4 leads patient cable	KPA 139 C	1

DEFIBRILLATOR, ECG & PACING

ITEM NUMBER: RT2-05-04

ITEM NAME: DEFIBRILLATOR, ECG & PACING

ACCESSORIES & OPTIONS:

RT2-05-04	Defibrillator Accessories and Consumables	Medhold Part Number	Quantity
KPA 137 C	5 lead patient cable	KPA 137 C	1
KPA 138 C	10 lead patient cable	KPA 138 C	1
DGS 600 A	Complete set of external adult and pediatric paddles	DGS 600 A	1
DGS 605 A	Complete set of external pediatric paddles (without external adult paddles)	DGS 605 A	1
DGS 630	Adult's external paddles	DGS 630	1
TKL 102 A	Adult disposable electrodes for the patient cable (bag of 50 units)	TKL 102 A	1
TKL 212 B	Pediatric disposable electrodes for the patient cable (bag of 50 units)	TKL 212 B	1
XEG 001 A	Thermosensible roll of paper, 50 mm and 30 m for the recorder x 10 rolls	XEG 001 A	1
XEH 001 A	Gel tube of 250ml	XEH 001 A	1
Test Load 50 Ohms	Test Load 50 Ohms	Test Load 50 Ohms	1

DEFIBRILLATOR, ECG & SPO2

ITEM NUMBER: RT2-05-05
ITEM NAME: DEFIBRILLATOR, ECG & SPO2

Make/Brand	Bexen Cardio
Model	Reanibex 700C
Medhold Part Number	12437/2020 700C with ECG, AED, SPO2
DOH SAP number	220486
Country of Origin	Spain



The REANIBEX 700 is a monitor defibrillator with AED, (Automated External Defibrillator) portable and lightweight ideal for rapid provision of the first treatment to patients suffering a sudden cardiac arrest.

Made with the latest technologies in the field of defibrillation (biphasic waveform up to 200J), incorporates a screen that allows not only the visualization of the user guidance messages but also of the ECG & SpO2 signal in real time.

Operates from the mains (AC), from a vehicle battery (DC) or from its own battery

FEATURING:

- Biphasic up to 200J
- Control of the defibrillator and the printer from the paddles
- Integrated paediatric paddles
- Autotest when switched on and during operation
- Simple and intuitive operation
- AED with ERC/AHA guidelines 2015 and its review of 2017
- Non Invasive pacing fixed and demand mode (optional)
- Pulse Oximetry SpO2 from Masimo SET
- Automatic daily selftests ensures the device continual readiness
- Optional 4, 5 or 10 Leads ECG
- Synchronised cardioversion
- Printer
- Optional "Reanibex Data Manager" software to manage and organise all the data gathered

UNIT INCLUDES:

Equipment supplied: Reanibex 700C with ECG, AED, and SPO2 P/n RT2-05-05

Manual and AED (adult and pediatric patients) biphasic defibrillator monitor REANIBEX Serie-700C with "MASIMO" PULSIOXIMETER (SPO2) module. with 50 mm 2 channel printer and with 5.7" TFT color screen 3 waves, ECG up to 12 leads, rechargeable battery, internal battery charger, external adult and pediatric paddles, 5 leads ECG patient cable, Includes adult finger transducer and PULSIOXIMETER transducer cable, compact flash memory card, a set of adult disposable cable-electrodes, a set of pediatric disposable cable-electrodes and standard accessories: main power and grounding cables, fuses, disposable ECG adult electrodes for the patient cable (1 bag of 50 units), 10 rolls of paper (30 m), 1 tube of gel (250 c.c.) and 1 Instruction Manual. Test Load 50 Ohms.

WARRANTY (2-YEAR): YEAR 1 AND YEAR 2

A 2-year warranty against poor workmanship and latent defects and parts. This must be all inclusive and include, BUT NOT LIMITED TO, amongst others, ALL PARTS (including Batteries), Labour, Traveling and Accommodation. The 2-year warranty must also include all quality check, quality assurance requirements, Preventative Maintenance/Calibrations, Software updates and upgrades to be included. This 2-year warranty will commence after formal Commissioning and handover of the equipment.

UNIT PRICE INCLUSIVE OF VAT

See attached pricing sheet

BEXEN CARDIO REANIBEX 700

(STANDARD ACCESSORIES & CONSUMABLES, AND 2-YEAR WARRANTY)
WEBSITE PRODUCT LINK:
[HTTPS://TINYURL.COM/HD5CY8J3](https://tinyurl.com/HD5CY8J3)

Scan me



ITEM NUMBER: RT2-05-05
ITEM NAME: DEFIBRILLATOR, ECG & SPO2
ACCESSORIES & OPTIONS:

RT2-05-05	Defibrillator, with ECG & SpO2 - - Accessories and Consumables	Medhold Part Number	Quantity
RT2-05-05a	3-Lead ECG patient cable (x1)	KPA 137 C	1
RT2-05-05b	ECG electrodes (25/Pack)	TKL 102 A	2 x 25
RT2-05-05c	AED pads (x2 pairs)	TKL 0311 A	2 of
RT2-05-05d	Cable for AED pads (x1)	KSA 0103 A	1
RT2-05-05e	Electrode gel (x1: 250ml)	XEH 001 A	1
RT2-05-05f	Printer paper (x 10 rolls)	XEG 001 A	10 rolls
RT2-05-05g	SpO2 transducer	TKC 0213 A	1
RT2-05-05h	SpO2 cable	TKC 0201 B	1
RT2-05-05i	Dedicated trolley	70126202DC	1
RT2-05-05j	Internal Defibrillation iD System (One set) ID Smart cable (Reusable) ID Electrode common for adult and pediatric patients (disposable) ID Paddle for adult patient (disposable) ID Paddle for pediatric patient (disposable) Total for one set	KSS 0730 TKC 0523 TKC 0520 1300	1 1 1 1 1
RT2-05-05k	Reusable external adult/paediatric defibrillator paddles	DGS 600 A	1

OPTIONAL FEATURES

PM	Pacemaker . Includes one set of disposable cable-electrodes.	PM	1
ACCESSORIES DGX 705 A	Transport bag with side pockets and shoulder strap.	DGX 705 A	1 1
R700+BAG - DGS 0850 B / R700ONLY - DGS 0830 B	Wall bracket for ambulance. 10 G certified (EN 1789)	R700+BAG - DGS 0850 B / R700ONLY - DGS 0830 B	1
DGS 5600 A	Rechargeable NiMH battery. Autonomy of more than 130 discharges at 200J. Charging time:3 hours	DGS 5600 A	1
DGS 0800 A	External battery charger	DGS 0800 A	1
DII 0020 A	PC application "REANIBEX- Data Manager" software : for transferring the ECG data and audio voices from the Defibrillator to the computer through a compact flash reader.	DII 0020 A	1
70126202DC	Trolley UniCart II with Roll 125mm (local) with basket HU30	70126202DC	1

DISPOSABLES AND SPARE PARTS

TKL 0311 A	Adult Disposable cable-electrodes, for defibrillation or for pacemaker option use.	TKL 0311 A	1
TKL 0313 A (TYCO)	Adult Preconnected disposable cable-electrodes.	TKL 0313 A (TYCO)	1

ITEM NUMBER: RT2-05-05
ITEM NAME: DEFIBRILLATOR, ECG & SPO2
ACCESSORIES & OPTIONS:

RT2-05-05	Defibrillator, with ECG & SpO2- - Accessories and Consumables	Medhold Part Number	Quantity
TKL 0312 B	Pediatric Disposable cable-electrodes (Not for AED use).	TKL 0312 B	1
KSA 0103 A	Reusable longer cable for ADULT disposable cable-electrodes for defibrillation (1,5 mts long).	KSA 0103 A	1
KAA 109 A	12VDC battery connection cable	KAA 109 A	1
TKC 0201 A	MASIMO Pulsioxymeter transducer longer cable (12 feet)	TKC 0201 A	1
TKC 0201 B	MASIMO Pulsioxymeter transducer longer cable (4 feet)	TKC 0201 B	1
TKC 0213 A	MASIMO Reusable Adult's finger Pulsioximeter transducer (>30 kg)	TKC 0213 A	1
TKC 0214 A	MASIMO Multisite reusable pulsioximeter transducer (> 1kg)	TKC 0214 A	1
TKC 0216 A	MASIMO Reusable pediatric pulsioximetre transducer (10-50 kg)	TKC 0216 A	1
KPA 139 C	4 leads patient cable	KPA 139 C	1
KPA 137 C	5 lead patient cable	KPA 137 C	1
KPA 138 C	10 lead patient cable	KPA 138 C	1
DGS 600 A	Complete set of external adult and pediatric paddles	DGS 600 A	1
DGS 605 A	Complete set of external pediatric paddles (without external adult paddles)	DGS 605 A	1
DGS 630	Adult's external paddles	DGS 630	1
TKL 102 A	Adult disposable electrodes for the patient cable (bag of 50 units)	TKL 102 A	1
TKL 212 B	Pediatric disposable electrodes for the patient cable (bag of 50 units)	TKL 212 B	1
XEG 001 A	Thermosensible roll of paper, 50 mm and 30 m for the recorder x 10 rolls	XEG 001 A	1
XEH 001 A	Gel tube of 250ml	XEH 001 A	1
	Test Load 50 Ohms		1

DEFIBRILLATOR, ECG, SPO2 & NIBP

ITEM NUMBER: RT2-05-06
ITEM NAME: DEFIBRILLATOR, ECG, SPO2 & NIBP

Make/Brand	Bexen Cardio
Model	Reanibex 700C
Medhold Part Number	12437/2020 700C with ECG, AED, SPO2 and NIBP
DOH SAP number	220487
Country of Origin	Spain



The REANIBEX 700 is a monitor defibrillator with AED, (Automated External Defibrillator) portable and lightweight ideal for rapid provision of the first treatment to patients suffering a sudden cardiac arrest.

Made with the latest technologies in the field of defibrillation (biphasic waveform up to 200J), incorporates a screen that allows not only the visualization of the user guidance messages but also of the ECG & optional SpO2 signal in real time and NIBP.

Operates from the mains (AC), from a vehicle battery (DC) or from its own battery

FEATURING:

- Biphasic up to 200J
- Control of the defibrillator and the printer from the paddles
- Integrated paediatric paddles
- Autotest when switched on and during operation
- Simple and intuitive operation
- AED with ERC/AHA guidelines 2015 and its review of 2017
- Non Invasive pacing fixed and demand mode (optional)
- Pulse Oximetry SpO2 from Masimo SET
- NIBP from Suntech
- Automatic daily selftests ensures the device continual readiness
- Optional 4, 5 or 10 Leads ECG
- Synchronised cardioversion
- Printer
- Optional "Reanibex Data Manager" software to manage and organise all the data gathered

UNIT INCLUDES:

Reanibex 700C with ECG, AED, SPO2 and NIBP P/n RT2-05-06

Manual and AED (adult and pediatric patients) biphasic defibrillator monitor REANIBEX Serie-700C Includes "MASIMO" PULSIOXIMETER (SPO2) module. Includes "SUNTECH" Non Invasive Blood Pressure - NIBP module. with 50 mm 2 channel printer and with 5.7" TFT color screen 3 waves, ECG up to 12 leads, rechargeable battery, internal battery charger, external adult and pediatric paddles, 5 leads ECG patient cable, One adult finger transducer and PULSIOXIMETER transducer cable. one reusable adult cuff (23 to 33 cm) and NIBP air hose. compact flash memory card, a set of adult disposable cable-electrodes, a set of pediatric disposable cable-electrodes and standard accessories: main power and grounding cables, fuses, disposable ECG adult electrodes for the patient cable (1 bag of 50 units), 10 rolls of paper (30 m), 1 tube of gel (250 c.c.) and 1 Instruction Manual. Test Load 50 Ohms

WARRANTY (2-YEAR): YEAR 1 AND YEAR 2

A 2-year warranty against poor workmanship and latent defects and parts. This must be all inclusive and include, BUT NOT LIMITED TO, amongst others, ALL PARTS (including Batteries), Labour, Traveling and Accommodation. The 2-year warranty must also include all quality check, quality assurance requirements, Preventative Maintenance/Calibrations, Software updates and upgrades to be included. This 2-year warranty will commence after formal Commissioning and handover of the equipment.

UNIT PRICE INCLUSIVE OF VAT

(STANDARD ACCESSORIES & CONSUMABLES, AND 2-YEAR WARRANTY)

WEBSITE PRODUCT LINK:

[HTTPS://TINYURL.COM/HD5CY8J3](https://tinyurl.com/HD5CY8J3)

See attached pricing sheet

BEXEN CARDIO REANIBEX 700

Scan me



ITEM NUMBER: RT2-05-06
ITEM NAME: DEFIBRILLATOR, ECG, SPO2 & NIBP
ACCESSORIES & OPTIONS:

RT2-05-06	Defibrillator, with ECG, SpO2, and NIBP Accessories	Medhold Part Number	Quantity
RT2-05-06a	3-Lead ECG patient cable (x1)	KPA 137 C	1
RT2-05-06b	ECG electrodes (25/Pack)	TKL 102 A	2 x 25
RT2-05-06c	AED pads (x2 pairs - Paeds & Adults)	TKL 0311 A & TKL 0312 B	2 of
RT2-05-06d	Cable for AED pads (x1)	KSA 0103 A	1
RT2-05-06e	Electrode gel (x1: 250ml)	XEH 001 A	1
RT2-05-06f	Printer paper (x 10 rolls)	XEG 001 A	10 rolls
RT2-05-06g	SpO2 sensor	TKC 0213 A	1
RT2-05-06h	SpO2 cable	TKC 0201 B	1
RT2-05-06i	BP hose	TKC 0710 A	1
RT2-05-06j	BP cuff (Adult)	TKC 0711 A	1
RT2-05-06k	Dedicated trolley	70126202DC	1
RT2-05-06l	Internal Defibrillation iD System		
	ID Smart cable (Reusable)	KSS 0730	1
	ID Electrode common for adult and pediatric patients (disposable)	TKC 0523	1
	ID Paddle for adult patient (disposable)	TKC 0520	1
	ID Paddle for pediatric patient (disposable)	1300	1
	Total for one set		1
RT2-05-06m	Reusable external adult/paediatric defibrillator paddles	DGS 600 A	1
OPTIONAL FEATURES			
PM	Pacemaker . Includes one set of disposable cable-electrodes. order upfront	PM	1
ACCESSORIES			
DGX 705 A	Transport bag with side pockets and shoulder strap.	DGX 705 A	1
R700+BAG - DGS 0850 B / R700ONLY - DGS 0830 B	Wall bracket for ambulance. 10 G certified (EN 1789)	R700+BAG - DGS 0850 B / R700ONLY - DGS 0830 B	1
DGS 5600 A	Rechargeable NiMH battery. Autonomy of more than 130 discharges at 200J. Charging time: 3 hours	DGS 5600 A	1
DGS 0800 A	External battery charger	DGS 0800 A	1
DII 0020 A	PC application "REANIBEX- Data Manager" software : for transferring the ECG data and audio voices from the Defibrillator to the computer through a compact flash reader.	DII 0020 A	1
70126202DC	Trolley UniCart II with Roll 125mm (local) with basket HU30	70126202DC	1

ITEM NUMBER: RT2-05-06
ITEM NAME: DEFIBRILLATOR, ECG, SPO2 & NIBP
ACCESSORIES & OPTIONS:

RT2-05-06	Defibrillator, with ECG, SpO2,	Medhold Part Number and NIBP Accessories	Quantity
DISPOSABLES AND SPARE PARTS			
TKL 0311 A	Adult Disposable cable-electrodes, for defibrillation or for pacemaker option use.	TKL 0311 A	1
TKL 0313 A (TYCO)	Adult Preconnected disposable cable-electrodes.	TKL 0313 A (TYCO)	1
TKL 0312 B	Pediatric Disposable cable-electrodes	TKL 0312 B	1
KSA 0103 A	Reusable longer cable for ADULT disposable cable-electrodes for defibrillation (1,5 mts long).	KSA 0103 A	1
KAA 109 A	12VDC battery connection cable	KAA 109 A	1
TKC 0201 A	MASIMO Pulsioxymeter transducer longer cable (12 feet)	TKC 0201 A	1
TKC 0201 B	MASIMO Pulsioxymeter transducer longer cable (4 feet)	TKC 0201 B	1
TKC 0213 A	MASIMO Reusable Adult's finger Pulsioximeter transducer (>30 kg)	TKC 0213 A	1
TKC 0214 A	MASIMO Multisite reusable pulsioximeter transducer (> 1kg)	TKC 0214 A	1
TKC 0216 A	MASIMO Reusable pediatric pulsioximetre transducer (10-50 kg)	TKC 0216 A	1
KPA 139 C	4 leads patient cable	KPA 139 C	1
KPA 137 C	5 lead patient cable	KPA 137 C	1
KPA 138 C	10 lead patient cable	KPA 138 C	1
DGS 600 A	Complete set of external adult and pediatric paddles	DGS 600 A	1
DGS 605 A	Complete set of external pediatric paddles (without external adult paddles)	DGS 605 A	1
DGS 630	Adult's external paddles	DGS 630	1
TKL 102 A	Adult disposable electrodes for the patient cable (bag of 50 units)	TKL 102 A	1
TKL 212 B	Pediatric disposable electrodes for the patient cable (bag of 50 units)	TKL 212 B	1
XEG 001 A	Thermosensible roll of paper, 50 mm and 30 m for the recorder x 10 rolls	XEG 001 A	1
XEH 001 A	Gel tube of 250ml	XEH 001 A	1
TKC 0710 A	3 mts NIBP air hose	TKC 0710 A	1
TKC 0712 A	infant reusable cuff (8 to 13 cm)	TKC 0712 A	1
TKC 0713 A	child reusable cuff (12 to 19 cm)	TKC 0713 A	1
TKC 0714 A	small adult reusable cuff (17 to 25 cm)	TKC 0714 A	1
TKC 0711 A	adult reusable cuff (23 to 33 cm)	TKC 0711 A	1
TKC 0715 A	large adult reusable cuff (31 to 40 cm)	TKC 0715 A	1
TKC 0716 A	Thigh reusable cuff (38 to 50 cm)	TKC 0716 A	1
	Test Load 50 Ohms		

ELECTROCARDIOGRAPH 12 CHANNEL WITH TROLLEY

ITEM NUMBER: RT2-05-07
ITEM NAME: ELECTROCARDIOGRAPH 12 CHANNEL WITH TROLLEY

Make/Brand	GE Healthcare
Model	MAC 2000 ECG recorder with trolley
Medhold Part Number	2063587-017
DOH SAP number	220490
Country of Origin	Germany/India


KEY FEATURES:

- Automatic capture of 10-second ECGs with one-touch operation
- 7" color display is configurable
- On-screen preview of 12-lead waveforms and ECG results helps streamline review and diagnosis
- Multiple connectivity and export options for easy data access, transfer and storage
- Marquette 12SL ECG analysis program for clinical decision support
- Supports paperless workflow with options for bidirectional communication, including Orders, ADT query features
- Seamless connectivity with the MUSE* cardiology information CardioSoft* program—or send reports to a shared directory on a PC
- Connectivity with LAN, Wi-Fi, serial port, modem, or SD card
- PDF and XML export capabilities
- Internal storage for up to 200 ECGs
- Marquette 12SL ECG analysis program with measurement and interpretation
- Arrhythmia detection
- ACI-TIPI algorithm for acute cardiac ischemia
- Gender-specific interpretation provides improved sensitivity to Acute Myocardial Infarction in women
- RR analysis for detection of patterns in heart rate variability

UNIT INCLUDES:

MAC 2000 ECG recorder with color graphics display, Resting ECG Measurements and Interpretation: Standard accessory kit, including IEC patient cable, with Leadwire Set/10, Grabber connector, 70 cm, IEC, and recording paper Power cord and user's manual. MAC 2000 RR Analysis Unicart II Trolley with basket for paper and accessories

WARRANTY (2-YEAR): YEAR 1 AND YEAR 2

A 2-year warranty against poor workmanship and latent defects and parts. This must be all inclusive and include, BUT NOT LIMITED TO, amongst others, ALL PARTS (including Batteries), Labour, Traveling and Accommodation. The 2-year warranty must also include all quality check, quality assurance requirements, Preventative Maintenance/Calibrations, Software updates and upgrades to be included. This 2-year warranty will commence after formal Commissioning and handover of the equipment.

UNIT PRICE INCLUSIVE OF VAT

See attached pricing sheet

(STANDARD ACCESSORIES & CONSUMABLES, AND 2-YEAR WARRANTY)
WEBSITE PRODUCT LINK:
[HTTPS://TINYURL.COM/5EPFHK6A](https://tinyurl.com/5EPFHK6A)

MAC 2000



ITEM NUMBER: RT2-05-07 ELECTROCARDIOGRAPH 12 CHANNEL WITH TROLLEY
ACCESSORIES & OPTIONS:

RT2-05-07	ELECTROCARDIOGRAPH (12 CHANNEL) ACCESSORIES and CONSUMABLES	Medhold Part Number	Quantity
RT2-05-07a	ECG cable with adaptor clips	2063587-713	1
RT2-05-07b	ECG electrodes (25/Pack) Box of 600 electrodes	1700-030	1
RT2-05-07c	Dedicated trolley	70126202DC	1
RT2-05-07d	Battery	2066261-085	1
STORAGE & COMMUNICATION OPTIONS			
2063587-104	Internal Storage 200 - SD Card XML/Hilltop	2063587-104	1
2063587-155	MUSE LAN Serial, Int Storage 200, SD Card PDF/XML/Hilltop	2063587-155	1
2063587-153	CardioSoft LAN Serial, Int Storage 200, SD Card PDF/XML/Hilltop	2063587-153	1
2063587-154	PC LAN, Int Storage 200, SD Card PDF/XML /Hilltop	2063587-154	1
2063587-713	STARTER KIT IEC TRUNK CABLE AND LEAD WIRES W/ ECG GRABBERS	2063587-713	1
REUSABLE ACCESSORIES VALUE STARTER SET			
2104726-001	ECG Value Cable 10 Lead Cable/Leadwire IEC	2104726-001	1
2104783-001	Reusable Value Bulb Electrode 6/set	2104783-001	1
2104785-001	Clamp Limb IEC Value 4/set	2104785-001	1
21708318	Electrode cream, 250-ml refill bottle	21708318	1
PATIENT CABLE AND LEADWIRE SET FOR DISPOSABLE ELECTRODES			
2104722-001	Trunk Cable 10-Lead IEC	2104722-001	1
2104748-001	Lead wire Set, Base 10, Grabber, IEC.	2104748-001	1
2104726-001	ECG Value Cable 10 Lead Cable/Leadwire IEC	2104726-001	1
2104729-001	Adapter Clip, Wide Mouth Clip, 14 pack	2104729-001	1
2001925-003	Leadwire, CAM/Universal, 40 in.	2001925-003	1
2001925-005	Leadwire, CAM/Universal, 26 in.	2001925-005	1
900178-102	Universal Leadwire Connector Set, Mactrode, IEC, 10/set	900178-102	1
2104724-002	Lead wire Set, Base 10, Banana, IEC.	2104724-002	1
2074294-001	MAC 2000 Computer based training CBT	2074294-001	1
DISPOSABLES			
1700-030	Cleartrace Clear Tape With Adhesive Gel Conmed Box of 600	1700-030	1
TKL 102 A	Adult disposable electrodes for the patient cable (bag of 50 units)	TKL 102 A	1
2104779-001	ECG Snap Electrode Foam TearDrop Wet Gel, 38.1mm/1.5", 600 Electrodes /Case. Adult-Packaging: 60 Electrodes/Pouches - 10 Pouches/Case. Pouches are re-sealable.	2104779-001	1
2104772-001	CONTRAST PAPER A4; 150 SHEETS	2104772-001	1
BATTERY			
2066261-085	14.54V 3.5AH LI-ION BATTERY PACK	2066261-085	1



RT2-2021

**SUPPLY, DELIVERY, INSTALLATION,
COMMISSIONING AND MAINTENANCE OF
MEDICAL EQUIPMENT**

**TO THE STATE FOR THE PERIOD
01 DECEMBER 2021
TO 01 DECEMBER 2024.**

SECTION C

DIGITAL SOLUTIONS

CONTACT DETAILS

Sibongile Leeu: s.leeu@medhold.co.za
Head Office: 011 966 0600

Or contact your local branch – details on the back page

CENTRICITY HIGH ACUITY ANAESTHESIA

ITEM NUMBER: RT2-01-060**ITEM NAME: AIMS (ANAESTHESIA INFORMATION MANAGEMENT SYSTEM)**

Centricity High Acuity Anaesthesia supports a compliant and high-quality patient record while helping staff focus on the patient during critical workflow transitions.

IT ENABLES:

- Automated device data documentation with clinically relevant visualization
- Risk assessment via tools such as the WHO Surgical Safety checklist
- Faster data entry with touchscreen and gesture inputs
- Digital patient archiving
- Consistent implementation of care protocols

UNIT INCLUDES:

- Centricity Anaesthesia 8 Theatre License
- Medical device interfacing*
- GCX Articulating Workstation Arms
- 23.8" Touch All-in-One Desktop Computer (x8)
- Production and Backup Server
- Basic Theatre interconnectivity
- Integration to HIS
- 5 Year Service and Maintenance Agreement
- Project Management (7 DAYS)
- Configuration Service (17 days)
- Interfacing to HIS (5 DAYS)
- Installation and Upgrade Service (20 Hours)

WARRANTY 5-YEAR**Warranty Exclusions:**

*Excluding negligence and Cybersecurity attacks + Load shedding electrical damage

*Excludes 3rd Party Integration Cost

*Excludes Anti-Virus

**UNIT PRICE INCLUSIVE OF VAT
(PLEASE REFER TO PRICING SHEET 5-YEAR WARRANTY)**

WEB SITE LINK

[HTTPS://TINYURL.COM/2XW3P679](https://tinyurl.com/2XW3P679)

CENTRICITY™ HIGH ACUITY ANAESTHESIA

Scan me



* Please confirm list of supported medical devices with your Medhold sales representative



Centricity High Acuity Anesthesia

HELPING TO ENHANCE CLINICAL DECISION MAKING
THROUGHOUT THE PERI-OPERATIVE WORKFLOW

CHALLENGE:

In the peri-operative environment, ensuring timely access to comprehensive patient data is a challenge. Disparate systems and fragmented data capture are the norm, frustrating the care team and impacting both the quality and efficiency of care delivery. Performance metrics suffer at all levels, from the clinical team to the unit to the overall health system.

SOLUTION:

Now, more than ever, immediate access to comprehensive patient data and clinical decision tools is key to providing a safe, effective clinical anesthesia environment for the patient.

Centricity High Acuity Anesthesia is a specialized clinical information management system that enhances clinical documentation and decision-making from pre-operative assessment throughout the peri-operative workflow.

Up to **84%**
reduction in
medication dose,
route errors ¹

Helps to
reduce delays
in procedure
due to incomplete
information ²

Helps to
**increase
confidence**
having access to
needed data ²

28 hrs
of clinical staff time
gained per week ²

1. <http://www.ahrq.gov/qual/aderia/aderia.htm#Costs>
2. Almere Case study. 8 ORs, 7 holding and 13 recovery beds
Performance may vary for different customers.



MEDHOLD

People. Products. Partnerships.

RT2-2021

**SUPPLY, DELIVERY, INSTALLATION,
COMMISSIONING AND MAINTENANCE OF
MEDICAL EQUIPMENT**

**TO THE STATE FOR THE PERIOD
01 DECEMBER 2021
TO 01 DECEMBER 2024.**

SECTION D

INFECTION CONTROL

CONTACT DETAILS

Sibongile Leeu: s.leeu@medhold.co.za
Head Office: 011 966 0600

Or contact your local branch – details on the back page

ITEM NUMBER: RT2-10-05
ITEM NAME: HEAT SEALING MACHINE
GETINGE PROSEAL PLUS OFFER 2

Make/Brand	Getinge
Model	Proseal Plus
Medhold Part Number	6003000194
DOH SAP number	220620
Country of Origin	Germany

HEAT SEALING MACHINE

GETINGE PROSEAL PLUS PRODUCT SPECIFICATION

Getinge Proseal Plus rotary heat sealing devices are designed with focus on user friendliness and can master high volume operations easily. Getinge Proseal Plus is a validateable rotary heat sealer which means that the critical parameters are controlled and monitored according to ISO 11607-2 standard.

If any of the pre-set parameters are out of range there will be an automatic clear and easy to understand indication, which ensures that each package gets an approved seal.

UNIT PRICE INCLUSIVE OF VAT AS PER RT2 CIRCULAR PRICELIST

Yr 1 - 2 Warranty and Preventative Maintenance.

WEB SITE LINK
HEAT SEALING MACHINE
[HTTPS://TINYURL.COM/5N838HK5](https://tinyurl.com/5N838HK5)
HEAT SEALING MACHINE


ITEM NUMBER: RT2-10-05
ITEM NAME: HEAT SEALING MACHINE
GETINGE PROSEAL PREMIUM OFFER 3

Make/Brand	Getinge
Model	Proseal Premium
Medhold Part Number	6003000195
DOH SAP number	220621
Country of Origin	Germany

HEAT SEALING MACHINE

GETINGE PROSEAL PREMIUM PRODUCT SPECIFICATION

Getinge Proseal Premium sealing devices are designed with focus on user friendliness and can master high volume operations easily. Proseal Premium devices are equipped with an integrated printer and automatic pouch width sensor for secure printing of parameters on the packaging. Getinge Proseal Premium with internal printer is our top version of rotary heat sealers and is validateable according to ISO 11607-2. This means that the critical parameters are controlled and monitored. If any of the pre-set parameters are out of range there will be an automatic clear and easy to understand indication, which ensures that each package gets an approved seal.

UNIT PRICE INCLUSIVE OF VAT AS PER RT2 CIRCULAR PRICELIST

Yr 1 - 2 Warranty and Preventative Maintenance.

WEB SITE LINK
HEAT SEALING MACHINE
[HTTPS://TINYURL.COM/2P8J3V6C](https://tinyurl.com/2P8J3V6C)
HEAT SEALING MACHINE


ITEM NUMBER: RT2-13-02
ITEM NAME: INSTRUMENT WASHER, SINGLE DOOR
GETINGE 8668 WD SD OFFER 4

Make/Brand	Getinge
Model	Getinge Instrument Washer 8668 SD
Medhold Part Number	S-86682001
DOH SAP Number	220646
Country of Origin	Sweden

GETINGE 8668 WD SD

**GETINGE INSTRUMENT
WASHER 8668 SD
PRODUCT SPECIFICATION**
Key features

- Drying is optimized with the use of variable speed fans (energy vs. flow) that maximizes the drying performance.
- The washer-disinfector uses wash carts with removable spray arms and shelves that are adaptable to different load heights. The wash carts have incorporated accessories for cleaning lumen/cannulated/laparoscopic instruments.
- The standard washer-disinfector is equipped with a full-view, power operated, vertical glass door. Double door, pass through models may be selected as well as Single door options.
- The washer-disinfector is controlled via a 7", color touchscreen. On pass-through models there is a touchscreen on both the clean and the soiled side.
- The washer's docking interface is optimized for use with the SMART 2.0 loading trolley that provides both ergonomic and efficient loading/unloading.
- 5 Level Wash Cart


UNIT PRICE INCLUSIVE OF VAT AS PER RT2 CIRCULAR PRICELIST

Yr 1 - 2 Warranty and Preventative Maintenance.

WEB SITE LINK
GETINGE 8668 SD WD
[HTTPS://TINYURL.COM/5N84BMSW](https://tinyurl.com/5N84BMSW)

GETINGE 8668 SD WD

Scan me



ITEM NUMBER: RT2-13-03
ITEM NAME: INSTRUMENT WASHER, DOUBLE DOOR
GETINGE 8668 WD DD OFFER 4

Make/Brand	Getinge
Model	Getinge Instrument Washer 8668 DD
Medhold Part Number	S-86682003
DOH SAP Number	220647
Country of Origin	Sweden

GETINGE 8668 WD DD

**GETINGE INSTRUMENT
WASHER 8668 DD
PRODUCT SPECIFICATION**
Key features

- Drying is optimized with the use of variable speed fans (energy vs. flow) that maximizes the drying performance.
- The washer-disinfector uses wash carts with removable spray arms and shelves that are adaptable to different load heights. The wash carts have incorporated accessories for cleaning lumen/cannulated/laparoscopic instruments.
- The standard washer-disinfector is equipped with a full-view, power operated, vertical glass door. Double door, pass through models may be selected as well as Single door options.
- The washer-disinfector is controlled via a 7", color touchscreen. On pass-through models there is a touchscreen on both the clean and the soiled side.
- The washer's docking interface is optimized for use with the SMART 2.0 loading trolley that provides both ergonomic and efficient loading/unloading.


UNIT PRICE INCLUSIVE OF VAT AS PER RT2 CIRCULAR PRICELIST

Yr 1 - 2 Warranty and Preventative Maintenance.

**WEB SITE LINK
GETINGE 8668 SD WD**
[HTTPS://TINYURL.COM/5N84BMSW](https://tinyurl.com/5N84BMSW)

GETINGE 8668 SD WD

Scan me



RT2-2021

**SUPPLY, DELIVERY, INSTALLATION,
COMMISSIONING AND MAINTENANCE OF
MEDICAL EQUIPMENT**

**TO THE STATE FOR THE PERIOD
01 DECEMBER 2021
TO 01 DECEMBER 2024.**

SECTION E

MATERNAL INFANT CARE & VENTILATION

CONTACT DETAILS

Sibongile Leeu: s.leeu@medhold.co.za
Head Office: 011 966 0600

Or contact your local branch – details on the back page

VENTILATION

ITEM NUMBER: RT2-02-02
ITEM NAME: VENTILATOR PAEDIATRIC & ADULT

Make/Brand	GE Healthcare
Model	Carescape R860
Medhold Part Number	1506-8600-000
DOH SAP number	167812
Country of Origin	USA (Various suppliers for accessories)

UNIT INCLUDES:

- Patient circuit
- Reusable, sterilisable Paeds Mask
- Reusable, sterilisable Adult Mask
- Flow sensor.
- Test lung
- Test circuit
- Bacteria filter (pk of 5)
- Air filters (pk of 5)


UNIT PRICE INCLUSIVE OF VAT

Offer price Includes: Start-up accessories & consumables as per tender spec, Plus Yr 1 - 2 Warranty and Preventative Maintenance.

Refer to pricing sheet

OPTIONAL ITEMS

Extended Warranty and Maintenance Contract - 5 years
 Optional Year 3 to 7 Warranty & Preventative Maintenance Including VAT

DOH ITEM NUMBER
RT2-02-02
DESCRIPTION
 The unit must be able to ventilate and monitor paediatric and adult patients.

MEDHOLD PART NUMBER
1506-8600-000
WEB SITE LINK
VENTILATOR PAEDIATRIC / ADULT
[HTTPS://TINYURL.COM/2P8ZKVMY](https://tinyurl.com/2P8ZKVMY)

CARESCAPE R860

Scan me



VENTILATION

ITEM NUMBER: RT2-02-02

ITEM NAME: VENTILATOR PAEDIATRIC & ADULT

RT2-2021 Item Number	Ventilator Paediatric Accessories	Medhold Part Number	Quantity
RT2-02-02a	Reusable patient circuit (Box of 10)	2096534-007	Box of 10
RT2-02-02b	Reusable, sterilizable Paeds Mask (Inflatable)	ANFM8848002	Each
RT2-02-02c	Reusable, sterilizable Adult Mask (Inflatable)	ANFM8848004	Each
RT2-02-02d	Flow sensor.	1505-3848-000	Each
RT2-02-02e	Test lung	LO600.6832D	Each
RT2-02-02f	Test circuit	M1019541	Each
RT2-02-02g	Bacteria filter (pk of 5)	557022500	Pack of 5
RT2-02-02h	Air filters (pk of 5)	2083208-001	Pack of 5
RT2-02-02i	Adult disposable patient circuit (Pk of 10)	M1019541	Pack of 10
RT2-02-02j	Paeds disposable patient circuits (Pk of 10)	M1019544	Pack of 10
RT2-02-02k	Temperature probes for humidifier	FIS900MR869	Each
Upgrades			
RT2-02-02l	Built-in nebulizer	Built into ventilator	Each
RT2-02-02m	Compressor air: 0.2 to 0.7 bar	Built into ventilator trolley	Each
RT2-02-02n	Humidifier The humidifier must be dual servo controlled with an alarm system All accessories for the humidifier must be included Humidifier must have a hose heating element, dual temperature probe and electrical adapter.	FISMR850	Each
RT2-02-02o	End Tidal CO2 curve	M1184206	Each

MATERNAL INFANT CARE

ITEM NUMBER: RT2-03-07
ITEM NAME: GIRAFFE WARMER

Make/Brand	GE Healthcare
Model	Giraffe Warmer
Medhold Part Number	M1118179
DOH SAP number	151321
Country of Origin	USA

UNIT INCLUDES:

- Main warmer unit with recessed head and fixed base
- Aimable procedure light
- Hands-free silence alarm
- Observation light
- Temperature probe panel with reusable temperature probe
- 1 x box of 50 x disposable heat reflecting patches
- Single side to side drawer
- Integrated T-Piece Resuscitation with blender
- 10 x disposable T-Piece Circuit kits
- IV Pole
- Integrated X-ray tray
- Monitor Tray
- E-cylinder holder
- Giraffe starter pack for developmental care
- User manual


UNIT PRICE INCLUSIVE OF VAT

Offer price Includes: Start-up accessories & consumables as per tender spec, Plus Yr 1 - 2 Warranty and Preventative Maintenance.

Refer to pricing sheet
OPTIONAL ITEMS

Extended Warranty and Maintenance Contract - 5 years
 Optional Year 3 to 7 Warranty & Preventative Maintenance Including VAT

DOH ITEM NUMBER
RT2-03-07
DESCRIPTION
ICU CRIB:
 A crib used to provide patient support and relief in the intensive care unit.

MEDHOLD PART NUMBER
M1118179
WEB SITE LINK
GIRAFFE WARMER
[HTTPS://TINYURL.COM/YD55C2Y3](https://tinyurl.com/YD55C2Y3)
GIRAFFE WARMER
Scan me


MATERNAL INFANT CARE

ITEM NUMBER: RT2-03-07

ITEM NAME: GIRAFFE WARMER

RT2-2021 Item Number	ICU Crib Warmer Accessories	Medhold Part Number	Quantity
RT2-03-07a	Temperature probe	2075796-001	Each
RT2-03-07c	IV pole	2043931-001	Each
RT2-03-07d	Shelf for monitor	2043933-001	Each
RT2-03-07e	Instrument tray	6600-0865-700	Each
RT2-03-07f	Battery	TBC	Each
RT2-03-07g	Air/oxygen blender (UPGRADE)	Built in	Each
RT2-03-07h	Built-in scale (UPGRADE)	M1125056	Each
RT2-03-07i	LED Phototherapy (UPGRADE)	2059118-001	Each
RT2-03-07j	Pulse Oxymeter (UPGRADE)	2064253-001	Each
RT2-03-07k	Disposable covers for skin temperature probe (box of 10)	2031980300	Each
RT2-03-07l	Integrated suction (UPGRADE)	Built in, integrated suction	Each

MATERNAL INFANT CARE

ITEM NUMBER: RT2-03-08
ITEM NAME: INCUBATOR ICU

Make/Brand	GE Healthcare
Model	Giraffe Incubator Carestation
Medhold Part Number	2082844-002
DOH SAP number	154907
Country of Origin	USA

UNIT INCLUDES:

- Closed incubator unit with elevating base
- Air boost protects the baby from heat loss
- Hands-free silence alarm
- Pressure diffusing mattress
- 1 x Temperature probe panel with reusable temperature probe
- 1 x Box of 50 x Disposable heat reflecting patches
- Single side to side drawer
- Patented Baby Susan rotating mattress
- Integrated X-ray tray
- Giraffe starter pack for developmental care
- User manual
- Internal battery for alarms
- Giraffe Instrument Shelf
- E-Cylinder Holder-for Careplus, Giraffe & IWS
- IV Pole
- Bacterial Filters (10)



UNIT PRICE INCLUSIVE OF VAT

Offer price Includes: Start-up accessories & consumables as per tender spec, Plus Yr 1 - 2 Warranty and Preventative Maintenance.

Refer to pricing sheet

OPTIONAL ITEMS

Extended Warranty and Maintenance Contract - 5 years
 Optional Year 3 to 7 Warranty & Preventative Maintenance Including VAT

DOH ITEM NUMBER
RT2-03-08

DESCRIPTION

Incubator ICU: A unit that provides a closed, controlled environment that warms an infant by circulating heated air around the infant's skin and is to be used in the neonatal intensive care unit.

MEDHOLD PART NUMBER
2082844-002

WEB SITE LINK
INCUBATOR ICU

[HTTPS://TINYURL.COM/3JTWH46S](https://tinyurl.com/3JTWH46S)

INCUBATOR ICU

Scan me



MATERNAL INFANT CARE

ITEM NUMBER: RT2-03-08

ITEM NAME: INCUBATOR ICU

RT2-2021 Item Number	Incubator ICU Accessories	Medhold Part Number	Quantity
RT2-03-08a	Bacterial filter (Pk of 4)	2107356-001	Each pack of 4
RT2-03-08b	Iris port sleeves (Pk of 10)	6600-0211-850	Each pack of 4
RT2-03-08d	Skin temperature probe (reusable)	2075796-001	Each
RT2-03-08e	Covers for temperature probe (Pk of 25)	2044385-001	Each pack of 25
RT2-03-08f	Internal battery for alarms	TBC	Each
RT2-03-08g	Scale (UPGRADE)	M1125056	Each
RT2-03-08h	Observation lamp (UPGRADE)	2111193-001	Each

MATERNAL INFANT CARE

ITEM NUMBER: RT2-03-09
ITEM NAME: INCUBATOR (STANDARD)

Make/Brand	GE Healthcare
Model	Giraffe Incubator Carestation
Medhold Part Number	2082844-002
DOH SAP number	155124
Country of Origin	USA

UNIT INCLUDES:

- Closed incubator unit with elevating base
- Air boost protects the baby from heat loss
- Hands-free silence alarm
- Pressure diffusing mattress
- 1 x Temperature probe panel with reusable temperature probe
- 1 x Box of 50 x Disposable heat reflecting patches
- Single side to side drawer
- Patented Baby Susan rotating mattress
- Integrated X-ray tray
- Giraffe starter pack for developmental care
- User manual
- Internal battery for alarms
- Giraffe Instrument Shelf
- E-Cylinder Holder-for Careplus, Giraffe & IWS
- IV Pole
- Bacterial Filters (10)


UNIT PRICE INCLUSIVE OF VAT

Offer price Includes: Start-up accessories & consumables as per tender spec, Plus Yr 1 - 2 Warranty and Preventative Maintenance.

Refer to pricing sheet
OPTIONAL ITEMS

Extended Warranty and Maintenance Contract - 5 years
 Optional Year 3 to 7 Warranty & Preventative Maintenance Including VAT

**DOH ITEM NUMBER
 RT2-03-09**
DESCRIPTION
 Incubator Standard: A unit that provides a closed, controlled environment that warms an infant by circulating heated air around the infant's skin.

**MEDHOLD PART NUMBER
 2082844-002**
**WEB SITE LINK
 INCUBATOR (STANDARD)**
[HTTPS://TINYURL.COM/3JTWH46S](https://tinyurl.com/3JTWH46S)
INCUBATOR (STANDARD)
Scan me


MATERNAL INFANT CARE

ITEM NUMBER: RT2-03-09

ITEM NAME: INCUBATOR (STANDARD)

RT2-2021 Item Number	Incubator Standard Accessories	Medhold Part Number	Quantity
RT2-03-09a	Bacterial filter (Pk of 4)	2107356-001	Each pack of 4
RT2-03-09b	Iris port sleeves (Pk of 10)	6600-0211-850	Each pack of 10
RT2-03-09d	Skin temperature probe (reusable)	2075796-001	Each
RT2-03-09e	Covers for temperature probe (Pk of 25)	2044385-001	Each pack of 25
RT2-03-09f	Internal battery for alarms	TBC	Each
RT2-03-09g	Scale (UPGRADE)	M1125056	Each
RT2-03-09h	Observation lamp (UPGRADE)	2111193-001	Each
RT2-03-09i	Blanket cover	6600-0846-800	Each
RT2-03-09j	X-ray tray	Built in	Each
RT2-03-09k	Accessory rail	Built in	Each

MATERNAL INFANT CARE

ITEM NUMBER: RT2-03-11
ITEM NAME: INFANT RADIANT WARMER

Make/Brand	GE Healthcare
Model	Lullaby Warmer
Medhold Part Number	2091100-001
DOH SAP number	149517
Country of Origin	India

UNIT INCLUDES:

- Main warmer unit with heater head and fixed base
- Lockable castors
- Heater head swivel 90 degree L or R
- Calrod heater element
- Examination lamp
- Integrated bed with firm mattress
- Continuous bed tilt with lock
- LED display control panel
- Manual and baby mode temperature measurement
- Temperature probe panel with reusable temperature probe
- 1 x box of 50 x disposable heat reflecting patches
- Dual front drawers
- Apgar timer
- Integrated X-ray tray
- IV Pole
- Utility pole
- Monitor Tray
- E-cylinder oxygen holder
- Giraffe starter pack for developmental care
- User Manual


UNIT PRICE INCLUSIVE OF VAT

Offer price Includes: Start-up accessories & consumables as per tender spec, Plus Yr 1 - 2 Warranty and Preventative Maintenance.

Refer to pricing sheet
OPTIONAL ITEMS

Extended Warranty and Maintenance Contract - 5 years
 Optional Year 3 to 7 Warranty & Preventative Maintenance Including VAT

DOH ITEM NUMBER
RT2-03-11
DESCRIPTION
 Infant Radiant Warmer: A warming device to provide uniform heat to the body using an overhead heating source.

MEDHOLD PART NUMBER
2091100-001
WEB SITE LINK
INFANT RADIANT WARMER
[HTTPS://TINYURL.COM/BDD59TBC](https://tinyurl.com/BDD59TBC)

INFANT RADIANT WARMER

Scan me



MATERNAL INFANT CARE

ITEM NUMBER: RT2-03-11

ITEM NAME: INFANT RADIANT WARMER

RT2-2021 Item Number	Infant Radiant Warmer Accessories	Medhold Part Number	Quantity
RT2-03-11a	Reusable Temperature probe	6600-0628-700	Each
RT2-03-11b	Mattress - Lullaby	0305-5060-300	Each
RT2-03-11b	Mattress - Giraffe	0305-5061-300	Each
RT2-03-11c	LED Phototherapy (UPGRADE)	2052663-001	Each
RT2-03-11d	IV pole	2043931-001	Each
RT2-03-11e	Disposable temperature probe (Pk of 10)	6600-0208-700	Each pack of 10
RT2-03-11f	Accessory Rail	Built in	Each
RT2-03-11g	Monitor stand	2043933-001	Each

MATERNAL INFANT CARE

ITEM NUMBER: RT2-03-11C
ITEM NAME: LULLABY LED PHOTOTHERAPY LIGHT

Make/Brand	GE Healthcare
Model	Lullaby LED Phototherapy Light
Medhold Part Number	2059118-001
DOH SAP number	149556
Country of Origin	India

UNIT INCLUDES:

- Lullaby LED Phototherapy
- Lullaby LED Phototherapy - Basic Unit (incl. Pole, bank light, service manual)
- Lullaby Phototherapy Power cord - India
- Lullaby LED Phototherapy - Operator Manual (English)

UNIT PRICE INCLUSIVE OF VAT

Offer price Includes: Start-up accessories & consumables as per tender spec,
 Plus Yr 1 - 2 Warranty and Preventative Maintenance.

Refer to pricing sheet

OPTIONAL ITEMS

Extended Warranty and Maintenance Contract - 5 years

Optional Year 3 to 7 Warranty & Preventative

Maintenance Including VAT


DOH ITEM NUMBER
RT2-03-11C
DESCRIPTION
 Phototherapy light (LED)

MEDHOLD PART NUMBER
2059118-001
WEB SITE LINK
PHOTOTHERAPY LIGHT (LED)
[HTTPS://TINYURL.COM/T37KNDMK](https://tinyurl.com/T37KNDMK)

PHOTOTHERAPY LIGHT (LED)

Scan me



MATERNAL INFANT CARE

ITEM NUMBER: RT2-03-13

ITEM NAME: BILISOFT 2.0

Make/Brand	GE Healthcare
Model	Bilisoft 2.0
Medhold Part Number	2108683-001
DOH SAP number	154866
Country of Origin	USA

UNIT INCLUDES:

- Bilisoft PT System English South Africa Plug Small Pad
Bilisoft Pad Cover, Disposable, Small (Box of 20)

UNIT PRICE INCLUSIVE OF VAT

Offer price Includes: Start-up accessories & consumables as per tender spec,
Plus Yr 1 - 2 Warranty and Preventative Maintenance.

Refer to pricing sheet

OPTIONAL ITEMS

Extended Warranty and Maintenance Contract - 5 years

Optional Year 3 to 7 Warranty & Preventative

Maintenance Including VAT



DOH ITEM NUMBER
RT2-03-13

DESCRIPTION
Phototherapy blanket LED:

MEDHOLD PART NUMBER
2108683-001

WEB SITE LINK
BILISOFT 2.0

[HTTPS://TINYURL.COM/4TKP9HWX](https://tinyurl.com/4TKP9HWX)

BILISOFT 2.0

Scan me



MATERNAL INFANT CARE

ITEM NUMBER: RT2-03-13

ITEM NAME: BILISOFT 2.0

RT2-2021 Item Number	Infant Radiant Warmer Accessories	Medhold Part Number	Quantity
RT2-03-13b	Effective treatment area: Large pad of 450-550cm2	2108752-001	Each



RT2-2021

**SUPPLY, DELIVERY, INSTALLATION,
COMMISSIONING AND MAINTENANCE OF
MEDICAL EQUIPMENT**

**TO THE STATE FOR THE PERIOD
01 DECEMBER 2021
TO 01 DECEMBER 2024.**

SECTION F

PATIENT MONITORING

CONTACT DETAILS

Sibongile Leeu: s.leeu@medhold.co.za
Head Office: 011 966 0600

Or contact your local branch – details on the back page

BASIC MONITOR

ITEM NUMBER: RT2-05-14_OPTION 6
ITEM NAME: BASIC MONITOR

Make/Brand	GE Healthcare
Model	CARESCAPE B650 Monitor with PDM Module
Medhold Part Number	5805688
DOH SAP number	154872
Country of Origin	Finland

PATIENT MONITOR INCLUDES:

B650 Monitor with PDM module

- Multi-link standard ECG cable
- 3 leadwire set with grabbers
- NIBP hose
- Reusable infant, child, small adult, adult and large adult cuffs
- SPO2 Interface cable
- Adult finger probe
- Multisite SPO2 sensor
- Temperature adaptor
- Oral and Skin temperature sensors
(No mounting bracket included)


WARRANTY (2-YEAR): YEAR 1 AND YEAR 2

Warranty Entails: A 2-year warranty against poor workmanship and latent defects and parts. This must be all inclusive and include, BUT NOT LIMITED TO, amongst others, ALL PARTS (including Batteries), Labour, Traveling and Accommodation. The 2-year warranty must also include all quality check, quality assurance requirements, Preventative Maintenance/Calibrations, Software updates and upgrades to be included. This 2-year warranty will commence after formal Commissioning and handover of the equipment.

UNIT PRICE INCLUSIVE OF VAT (STANDARD ACCESSORIES & CONSUMABLES, AND 2-YEAR WARRANTY)

See attached pricing sheet

WEB SITE LINK <https://tinyurl.com/YRXXX269>

Scan me


OPTIONAL ACCESSORIES:

DOH Item Number	Description	Medhold Part Number	Quantity
RT2-05-14a	Paediatric 3/5 lead ECG cables	412682-003	Each
RT2-05-14b	Adult 3/5 lead ECG cables	412682-004	Each
RT2-05-14c	ECG interface cable	2017003-003	Each
RT2-05-14d	Reusable latex-free paediatric cuff	DUR-P2-2A-E	Each
RT2-05-14e	Reusable latex-free adult cuff	DUR-A2-2A-E	Each
RT2-05-14f	Reusable latex-free neonate cuff	DUR-P1-2A-E	Each
RT2-05-14g	Reusable latex-free obese cuff	DUR-T1-2A-E	Each
RT2-05-14h	NIBP hose	2058203-002	Each
RT2-05-14i	Finger sensor for Paeds (Nellcor)	D-YS	Each
RT2-05-14i-1	Finger sensor for Paeds (Masimo)	4051	Each
RT2-05-14j	Finger sensor for adults (Nellcor)	FLEXMAX	Each
RT2-05-14j-1	Finger sensor for adults (Masimo)	4050	Each
RT2-05-14k	Finger wrap-around probe for neonates (Nellcor)	D-YS	Each
RT2-05-14k-1	Finger wrap-around probe for neonates (Masimo)	4054	Each
RT2-05-14l	Spo2 interface cable (Nellcor)	2021406-001	Each
RT2-05-14l-1	Spo2 interface cable (Masimo)	MASIMO-RD-4085	Each
RT2-05-14m	Two reusable temperature probes	M1024254 / M1024247	Each
RT2-05-14	Monitor Front bezel	2109443-001	Each
RT2-05-14o	Battery	M1168356	Each
RT2-05-14p	NIBP pump module	2031069-007	Each
RT2-05-14q	Wall bracket, height adjustable on which the monitor can be mounted.	WC-0002-05 / WMM-0002-01C/02C/03D/ WMM-0006-04D	Each
RT2-05-14r	Roll stand	ME-0030-63-DCK	Each

CENTRAL MONITORING

ITEM NUMBER: RT2-05-17_OPTION 6
ITEM NAME: CENTRAL MONITORING, 8-BED PACKAGE

Make/Brand	GE Healthcare
Model	CARESCAPE B650 Monitor with PDM Module and 8 Bed CARESCAPE Central Station
Medhold Part Number	5805688 & 5514385
DOH SAP number	220550
Country of Origin	Finland and USA

PATIENT MONITOR INCLUDES:

B650 Monitor with PDM module

- Multi-link standard ECG cable
- 3 leadwire set with grabbers
- NIBP hose and reusable adult cuff
- SPO2 Interface cable
- Adult finger probe
- Temperature adaptor
- Skin temperature sensors
- Dual invasive pressure adaptor
- 3 invasive pressure cables (no consumables)
- Wall mount bracket with channel and cable hooks


CENTRAL STATION INCLUDES:

CARESCAPE Central Station MAI800 Integrated System

- Keyboard and Mouse
- Dual 54.61 cm (21, 5 in) medical grade flat panel displays
- Laser printer
- UPS (Uninterrupted Power Supply)
- Networking
- 8 bed licenses
- 144 hours Full disclosure per bed


WARRANTY (2-YEAR): YEAR 1 AND YEAR 2

Warranty Entails: A 2-year warranty against poor workmanship and latent defects and parts. This must be all inclusive and include, BUT NOT LIMITED TO, amongst others, ALL PARTS (including Batteries), Labour, Traveling and Accommodation. The 2-year warranty must also include all quality check, quality assurance requirements, Preventative Maintenance/Calibrations, Software updates and upgrades to be included. This 2-year warranty will commence after formal Commissioning and handover of the equipment.

UNIT PRICE INCLUSIVE OF VAT
(STANDARD ACCESSORIES & CONSUMABLES, AND 2-YEAR WARRANTY)

See attached pricing sheet

WEB SITE LINK
B 650
[HTTPS://TINYURL.COM/YRXXX269](https://tinyurl.com/YRXXX269)
CENTRAL STATION
[HTTPS://TINYURL.COM/22N7RPHR](https://tinyurl.com/22N7RPHR)

B650

Scan me



CENTRAL STATION

Scan me



CENTRAL MONITORING

ITEM NUMBER: RT2-05-17_OPTION 9
ITEM NAME: CENTRAL MONITORING, 8-BED PACKAGE

Make/Brand	GE Healthcare
Model	CARESCAPE B650 Monitor with CS1 Module and 8 Bed CARESCAPE Central Station
Medhold Part Number	5805688 & 5514385
DOH SAP number	220551
Country of Origin	Finland and USA

PATIENT MONITOR INCLUDES:

B650 Monitor with CS1 module

- 3 leadwire set with grabbers
- NIBP hose and reusable adult cuff
- Adult finger probe
- Skin temperature sensors
- 2 invasive pressure cables (no consumables)
- Wall mount bracket with channel and cable hooks


CENTRAL STATION INCLUDES:

CARESCAPE Central Station MAI800 Integrated System

- Keyboard and Mouse
- Dual 54.61 cm (21, 5 in) medical grade flat panel displays
- Laser printer
- UPS (Uninterrupted Power Supply)
- Networking
- 8 bed licenses
- 144 hours Full disclosure per bed


WARRANTY (2-YEAR): YEAR 1 AND YEAR 2

Warranty Entails: A 2-year warranty against poor workmanship and latent defects and parts. This must be all inclusive and include, BUT NOT LIMITED TO, amongst others, ALL PARTS (including Batteries), Labour, Traveling and Accommodation.

The 2-year warranty must also include all quality check, quality assurance requirements, Preventative Maintenance/Calibrations, Software updates and upgrades to be included. This 2-year warranty will commence after formal Commissioning and handover of the equipment.

UNIT PRICE INCLUSIVE OF VAT
(STANDARD ACCESSORIES & CONSUMABLES, AND 2-YEAR WARRANTY)


See attached pricing sheet

WEB SITE LINK
B650
[HTTPS://TINYURL.COM/YRXXX269](https://tinyurl.com/YRXXX269)
CENTRAL STATION
[HTTPS://TINYURL.COM/22N7RPHR](https://tinyurl.com/22N7RPHR)

B650



CENTRAL STATION



CENTRAL MONITORING

ITEM NUMBER: RT2-05-17a_OPTION 6
ITEM NAME: CENTRAL MONITORING, 12-BED PACKAGE

Make/Brand	GE Healthcare
Model	CARESCAPE B650 Monitor with PDM Module and 12 Bed CARESCAPE Central Station
Medhold Part Number	5805688 & 5514385
DOH SAP number	220564
Country of Origin	Finland and USA

PATIENT MONITOR INCLUDES:

B650 Monitor with PDM module

- Multi-link standard ECG cable
- 3 leadwire set with grabbers
- NIBP hose and reusable adult cuff
- SPO2 Interface cable
- Adult finger probe
- Temperature adaptor
- Skin temperature sensors
- Dual invasive pressure adaptor
- 3 invasive pressure cables (no consumables)
- Wall mount bracket with channel and cable hooks


CENTRAL STATION INCLUDES:

CARESCAPE Central Station MAI800 Integrated System

- Keyboard and Mouse
- Dual 54.61 cm (21, 5 in) medical grade flat panel displays
- Laser printer
- UPS (Uninterrupted Power Supply)
- Networking
- 12 bed licenses
- 144 hours Full disclosure per bed


WARRANTY (2-YEAR): YEAR 1 AND YEAR 2

Warranty Entails: A 2-year warranty against poor workmanship and latent defects and parts. This must be all inclusive and include, BUT NOT LIMITED TO, amongst others, ALL PARTS (including Batteries), Labour, Traveling and Accommodation. The 2-year warranty must also include all quality check, quality assurance requirements, Preventative Maintenance/Calibrations, Software updates and upgrades to be included. This 2-year warranty will commence after formal Commissioning and handover of the equipment.

UNIT PRICE INCLUSIVE OF VAT
(STANDARD ACCESSORIES & CONSUMABLES, AND 2-YEAR WARRANTY)

See attached pricing sheet

WEB SITE LINK
B650
[HTTPS://TINYURL.COM/YRXXX269](https://tinyurl.com/YRXXX269)
CENTRAL STATION
[HTTPS://TINYURL.COM/22N7RPHR](https://tinyurl.com/22N7RPHR)

B650

Scan me



CENTRAL STATION

Scan me



CENTRAL MONITORING

ITEM NUMBER: RT2-05-17a_OPTION 9
ITEM NAME: CENTRAL MONITORING, 12-BED PACKAGE

Make/Brand	GE Healthcare
Model	CARESCAPE B650 Monitor with CS1 Module and 12 Bed CARESCAPE Central Station
Medhold Part Number	5805688 & 5514385
DOH SAP number	220563
Country of Origin	Finland and USA

PATIENT MONITOR INCLUDES:

B650 Monitor with CS1 module

- 3 leadwire set with grabbers
- NIBP hose and reusable adult cuff
- Adult finger probe
- Skin temperature sensors
- 2 invasive pressure cables (no consumables)
- Wall mount bracket with channel and cable hooks


CENTRAL STATION INCLUDES:

CARESCAPE Central Station MAI800 Integrated System

- Keyboard and Mouse
- Dual 54.61 cm (21, 5 in) medical grade flat panel displays
- Laser printer
- UPS (Uninterrupted Power Supply)
- Networking
- 12 bed licenses
- 144 hours Full disclosure per bed


WARRANTY (2-YEAR): YEAR 1 AND YEAR 2

Warranty Entails: A 2-year warranty against poor workmanship and latent defects and parts. This must be all inclusive and include, BUT NOT LIMITED TO, amongst others, ALL PARTS (including Batteries), Labour, Traveling and Accommodation. The 2-year warranty must also include all quality check, quality assurance requirements, Preventative Maintenance/Calibrations, Software updates and upgrades to be included. This 2-year warranty will commence after formal Commissioning and handover of the equipment.

UNIT PRICE INCLUSIVE OF VAT
(STANDARD ACCESSORIES & CONSUMABLES, AND 2-YEAR WARRANTY)

See attached pricing sheet

WEB SITE LINK
B 650
[HTTPS://TINYURL.COM/YRXXX269](https://tinyurl.com/YRXXX269)
CENTRAL STATION
[HTTPS://TINYURL.COM/22N7RPHR](https://tinyurl.com/22N7RPHR)
B650

Scan me


CENTRAL STATION

Scan me



CENTRAL MONITORING


ITEM NUMBER: RT2-05-17b_OPTION 6
ITEM NAME: CENTRAL MONITORING, 16-BED PACKAGE

Make/Brand	GE Healthcare
Model	CARESCAPE B650 Monitor with PDM Module and 16 Bed CARESCAPE Central Station
Medhold Part Number	5805688 & 5514385
DOH SAP number	220567
Country of Origin	Finland and USA

PATIENT MONITOR INCLUDES:

B650 Monitor with PDM module

- Multi-link standard ECG cable
- 3 leadwire set with grabbers
- NIBP hose and reusable adult cuff
- SPO2 Interface cable
- Adult finger probe
- Temperature adaptor
- Skin temperature sensors
- Dual invasive pressure adaptor
- 3 invasive pressure cables (no consumables)
- Wall mount bracket with channel and cable hooks


CENTRAL STATION INCLUDES:

CARESCAPE Central Station MAI800 Integrated System

- Keyboard and Mouse
- Dual 54.61 cm (21, 5 in) medical grade flat panel displays
- Laser printer
- UPS (Uninterrupted Power Supply)
- Networking
- 16 bed licenses
- 144 hours Full disclosure per bed


WARRANTY (2-YEAR): YEAR 1 AND YEAR 2

Warranty Entails: A 2-year warranty against poor workmanship and latent defects and parts. This must be all inclusive and include, BUT NOT LIMITED TO, amongst others, ALL PARTS (including Batteries), Labour, Traveling and Accommodation. The 2-year warranty must also include all quality check, quality assurance requirements, Preventative Maintenance/Calibrations, Software updates and upgrades to be included. This 2-year warranty will commence after formal Commissioning and handover of the equipment.

UNIT PRICE INCLUSIVE OF VAT
(STANDARD ACCESSORIES & CONSUMABLES, AND 2-YEAR WARRANTY)

See attached pricing sheet

WEB SITE LINK
B 650
[HTTPS://TINYURL.COM/YRXXX269](https://tinyurl.com/YRXXX269)
CENTRAL STATION
[HTTPS://TINYURL.COM/22N7RPHR](https://tinyurl.com/22N7RPHR)

B650

Scan me



CENTRAL STATION

Scan me



CENTRAL MONITORING

ITEM NUMBER: RT2-05-17b_OPTION 9
ITEM NAME: CENTRAL MONITORING, 16-BED PACKAGE

Make/Brand	GE Healthcare
Model	CARESCAPE B650 Monitor with CS1 Module and 16 Bed CARESCAPE Central Station
Medhold Part Number	5805688 & 5514385
DOH SAP number	220569
Country of Origin	Finland and USA

PATIENT MONITOR INCLUDES:

B650 Monitor with CS1 module

- 3 leadwire set with grabbers
- NIBP hose and reusable adult cuff
- Adult finger probe
- Skin temperature sensors
- 2 invasive pressure cables (no consumables)
- Wall mount bracket with channel and cable hooks


CENTRAL STATION INCLUDES:

CARESCAPE Central Station MAI800 Integrated System

- Keyboard and Mouse
- Dual 54.61 cm (21, 5 in) medical grade flat panel displays
- Laser printer
- UPS (Uninterrupted Power Supply)
- Networking
- 16 bed licenses
- 144 hours Full disclosure per bed


WARRANTY (2-YEAR): YEAR 1 AND YEAR 2

Warranty Entails: A 2-year warranty against poor workmanship and latent defects and parts. This must be all inclusive and include, BUT NOT LIMITED TO, amongst others, ALL PARTS (including Batteries), Labour, Traveling and Accommodation. The 2-year warranty must also include all quality check, quality assurance requirements, Preventative Maintenance/Calibrations, Software updates and upgrades to be included. This 2-year warranty will commence after formal Commissioning and handover of the equipment.

UNIT PRICE INCLUSIVE OF VAT
(STANDARD ACCESSORIES & CONSUMABLES, AND 2-YEAR WARRANTY)

See attached pricing sheet

WEB SITE LINK
B650
[HTTPS://TINYURL.COM/YRXXX269](https://tinyurl.com/YRXXX269)
CENTRAL STATION
[HTTPS://TINYURL.COM/22N7RPHR](https://tinyurl.com/22N7RPHR)

B650



CENTRAL STATION



VITAL SIGNS MONITORING

ITEM NUMBER: RT2-05-19_OPTION 4
ITEM NAME: VITAL SIGNS MONITOR WITH OXIMETER AND TEMP

Make/Brand	GE Healthcare
Model	CARESCAPE V100 (TruSignal SPO2)
Medhold Part Number	2038172-521
DOH SAP number	220581
Country of Origin	USA (Mexico)

UNIT INCLUDES:

CARESCAPE V100 monitor

- NIBP hose, 3.65 metres long
- Reusable infant, child, small adult, adult and large adult Dura-Cuffs
- Ohmeda pulse oximetry cable
- Ohmeda reusable finger probe
- Ohmeda multisite sensor
- Reusable Exergen Temporal scanner
- Power supply with bracket and power cord
- User and technical manuals
- Roll stand with accessories basket


WARRANTY (2-YEAR): YEAR 1 AND YEAR 2

Warranty Entails: A 2-year warranty against poor workmanship and latent defects and parts. This must be all inclusive and include, BUT NOT LIMITED TO, amongst others, ALL PARTS (including Batteries), Labour, Traveling and Accommodation. The 2-year warranty must also include all quality check, quality assurance requirements, Preventative Maintenance/Calibrations, Software updates and upgrades to be included. This 2-year warranty will commence after formal Commissioning and handover of the equipment.

UNIT PRICE INCLUSIVE OF VAT
(STANDARD ACCESSORIES & CONSUMABLES, AND 2-YEAR WARRANTY)

See attached pricing sheet

WEB SITE LINK
[HTTPS://TINYURL.COM/4BJJW9AP](https://tinyurl.com/4BJJW9AP)

V100

Scan me


OPTIONAL ACCESSORIES:

DOH Item Number	Description	Medhold Part Number	Quantity
RT2-05-19a	Reusable latex-free paediatric cuff	DUR-P2-2A-E	Each
RT2-05-19b	Reusable latex-free adult cuff	DUR-A2-2A-E	Each
RT2-05-19c	Reusable latex-free neonate cuff	DUR-P1-2A-E	Each
RT2-05-19d	Reusable latex-free obese cuff	DUR-T1-2A-E	Each
RT2-05-19e	NIBP hose (adult/Paeds)	2058205-001	Each
RT2-05-19f	Finger sensor for Paeds	TS-W-D	Each
RT2-05-19g	Finger sensor for adults	TS-F-D	Each
RT2-05-19h	Finger wrap-around probe for neonates	TS-SE-3	Each
RT2-05-19i	Spo2 interface cable	TS-G3	Each
RT2-05-19j	Two reusable temperature probes	Exergen Temporal Scanner	Each
RT2-05-19k	Front bezel	2037103-062	Each
RT2-05-19l	Back cover	2037103-002	Each
RT2-05-19m	Battery	PU044	Each
RT2-05-19n	NIBP pump/module	2037103-015	Each

VITAL SIGNS MONITORING

ITEM NUMBER: RT2-05-19_OPTION 5
ITEM NAME: VITAL SIGNS MONITOR WITH OXIMETER AND TEMP

Make/Brand	GE Healthcare
Model	CARESCAPE V100 (Nellcor SPO2)
Medhold Part Number	2038172-529
DOH SAP number	220585
Country of Origin	USA (Mexico)

UNIT INCLUDES:

CARESCAPE V100 monitor

- NIBP hose, 3.65 metres long
- Reusable infant, child, small adult, adult and large adult Dura-Cuffs
- Nellcor pulse oximetry cable
- Nellcor reusable finger probe
- Nellcor multisite sensor
- Reusable Exergen Temporal scanner
- Power supply with bracket and power cord
- User and technical manuals
- Roll stand with accessories basket


WARRANTY (2-YEAR): YEAR 1 AND YEAR 2

Warranty Entails: A 2-year warranty against poor workmanship and latent defects and parts. This must be all inclusive and include, BUT NOT LIMITED TO, amongst others, ALL PARTS (including Batteries), Labour, Traveling and Accommodation. The 2-year warranty must also include all quality check, quality assurance requirements, Preventative Maintenance/Calibrations, Software updates and upgrades to be included. This 2-year warranty will commence after formal Commissioning and handover of the equipment.

UNIT PRICE INCLUSIVE OF VAT
(STANDARD ACCESSORIES & CONSUMABLES, AND 2-YEAR WARRANTY) See attached pricing sheet

WEB SITE LINK
[HTTPS://TINYURL.COM/4BJJW9AP](https://tinyurl.com/4BJJW9AP)

V100

Scan me


OPTIONAL ACCESSORIES:

DOH Item Number	Description	Medhold Part Number	Quantity
RT2-05-19a	Reusable latex-free paediatric cuff	DUR-P2-2A-E	Each
RT2-05-19b	Reusable latex-free adult cuff	DUR-A2-2A-E	Each
RT2-05-19c	Reusable latex-free neonate cuff	DUR-P1-2A-E	Each
RT2-05-19d	Reusable latex-free obese cuff	DUR-T1-2A-E	Each
RT2-05-19e	NIBP hose (adult/Paeds)	2058205-001	Each
RT2-05-19f	Finger sensor for Paeds	D-YS	Each
RT2-05-19g	Finger sensor for adults	DS-100A	Each
RT2-05-19h	Finger wrap-around probe for neonates	D-YS	Each
RT2-05-19i	Spo2 interface cable	2021406-001	Each
RT2-05-19j	Two reusable temperature probes	Exergen Temporal Scanner	Each
RT2-05-19k	Front bezel	2037103-062	Each
RT2-05-19l	Back cover	2037103-002	Each
RT2-05-19m	Battery	PU044	Each
RT2-05-19n	NIBP pump/module	2037103-015	Each

VITAL SIGNS MONITORING

ITEM NUMBER: RT2-05-19_OPTION 6

ITEM NAME: VITAL SIGNS MONITOR WITH OXIMETER AND TEMP

Make/Brand	GE Healthcare
Model	CARESCAPE V100 (Masimo SPO2)
Medhold Part Number	2038172-525
DOH SAP number	220599
Country of Origin	USA (Mexico)

UNIT INCLUDES:

CARESCAPE V100 monitor

- NIBP hose, 3.65 metres long
- Reusable infant, child, small adult, adult and large adult Dura-Cuffs
- Masimo pulse oximetry cable
- Masimo reusable finger probe
- Masimo multisite sensor
- Reusable Exergen Temporal scanner
- Power supply with bracket and power cord
- User and technical manuals
- Roll stand with accessories basket



WARRANTY (2-YEAR): YEAR 1 AND YEAR 2

Warranty Entails: A 2-year warranty against poor workmanship and latent defects and parts. This must be all inclusive and include, BUT NOT LIMITED TO, amongst others, ALL PARTS (including Batteries), Labour, Traveling and Accommodation. The 2-year warranty must also include all quality check, quality assurance requirements, Preventative Maintenance/Calibrations, Software updates and upgrades to be included. This 2-year warranty will commence after formal Commissioning and handover of the equipment.

UNIT PRICE INCLUSIVE OF VAT

(STANDARD ACCESSORIES & CONSUMABLES, AND 2-YEAR WARRANTY) See attached pricing sheet

WEB SITE LINK

[HTTPS://TINYURL.COM/4BJJW9AP](https://tinyurl.com/4BJJW9AP)

V100

Scan me



OPTIONAL ACCESSORIES:

DOH Item Number	Description	Medhold Part Number	Quantity
RT2-05-19a	Reusable latex-free paediatric cuff	DUR-P2-2A-E	Each
RT2-05-19b	Reusable latex-free adult cuff	DUR-A2-2A-E	Each
RT2-05-19c	Reusable latex-free neonate cuff	DUR-P1-2A-E	Each
RT2-05-19d	Reusable latex-free obese cuff	DUR-T1-2A-E	Each
RT2-05-19e	NIBP hose (adult/Paeds)	4051	Each
RT2-05-19f	Finger sensor for Paeds	4050	Each
RT2-05-19g	Finger sensor for adults	4054	Each
RT2-05-19h	Finger wrap-around probe for neonates	MASIMO-RD-4085	Each
RT2-05-19i	Spo2 interface cable	2021406-001	Each
RT2-05-19j	Two reusable temperature probes	Exergen Temporal Scanner	Each
RT2-05-19k	Front bezel	2037103-062	Each
RT2-05-19l	Back cover	2037103-002	Each
RT2-05-19m	Battery	PU044	Each
RT2-05-19n	NIBP pump/module	2037103-015	Each

RT2-2021

**SUPPLY, DELIVERY, INSTALLATION,
COMMISSIONING AND MAINTENANCE OF
MEDICAL EQUIPMENT**

**TO THE STATE FOR THE PERIOD
01 DECEMBER 2021
TO 01 DECEMBER 2024.**

SECTION G

THEATRE TABLES & LIGHTS

CONTACT DETAILS

Sibongile Leeu: s.leeu@medhold.co.za
Head Office: 011 966 0600

Or contact your local branch – details on the back page

THEATRE TABLE

ITEM NUMBER: RT2--14-02 OFFER 2
ITEM NAME: MAQUET MEERA BARIATRIC THEATRE TABLE

Make/Brand	Getinge Maquet
Model	MEERA Bariatric Theatre Table
Medhold Part Number	7200.01B0
DOH SAP number	
Country of Origin	Germany

UNIT INCLUDES:

- MEERA electrohydraulic, without cable- connected hand control, head rest or extension segment, with SFC padding”
- Cable-connected hand control
- Adaptors for head rests, pair, for 1131, 1132, 1133, 1433
- Extension segment, with transponder, coded
Suitable for: 1118.06, 1131, 1132, 1133, 1433”
- Head rest, inclinable, SFC padding
- Leg plates, pair, Gas-spring assisted, two pieces, hinged and abductable, SFC padding Suitable for: 1118, 1131, 1132, 1133, 1434”
- Swivel arm board 60mm - pair
- Goepel knee crutch, adjustable via ball- and-socket joint, with fastening strap, SFC pad (pair.) Included 2 x 1003.23C0 radial setting clamps
- Anaesthesia screen without clamp Suitable for: side rail
- Accessories stand


OPTIONAL ITEMS
Extended Warranty and Maintenance Contract - 5 years

Optional Year 3 to 7 Warranty & Preventative

Maintenance Including VAT AS per RT2 Circular and attached pricelist

UNIT PRICE INCLUSIVE OF VAT
AS PER RT2 CIRCULAR

Offer price Includes: Start-up accessories & consumables as per tender spec, Plus Yr 1 - 2 Warranty and Preventative Maintenance.

WEB SITE LINK
GETINGE MAQUET
[TINYURL.COM/A79BVXXB](https://tinyurl.com/A79BVXXB)

GETINGE MAQUET

Scan me



THEATRE TABLE ACCESSORIES

ITEM NUMBER: RT2--14-02 OFFER 2

ITEM NAME: MAQUET MEERA BARIATRIC THEATRE TABLE

Item Number	Spec Number	Description	Medhold Part Number	Quantity
RT2-14-02	1.1	MEERA electrohydraulic, without cable- connected hand control, head rest or extension segment, with SFC padding	7200.01B0	
	1.1	Cable-connected hand control	7200.90A0	
	1.1	Adaptors for head rests, pair, for 1131, 1132, 1133, 1433	1130.81A0	
	1.1	Extension segment, with transponder, coded Suitable for: 1118.06, 1131, 1132, 1133, 1433	1131.31BC	
	1.1	Head rest, inclinable, SFC padding	1130.64C0	
	1.1	Leg plates, pair, Gas-spring assisted, two pieces, hinged and abductable, SFC padding Suitable for: 1118, 1131, 1132, 1133, 1434	1118.60B0	
	4.1	Swivel arm board 60mm - pair	SU01.04990002	
	4.2	Goepel knee crutch, adjustable via ball- and-socket joint, with fastening strap, SFC pad (pair.) Included 2 x 1003.23C0 radial setting clamps	1001.65A0	
	4.4	Anaesthesia screen without clamp Suitable for: side rail	1002.57A0	
	4.5	Accessories stand	1132.40A0	

ITEM NUMBER: RT2--14-02

ITEM NAME: BARIATRIC ACCESSORIES

Code	Spec Number	Description	Quantity
1001.86B0	RT2-14-02-c	Foot plates, with radial setting	1
1133.58BC	RT2-14-02d	Leg plate, one-section, incl. SFC padding Suitable for: 1118.01, 1118.06, 1131, 1132, 1133, 1433	1
1001.73A0		Leg rest for obese patients, without clamp, Pair included 2 x Radial setting clamps 1003.23C0	1
1001.4700	RT2-14-02-b	Leg restraint cuff - pair	1
1101.6485	RT2-14-02-b	Leg restraint cuff, pair	1
1001.75A0	RT2-14-02-a	Table width extender, long, pair, 300 x 100 mm Suitable for: side rail	1
1001.76A0	RT2-14-02-a	Table width extender, short, pair, 200 x 100 mm Suitable for: side rail	2

THEATRE TABLE ACCESSORIES

ITEM NUMBER: RT2-14-02A

ITEM NAME: ORTHOPAEDIC SET ACCESSORIES SET

Code	Spec Number	Description	Quantity
		ASSIST SYSTEM Extension/Traction Device	*
1419.01JC	1.2	Extension device/Traction device, CrNi steel design, incl. cart for transport and mounting, accessory basket and basic equipment, with seat plate and buttock support sections, SFC padding, without traction stirrup clamp with rotation or rotating and tilting clamp, without foot plates or leg plates, without special accessories for MEERA/Alphastar PRO (1132.21/22) incl. basic equipment: 1 pair of adaptors 1 transport cart with basket 3114.6023 1 traction bar, short 3114.6723 1 traction bar, long 1003.3700 1 screw tension device 1003.49B0 1 foot plate mount 1002.18A0 2 supporting bars 1003.23C0 2 radial setting clamps 3112.3134 1 countertraction post for neck of femur 9080.5544 1 padded roll for 3112.3134 1004.91A0 1 side rail extension"	1
1003.34A0	1.2	Rotation and tilting clamp	2
1001.97A0	1.8	Foot plates, pair	1
1150.64B0	1.2	Leg plates, pair	1
1003.54A0	1.2	Bar elongation	2
1130.54B0	1.1	Connecting bracket for 1002.65A0/.07A0/.10A0 and skull clamp accessories Suitable for: head rest mounting point"	1

ORTHOPAEDIC HAND AND ARM ACCESSORIES OPTIONS

1001.63A0	1.4	Hand operating table, 845 x 380 mm, metal- free, radiolucent, with side rails, without clamp, with SFC padding need to add 1003.23C0 radial setting clamp	1
1003.23C0	1.4	Radial setting clamp Suitable for: side rail"	1

ACCESSORIES

1132.40A0	1.6	Accessories stand	1
1001.59A0	1.7	Body strap New type Metal clamp	1
1002.65A0	1.1	Connecting fixture, with rail for fixation of head rests	1
1002.72A0	1.1	Horseshoe-shaped head rest, two-part	1
1131.55BC	1.3	Seat plate extension, with transponder, coded Suitable for: 1118.06, 1131, 1132, 1133, 1433"	1

THEATRE TABLE ACCESSORIES

ITEM NUMBER: RT2-14-01B

ITEM NAME: NEUROSURGERY ACCESSORIES SET

Code	Spec Number	Description	Quantity
		Neuro ACCESSORIES	*
1130.54B0	1.8	Connecting bracket for 1002.65A0/07A0/10A0 and skull clamp accessories Suitable for: head rest mounting point	1
1005.51A0	1.5	Connection clamp for skull clamp	1
1005.2300	1.6	Cross-bar attachment for positioning the sitting patient, without clamps	1
1005.56A0	1.1	Horseshoe-shaped head rest with gel pad	1
1005.53B0	1.2	Pins, adults, scope of delivery includes 3 pcs, to be used with skull clamp, reusable product	1
1005.52C0	1.2	Skull clamp Doro	1
1005.50A0	1.4	Standard holder for skull clamp, width- adjustable	1
1132.40A0	1.7	Accessories stand	1
1002.65A0	1.3	Connecting fixture, with rail for fixation of head rests	1
1002.72A0	1.3	Horseshoe-shaped head rest, two-part	1

ITEM NUMBER: RT2-14-01C

ITEM NAME: SPINAL ACCESSORIES SET

Code	Spec Number	Description	Quantity
		ASSIST SYSTEMS - Universal Frame	
1007.24A1	1.9	Universal Frame for patient positioning, 2 carbon fibre bars, 360° metalfree translucent, mirror to control position of tracheal tube, 2 arm supports, 1 perspex head rest, 1 accessory clamp, 1 head rest adapter, 1 trolley, 1 universal mounting point. L x W x H: 1600 x 630 x 600/900 mm"	1
1007.25A0	1.9	Side rail clamp. side rail format for EU, for universal frame	1
1007.26A0	1.9	Short positioning plate. for Universal Frame L x W x H: 850 x 460 x 130 mm	1
1007.26B0	1.9	Long positioning plate. for Universal Frame L x W x H: 1350 x 460 x 130 mm	1
1007.28B0	1.9	Adjustable body supports. with SFC pad for universal frame"	1
1007.29A0	1.9	Single-use covers. for adjustable body supports, 5 sets with 4 pieces each"	1
1007.30A0	1.9	Adapter. for Radiolucent Spine Surgery Frame	1
1007.31A0	1.9	Radiolucent Spine Surgery Frame	1
1002.72A0	1.9 and 1.1	Horseshoe-shaped head rest, two-part	1

THEATRE TABLE ACCESSORIES

ITEM NUMBER: RT2-14-01C

ITEM NAME: SPINAL ACCESSORIES SET

Code	Spec Number	Description	Quantity
		Neuro ACCESSORIES	
1130.54B0	1.1	Connecting bracket for 1002.65A0/.07A0/.10A0 and skull clamp accessories Suitable for: head rest mounting point"	1
1005.51A0	1.5	Connection clamp for skull clamp	1
1005.2300	1.6	Cross-bar attachment for positioning the sitting patient, without clamps	1
1005.53B0	1.2	Pins, adults, scope of delivery includes 3 pcs, to be used with skull clamp, reusable product	1
1005.52C0	1.2	Skull clamp Doro	1
1005.50A0	1.4	Standard holder for skull clamp, width- adjustable	1
		ACCESSORIES	
1132.40A0	1.7	Accessories stand	1
1002.65A0	1.1	Connecting fixture, with rail for fixation of head rests	1

ITEM NUMBER: RT2-14-01C


ITEM NAME: MAQUET REPLACEMENT ACCESSORIES

Code	Spec Number	Description	Quantity
SU01.04990002	RT2-14-01d	Swivel arm board 60mm	2
1001.65A0	RT2-14-01e	Goepel knee crutch, adjustable via ball-and-socket joint, with fastening strap, SFC pad, without clamp, 1 piece (please order 2 pcs.)including 2 x 1003.23C0 radial setting clamp"	2
1003.23C0	RT2-14-01f	Radial setting clamp Suitable for: side rail	2
1002.57A0	RT2-14-01g	Anaesthesia screen With Radial setting clamp 1003.23C0	1
1132.40A0	RT2-14-01h	Accessories stand	1
1003.45D0	RT2-14-01i	TUR Drainage set, drainage bowl with frame, sieve and discharge hose	1
1003.46A0	RT2-14-01i	Wheeled stand with bucket for TUR drainage set	1
1131.55BC	RT2-14-01j	Seat plate extension, with transponder, coded Suitable for: 1118.06, 1131, 1132, 1133, 1433	1
1130.64C0	RT2-14-01k	Head rest, inclinable, SFC padding	1
1005.53B0	RT2-14-01l	Pins, adults, scope of delivery includes 3 pcs, to be used with skull clamp, reusable product	1
1005.52C0	RT2-14-01l	Skull clamp Doro	1
1002.72A0	RT2-14-01m	Horseshoe-shaped head rest, two- part	1
1002.65A0	RT2-14-01m	Connecting fixture, with rail for fixation of head rests	1
1130.54B0	RT2-14-01m	Connecting bracket for 1002.65A0/.07A0/.10A0 and skull clamp accessories Suitable for: head rest mounting point	1
1005.50A0	RT2-14-01n	Standard holder for skull clamp, width-adjustable	1
1005.51A0	RT2-14-01o	Connection clamp for skull clamp	1
1005.2300	RT2-14-01p	Cross-bar attachment for positioning the sitting patient, without clamps	1
1132.40A0	RT2-14-01q	Accessories stand	1
1130.54B0	RT2-14-01r	Connecting bracket for 1002.65A0/.07A0/.10A0 and skull clamp accessories Suitable for: head rest mounting point	1

EXAMINATION LIGHT

ITEM NUMBER: RT2-05-10 OFFER 1

ITEM NAME: MAQUET LUCEA 50 MOBILE WITH BATTERY BACK-UP

Make/Brand	Getinge Maquet	
Model	Lucea 50 M	
Medhold Part Number	ARDLCA309005A	
DOH SAP number		
Country of Origin	France	
UNIT INCLUDES:		
• Illumination 1m(Lux)	60 000	
• Field Diameter (cm/inch)*	22/8.6	
• Depth of illumination at 20%(cm/inch)	180/70.9	
• Depth of illumination at 60%(cm/inch)	120/47.2	
• Colour Temperature	4300 K	
• Colour Rendering index (Index/CRI)(Ra)	96	
• Irradiance at nominal illuminance (W/m2)	234	
• Service Life	60 000	
OPTIONAL ITEMS		
Extended Warranty and Maintenance Contract - 5 years		
Optional Year 3 to 7 Warranty & Preventative Maitenance Including VAT		
UNIT PRICE INCLUSIVE OF VAT		
AS PER RT2 CIRCULAR		
		Offer price Includes: Start-up accessories & consumables as per tender spec, Plus Yr 1 - 2 Warranty and Preventa- tive Maintenance.

WEB SITE LINK
VOLISTA STANDOP 64
[HTTPS://TINYURL.COM/4AZS3PPA](https://tinyurl.com/4AZS3PPA)


LUCEA 50

Scan me



EXAMINATION LIGHT

ITEM NUMBER: RT2-05-11 OFFER 2
ITEM NAME: MAQUET LUCEA 50 WALL MOUNT

Make/Brand	Getinge Maquet	
Model	Lucea 50 WM	
Medhold Part Number	ARDLCA109004A-C	
DOH SAP number		
Country of Origin	France	
UNIT INCLUDES: <ul style="list-style-type: none"> • Illumination 1m(Lux) 60 000 • Field Diameter (cm/inch)* 22/8.6 • Depth of illumination at 20%(cm/inch) 180/70.9 • Depth of illumination at 60%(cm/inch) 120/47.2 • Colour Temperature 4300 K • Colour Rendering index (Index/CRI)(Ra) 96 • Irradiance at nominal illuminance (W/m2) 234 • Service Life 60 000 		
OPTIONAL ITEMS Extended Warranty and Maintenance Contract - 5 years Optional Year 3 to 7 Warranty & Preventative Maintenance Including VAT		UNIT PRICE INCLUSIVE OF VAT AS PER RT2 CIRCULAR Offer price Includes: Start-up accessories & consumables as per tender spec, Plus Yr 1 - 2 Warranty and Preventa- tive Maintenance.

WEB SITE LINK
GETINGE MAQUET
[HTTPS://TINYURL.COM/4AZS3PPA](https://tinyurl.com/4AZS3PPA)


LUCEA 50

Scan me



CEILING MOUNTED THEATRE LIGHTS

ITEM NUMBER: RT2--13-04 OFFER 1
ITEM NAME: MAQUET VOLISTA ACCESS 44-THEATER LIGHT, CEILING-MOUNT

Make/Brand	Getinge Maquet	
Model	Volista Access 44	
Medhold Part Number	ARDVCS209005A-C	
DOH SAP number		
Country of Origin	France	
UNIT INCLUDES:		
• Nominal Luminance(Lux)	130 000	
• Nominal Luminance with Boost Mode	160 000	
• Dimming Range	10-100%	
• Light Patch Diameter (cm/inch)	20-25 cm / 7.9-9.8 inches	
• Colour Temperature	4300 K	
• Colour Rendering index (Index/CRI)(Ra)	95	
• Shadow Dilution with one offset mask(%)	77	
• Shadow Dilution with one offset 2 Masks(%)	50	
• Irradiance at nominal illuminance (W/m2)	<500	
• Standard Features	Ambient Light	
OPTIONAL ITEMS		
Extended Warranty and Maintenance Contract - 5 years		
Optional Year 3 to 7 Warranty & Preventative Maintenance Including VAT		
UNIT PRICE INCLUSIVE OF VAT		
AS PER RT2 CIRCULAR		
Offer price Includes: Start-up accessories & consumables as per tender spec, Plus Yr 1 - 2 Warranty and Preventative Maintenance.		


WEB SITE LINK
VOLISTA ACCESS 44
[HTTPS://TINYURL.COM/2JRDVTP3](https://tinyurl.com/2JRDVTP3)
VOLISTA ACCESS 44

Scan me



CEILING MOUNTED THEATRE LIGHTS

ITEM NUMBER: RT2--13-04 OFFER 2
ITEM NAME: MAQUET VOLISTA ACCESS 64-THEATER LIGHT, CEILING-MOUNT

Make/Brand	Getinge Maquet	
Model	Volista Access 64	
Medhold Part Number	ARDVCS209007A-C	
DOH SAP number		
Country of Origin	France	
UNIT INCLUDES:		
<ul style="list-style-type: none">• Nominal Luminance(Lux)130 000• Nominal Luminance with Boost Mode160 000• Dimming Range10-100%• Light Patch Diameter (cm/inch)20-25 cm / 7.9-9.8 inches• Colour Tembperature4300 K• Colour Rendering index (Index/CRI)(Ra)95• Shadow Dilution with one offset mask(%)77• Shadow Dilution with one offset 2 Masks(%)50• Irradiance at nominal illuminance (W/m2)<500• Standard FeaturesAmbient Light		
OPTIONAL ITEMS		
Extended Warranty and Maintenance Contract - 5 years		
Optional Year 3 to 7 Warranty & Preventative Maitenance Including VAT		
UNIT PRICE INCLUSIVE OF VAT		AS PER RT2 CIRCULAR
Offer price Includes: Start-up accessories & consumables as per tender spec, Plus Yr 1 - 2 Warranty and Preventa-tive Maintenance.		


WEB SITE LINK
VOLISTA ACCESS 64
[HTTPS://TINYURL.COM/2JRDVTP3](https://tinyurl.com/2JRDVTP3)
VOLISTA ACCESS 64

Scan me



CEILING MOUNTED THEATRE LIGHTS

ITEM NUMBER: RT2--13-04 OFFER 2
ITEM NAME: MAQUET VOLISTA STANDOP 64-THEATER LIGHT, CEILING-MOUNT

Make/Brand	Getinge Maquet	
Model	Volista StandOP 64	
Medhold Part Number	ARDVST229009A-C	
DOH SAP number		
Country of Origin	France	
UNIT INCLUDES:		
<ul style="list-style-type: none">Nominal Luminance(Lux)130 000Nominal Luminance with Boost Mode160 000Dimming Range10-100%Light Patch Diameter (cm/inch)20-25 cm / 7.9-9.8 inchesColour Tembperature3900 KColour Rendering index (Index/CRI)(Ra)95Shadow Dilution with one offset mask(%)77Shadow Dilution with one offset 2 Masks(%)50Irradiance at nominal illuminance (W/m2)<500Standard Features : Automatic Illumination Management, Ambient Light		
OPTIONAL ITEMS		
Extended Warranty and Maintenance Contract - 5 years		
Optional Year 3 to 7 Warranty & Preventative Maitenance Including VAT		
UNIT PRICE INCLUSIVE OF VAT		
AS PER RT2 CIRCULAR		
Offer price Includes: Start-up accessories & consumables as per tender spec, Plus Yr 1 - 2 Warranty and Preventa-tive Maintenance.		

WEB SITE LINK
VOLISTA STANDOP 64
[HTTPS://TINYURL.COM/4UPMVMtv](https://tinyurl.com/4UPMVMtv)
VOLISTA STANDOP 64

Scan me



CEILING MOUNTED THEATRE LIGHTS

ITEM NUMBER: RT2-13-04-OFFER 4
ITEM NAME: MAQUET POWERLED II THEATRE LIGHT - CEILING MOUNTED

Make/Brand	Getinge Maquet	
Model	POWERLED II SB 55	
Medhold Part Number	ARDPWT209038A-C	
DOH SAP number		
Country of Origin	France	
UNIT INCLUDES:		
• Nominal Luminance(Lux)	140 000	
• Nominal Luminance with Boost Mode	160 000	
• Dimming Range	10-100% with sensitivity response	
• Light Patch Diameter (cm/inch)	13-20 / 5.1"-7.9"	
• Colour Tembperature	4300 K	
• Colour Rendering index (Index/CRI)(Ra)	96	
• Shadow Dilution with one offset mask(%)	100	
• Shadow Dilution with one offset 2 Masks(%)	80	
• Irradiance at nominal illuminance (W/m2)	<500	
• Standard Features	: AIM, FSP, integrated antibacterial coating, laser positioning	
OPTIONAL ITEMS		
Extended Warranty and Maintenance Contract - 5 years		
Optional Year 3 to 7 Warranty & Preventative Maitenance Including VAT		
UNIT PRICE INCLUSIVE OF VAT		
AS PER RT2 CIRCULAR		
Offer price Includes: Start-up accessories & consumables as per tender spec, Plus Yr 1 - 2 Warranty and Preventa-tive Maintenance.		

WEB SITE LINK
POWERLED II SB 55
[HTTPS://TINYURL.COM/BDFZRPXF](https://tinyurl.com/BDFZRPXF)
POWERLED II SB 55

Scan me





CONTACTS



WWW.MEDHOLD.CO.ZA

INFO@MEDHOLD.CO.ZA

HEAD OFFICE GAUTENG

Telephone: +27 (11) 966 0600
General Fax: +27 (11) 975 3870
MSI Business Park,
68 Rigger Rd, Spartan,
Kempton Park,
1619

EAST LONDON

Branch Manager: Francois Fourie
Telephone: +27 (043) 748-4986
Fax: +27 (043) 748-4981
Suite No. 18 Beacon Park
89 Beacons Road
East London
5252

GQEBERHA

Branch Manager: Francois Fourie
Telephone: +27 (041) 364-0168
Fax: +27 (041) 365-2591
27 Shirley Street
Newton Park
Gqeberha
6025

FREE STATE

Branch Manager: Dawid Venter
Telephone: +27 (051) 406-2400
Fax: +27 (051) 448-9319
Palm Park Suite 4&5
94 Kellner Street
Westdene
Bloemfontein
9301

KWA-ZULU NATAL

Branch Manager: Wayne Gimingham
Telephone: +27 (031) 764-3638
Fax: +27 (031) 764-4891
78 Old Main Road,
Kloof
Durban
3610

WESTERN CAPE

Branch Manager: Hayley Bloomer
Telephone: +27 (021) 976-7499
Fax2mail: 087 231 4668
8 Langenhoven Street
Durbanville
7550