

Compact Trolley

Compact trolley designed to provide mobility for MAC™ 5 Resting ECG systems



Improved cable management



Dimensions (without resting ECG system)

A4 trolley (without resting ECG system)

Height	37.8 in (96 cm) 56.4 in (143 cm) with acquisition module and holder
Width	23.6 in (60 cm)
Depth	23.6 in (60 cm)
Weight	37 lbs (16.9 kg) — basic trolley 2.4 lbs (1.1 kg) — acquisition simply arm and barcode holder

A5 and Lite trolley (without resting ECG system)


Height	37.8 in (96 cm) 56.4 in (143 cm) with acquisition module and holder
Width	23.6 in (60 cm)
Depth	23.6 in (60 cm)
Weight	35 lbs (15.9 kg) — basic trolley 2.4 lbs (1.1 kg) — acquisition simply arm and barcode holder

Supported ECG Systems

MAC 5

Optional accessories

- **KISS suction system (option where available)**
- **Second tray**
 - Holds A4, A5 and letter size paper
 - Inside dimensions of tray: 13 in (33 cm) x 4.3 in (110 cm) x 5.7 in (144 cm)



©2021 General Electric Company – All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Not 510k cleared. GE, the GE Monogram, and MAC are registered trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. Compact trolley is produced by SUZHOU TIANYE MECHANICAL INDUSTRY CO. GE Medical Systems, Inc., doing business as GE Healthcare.

JB18569XX

