



CrystalView™ Pro Fluid Management System

CrystalView™ Pro is designed to be simple and intuitive for the user, while providing reliable fluid management in arthroscopic procedures. With CrystalView™ Pro, you can **focus on the patient, not the pump.**

NEW




Hall® MicroFree™
Small Bone Power System



TruShot™ with Y-Knot®
Tissue Fixation System



Lanza™
Tapered Specialty Shaver



To learn more about these
and other products,
call 800-448-6506 or
visit CONMED.com.

FLUID MANAGEMENT SYSTEM

A Cut Above

CrystalView™ Pro Fluid Management System combines high performance inflow and outflow for optimal distension and clear visualization with intuitive features and customizable surgeon presets to help simplify the OR.

Simple

- The CrystalView™ Pro was designed so your fluid pump is one less thing you have to worry about. In fact, it only takes three clicks for it to start pumping fluid.

Intuitive

- For maximum ease-of-use, the user-friendly touchscreen works like an app on your phone and is designed to simplify pump setup.

Reliable

- New direct drive motor design is intended to eliminate gear seizing, reduce noise, and increase reliability.

Customizable

- Customizable surgeon presets allow for quick set-up. Simply select the surgeon and procedure type and that's it. These settings can be easily exported to other CrystalView™ Pro units, saving installation time for multi-room facilities.

Schedule an evaluation today to simplify your OR.

CONMED.com/CrystalView

CrystalView™ Pro is compatible with CONMED's cost-effective Eco-Flow™ Day-Use Two-Part Tubing that allows only the patient portion of the tubing to be replaced after each case so the pump tubing can stay in place throughout the day.



Additional Features

- Action-oriented alerts for easy, intuitive troubleshooting
- Two unique rinse modes
- Wave mode available in all procedure types
- Wider distance between rotors for easier cassette installation
- Easily toggle between set and actual pressure as well as mmHg or H₂O
- Autoclavable remote control



Ordering Information

CrystalView™ Pro Fluid Management System

Description	Product Code	Quantity
CrystalView™ Pro Irrigation Console	CV4000	1 each
CrystalView™ Pro Remote	CV4040	1 each
Pump Interface Cable	D4024	1 each
Shaver Sensor	24k900	1 each
Shaver Sensor Extension	24k191	1 each
Outflow/Shaver Suction/Waste Tubing	24k100	10 per box
Inflow Arthroscopy Tube Set	10k100	10 per box
10K Arthroscopy Tube Set with Suction Line	10k150	10 per box
Day-Use Pump Tubing	10k600	4 per box
Single-Use Patient Tubing	10k605	12 per box

CONMED Corporation
525 French Road, Utica, NY 13502

Local 727-392-6464
Toll Free 1-800-237-0169

CONMED.com
customer_service@CONMED.com