


Unter Haßlen 20 · D-78532 Tuttlingen  
Phone: +49 (0)7462/9470-0  
Fax: +49 (0)7462/9470-50  
e-mail: [sales@rz-medizintechnik.com](mailto:sales@rz-medizintechnik.com)  
<http://www.rz-medizintechnik.com>



RZ Medizintechnik  
GmbH *Progress in Surgery*

## RZ MAMMAPLAST



## RZ MAMMAPLAST

In order to proceed a modern breast reduction surgery, the deepithelization technique is still the first step to start with. The RZ MAMMAPLAST self-retaining retractor helps to provide constant pressure on the breast so that deepithelization can be done without further assistance: the breast will be kept tightly and safely in its position.

I 260-500-000  
Breast Elevator for Mamma Plastic,  
detachable

