



progress in surgery

360° rotatable
BIPOLAR
FORCEPS



ENT

BIPOLAR MICRO-COAGULATION FORCEPS

RZ - Medizintechnik adds the state-of-the-art 360° Bipolar Micro-Coagulation Forceps to comply with the demand of the surgeons for a bipolar Micro Forceps with a minimal diameter of 3mm and using bipolar electricity.

This bipolar forceps can be used in the field of Ear Nose and Throat or in any other operations with very limited access for larger instruments. The large variety of jaws allows the surgeon to use his preferred jaw design.


To comply with the demand of hygienists for easy cleaning and sterilisation we added a disposable single use cleaning and flushing tube, which allows the surgeon a better view during surgeries.

Advantages:

- For operations with difficult accessible operation areas
- Bipolar technology for a maximum of patient safety
- Various popular electrode inserts, available
- Ergonomic, surgeon friendly handle design
- 3 piece take apart system with an additional single use cleaning and flushing tube




1:1




360° rotatable shaft



Detachable system



Ø 3mm shaft



Additional Suction Tube



SHAFT AND HANDLE



1

1 252-808-001

Shaft, 150mm, Ø 3mm



2

2 252-808-000

Ring Handle

ACCESSORIES



252-808-007

Disposable Suction Tube

INSERTS

150mm working length



252-808-002
Grasping Forceps, fenestrated



252-808-003
Spoon Forceps Ø2mm, curved



252-808-004
Spoon Forceps Ø2.7mm, curved



252-808-005
Maryland dissector, curved



252-808-006
Dissector, angled

HF CABLES



Generator end



2-Banana *

3m 700-350-030
5m 700-350-050



Erbe

3m 700-738-003
5m 700-738-005



Valleylab
Lamidey, EMC

3m 700-736-003
5m 700-736-005



Martin, Aesculap
Berchtold

3m 700-737-003
5m 700-737-005



U20

3m 700-252-818

* Not available with CE marking within Europe



progress in surgery

RZ Medizintechnik GmbH

Unter Hasslen 20
D-78532 Tuttlingen

Phone: +49 (0)7462/9470-0
Fax: +49 (0)7462/9470-50

sales@rz-medizintechnik.com
www.rz-medizintechnik.com

Media Code Number: 021550