


RZ Medizintechnik GmbH
Unter Hasslen 20 | D-78532 Tuttlingen
Phone: +49 (0)7462/9470-0
Fax: +49 (0)7462/9470-50
sales@rz-medizintechnik.com
www.rz-medizintechnik.com




BIPOLAR Tonsillectomy Forceps



news
07.2015

BIPOLAR TONSILLECTOMY FORCEPS



Ordinary-shaped, healthy tonsil

Tonsil, indicated to surgical extraction

Conventional COLD STEEL TONSILLECTOMY



Disadvantages

- use of many surgical instruments
- limited visibility due to heavy bleedings
- vessel sealing only with tonsil snare
- frequent changes between instruments
- sometimes stitches are needed

New procedure BIPOLAR TONSILLECTOMY



Advantages

- use of only 1 instrument
- significantly less bleeding
- better view
- almost bloodless operation
- immediate atrophy
- patient-friendly
- shorter operation time
- less post-operative pain
- no stitches needed

430-701-000
Bipolar Tonsillectomy Forceps
82mm working length, **with** suction



Suction on/off with finger rest



Suction directly on tip







Suction plug

430-701-001
Bipolar Tonsillectomy Forceps
82mm working length, **without** suction

BIPOLAR CABLES



Bipolar Cables for RZ Tonsillar Forceps
4.5m cable length

-  700-740-001
2-Banana
-  700-740-002
Erbe/Storz/Wolf
-  700-740-003
Martin
-  700-740-004
Valleylab

