

RZ Medizintechnik GmbH
Unter Hasslen 20 | D-78532 Tuttlingen

Phone: +49 (0)7462/9470-0
Fax: +49 (0)7462/9470-50

sales@rz-medizintechnik.com
www.rz-medizintechnik.com



RZ
Medizintechnik
GmbH

*over 20 years of
progress in surgery*


Non-Stick monopolar SUCTION COAGULATOR



news
04.2016

ENT NON-STICK SUCTION COAGULATOR


Stop with the new Suction Coagulator bleedings in an effective, very comfortable and safety way. Parallel coagulation and suction without interruptions are integrated into a comfortable and well designed handle which improves every surgery. The superior non-stick tip ensures successive coagulation without complications.



Universal suction tubing connection
Standard international
4mm electrical connection



Insulation coating for maximum safety



Suction thumb control




Non-stick and non-adherencing safety valve




Ergonomic handle for effortless working

SUCTION COAGULATOR

Order Number	Working length	Feature
700-420-101	108mm	2x angled











Order Number	Working length	Feature
700-420-102	175mm	angled



Order Number	Working length	Feature
700-420-103	176mm	straight, malleable

MONOPOLAR CABLES

			
		4mm socket	4mm socket, ang.
Generator end			
	3m	700-501-303	
	5m	700-501-503	
RZ HF Unit Standard Plug			
	3m	700-501-302	
	5m	700-501-502	
RZ HF Unit SDS Plug, programmable			
	3m	700-280-030	700-290-030
	5m	700-280-050	700-290-050
Erbe, Martin, Berchtold			
	3m	700-360-031	700-361-030
	5m	700-360-051	700-361-050
Bovie, Valleylab, Conmed			
	3m	700-285-030	
	5m	700-285-050	
Martin, Berchtold			
	3m	700-101-030	700-102-030
	5m	700-101-050	
Erbe, ACC/ICC			



RZ HIGH FREQUENCY CABLES

RZ Medizintechnik GmbH offers an extensive range of high-quality cables for use with established devices in bipolar and monopolar HF surgery. All cables meet highest safety standards and provide best performance at persistent reliability.

Advantages:

- FEP Insulation
- Poreless silver plated wires
- Silicone jacket
- Reliable and strong sheathing
- Latex-free

700-350-400
Electrosurgical Generator G400



Media Code Number: 021310

© 04/2016 Copyright by RZ Medizintechnik GmbH · Unter Hasslen 20 · D-78532 Tuttlingen. All rights reserved.
Modifications of technical specification or design are subject to change without prior notice.