


RZ Medizintechnik GmbH
Unter Hasslen 20 | D-78532 Tuttlingen
Phone: +49 (0)7462/9470-0
Fax: +49 (0)7462/9470-50
sales@rz-medizintechnik.com
www.rz-medizintechnik.com



RZ
Medizintechnik
GmbH

Progress in Surgery

Bone Fragment DISTRACTORS



news
03.2014


RZ BONE FRAGMENT DISTRACTORS

Flexible positionable forceps with closed and outspread arms and two different pin holes for K-Wire pins to spread bone fragments.

- two different K-Wire hole sizes
- two arm designs
- 150mm and 210mm sizes
- variety of indications
- flexible positioning



2.5mm and 1.6mm pin holes



Two different models available - outspread model for maximum spreading width



CLOSED

210-974-301	150mm
210-974-311	210mm



OUTSPREAD

210-974-300	150mm
210-974-310	210mm