





progress in surgery

ArthVANCE[®]


SHOULDER
INSTRUMENTS



ARTHROSCOPY

ArthVANCE® SHOULDER INSTRUMENTS

SUTURE PASSER



Suture Passer

170mm working length

225-820-217



Needle with lateral opening, 0.3 x 170mm

for Suture Passer

170mm working length

225-820-219


KNOT PUSHER



Knot Pusher, single-hole, for sliding knots and half-stitches,
modified distal tip for easy slipknots

225-259-000

SUTURE CUTTER / KNOT PUSHER



- special locking mechanism
- secure knot pushing
- easy cutting behind the knot
- ergonomic axial handle

Suture Cutter, axial handle, \varnothing 4.2 mm, v-notch,
left opened, 130 mm working length

225-820-100



Suture Cutter, ring handle, \varnothing 4.2 mm, v-notch,
left opened, 130 mm working length

225-820-101

with TC inserts, golden rings

225-820-103

ArthVANCE® SHOULDER INSTRUMENTS


SUTURE RETRIEVER



Fiber Retriever, straight, 155mm working length,
retriever \varnothing 1.5mm pin

225-820-259


DEBRIDEMENT RASP



Debridement Rasp, razor sharp, ergonomic PEEK handle


225-600-630

TISSUE ELEVATOR



Tissue Elevator, precise, 15° angled,
ergonomic PEEK handle

225-600-650



SMART HOOK

Best support for rotator cuff repair.




Left curved jaw
225-820-307



Right curved jaw
225-820-309

SIXER

Best support for labral repair and capsule shift.




Twisted left
225-820-306



Twisted right
225-820-308

ArthVANCE® SHOULDER INSTRUMENTS



BIRDY

Best support for slap repair, labral repair and rotator cuff repair.



Straight jaw, upcurved shaft

A 225-820-400


B 225-840-400



20° upwards angled

A 225-820-402

B 225-840-402



40° upwards angled

A 225-820-404

B 225-840-404



45° right angled

A 225-820-406

B 225-840-406



45° left angled

A 225-820-408

B 225-840-408

ROTATOR CUFF TESTER





Rotator Cuff Tester

(available on request)

225-225-225

Please do not test the Suture Passer without adequate material in between the jaw. The needle may break due to the friction. We offer and recommend you this Rotator Cuff Tester 225-225-225.





progress in surgery

RZ Medizintechnik GmbH

Unter Hasslen 20
78532 Tuttlingen, Germany

Phone: +49 (0)7462/9470-0
Fax: +49 (0)7462/9470-50

sales@rz-medizintechnik.com
www.rz-medizintechnik.com

Media Code Number: 021660

© 11/2018 Copyright by RZ Medizintechnik GmbH - Unter Hasslen 20 - D-78532 Tuttlingen. All rights reserved.
Modifications of technical specification or design are subject to change without prior notice.