

Neuromuscular Transmission Module, E-NMT

For integrated measurement of the level of neuromuscular block



The neuromuscular transmission module, E-NMT, integrates the measurement of the level of neuromuscular block into a single-width, plug-in module. The module complies with IEC 60601-1 3rd edition.

Features

- Provides a quantitative, automatic measurement of muscle response to an electrical stimulus
- Offers all common stimulation modes: train of four (TOF), single twitch (ST), double burst stimulation (DBS) and tetanic stimulation
- Covers the whole range of neuromuscular blockade
- Two measurement sensor options for signal acquisition: MechanoSensor for kinemyography (KMG) and ElectroSensor for electromyography (EMG)
- Offers two MechanoSensor options: for adults and for pediatric patients
- Provides automatic setting of supramaximal current
- Provides automatic measurement at user-selectable interval
- Recovery Note can be set to notify when the block is wearing off
- Recall function enables moving the module with the patient without losing the reference values and the measured supramaximal current
- In addition to the integrated neuromuscular block measurement, the module enables nerve location for regional block by utilizing a comfortable 40 μ s pulse

Technical specifications

Direct function keys

Start-up	Automatically initiates the measurement by setting supramaximal current, reference value and starting cycle
Stop/Continue	Stops / continues measurement with same patient

NMT

Stimulation modes	Train of four (TOF) Double burst (3,3), (DBS) Single twitch (ST) 50 Hz tetanic (post tetanic count)
Numeric display	TOF%/DBS%, Count, T1%, PTC

Measurement intervals for TOF/DBS

CARESCAPE™ modular monitors with OR and PACU software packages: Manual, 10 s, 12 s, 15 s, 20 s, 1 min, 5 min or 15 min.

CARESCAPE modular monitors with ICU and ED software packages: Manual, 20 s, 1 min, 5 min, 15 min, 30 min, 60 min or 120 min.

Measurement intervals for ST

Manual, 1 s, 10 s, 20 s

Stimulus pulse	Square wave, constant current
Pulse width	100, 200 or 300 μ s
Stimulus current range	Supramax: 1 to 70 mA Manual: 1 to 70 mA, in 1 mA steps
Stimulus current accuracy	20% or \pm 3 mA, whichever is greater
Load for 70mA	3 k Ω (max)
Voltage	300 V (max)

Regional Block mode

Stimulation mode	Single twitch
Stimulation intervals	1, 2, 3 s
Stimulus pulse	Square wave, constant current
Pulse width	40 μ s
Stimulus current range	0 to 5.0 mA with 0.1 mA steps
Stimulus current accuracy	20% or \pm 0.3 mA, whichever is greater

Monitor compatibility

CARESCAPE modular monitors with OR, PACU, ED and/or ICU software

Environmental specifications

Operating conditions

Temperature	10 to 40°C (50 to 104°F)
Relative humidity	10 to 90% non-condensing

Storage conditions

Temperature	-20 to 60°C (-4 to 140°F)
Relative humidity	10 to 90% non-condensing

Physical specifications

Dimensions (H x W x D)	11.2 x 3.7 x 18.6 cm (4.4 x 1.5 x 7.3 in)
Weight	0.35 kg (0.77 lb)



Imagination at work

Product may not be available in all countries and regions. Full product technical specification is available upon request. Contact a GE Healthcare Representative for more information. Please visit www.gehealthcare.com/promotional-locations.

Data subject to change.

© 2013 - 2020 General Electric Company.

GE, GE Monogram, Imagination at work and CARESCAPE are trademarks of General Electric Company.

All other third-party trademarks are the property of their respective owners.

Reproduction in any form is forbidden without prior written permission from GE. Nothing in this material should be used to diagnose or treat any disease or condition. Readers must consult a healthcare professional.

DOC2057228 Rev 2 2/20